



# Coloplast Wound Care

Product overview



*"I remember I was hospitalized on the 3<sup>rd</sup> floor. The doctor was planning to amputate my leg to the knee. Monica changed the whole treatment and told the management that she would take care of my foot for one month, and then they could reassess the situation. She brought a light at the end of the tunnel. Her dedication saved me"*

Celso, 55, Brazil, Diabetic foot ulcer

<b>Treatment of non-infected wounds</b>	<b>4</b>
Biatain® Silicone with 3DFit Technology™	5
Biatain® Silicone Multishape	8
Biatain® Silicone Lite	10
Biatain® Non-Adhesive	11
Biatain® Adhesive	12
Biatain® Soft-Hold	13
Biatain® Super	14
Biatain® Alginate	16
Triad™	18

<b>Treatment of infected wounds</b>	<b>19</b>
Biatain® Ag	20
Biatain® Alginate Ag	22

<b>Protection of wound and skin</b>	<b>25</b>
Comfeel® Plus Transparent	26
Comfeel® Plus	27
Physiotulle®	28

<b>Wound bed preparation: cleansing and debridement</b>	<b>29</b>
Sea-Clens®	30
Purilon®	31
Woun'Dres®	32

<b>Brief statements</b>	<b>33</b>
-------------------------	-----------

# Treatment of non-infected wounds



## **Biatain® Silicone**

Foam Dressings

Biatain® Silicone with 3DFit Technology™ conforms to the wound bed to fill the gap and reduce exudate pooling for optimal healing conditions.



### **Conforms to the wound bed**

Upon contact with exudate, the unique foam with 3DFit Technology conforms to the wound bed and reduces exudate pooling.



### **Absorbs vertically**

Microcapillaries go into action to absorb exudate vertically and create optimal healing conditions.



### **Retains exudate**

Exudate is locked away and retained even under compression, reducing risk of leakage and maceration of the wound edges and periwound skin.

# Biatain® Silicone

## A perfect fit for treatment + prevention

Biatain® Silicone is a soft and conformable polyurethane foam dressing with an advanced 5-layer construction designed to address exudate management and the extrinsic risk factors that can lead to skin breakdown. It can be used on a broad range of exuding chronic and acute wounds, and the Sacral and Multishape dressings may be used as part of a protocol to prevent pressure injuries.

### Promotes optimal moist wound healing

- Conforms to the wound bed
- Absorbs vertically
- Retains exudate

### Protects against pressure injury risk factors

- Redistributes pressure
- Reduces friction and shear
- Manages microclimate

### Usage

Biatain Silicone is indicated for a wide range of exuding chronic and acute wounds. They include leg ulcers, pressure ulcers, non-infected diabetic foot ulcers, donor sites, postoperative wounds and traumatic wounds. It can be used in combination with compression therapy and may be left in place for up to 7 days depending on the amount of exudate, dressing conditions and type of wound.



reddot award 2017  
winner



Due to its outstanding design, Biatain Silicone has been awarded the Red Dot Award, an internationally recognized consumer award that focuses on user experience and aesthetics of products.

### Biatain® Silicone Foam Dressing

Code	Size	Units	HCPCS <sup>1</sup>
33434	3 x 3" (7.5 x 7.5 cm)	10	A6212
33435	4 x 4" (10 x 10 cm)	10	A6212
33436	5 x 5" (12.5 x 12.5 cm)	10	A6212
33437	6 x 6" (15 x 15 cm)	5	A6213
33438	7 x 7" (17.5 x 17.5 cm)	5	A6213
33400	4 x 8" (10 x 20 cm)	5	A6245
33401	4 x 12" (10 x 30 cm)	5	A6246
33406	Heel 7 x 7" (18 x 18 cm)	5	A6245
39000*	Sacral 9.8 x 9.8" (25 x 25 cm)	5	A6213
39001*	Sacral 6 x 7.5" (15 x 19 cm)	5	A6210
39010*	Multishape 5.5 x 7.6" (14 x 19.5 cm)	5	A6212

 Exudate Level

\*Indicated for pressure injury prevention



# Biatain® Silicone Multishape

Award-winning, user-friendly design  
for difficult to dress areas



Scan for application video

## Preparation

- Cleanse the wound and periwound skin in accordance with local guidelines, e.g. lukewarm water or physiological saline solution
- Gently dry the periwound skin
- If any film, cream, ointment or similar product is used, allow the periwound skin to dry before applying the product

## Application



Use protective films to avoid touching adhesive side and ensure aseptic application.  
Remove center protective film.



Apply adhesive side towards wound.



Remove remaining protective films, one at a time, gently securing each side to ensure an even and smooth fit on skin.



Gently run your fingers around border to confirm securement.



Please refer to the Instructions for Use (IFU) packaged with the product for further information

## Biatain® Silicone Lite

Biatain® Silicone Lite is a thin, ultra flexible and absorbent polyurethane foam dressing with a waterproof topfilm and a gentle silicone adhesive layer. The soft and thin foam pad is designed with a low profile to offer patients increased mobility.



- Conforms to the wound bed
- Absorbs vertically
- Retains exudate

### Usage

Biatain Silicone Lite can be used for a wide range of low exuding chronic and acute wounds. They include leg ulcers, pressure ulcers, non-infected diabetic foot ulcers, donor sites, postoperative wounds and skin abrasions. It can be used in combination with compression therapy and may be left in place for up to 7 days depending on the amount of exudate, dressing conditions and type of wound.

### Biatain® Silicone Lite Foam Dressing

Code	Size	Units	HCPCS <sup>1</sup>
33452	2 x 2" (5 x 5 cm)	5	A6413*
33453	2 x 5" (5 x 12.5 cm)	5	A6212
33444	3 x 3" (7.5 x 7.5 cm)	10	A6212
33445	4 x 4" (10 x 10 cm)	10	A6212
33446	5 x 5" (12.5 x 12.5 cm)	10	A6212

 —  Exudate Level

\* Not covered by Medicare

## Biatain® Non-Adhesive

Biatain® Non-Adhesive is a soft polyurethane foam dressing with beveled edges and a waterproof backing designed for extra fragile skin. Biatain Non-Adhesive absorbs and retains exudate to promote an optimal moist healing environment and is suitable for use under compression.



- Conforms to the wound bed
- Absorbs vertically
- Retains exudate

### Usage

Biatain Non-Adhesive can be used for a wide range of exuding wounds, including leg ulcers, pressure ulcers, non-infected diabetic foot ulcers, second degree burns, donor sites, postoperative wounds and skin abrasions. It can be used in combination with compression therapy and may be left in place for up to 7 days depending on the amount of exudate, dressing conditions and type of wound.


### Biatain® Non-Adhesive Foam Dressing

Code	Size	Units	HCPCS <sup>1</sup>
6105	2 x 2.75" (5 x 7 cm)	10	A6209
3410	4 x 4" (10 x 10 cm)	10	A6209
3413	6 x 6" (15 x 15 cm)	5	A6210
3416	8 x 8" (20 x 20 cm)	5	A6211

 —  Exudate Level

## Biatain® Adhesive

Biatain® Adhesive is a soft and conformable polyurethane foam dressing with a hydrocolloid-based adhesive border designed for wounds with healthy periwound skin that need extra adhesion. Biatain Adhesive absorbs and retains exudate to promote an optimal moist wound healing environment.



- Conforms to the wound bed
- Absorbs vertically
- Retains exudate

### Usage

Biatain Adhesive can be used for a wide range of exuding wounds, including leg ulcers, pressure ulcers, second degree burns, donor sites, postoperative wounds and skin abrasions. Biatain Adhesive can be used in combination with compression therapy and may be left in place for up to 7 days depending on the amount of exudate, dressing conditions and type of wound.

### Biatain® Adhesive Foam Dressing

Code	Size	Units	HCPCS <sup>1</sup>
3430	4 x 4" (10 x 10 cm) circle pad	10	A6212
3420	5 x 5" (12.5 x 12.5 cm)	10	A6212
3423	7 x 7" (18 x 18 cm)	5	A6213
3485	Sacral 9 x 9" (23 x 23 cm)	5	A6213
3488	Heel 7.5 x 8" (19 x 20 cm)	5	A6212

 Exudate Level

## Biatain® Soft-Hold

Biatain® Soft-Hold is a soft and conformable polyurethane foam dressing with an adherent layer covering less than 50% of the foam surface designed to hold the dressing in place during application and removal. Biatain Soft-Hold absorbs and retains exudate to promote an optimal moist wound healing environment and is suitable for use under compression.



- Conforms to the wound bed
- Absorbs vertically
- Retains exudate

### Usage

Biatain Soft-Hold can be used for a wide range of exuding wounds, including leg ulcers, pressure ulcers, second-degree burns, postoperative wounds and skin abrasions. It can be used in combination with compression therapy and may be left in place for up to 7 days depending on the amount of exudate, dressing conditions and type of wound. Biatain Soft-Hold should not be used for donor sites.

### Biatain® Soft-Hold Foam Dressing

Code	Size	Units	HCPCS <sup>1</sup>
3473	2 x 2.75" (5 x 7 cm)	5	A6209
3470	4 x 4" (10 x 10 cm)	5	A6209
3475	6 x 6" (15 x 15 cm)	5	A6210

 Exudate Level

# Biatain® Super

Biatain® Super is a highly-absorbent hydrocapillary dressing available in adhesive and non-adhesive versions. Biatain Super consists of a wound contact layer, superabsorbent pad and waterproof topfilm to promote an optimal healing environment for highly exuding wounds.



- Fast absorption, distribution and retention of exudate
- Biatain Super Non-Adhesive is suitable for highly exuding wounds surrounded by fragile skin

## Usage

Biatain Super can be used for the treatment of a wide range of exuding wounds, including leg ulcers, pressure ulcers, diabetic foot ulcers, second degree burns, surgical wounds and skin abrasions. Biatain Super is suitable for use under compression therapy and may be left in place for up to 7 days depending on amount of exudate, dressing conditions and type of wound.

## Biatain® Super, Adhesive

Code	Size	Units	HCPCS <sup>1</sup>
46100	4 x 4" (10 x 10 cm)	10	A6203
46120	5 x 5" (12.5 x 12.5 cm)	10	A6203
46250	5 x 8" (12.5 x 20 cm)	10	A6204
46150	6 x 6" (15 x 15 cm)	10	A6203
46200	8 x 8" (20 x 20 cm)	10	A6204

💧💧💧 – 💧💧💧 Exudate Level

## Biatain® Super, Non-Adhesive

Code	Size	Units	HCPCS <sup>1</sup>
46300	4 x 4" (10 x 10 cm)	10	A6251
46320	5 x 5" (12.5 x 12.5 cm)	10	A6251
46450	5 x 8" (12.5 x 20 cm)	10	A6252
46350	6 x 6" (15 x 15 cm)	10	A6251
46390	8 x 8" (20 x 20 cm)	10	A6252


💧💧💧 – 💧💧💧 Exudate Level



# Biatain® Alginate

Biatain® Alginate is a highly-absorbent alginate dressing for moderately to heavily exuding wounds of many sizes, shapes and depths. Biatain Alginate consists of calcium alginate and carboxymethylcellulose (CMC) and converts to a soft, cohesive gel upon contact with exudate.

- Recommended for cavity filling
- High dressing integrity for one-piece removal
- Documented haemostatic effect<sup>2</sup>



## Usage

Biatain Alginate can be used on moderate to heavily exuding partial to full thickness wounds.

Biatain Alginate may assist in supporting the control of minor bleeding in superficial wounds.

Biatain Alginate can be used with compression therapy.

Biatain Alginate should not be used on dry or lightly exuding wounds or to control heavy bleeding. Biatain Alginate is indicated for pressure ulcers, venous and arterial leg ulcers, diabetic ulcers, cavity wounds, lacerations, skin abrasions, graft and donor sites, post operative surgical wounds, superficial and partial thickness burns. It can remain in situ up to 7 days depending on wound condition and level of exudate.



## Biatain® Alginate Dressing

Code	Size	Units	HCPCS
3705	2 x 2" (5 x 5 cm)	30	A6196
3710	4 x 4" (10 x 10 cm)	10	A6196
3715	6 x 6" (15 x 15 cm)	10	A6197

💧—💧💧 Exudate Level



## Biatain® Alginate Rope

Code	Size	Units	HCPCS
3740	1 x 17.5" rope (3 x 44 cm)	6	A6199

💧—💧💧 Exudate Level

## Triad®

Triad® is a zinc-oxide based hydrophilic wound dressing designed to manage low to moderate levels of exudate, while promoting a moist wound environment to facilitate autolytic debridement.



- Ideal alternative for difficult-to-dress wounds
- Can be applied directly to wound or impregnated into gauze
- Unique preservative-free blend of four ingredients adheres to wet, eroded skin

### Usage

Triad is indicated for the local management of pressure and venous stasis ulcers, dermal lesions/injuries, superficial wounds, scrapes, first- and second-degree burns, partial- and full-thickness wounds. Contraindicated for third-degree burns or infected wounds. May be left in place up to 7 days.

### Triad™ Hydrophilic Wound Dressing

Code	Size	Units	HCPCS <sup>1</sup>
1964	2.5 oz. tube (71 g)	12	A6240
1967	6 oz. tube (170 g)	12	A6240

💧—💧 Exudate Level

## Treatment of infected wounds

# Biatain® Ag

## Combats infection and biofilms where it matters

Biatain® Ag manages moderate to high amounts of exudate and conforms to the wound bed to reduce exudate pooling. Biatain Ag has been shown to kill 99.99% of mature biofilms\* (*P. aeruginosa*) and to prevent biofilm formation (shown in vitro)<sup>3</sup>. The dressing is also effective against a broad spectrum of bacteria and fungi for up to 7 days.<sup>3</sup>

- Conforms to the wound bed
- Absorbs vertically
- Retains exudate
- Available in non-adhesive and adhesive version with hydrocolloid border

### Usage

Biatain Ag foam dressings are indicated for a wide range of exuding wounds including leg ulcers, pressure ulcers, second-degree burns, donor sites, postoperative wounds and traumatic wounds with delayed healing due to bacteria, or where there is a risk of infection. Biatain Ag Non-Adhesive foam dressings are additionally indicated for diabetic foot ulcers. They can be used in combination with compression therapy and may be left in place for up to 7 days depending on the amount of exudate, dressing conditions and type of wound.

\* Mature biofilms are defined as minimum 24 hours old.



### Biatain® Ag Non-Adhesive Foam Dressing

Code	Size	Units	HCPCS
9622	4 x 4" (10 x 10 cm)	5	A6209
9625	6 x 6" (15 x 15 cm)	5	A6210

💧—💧💧 Exudate Level



### Biatain® Ag Adhesive Foam Dressing

Code	Size	Units	HCPCS
9632	5 x 5" (12.5 x 12.5 cm)	5	A6212
9635	7 x 7" (18 x 18 cm)	5	A6213
9641	Sacral 9 x 9" (23 x 23 cm)	5	A6213
9643	Heel 7.5 x 8" (19 x 20 cm)	5	A6212

💧—💧💧 Exudate Level

# Biatain® Alginate Ag

Biatain® Alginate Ag is a highly-absorbent silver alginate dressing for moderately to heavily exuding infected wounds of many sizes, shapes and depths. Biatain Alginate Ag may be used on infected wounds under the discretion of a healthcare professional.

- Recommended for cavity filling
- Contains an ionic silver complex with a broad antimicrobial profile
- High dressing integrity for one-piece removal
- Documented haemostatic effect<sup>2</sup>

## Usage

Biatain Alginate Ag can be used on moderate to heavily exuding partial to full thickness wounds. Biatain Alginate Ag may assist in supporting the control of minor bleeding in superficial wounds. Biatain Alginate Ag may be used on infected wounds or wounds at risk of infection. Biatain Alginate Ag can be used under compression bandages. Biatain Alginate Ag should not be used on dry or lightly exuding wounds or to control heavy bleeding. Biatain Alginate Ag is indicated for post-operative wounds, trauma wounds (dermal lesions, trauma injuries or incisions), leg ulcers, pressure ulcers, diabetic ulcers, graft and donor sites, cavity wounds, superficial and partial thickness burns. Dressing change frequency will depend on wound condition and the level of exudate. Initially it may be necessary to change the dressing every 24 hours.



## Biatain® Alginate Ag Dressing

Code	Size	Units	HCPCS
3755	2 x 2" (5 x 5 cm)	30	A6196
3760	4 x 4" (10 x 10 cm)	10	A6196
3765	6 x 6" (15 x 15 cm)	10	A6197

💧—💧💧 Exudate Level



## Biatain® Alginate Ag Rope

Code	Size	Units	HCPCS
3780	1 x 17.5" rope (3 x 44 cm)	10	A6199

💧—💧💧 Exudate Level






*"These dressings have a dual purpose for me. It is full time protection and full time treatment when I need it – and I can manage it totally independently. So it means everything to me and makes me live a life with my leg despite the injury. That's the bottom line."*

Janice, 63, United Kingdom

## Protection of wound and skin



## Comfeel® Plus Transparent

Comfeel® Plus Transparent is a thin and flexible semi-permeable hydrocolloid dressing designed for visual inspection without disturbing the wound. May be used as primary or secondary dressing.

- Protection of wound and skin
- Visual wound monitoring with mapping grid
- User-friendly design



### Usage

Comfeel Plus Transparent can be used for treatment of non to low exuding chronic wounds and superficial acute wounds in the final stages of wound healing. Comfeel Plus Transparent is indicated for pressure ulcers, leg ulcers, superficial burns, donor sites, traumatic wounds, postoperative wounds, and skin abrasions. The dressing can stay on for up to 7 days depending on the specific wound characteristics.

### Comfeel® Plus Transparent Dressing (Thin)


Code	Size	Units	HCPCS <sup>1</sup>
33530	2 x 2.75" (5 x 7 cm)	10	A6234
33533	4 x 4" (10 x 10 cm)	10	A6234
33536	3.5 x 5.5" (9 x 14 cm)	10	A6235
33539	6 x 6" (15 x 15 cm)	5	A6235

 Exudate Level

## Comfeel® Plus

Comfeel® Plus is a semi-permeable hydrocolloid dressing containing calcium alginate for enhanced absorption capacity. May be used as primary or secondary dressing.

- Protection of wound and skin
- Visual wound monitoring with mapping grid
- User-friendly design



### Usage

Comfeel Plus can be used for treatment of low to moderately exuding chronic wounds and superficial acute wounds in the final stages of wound healing. Comfeel Plus is indicated for pressure ulcers, leg ulcers, superficial burns, donor sites, post-operative wounds, and skin abrasions. The dressing can stay on for up to 7 days depending on the specific wound characteristics.

### Comfeel® Plus Ulcer Dressing

Code	Size	Units	HCPCS <sup>1</sup>
33146	1.5 x 2.5" (4 x 6 cm)	30	A6234
33110	4 x 4" (10 x 10 cm)	10	A6234
33115	6 x 6" (15 x 15 cm)	5	A6235
33120	8 x 8" (20 x 20 cm)	5	A6236
33285	Sacral 7 x 7" (17 x 17 cm)	5	A6235

  Exudate Level

### Comfeel® Plus Contour Dressing

Code	Size	Units	HCPCS <sup>1</sup>
33280	24 sq. in. (6 x 8 cm)	5	A6235
33283	42 sq. in. (9 x 11 cm)	5	A6235

  Exudate Level

# Physiotulle®

Physiotulle® is a non-adherent wound contact layer comprised of a knitted polyester fabric impregnated by hydrocolloid particles (CMC) suspended in petrolatum.

- Easy to handle and easy to use
- Offers comfort and protection
- Does not stick to wound bed for atraumatic removal



## Usage

Physiotulle is indicated for low to highly exuding chronic and acute wounds; superficial burns, donor sites, postoperative wounds, skin abrasions, pressure ulcers and leg ulcers. Can be used with a broad range of absorbent secondary dressings and also under compression therapy. The recommended dressing change interval is approximately every 2 days, depending on the wound type, exudate level and progression.

## Physiotulle® Wound Contact Layer

Code	Size	Units	HCPCS <sup>1</sup>
3910	4 x 4" (10 x 10 cm)	10	A6206
3912	4 x 8" (10 x 20 cm)	3	A6207
3920	6 x 6" (15 x 15 cm)	10	A6207
3915	6 x 8" (15 x 20 cm)	10	A6207

💧—💧💧💧 Exudate Level

# Wound bed preparation: cleansing and debridement



## Sea-Cleans®

Sea-Cleans® saline-based wound cleanser softens and removes necrotic tissue, assisting in maintaining a moist wound environment.

- Utilizes a nonionic surfactant in a non-toxic, non-irritating formula
- Easy-to-use spray bottle with SPRAY and STREAM setting with irrigation pressure
- No-rinse, isotonic, pH-balanced formula

### Usage

Sea-Cleans is indicated to soften and remove necrotic tissue and assist in maintaining a moist wound environment.



### Sea-Cleans® Wound Cleanser

Code	Size	Units	HCPCS <sup>1</sup>
1063	6 fl. oz. (178 mL)	12	A6260
1061	12 fl. oz. (355 mL)	12	A6260

💧—💧💧 Exudate Level

## Purilon®

Purilon® is a gentle hydrogel consisting of natural ingredients without any additives. Purilon provides a unique match of hydrating and absorbing properties to promote a moist wound healing environment that can facilitate autolytic debridement.

- Cohesive and viscous gel to minimize risk of leakage and maceration
- Easy to apply and remove
- Sterile
- Contains alginate



### Usage

Purilon gel is indicated for dry and sloughy necrotic wounds as well as wounds with a mix of necrotic and granulated tissue such as leg ulcers, pressure ulcers, non-infected diabetic foot ulcers, and it may be used on first- and second-degree burns. The gel may be used throughout the healing process to provide a moist healing environment. Purilon gel should be used in conjunction with a secondary dressing and changed at least every 3 days. Do not use with third-degree burns.

### Purilon® Gel

Code	Size	Units	HCPCS <sup>1</sup>
3906	0.28 oz. (8 g)	10	A6248
3900	0.5 oz. (15 g)	10	A6248
3903	0.88 oz. (25 g)	10	A6248

💧—💧 Exudate Level



# Woun'Dres®

Woun'Dres® Collagen Hydrogel promotes a moist wound healing environment to facilitate autolytic debridement and wound healing.

- Hydrates dry wounds and eschar
- Resealable cap and easy application
- Multi-use, cohesive gel



## Usage

Woun'Dres is indicated for superficial wounds and abrasions and minor burns.

Use under the supervision of a healthcare professional for the local management of partial- and full-thickness wounds including pressure and diabetic ulcers; lower extremity ulcers including those of venous, arterial and mixed etiology; surgical wounds; first- and second-degree burns.

## Woun'Dres® Collagen Hydrogel

Code	Size	Units	HCPCS <sup>1</sup>
1166	1 oz. tube (28 g)	36	A6248
7690	3 oz. tube (84 g)	12	A6248

 Exudate Level

## Brief Statements

### BIATAIN® AG FOAM ADHESIVE AND NON-ADHESIVE DRESSINGS

Biatain® Ag Foam Adhesive and Non Adhesive Dressings are indicated for use in patients with moderately to highly exuding leg ulcers and pressure sores, second degree burns, donor sites, post-operative wounds and skin abrasions. Biatain® Ag Foam Non-Adhesive Dressings are also indicated for diabetic foot ulcers. Patients being treated for wounds should follow the advice of, and direct questions about use of the product to, their medical professional. Before using the device, carefully read the product labels and information accompanying the device including the instructions for use which contain additional safety information. Coloplast's wound dressings are for single-use only; discard it after use. The use of Biatain Ag dressings during pregnancy, lactation and on children has not been demonstrated. Biatain Ag foam dressings should not be used on patients with a known sensitivity to silver. If you experience symptoms of an allergic reaction, contact your healthcare professional. The risk information provided here is not comprehensive. To learn more, talk to your healthcare provider. For further information, call Coloplast Corp. at 1-800-533-0464 and/or consult the company website at [www.coloplast.us](http://www.coloplast.us).

### BIATAIN® ALGINATE AG DRESSINGS

Biatain® Alginate Ag Dressings are indicated for wounds with moderately to heavily exuding partial to full thickness wounds such as post-operative wounds, trauma wounds, leg ulcers, pressure ulcers, diabetic ulcers, post-operative wounds and skin abrasions with delayed healing due to bacteria, or where there is a risk of infection. The dressings may be used to support moist wound healing in patients who are in treatment for a local or systemic infection under the discretion of a healthcare professional; are suitable for use in combination with compression. Patients being treated for wounds should follow the advice of, and direct questions about use of the product to, their medical professional. Before using the device, carefully read the product labels and information accompanying the device including the instructions for use which contain additional safety information. Coloplast's wound dressings are for single-use only; discard it after use. The use of Biatain Ag dressings during pregnancy, lactation and on children has not been demonstrated. Biatain Ag foam dressings should not be used on patients with a known sensitivity to silver. If you experience symptoms of an allergic reaction, contact your healthcare professional. The risk information provided here is not comprehensive. To learn more, talk to your healthcare provider. For further information, call Coloplast Corp. at 1-800-533-0464 and/or consult the company website at [www.coloplast.us](http://www.coloplast.us).

### TRIAD® HYDROPHILIC WOUND DRESSING

Triad™ hydrophilic wound dressing is indicated for the local management of pressure and venous stasis ulcers, dermal lesions/injuries, superficial wounds, scrapes, first- and second-degree burns, partial- and full-thickness wounds. Suitable for use on neonates. Triad™ is contraindicated for third degree burns or infected wounds. See device label for complete information including warnings and precautions. For further information call Coloplast Corp. at 1-800-533-0464 and/or consult the company website at [www.coloplast.us](http://www.coloplast.us).

Notes:

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

Notes:

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

Coloplast develops products and services that make life easier for people with very personal and private medical conditions. Working closely with the people who use our products, we create solutions that are sensitive to their special needs. We call this intimate healthcare.

Our business includes Ostomy Care, Continence Care, Wound and Skin Care and Urology Care. We operate globally and employ about 11,000 employees.

#### References

1. PDAC verification as of material print date of 11.18. Refer to [www.dmepdac.com](http://www.dmepdac.com) for most up to date information.
2. Segal et al. Journal of Biomaterials Applications 1998;249.
3. Christiansen et al. Poster presented at EWMA 2018.

Ostomy Care / Continence Care / Wound & Skin Care /  
Interventional Urology

Coloplast Corp. Minneapolis, MN 55411 / 1-800-533-0464

[www.coloplast.us](http://www.coloplast.us) The Coloplast logo is a registered trademark of Coloplast A/S.

© 2018 Coloplast Corp. All rights reserved.



PM-06123 11.18