

Ostomy Product Guide

Effective December 2022



Stina

Making life easier_

Welcome to Coloplast

Coloplast is a leader in the market within intimate healthcare. For more than sixty years, we've worked closely with healthcare professionals and users to develop products and services that make life easier for people with intimate healthcare needs.

Our Mission

***Making life easier for
people with intimate
healthcare needs***

Table of Contents

Accessories.....	7-16
Irrigation.....	17-18
2-Piece Click.....	19-40
2-Piece Flex.....	41-50
1-Piece.....	51-80
Hospital Assortment / Wound Management	81-102
Pediatric	103-114
Reference Tools.....	115-124
Education	
Ostomy Care	115-118
Ostomy Education.....	119
Medicare Utilization	
Ostomy.....	120-121
HCPCS Codes	122-124
Index.....	125-132

Through BodyFit Technology, we offer a secure fit to real bodies *to help prevent leakage and skin issues*

Leakage and skin issues are closely related challenges. The **right fit** helps **prevent** both.



Leakage

76%

of people with a stoma have experienced leakage within the last six months.¹

Leakage

is often caused by a poor fit to the body.²



Skin issues

Skin issues

There is a significant correlation between leakage frequency and the severity of the peristomal skin condition.²

61%

of users attributed their most recent skin issue to leakage under the barrier.³

1. Results from the 2014 Ostomy Life Study - an international survey aimed to better understand the challenges people living with an ostomy face (n=4138 ostomates). Published: Claessens I, Probert R, et. al. The Ostomy Life Study: the everyday challenges faced by people living with a stoma in a snapshot. *Gastrointestinal Nursing*. 2015;13(5), 18-25.
2. Results from a multinational, randomized controlled cross-over trial comparing SenSura Mio Convex to the study participants' own appliance (n=129). Published: Walker H, Hopkins G, et. al. Raising the bar: New flexible convex ostomy appliance - a randomised controlled trial. *World Council of Enterostomal Therapists Journal*. 2016;36(1) S6-S11.
3. Results from the 2018-2019 Ostomy Life Study - an international survey (n=3112 ostomates); Coloplast market research, Data on file.

The #1 ostomy brand in the world, **SenSura® Mio**, is now **Extended Wear***!



To learn more, visit us at www.sensuramio.us or call 1-888-726-7872.

SenSura Mio was designed to fit individual bodies – and now, a recent 2,359-patient survey demonstrated that SenSura Mio users report less leakage onto their clothes than other extended wear products, while experiencing similar or even longer wear time using SenSura Mio¹.

While SenSura Mio was initially launched into the standard wear category, recent user experience¹, clinician evaluation², and lab-based data³ reinforces SenSura Mio's performance as an extended wear product.

With this change, users who need or would prefer an extended wear product will now have access to SenSura Mio's unique technology, designed to reduce the worry about leakage & provide a comfortable fit:

- Full-circle filter to reduce ballooning
- Discreet neutral-gray, soft woven textile
- Wide outlet
- A flexible fit, thanks to BodyFit Technology, an elastic adhesive that is unique in all SenSura Mio barriers.

*Excluding 1-piece closed-end pouches

1. Ostomy pouch user survey, 2021, n=2,359. Dataset collected by Coloplast but includes all brands. Data on file.

2. Product evaluation of SenSura Mio Convex, 2016, n=98 clinicians.

3. Data from 48h bench testing performed by Coloplast; may not be indicative of clinical performance. Data on file.

Products of Distinction



SenSura® Mio Convex Flip

Designed to fit curves, bulges and hernias with less creasing and folding



SenSura® Mio Flex

- Floating flange
- Easy to apply
- Flexible coupling



Soft
6mm

Light
6mm

Deep
9mm

SenSura® Mio Convex

Three types and two depths to fit different bodies



SenSura® Mio High Output

The only 1-piece convex option with a spout outlet for fecal output. Also available in 1-piece flat, 2-piece Click, and 2-piece Flex

Brava[®]

Accessories

Security to live *everyday life*



“I was determined to play sports again. The Brava[®] Protective Seal gave me confidence.”

– Kevin*

*Kevin is a Brava user who has received compensation from Coloplast to provide this information.

Brava® supporting products

For customizing pouching system fit and for issues such as leakage and skin irritation.

NEW! Brava® Protective Seal Convex

The new Brava Protective Seal Convex is designed to protect the skin against leakage for flush or retracted stomas and patients with folds, creases or scarring around their stoma.



Brava® Protective Seal

- Easily shaped
- Leaves minimal residue behind
- Follows body movements
- Significantly reduces leakage compared to competitor rings*

The **Brava® Ostomy Support Belt** brings relief by relieving the heavy sensation of hernias.

- Comfortable and breathable fabric
- Stretchable to fit body contours
- Easy-to-use pocket closure
- Silicone grip to prevent sliding and rolling



*Based on bench test data. May not be indicative of clinical performance.

Think again. Think Brava®

Designed to address leakage and care for your skin.

Brava® Elastic Barrier Strips

Skin-friendly strips are used to “picture frame” around an ostomy barrier, designed for extra security.

- Elastic (follow the contours of the body)
- Ensure the position of the barrier and designed to minimize roll-up
- Skin-friendly and secure alternative to tape
- **Not made with natural rubber latex**



Brava® Elastic Barrier Strips- Curved

For round or oval barrier shapes

Code	Size	Units	HCPCS*
120700**	5½" strips	20 strips	A4362

**2 strips = 1 billable unit of A4362

Brava® Elastic Barrier Strips- Y shape

Designed for users with challenging body shapes, e.g. hernias or skin folds



Code	Size	Units	HCPCS*
120721**	5¾" strips	30 strips	A5121

**3 strips = 1 billable unit of A5121

Brava® Elastic Barrier Strips- Straight

For square or rectangular barrier shapes



Code	Size	Units	HCPCS*
120740**	5¼" strips	20 strips	A4362

**2 strips = 1 billable unit of A4362

Brava® Elastic Barrier Strip XL

Designed for increased security & support around the barrier. Wider shape accommodates the petals of Convex Flip.



Code	Size	Units	HCPCS*
12076**	6¾" strips	30 strips	A5121

**3 strips = 1 billable unit of A5121

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

Brava® Protective Seal

Moldable barrier ring designed to adhere securely and protect the skin around the stoma.

- Soft, flexible material adapts to body shape as the wearer bends and stretches
- Can be easily and cleanly removed
- Alcohol-free and sting-free
- Multiple sizes for a custom fit
- **Not made with natural rubber latex**



Code	Stoma size	Adhesive diameter	Thickness	Seal width	Units	HCPCS*
12035	3/4"	1 7/8"	2.5 mm	Standard (5/8")	10	A4385
12036	3/4"	2 1/4"	2.5 mm	Wide (3/4")	10	A4385
12038	3/4"	2 1/2"	2.5 mm	Wide (7/8")	10	A4385
12037	1 1/8"	2 1/4"	2.5 mm	Standard (5/8")	10	A4385
12039	1 3/8"	2 1/2"	2.5 mm	Standard (5/8")	10	A4385
12045	3/4"	1 7/8"	4.2 mm	Standard (5/8")	10	A4385
12046	3/4"	2 1/4"	4.2 mm	Wide (3/4")	10	A4385
12048	3/4"	2 1/2"	4.2 mm	Wide (7/8")	10	A4385
12047	1 1/8"	2 1/4"	4.2 mm	Standard (5/8")	10	A4385
12049	1 3/8"	2 1/2"	4.2 mm	Standard (5/8")	10	A4385
12032	3/4"	3"	2.5 mm	Extra wide (1 1/8")	10	A4385
12033	1 1/8"	3"	2.5 mm	Extra wide (1 5/16")	10	A4385

Brava® Protective Seal Convex

Moldable barrier ring designed to provide convexity for flush or retracted stomas and patients with folds, creases or scarring around their stoma.

- High-tack adhesive that flows quickly into skin, filling small creases
- Erosion resistant composition
- High flexibility to maintain seal during body movement
- Steep sides designed to provide convexity to peristomal area
- Thin edges fit snugly around peristomal area
- **Not made with natural rubber latex**



Code	Stoma size	Thickness	Units	HCPCS*
12090	3/4"	8.2mm	10	A4411
12091	1"	8.2mm	10	A4411
12092	1 3/16"	8.2mm	10	A4411
12093	1 3/8"	8.2mm	10	A4411
12094	1 9/16"	8.2mm	10	A4411
12095	1 3/4"	8.2mm	10	A4411

Brava® Moldable Ring

Designed for dual-action leakage reduction - used to fill in uneven skin surfaces and protect the skin.

- Molds and adheres securely to dips, creases, and folds to help create a tight seal around the stoma
- Helps keep skin healthy by protecting it from output and absorbing moisture from the skin
- Easy to handle, mold and apply
- Alcohol-free and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
120307	2 mm (thin)	10	A4385
120427	4.2 mm (thick)	10	A4385

Brava® Strip Paste

Pectin-based product that protects peristomal skin.

- Protects peristomal skin from stoma output
- Used to fill uneven skin surfaces to create a flat pouching surface
- Alcohol-free and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
26555	2 oz. total package (.2 oz. per strip)	10 strips	A4406

Brava® Paste

A sting-free alternative to conventional alcohol-based paste.

- Creates a tight, effective seal between peristomal skin and the ostomy/fistula appliance
- **No alcohol** = no need to "set"; **contains pectin**
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
12050	2 oz.	1	A4406

Coloplast® Paste

Used to fill in uneven skin surfaces to create a flat pouching surface.

- Provides tight sealing around the stoma
- Protects against stoma output, minimizing risk of maceration
- **Low alcohol content; no pectin**
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
2650	2 oz.	12	A4405

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

Brava® Skin Cleanser Wipes

The wipe that cleans and cares. A convenient way to cleanse and maintain healthy skin when changing the ostomy pouching system.

- Contains ProVitamin B5 to provide a soothing and moisturizing effect
- Does not contain alcohol or perfumes
- No residue left on the skin that could affect barrier adhesion
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
12080	1 wipe	15 wipes/pack	A9286**

Brava® Lubricating Deodorant

A pouch lubricant + deodorizer in one. Powder-fresh liquid neutralizes odor and allows output to slide to the bottom of the pouch for easier emptying.

- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
12060	.25 fl. oz. (7.5 mL) sachet	20 sachets	A4394
12061	8 oz. bottle	16	A4394

Brava® Adhesive Remover Spray

- Designed to lift the edges of the barrier for ease in removal
- Dries within seconds without leaving an oily residue behind
- Alcohol-free, silicone-based, and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
120105	1.7 oz. (50 mL)	1	A4455**

Brava® Adhesive Remover Wipes

- Wipe away adhesive residue
- Dries within seconds without leaving an oily residue behind
- Alcohol-free, silicone-based, and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
120115	1 wipe	30 wipes/box	A4456**

Brava® Skin Barrier Spray

- Protects against stoma output and adhesives
- Dry within seconds without building a thick or tacky layer
- Alcohol-free, silicone-based, and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
120205	1.7 oz. (50 mL)	1	A4369**

Brava® Skin Barrier Wipes

- Protects against stoma output and adhesives
- Dry within seconds without building a thick or tacky layer
- **Tip:** wring the wipe like a wet rag to expel the liquid, then wipe to disperse on skin
- Alcohol-free, silicone-based, and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
120215	1 wipe	30 wipes/box	A5120**

Brava® Protective Sheet

Hydrocolloid skin barrier that protects the skin from stoma output.

- Offers effective protection for peristomal skin and "3-way stretch"
- Provides a ready surface for good adhesion
- Adheres to moist skin
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
32105	4 x 4" (10 x 10 cm)	10	A4362
32155	6 x 6" (15 x 15 cm)	5	A5121
32205	8 x 8" (20 x 20 cm)	5	A5122

Brava® Powder

Non-sterile protective powder for moisture absorption.

- May use in conjunction with Brava® Skin Barrier Spray or PREP™ to "crust", thereby protecting moist skin
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
19075	1 oz. (25 g)	16	A4371

Hex-On®

Odor Antagonist

- Helps control airborne odors
- Non-aerosol spray
- **Not made with natural rubber latex**



Product	Code	Size	Units
Fresh Linen	7583	2 fl. oz./59 mL	12
Fragrance Free	67589	2 fl. oz./59 mL	12

Ostomy Scissors

Curved ostomy scissors.

- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
95050	One size	30/bx	N/A

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

**Product has not been through PDAC verification

Coloplast® Skin Barrier Rings

Soft and flexible barrier rings provide skin protection for stomas, fistulas, and draining wounds.

- Non-pectin based
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
2310	3/8" (10 mm)	30	A4404
2315	3/5" (15 mm)	30	A4404
2320	3/4" (20 mm)	30	A4404
2325	1" (25 mm)	30	A4404
2330	1 1/8" (30 mm)	30	A4404
2340	1 3/5" (40 mm)	30	A4404
2350	2" (50 mm)	30	A4404

Prep™ Protective Skin Barrier

Clear film forms a breathable barrier protecting against skin irritants such as adhesives, urine and feces. No need to remove between applications.

- **Not made with natural rubber latex**



Wipes

Code	Size	Units	HCPCS*
2041	54 wipes/pack	12 packs/box	A5120

Dabber

Code	Size	Units	HCPCS*
0925	2 fl. oz. (59 mL)	12	A4369

Pouch Filters (Filtrodor®)

Charcoal-activated, deodorizing filter with self-sealing foam center; adheres to pouch.

- **Not made with natural rubber latex**



Code	Units	HCPCS*
0509	50 <i>each box contains 5 packs of 10</i>	A4368

Drainable Pouch Clamp

Simple closure clamp is easy to manipulate, even for those with low dexterity.

- **Not made with natural rubber latex**



Code	Units	HCPCS*
9500	20	A4363

Brava® Ostomy Support Belt

This belt relieves the heavy sensation of hernias and may help prevent hernias post-operatively. Comfortable and breathable fabric, stretchable to fit body contours, intuitive pocket closure, and includes a silicone grip to prevent sliding and rolling.



- Not made with natural rubber latex
- All sizes are an 8" width

Code	Size	Units	HCPCS*
12003	S (29-32")	1	A4396
12004	M (31-35")	1	A4396
12005	L (34-39")	1	A4396
12006	XL (37-42")	1	A4396
12007	2XL (40-46")	1	A4396
12008	3XL (44-51")	1	A4396
12009	4XL (50-58")	1	A4396

Brava® Belt for SenSura® Mio

Neutral gray color was designed to blend in under clothing. Four-tab connection provides extra stability across the entire peristomal area. Designed specifically for use with SenSura® Mio products with belt tabs.



- Not made with natural rubber latex

Code	Size	Units	HCPCS*
4237	40" length	15	A4367
4247	61" length	15	A4367

Brava® Belt

Neutral gray color was designed to blend in under clothing. Comfortable, adjustable belt offers extra security for people with a flush or retracted stoma by supporting a convex barrier and allowing the stoma to protrude.



- Single tab attachment is designed to be compatible with most conventional ostomy products with belt tabs.
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
4215	43 $\frac{2}{3}$ " length	15	A4367
4220	49" length	15	A4367

Urostomy Night Bag

Includes a fixed connector that fits directly to SenSura® Mio, SenSura® and Assura® urostomy pouches. The bag holds up to 2 liters and has a 120 cm long flexible tube with anti-reflux valve. Bag hanger (#5070) not included. Order separately.



- Not made with natural rubber latex

Code	Units	HCPCS*
21365	10	A4357

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

Bed Drainage Bag

Wide-lumen bag tubing specifically designed for fecal output. Holds up to 2000 mL. Long tubing can be cut easily and allows nurses to find the right length for bed drainage. Can be used with fistula / wound management systems or high-output pouches. Bag hanger (#5070) not included. Order separately.



- **Not made with natural rubber latex**

Code	Units	HCPCS*
14010	6	A4357

Bag Hanger

Reusable, sturdy, plastic drainage bag hanger secures tubing and ensures bag stability.



- **Not made with natural rubber latex**

Code	Size	Units	HCPCS*
5070	One Size	1	Medicare N/A Other A4421

Extra Drain Port

For use with Fistula Wound Management Systems (page 89). Drain port can be easily applied to the lid to allow for catheter insertion. For catheters ranging from 18FR - 22FR, the drain port is ready for use. For catheters ranging from 24FR-30FR, pull off the top of the drain port - DO NOT CUT.



- **Not made with natural rubber latex**

Code	Units
14015	12

Coloplast® Transparent Fistula Lids

Replacement lids match the transparent lid included with each fistula pouch.



- **Not made with natural rubber latex**

Code	Size	Units
14021	Mini (matches 14050)	10
14031	Midi (matches 14060)	10
14041	Maxi (matches 14070)	10

Irrigation Sets

Coloplast® Economy Version

Durable water bag with integrated thermometer, bold print, large opening. Can be used with 1- or 2-piece ostomy systems.

- **Not made with natural rubber latex**

Code	Units	Includes	HCPCS*
1500	1	1- Water Bag/Regulator	A4398
		1- Cone (#1110)	A4399

Assura® Hospital Version

Durable water bag with integrated thermometer, bold print, large opening. Can be used with 1- or 2-piece ostomy systems.

- **Not made with natural rubber latex**

Code	Units	Includes	HCPCS*
1500H	1	1- Water Bag/Regulator	A4398
		1- Cone (#1110)	A4399
		4- Adhesive Irrigation Sleeves (#1003)	A4437
		2- Closure Clips (#9500)	



SenSura® or Assura® Deluxe Version

Components supplied in a travel bag.

- **Not made with natural rubber latex**

Code	Units	Includes	HCPCS*
12830	1	1- Water Bag/Regulator	A4398
		1- Cone (#1110)	A4399
		1- Irrigation Faceplate (#12820)	A4361
		1- Belt (#4215)	A4367
		2- Irrigation Sleeves (#12836)	A4436

Irrigation Sleeves

Irrigation Sleeves (with adhesive)

Disposable irrigation sleeve (adhesive-backed).

- **Not made with natural rubber latex**

Transparent

Code	Units	HCPCS*
1003	10	A4437






*For important information regarding the use of HCPCS codes listed and code definitions see page 120




**Irrigation Sleeves
(without adhesive)**



Irrigation sleeves for use with SenSura® Click or Assura® barriers or faceplate. Consult chart below to order the correct transparent sleeves.

• **Not made with natural rubber latex**

Code	Size	Units	HCPCS*
12834		5	A4436
12835		5	A4436
12836		5	A4436

Irrigation Sleeve #	Use with SenSura® Click Product Code	Units
12834	 - 10011, 10012, 10013, 10015-10017, 11011, 11012, 11015, 11016	5
12835	 - 10021, 10022, 10023, 10025-10027, 11021-11025	5
12836	 - 10031, 10032, 10035, 10036, 11031, 11032, 11035, 11036	5

Irrigation Accessories

Assura® Irrigation Faceplate

For use with ostomy belt 4215 or 4220.
(Use with #12836 Irrigation Sleeve) Note: not compatible with Brava® Belt for SenSura® Mio



• **Not made with natural rubber latex**

Code	Units	HCPCS*
12820	1	A4361

Coloplast Stoma Cone

Soft and flexible stoma cone with tubing.



• **Not made with natural rubber latex**

Code	Units	HCPCS*
1110	1	A4399



2-piece Click

(Mechanical coupling)



Pouch can be repositioned while it's still on the body



Tactile experience of pressing two pieces together can be reassuring



Lock ring gives an audible cue that it's securely locked

SenSura[®] Mio

2-Piece Click

The fit is
the difference
between sitting on the
edge and diving right in.

Mike is a SenSura[®] Mio Click user and has received compensation from Coloplast.

SenSura[®] Mio Click

BodyFit Technology is an elastic adhesive barrier found in all SenSura[®] Mio Barriers. This unique barrier quickly adheres to the skin, flows into creases and folds, is highly stretchable, and flexes with the body. This may help prevent channels and gaps and may result in reducing leakage and skin issues.



NOTE: SenSura® Mio Click is NOT compatible with SenSura® Click or Assura® systems.

SenSura® Mio Click Barrier

- Extended wear
- Flat
- Belt tabs for use with 4237, 4247
- **Not made with natural rubber latex**



Code	Stoma Size	Color	Units: 5 HCPCS*
Cut-to-fit			
10502	3/8-1 3/8" (10-35 mm)		A4409
10512	3/8-1 3/4" (10-45 mm)		A4409
10522	3/8-2 1/8" (10-55 mm)		A4410
10531	3/8-2 9/16" (10-65 mm)		A4410
Pre-cut			
10503	3/4" (20 mm)		A4409
10506	7/8" (22 mm)		A4409
10504	1" (25 mm)		A4409
10505	1 1/8" (28 mm)		A4409
10513	1 3/16" (30 mm)		A4409
10514	1 1/4" (32 mm)		A4409
10515	1 3/8" (35 mm)		A4409
10516	1 1/2" (38 mm)		A4409
10523	1 9/16" (40 mm)		A4409
10524	1 3/4" (45 mm)		A4409
10525	2" (50 mm)		A4409

SenSura® Mio Click Convex Light Barrier

- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Not made with natural rubber latex**



Code	Stoma Size	Color	Units: 5 HCPCS*
Cut-to-fit			
16901	3/8 - 3/4" (10-20 mm)		A4407
16911	5/8 - 1 3/16" (15-30 mm)		A4407
16921	5/8 - 1 9/16" (15-40 mm)		A4407
16931	5/8 - 2 1/16" (15-53 mm)		A4408
Pre-cut			
16903	1 1/16" (18 mm)		A4407
16913	1 3/16" (21 mm)		A4407
16914	1" (25 mm)		A4407
16915	1 1/8" (28 mm)		A4407
16923	1 1/4" (31 mm)		A4407
16924	1 3/8" (35 mm)		A4407
16933	1 1/2" (38 mm)		A4408
16934	1 5/8" (41 mm)		A4408













Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

SenSura® Mio Click Deep Convex Barrier 

- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Not made with natural rubber latex**

















Code	Stoma Size	Color	Units: 5 HCPCS*
Cut-to-fit			
16941	3/8 - 3/4" (10-20 mm)		A4407
16951	5/8 - 1 1/16" (15-30 mm)		A4407
16961	5/8 - 1 9/16" (15-40 mm)		A4407
16971	5/8 - 2 1/16" (15-53 mm)		A4408
Pre-cut			
16942	1 1/16" (18 mm)		A4407
16953	1 3/16" (21 mm)		A4407
16954	1" (25 mm)		A4407
16955	1 1/8" (28 mm)		A4407
16963	1 1/4" (31 mm)		A4407
16964	1 3/8" (35 mm)		A4407
16973	1 1/2" (38 mm)		A4408
16974	1 5/8" (41 mm)		A4408

SenSura® Mio Click Convex Flip Barrier 

- Extended wear
- Designed for curves, bulges and hernias
- Curved, star-shaped barrier conforms to outward areas with less creasing and folding
- "Flip to fit" design allows for easy application
- Belt tabs for use with 4237, 4247
- **Not made with natural rubber latex**



Code	Stoma Size	Color	Units: 5 HCPCS*	
Cut-to-fit				
18500	3/8 - 1 3/8" (10-35 mm)		A4407	
18510	3/8 - 1 3/4" (10-45 mm)		A4408	
18520	3/8 - 2 1/8" (10-55 mm)		A4408	
18530	3/8 - 2 9/16" (10-66 mm)		A4408	
Code	Stoma Size	Max Cut	Color	HCPCS*
Pre-cut				
18501	3/4" (20 mm)	1 3/8" (35 mm)		A4407
18511	1" (25 mm)	1 3/4" (45 mm)		A4407
18512	1 3/16" (30 mm)	1 3/4" (45 mm)		A4407
18513	1 3/8" (35 mm)	1 3/4" (45 mm)		A4407
18521	1 3/16" (30 mm)	2 1/8" (55mm)		A4407
18522	1 3/8" (35 mm)	2 1/8" (55 mm)		A4407
18523	1 9/16" (40 mm)	2 1/8" (55 mm)		A4407
18524	1 3/4" (45 mm)	2 1/8" (55 mm)		A4407
18531	1 3/4" (45 mm)	2 9/16" (66 mm)		A4408
18532	2" (50 mm)	2 9/16" (66 mm)		A4408

NOTE: SenSura® Mio Click is NOT compatible with SenSura® Click or Assura® systems.

SenSura® Mio Click MAXI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo inspection window for easy application
- **Not made with natural rubber latex**



11" (28 cm) / 655 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent			
Opaque			
11463	11462		A4427
11473	11472		A4427
11483	11482		A4427
11487	11486		A4427

SenSura® Mio Click MAXI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Without filter**
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**



11" (28 cm) / 655 mL

Units: 10

Code	Color	HCPCS*
Transparent		
11474		A4426
11484		A4426
11488		A4426

SenSura® Mio Click MIDI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window for easy application
- **Not made with natural rubber latex**



10¾" (26 cm) / 485 mL

Units: 10

Code	Color	HCPCS*
Opaque		
11461		A4427
11471		A4427



Indicates convexity



Indicates Convex Flip

*For important information regarding the use of HCPCS codes listed and code definitions see page 120





SenSura® Mio Click High Output Pouch with Tap Outlet

- Special full-circle filter designed specifically for ileostomy output
- Neutral gray, water-resistant textile
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- **Not made with natural rubber latex**



12½" (32 cm) / 1080 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
18610	18600		A4413
18611	18601		A4413
18612	18602		A4413
18613	18603		A4413




SenSura® Mio Click High Output Pouch with Soft Outlet

- Special full-circle filter designed specifically for ileostomy output
- Neutral gray, water-resistant textile
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Soft outlet provides effective drainage; connects to a standard drain bag or Coloplast fistula drain bag #14010
- **Not made with natural rubber latex**



12 ½" (32 cm) / 1080 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
18640	18630		A4413
18641	18631		A4413
18642	18632		A4413




SenSura® Mio Click Drainable Pouch with Tap Outlet

- Special full-circle filter designed specifically for ileostomy output
- Neutral gray, water-resistant textile
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- **Not made with natural rubber latex**



9" (23 cm) / 650 mL

Units: 10

Code	Color	HCPCS*
Opaque		
18620		A4427
18621		A4427
18622		A4427





SenSura® Mio Click MAXI Closed Pouch

- Full-circle filter
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo inspection window for easy application
- **Not made with natural rubber latex**



8¼" (20½ cm) / 570 mL

Units: 30

Code	Code	Color	HCPCS*
Transparent			
Opaque			
N/A	11402		A4423
N/A	11413		A4423
N/A	11422		A4423
11426	11425		A4423




SenSura® Mio Click MIDI Closed Pouch

- Full-circle filter
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window for easy application
- **Not made with natural rubber latex**



7¼" (18½ cm) / 275 mL

Units: 30

Code	Color	HCPCS*
Opaque		
11401		A4423
11412		A4423
11421		A4423


*For important information regarding the use of HCPCS codes listed and code definitions see page 120

Urostomy

SenSura® Mio Uro

Wide anti-reflux valves support varying ways of managing stents:

- Stents threaded into a pouch connected to a night drainage bag
- Rolling the stents in the top of the pouch, making a loop
- Stents that have been cut



Welded "multi-chambers" minimize ballooning




SenSura® Mio Click MAXI Urostomy Pouch

- Multi-chambers minimize sloshing
- Wider anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**



10½" (26½ cm) / 480 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
11492	11491		A4433
11498	11497		A4433
11500	11499		A4433

NOTE: SenSura® Mio Click is NOT compatible with SenSura® Click or Assura® systems. SenSura® Click and Assura® are compatible.

SenSura® Click Flat Barrier



- Standard wear
- Belt tabs for use with 4215, 4220
- **Not made with natural rubber latex**








Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit

10011	3/8 - 1 3/8" (10-35 mm)		A4414
10021	3/8 - 1 3/4" (10-45 mm)		A4414
10031	3/8 - 2 1/8" (10-55 mm)		A4414
10041	3/8 - 2 1/2" (10-65 mm)		A4415**

Pre-cut

10012	1 3/16" (20 mm)		A4414
10013	1" (25 mm)		A4414
10022	1 3/16" (30 mm)		A4414
10023	1 3/8" (35 mm)		A4414
10032	1 9/16" (40 mm)		A4414

SenSura® Click Xpro Flat Barrier



- Extended wear
- Belt tabs
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit

10015	3/8 - 1 3/8" (10-35 mm)		A4409
10025	3/8 - 1 3/4" (10-45 mm)		A4409
10035	3/8 - 2 1/4" (10-55 mm)		A4409
10045	3/8 - 2 1/2" (10-65 mm)		A4410

Pre-cut

10016	1 3/16" (20 mm)		A4409
10017	1" (25 mm)		A4409
10026	1 3/16" (30 mm)		A4409
10027	1 3/8" (35 mm)		A4409
10036	1 9/16" (40 mm)		A4409










*For important information regarding the use of HCPCS codes listed and code definitions see page 120

**Product has not been through PDAC verification

SenSura® Click Convex Light Barrier 

- Standard wear
- Belt tabs
- Not made with natural rubber latex












Code	Stoma Size	Color	Units: 5 HCPCS*
Cut-to-fit			
11011	5/8 - 7/8" (15-22 mm)		A4373
11021	5/8 - 1 5/16" (15-33 mm)		A4373
11031	5/8 - 1 11/16" (15-43 mm)		A4373
11041	5/8 - 2 1/16" (15-53 mm)		A4373**
Pre-cut			
11012	7/8" (22 mm)		A4373
11022	1" (25 mm)		A4373
11023	1 1/8" (28 mm)		A4373
11024	1 1/4" (31 mm)		A4373
11032	1 3/8" (35 mm)		A4373

SenSura® Click Xpro Convex Light Barrier 

- Extended wear
- Belt tabs
- Not made with natural rubber latex



Code	Stoma Size	Color	Units: 5 HCPCS*
Cut-to-fit			
11015	5/8 - 7/8" (15-22 mm)		A4407
11025	5/8 - 1 5/16" (15-33 mm)		A4407
11035	5/8 - 1 11/16" (15-43 mm)		A4407
11045	5/8 - 2 1/16" (15-53 mm)		A4408
Pre-cut			
11016	7/8" (22 mm)		A4407
11026	1" (25 mm)		A4407
11027	1 1/8" (28 mm)		A4407
11028	1 1/4" (31 mm)		A4407
11036	1 3/8" (35 mm)		A4407

For Assura® CONVEX barriers, see pages 33-34.





SenSura® Click MAXI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11½" (29½ cm) / 655 mL

Units: 20

Code	Code	Color	HCPCS*
Transparent	Opaque		
11194	11124		A4427
11195	11125		A4427
11196	11126		A4427
11197	11127		A4427**



SenSura® Click MIDI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**



10½" (26½ cm) / 485 mL

Units: 20

Code	Color	HCPCS*
Opaque		
11114		A4427
11115		A4427



SenSura® Click MAXI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Without filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11½" (29½ cm) / 655 mL

Units: 20

Code	Code	Color	HCPCS*
Transparent	Opaque		
11185	11135		A4426
11186	11136		A4426

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

**Product has not been through PDAC verification

NOTE: SenSura® Click is compatible with Assura® but NOT with SenSura® Mio Click.




SenSura® Click Magnum High Output Drainable Pouch

- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



12½" (32 cm) / 1160 mL

Units: 10

Code	Code	Color	HPCPS*
Transparent	Opaque		
19050	19040		A4413**
19051	19041		A4413**
19052	19042		A4413**

For Assura® Ileo Night Bag, see page 38.





SenSura® Click MAXI Closed Pouch

- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



8½" (21 cm) / 570 mL

Units: 30

Code	Code	Color	HPCPS*
Transparent	Opaque		
NA	10164		A4423
NA	10165		A4423
NA	10166		A4423
10187	10167		A4423**




SenSura® Click MIDI Closed Pouch

- **Filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**



7" (18 cm) / 275 mL

Units: 30

Code	Color	HPCPS*
Opaque		
10154		A4423
10155		A4423
10156		A4423

For Assura® stoma caps, see page 39.

Urostomy

SenSura® Click Uro MAXI Urostomy Pouch

- Multi-chambers minimize sloshing
- Wider anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



10³/₈" (26 cm) / 480 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
11854	11844		A4433
11855	11845		A4433
11856	11846		A4433

SenSura® Click Uro MIDI Urostomy Pouch

- Multi-chambers minimize sloshing
- Wider anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



9¹/₂" (24 cm) / 430 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
11851	11841		A4433

For Assura® urostomy micro-pouches and minicaps, see page 39.

Urostomy Night Bag and Bag Hanger

See Accessories page 16 for more information.



NOTE: Assura® is compatible with SenSura® Click , but **NOT** with SenSura® Mio Click.










*For important information regarding the use of HCPCS codes listed and code definitions see page 120

**Product has not been through PDAC verification

Assura® Flat Barrier

- Standard wear
- Belt tabs for use with 4215, 4220
- **Not made with natural rubber latex**













Code	Stoma Size	Color	Units: 5 HCPCS*
Cut-to-fit			
2881	3/8 - 1 3/8" (10-35 mm)		A4414
2882	3/8 - 1 3/4" (10-45 mm)		A4414
2883	3/8 - 2 1/4" (10-55 mm)		A4414
Pre-cut			
12841	1 3/16" (20 mm)		A4414
12842	1" (25 mm)		A4414
12844	1 3/16" (30 mm)		A4414
12845	1 3/8" (35 mm)		A4414
12847	1 9/16" (40 mm)		A4414
12848	1 3/4" (45 mm)		A4414

Assura® Flat Barrier

- Extended wear
- Belt tabs for use with 4215, 4220
- **Not made with natural rubber latex**



Code	Stoma Size	Color	Units: 5 HCPCS*
Cut-to-fit			
2831	3/8 - 1 3/8" (10-35 mm)		A4409
2832	3/8 - 1 3/4" (10-45 mm)		A4409
2833	3/8 - 2 1/4" (10-55 mm)		A4409
Pre-cut			
14250	5/8" (15 mm)		A4409
14251	1 3/16" (20 mm)		A4409
14252	1" (25 mm)		A4409
14254	1 3/16" (30 mm)		A4409
14255	1 3/8" (35 mm)		A4409
14257	1 9/16" (40 mm)		A4409
14258	1 3/4" (45 mm)		A4409












NOTE: Assura® is compatible with SenSura® Click, but NOT with SenSura® Mio Click.

Assura® Convex Light Barrier 

- Standard wear
- Belt tabs for use with 4215, 4220
- **Not made with natural rubber latex**



Units: 5












Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
14261	5/8 - 15/16" (15-23 mm)		A4373
14262	5/8 - 15/16" (15-33 mm)		A4373
14263	5/8 - 111/16" (15-43 mm)		A4373
Pre-cut			
14271	3/4" (18 mm)		A4373
14272	7/8" (22 mm)		A4373
14273	1" (25 mm)		A4373
14274	1 1/8" (28 mm)		A4373
14275	1 1/4" (31 mm)		A4373
14276	1 3/8" (35 mm)		A4373
14277	1 1/2" (38 mm)		A4373
14278	1 5/8" (41 mm)		A4373

Assura® Convex Light Barrier 

- Extended wear
- Belt tabs for use with 4215, 4220
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
14281	5/8 - 15/16" (15-23 mm)		A4407
14282	5/8 - 15/16" (15-33 mm)		A4407
14283	5/8 - 111/16" (15-43 mm)		A4407
Pre-cut			
14291	3/4" (18 mm)		A4407
14292	7/8" (22 mm)		A4407
14293	1" (25 mm)		A4407
14294	1 1/8" (28 mm)		A4407
14295	1 1/4" (31 mm)		A4407
14296	1 3/8" (35 mm)		A4407
14297	1 1/2" (38 mm)		A4407
14298	1 5/8" (41 mm)		A4407

 Indicates convexity












*For important information regarding the use of HCPCS codes listed and code definitions see page 120

Assura® Deep Convex Barrier 

- Standard wear
- Belt tabs for use with 4215, 4220
- **Not made with natural rubber latex**



Units: 5












Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
12716	5/8 - 1 ⁵ / ₁₆ " (15-33 mm)		A4373
12719	5/8 - 1 ¹¹ / ₁₆ " (15-43 mm)		A4373
Pre-cut			
12701	5/8" (15 mm)		A4373
12702	3/4" (18 mm)		A4373
12703	7/8" (22 mm)		A4373
12704	1" (25 mm)		A4373
12705	1 ¹ / ₈ " (28 mm)		A4373
12706	1 ¹ / ₄ " (31 mm)		A4373
12707	1 ³ / ₈ " (35 mm)		A4373
12708	1 ¹ / ₂ " (38 mm)		A4373
12709	1 ⁵ / ₈ " (41 mm)		A4373

Assura® Deep Convex Barrier 

- Extended wear
- Belt tabs for use with 4215, 4220
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
14243	5/8 - 1 ⁵ / ₁₆ " (15-23 mm)		A4407
14246	5/8 - 1 ⁵ / ₁₆ " (15-33 mm)		A4407
14249	5/8 - 1 ¹¹ / ₁₆ " (15-43 mm)		A4407
Pre-cut			
14232	3/4" (18 mm)		A4407
14233	7/8" (22 mm)		A4407
14234	1" (25 mm)		A4407
14235	1 ¹ / ₈ " (28 mm)		A4407
14236	1 ¹ / ₄ " (31 mm)		A4407
14237	1 ³ / ₈ " (35 mm)		A4407
14238	1 ¹ / ₂ " (38 mm)		A4407
14239	1 ⁵ / ₈ " (41 mm)		A4407

 Indicates convexity




NOTE: Assura® is compatible with SenSura® Click, but NOT with SenSura® Mio Click.

Assura® MAXI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11½" (29½ cm) / 655 mL Units: 10



Code	Color	HCPCS*
Transparent	Opaque	
15974	15984 	A4427
15975	15985 	A4427
15976	15986 	A4427

Assura® MAXI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Without filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**



11½" (29½ cm) / 655 mL Units: 10




Code	Color	HCPCS*
Opaque		
15982		A4426
15983		A4426

Assura® MIDI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**



10½" (26½ cm) / 485 mL Units: 10

Code	Color	HCPCS*
Opaque		
15954		A4427
15955		A4427
15956		A4427

*For important information regarding the use of HCPCS codes listed and code definitions see page 120




NOTE: Assura® is compatible with SenSura® Click, but NOT with SenSura® Mio Click.

Assura® MAXI Drainable Pouch

- EasiClose™ Outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11½" (29½ cm) / 655 mL Units: 10



Code	Code	Color	HCPCS*
Transparent			
13974	Opaque		A4427
13975	13985		A4427
13976	13986		A4427

Assura® MIDI Drainable Pouch

- EasiClose™ Outlet
- **Filter**
- Soft cloth on front and back
- Graphic pattern
- **Not made with natural rubber latex**



10½" (26½ cm) / 485 mL Units: 10



Code	Color	HCPCS*
Opaque		
13964		A4427
13965		A4427
13966		A4427

Assura® MINI Drainable Pouch

- EasiClose™ Outlet
- **Filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**



9¾" (23½ cm) / 280 mL Units: 10

Code	Color	HCPCS*
Opaque		
13924		A4427
13925		A4427

NOTE: Assura® is compatible with SenSura® Click, but NOT with SenSura® Mio Click.




Assura® MAXI Drainable Pouch

- Clamp closure
- **Without filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**



11¼" (28½ cm) / 760 mL

Units: 10

Code	Color	HCPCS*
Opaque		
12577		A4426
12578		A4426
12579		A4426




Assura® MIDI Drainable Pouch

- Clamp closure
- **Without filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**



9¾" (24½ cm) / 400 mL

Units: 10

Code	Color	HCPCS*
Opaque		
12527		A4426
12528		A4426
12529		A4426




Assura® High Output Drainable Pouch

- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



12½" (32 cm) / 750 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent			
Opaque			
8114	NA		A4413
8115	2846		A4413
8116	2847		A4413



For bed drainage bag and bag hanger, see page 16.

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

Assura® Ileo Night Bag

- 2 liter capacity for very high output stomas
- Clamp closure
- Comes with bag hanger
- **Filter**
- **Not made with natural rubber latex**






2 L		Units: 5	
Code	Color	HCPCS*	
Opaque			
2836		A4413	
2837		A4413	

Assura® MAXI Closed Pouch

- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**






8½" (21 cm) / 570 mL			Units: 30
Code	Code	Color	HCPCS*
Transparent		Opaque	
12374	12384		A4423
12375	12385		A4423
12376	12386		A4423

Assura® MIDI Closed Pouch

- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



7" (18 cm) / 275 mL			Units: 30
Code	Code	Color	HCPCS*
Transparent		Opaque	
12344	12354		A4423
12345	12355		A4423
12346	12356		A4423

NOTE: Assura® is compatible with SenSura® Click, but NOT with SenSura® Mio Click.




Ileostomy / Colostomy

Assura® Stoma Cap

- **Filter**
- Small size is optimal for intimate moments or colostomates who irrigate
- **Not made with natural rubber latex**



Units: 30

Code	Color	HCPCS*
Opaque		
2801		A5055
2802		A5055
2803		A5055**

Urostomy




Assura® Urostomy Micro-Pouch

- Miniature size prevents twisting
- Soft outlet can be connected to a drain bag for nighttime use, or a leg bag for patients in a wheelchair
- Soft cloth on front
- **Not made with natural rubber latex**



6" (15 cm)

Units: 10

Code	Color	HCPCS*
Transparent		
14204		A4433
14205		A4433
14206		A4433

For nighttime, use a urostomy night bag found on page 15.



Assura® Urostomy Minicap

- Ultra-absorbent material inside
- Small size is optimal for intimate use
- **Not made with natural rubber latex**



100 mL

Units: 30

Code	Color	HCPCS*
Opaque		
2807		A5055
2808		A5055

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

**Product has not been through PDAC verification

NOTE: Assura® is compatible with SenSura® Click, but NOT with SenSura® Mio Click.




Assura® MAXI Urostomy Pouch

- Multi-chambers minimize sloshing
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



10½" (26½ cm) / 488 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent			
Opaque			
14227	14224		A4433
14228	14225		A4433
14229	14226		A4433




Assura® MAXI Urostomy Pouch

- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



10¾" (27¼ cm) / 375 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent			
Opaque			
1757	1754		A4433
1758	1755		A4433
1759	1756		A4433



Assura® MIDI Urostomy Pouch

- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back
- **Not made with natural rubber latex**



9½" (24 cm) / 150 mL

Units: 10

Code	Color	HCPCS*
Transparent		
1751		A4433
1752		A4433

Urostomy Night Bag and Bag Hanger

See Accessories page 16 for more information.





2-piece Flex

(Adhesive coupling)

2-Piece Flex



Flexible system bends with body movements



The low profile of a 1-piece – the convenience of a 2-piece



Floating flange = minimal pressure to apply



Built-in security channel + strong adhesive = strong leakage reduction

SenSura® Mio

2-Piece Flex

The fit is
the difference
between staying in
and working out

Deirdre is a SenSura® Mio user.

SenSura® Mio Flex

The BodyFit Technology in all SenSura Mio barriers combines an elastic adhesive, that bends and stretches with the body, with an innovative barrier designed to provide a secure and comfortable fit.



NOTE: SenSura® Mio Flex and SenSura® Flex are compatible with one another

SenSura® Mio Flex Flat Barrier

- Extended wear
- **Not made with natural rubber latex**



Units: 5
18390 - 10 units

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
10551	3/8 - 1 5/16" (10-33 mm)		A4409
10561	3/8 - 1 7/8" (10-48 mm)		A4410
10571	3/8 - 2 11/16" (10-68 mm)		A4410
18390	3/8 - 3 1/2" (10-88mm)		A4410

Pre-cut

10562	1" (25 mm)		A4409
10566	1 1/8" (28 mm)		A4409
10563	1 3/16" (30 mm)		A4409
10567	1 1/4" (32 mm)		A4409
10564	1 3/8" (35 mm)		A4409
10568	1 1/2" (38 mm)		A4409
10565	1 9/16" (40 mm)		A4409
10572	1 3/4" (45 mm)		A4410
10573	2" (50 mm)		A4410

SenSura® Mio Flex Convex Light Barrier

- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
16471	3/8 - 3/4" (10-20 mm)		A4407
16481	5/8 - 1 3/16" (15-30 mm)		A4407
16483	5/8 - 1 9/16" (15-40 mm)		A4407
16491	5/8 - 2 1/16" (15-53 mm)		A4408

Pre-cut

16472	1 1/16" (18 mm)		A4407
16484	1 3/16" (21 mm)		A4407
16485	1" (25 mm)		A4407
16486	1 1/8" (28 mm)		A4407
16487	1 1/4" (31 mm)		A4407
16488	1 3/8" (35 mm)		A4407
16493	1 1/2" (38 mm)		A4408
16494	1 5/8" (41 mm)		A4408

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

NOTE: SenSura® Mio pouches **are compatible** with SenSura Xpro extended wear barriers





SenSura® Mio Flex MAXI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo inspection window for easy application
- **Not made with natural rubber latex**



11" (28 cm) / 655 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent			
Opaque			
NA	12262		A4425
12273	12272		A4425
12283	12282		A4425
18389	NA		A4425



SenSura® Mio Flex MAXI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Without filter**
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**



11" (28 cm) / 655 mL

Units: 10

Code	Color	HCPCS*
Transparent		
12274		A5063
12284		A5063


SenSura® Mio Flex MIDI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window for easy application
- **Not made with natural rubber latex**



10¼" (26 cm) / 485 mL

Units: 10

Code	Color	HCPCS*
Opaque		
12281		A4425



SenSura® Mio Flex High Output Pouch with Tap Outlet

- Special full-circle filter designed specifically for ileostomy output
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- **Not made with natural rubber latex**



12½" (32 cm) / 1080 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
18655	18650		A4413
18656	18651		A4413



SenSura® Mio Flex High Output Pouch with Soft Outlet

- Special full-circle filter designed specifically for ileostomy output
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- **Not made with natural rubber latex**



12½" (32 cm) / 1080 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
18657	18652		A4413
18658	18653		A4413

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

NOTE: SenSura® Mio Flex and SenSura® Flex are compatible with one another




SenSura® Mio Flex MAXI Closed Pouch

- Full-circle filter
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo inspection window for easy application
- Not made with natural rubber latex



8¼" (20½ cm) / 570 mL

Units: 30

Code	Code	Color	HCPCS*
Transparent			
Opaque			
12203	12202		A4419
12216	12213		A4419
12223	12222		A4419



SenSura® Mio Flex MAXI Urostomy Pouch

- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- Not made with natural rubber latex



10½" (27 cm) / 400 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent			
Opaque			
12292	12291		A5073**
12298	12297		A5073**

SenSura® Mio Flex MIDI Urostomy Pouch

- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- Not made with natural rubber latex



9½" (24 cm) / 330 mL

Units: 10

Code	Color	HCPCS*
Opaque		
12296		A5073**





NOTE: SenSura® Mio Flex and SenSura® Flex are compatible with one another

SenSura® Flex Flat Barrier

- Standard wear
- **Not made with natural rubber latex**



Units: 10
10108 - 5 units




Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
10101	3/8 - 1 5/16" (10-33 mm)		A4414
10102	3/8 - 1 7/8" (10-48 mm)		A4414
10103	3/8 - 2 3/4" (10-68 mm)		A4415
10108	3/8 - 3 1/2" (10-88 mm)		A4415

SenSura® Flex Xpro Flat Barrier

- Extended wear
- Compatible with SenSura Mio pouches
- **Not made with natural rubber latex**



Units: 10





Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
10104	3/8 - 1 5/16" (10-33 mm)		A4409
10105	3/8 - 1 7/8" (10-48 mm)		A4409
10106	3/8 - 2 3/4" (10-68 mm)		A4410

SenSura® Flex Convex Light Barrier

- Standard wear
- Belt tabs for use with 4215, 4220
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
11301	5/8 - 1 5/16" (15-23 mm)		A4373
11302	5/8 - 1 5/16" (15-33 mm)		A4373
11303	5/8 - 1 11/16" (15-43 mm)		A4373
11307	5/8 - 2 1/16" (15-53 mm)		A4373

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 120





**Product has not been through PDAC verification

SenSura® Flex Convex Light Barrier

- Standard wear
- Belt tabs for use with 4215, 4220
- **Not made with natural rubber latex**



Units: 5





Code	Stoma Size	Color	HCPCS*
Pre-cut			
11312	7/8" (22 mm)		A4373
11313	1" (25 mm)		A4373
11314	1 1/8" (28 mm)		A4373
11315	1 1/4" (31 mm)		A4373

SenSura® Flex Xpro Convex Light Barrier

- Extended wear
- Compatible with SenSura Mio pouches
- Belt tabs for use with 4215, 4220
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
11304	5/8 - 15/16" (15-23 mm)		A4407
11305	5/8 - 1 1/16" (15-33 mm)		A4407
11306	5/8 - 1 11/16" (15-43 mm)		A4407
11308	5/8 - 2 1/16" (15-53 mm)		A4408





SenSura® Flex MAXI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11 1/2" (29 1/2 cm) / 655 mL

Units: 20

Code	Code	Color	HCPCS*
Transparent			
NA	11514		A4425
11512	11515		A4425
11517	11518		A4425
11 1/2" (29 1/2 cm) / 750 mL XXL pouch			
11590	11591		A4425

 Indicates convexity

NOTE: SenSura® Mio Flex and SenSura® Flex are compatible with one another

SenSura® Flex MAXI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Without filter**
- Soft cloth on back
- **Not made with natural rubber latex**




11½" (29½ cm) / 655 mL

Units: 20

Code	Color	HCPCS*
------	-------	--------

Transparent

11612		A5063
-------	---	-------

11½" (29½cm) / 750 mL XXL pouch

11690		A5063
-------	---	-------

SenSura® Flex MAXI Closed Pouch

- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



8½" (21 cm) / 570 mL

Units: 30

Code	Code	Color	HCPCS*
------	------	-------	--------

Transparent Opaque

10911	10914		A4419
-------	-------	--	-------

10912	10915		A4419
-------	-------	---	-------

10917	10918		A4419
-------	-------	---	-------

*For important information regarding the use of HCPCS codes listed and code definitions see page 120



1-piece



The choice for flexibility:
bends with body movements



No coupling = pouch and
barrier are welded together
and apply as one piece



Low-profile

SenSura[®] Mio

The fit is
the difference
between asking for help
and doing it yourself

Walton is a SenSura[®] Mio Convex Flip user and has received compensation from Coloplast.

SenSura[®] Mio 1-piece Pouching System

The BodyFit Technology in all SenSura Mio barriers combines an elastic adhesive, that bends and stretches with the body, with an innovative barrier designed to provide a secure and comfortable fit.



1-piece

SenSura® Mio Flat MAXI Drainable Pouch

- Extended wear
- EasiClose™ WIDE Outlet
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo inspection window for easy application
- **Not made with natural rubber latex**



11" (28 cm) / 655 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

10481	10471	Cut-to-fit	$\frac{3}{8}$ - $2\frac{1}{8}$ " (10-55 mm)	A5056
10482	10472	Pre-cut	1" (25 mm)	A5056
10488	10473	Pre-cut	$1\frac{1}{8}$ " (28 mm)	A5056
10483	10474	Pre-cut	$1\frac{3}{16}$ " (30 mm)	A5056
10490	10475	Pre-cut	$1\frac{1}{4}$ " (32 mm)	A5056
10484	10476	Pre-cut	$1\frac{3}{8}$ " (35 mm)	A5056
10491	10480	Pre-cut	$1\frac{1}{2}$ " (38 mm)	A5056
10487	10479	Pre-cut	$1\frac{9}{16}$ " (40 mm)	A5056
10485	10477	Pre-cut	$1\frac{3}{4}$ " (45 mm)	A5056
10486	10478	Pre-cut	2" (50 mm)	A5056

SenSura® Mio Flat MIDI Drainable Pouch

- Extended wear
- EasiClose™ WIDE Outlet
- **Full-circle filter**
- Peek-a-boo inspection window allows easy application
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**

**Opaque** 10 $\frac{1}{4}$ " (26 cm) / 485 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
10461	Cut-to-fit	$\frac{3}{8}$ - $1\frac{3}{4}$ " (10-45 mm)	A5056

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

SenSura® Mio Flat MAXI Drainable Pouch

- Extended wear
- EasiClose™ WIDE Outlet
- **Without filter**
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**



Transparent 11" (28 cm) / 655 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
10489	Cut-to-fit	3/8" - 2 1/8" (10-55 mm)	A4388
10492	Pre-cut	1" (25 mm)	A4388
10493	Pre-cut	1 1/8" (28 mm)	A4388
10494	Pre-cut	1 1/4" (32 mm)	A4388
10495	Pre-cut	1 3/8" (35 mm)	A4388
10496	Pre-cut	1 1/2" (38 mm)	A4388

SenSura® Mio Soft Convex MAXI Drainable Pouch

- Extended wear
- Belt tabs for use with 4237, 4247
- EasiClose™ WIDE Outlet
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**



11" (28 cm) / 655 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent				
16715	16705	Cut-to-fit	5/8" - 1 5/16" (15-33 mm)	A5057
16716	16706	Cut-to-fit	3/8" - 2" (10-50 mm)	A5057
Opaque				
NA	13671	Pre-cut	7/8" (21 mm)	A5057
NA	13672	Pre-cut	1" (25 mm)	A5057
NA	13673	Pre-cut	1 1/8" (28 mm)	A5057
NA	13674	Pre-cut	1 1/4" (31 mm)	A5057
NA	13675	Pre-cut	1 3/8" (35 mm)	A5057
NA	13676	Pre-cut	1 1/2" (38 mm)	A5057

 Indicates convexity

SenSura® Mio Convex Light MAXI Drainable Pouch

- Extended wear
- Belt tabs for use with 4237, 4247
- EasiClose™ WIDE Outlet
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**



11" (28 cm) / 655 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
16737	16725	Cut-to-fit	$\frac{3}{8}$ - $\frac{7}{8}$ " (10-22 mm)	A5057
16738	16726	Cut-to-fit	$\frac{5}{8}$ - $1\frac{5}{16}$ " (15-33 mm)	A5057
16739	16733	Cut-to-fit	$\frac{3}{8}$ - $1\frac{11}{16}$ " (10-43 mm)	A5057
16740	16730	Pre-cut	$1\frac{3}{16}$ " (21 mm)	A5057
16741	16731	Pre-cut	1" (25 mm)	A5057
16742	16732	Pre-cut	$1\frac{1}{8}$ " (28 mm)	A5057
NA	16727	Pre-cut	$1\frac{1}{4}$ " (31 mm)	A5057
NA	16734	Pre-cut	$1\frac{3}{8}$ " (35 mm)	A5057
NA	16735	Pre-cut	$1\frac{1}{2}$ " (38 mm)	A5057

SenSura® Mio Convex Light MAXI Drainable Pouch

- Extended wear
- Belt tabs for use with 4237, 4247
- EasiClose™ WIDE Outlet
- FlexShell is designed for stability and flexibility
- **Without filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**



Transparent 11" (28 cm) / 655 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
16748	Cut-to-fit	$\frac{3}{8}$ - $1\frac{11}{16}$ " (10-43 mm)	A4390

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

SenSura® Mio Deep Convex MAXI Drainable Pouch

- Extended wear
- Belt tabs for use with 4237, 4247
- EasiClose™ WIDE Outlet
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**



11" (28cm) / 655 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
16765	16755	Cut-to-fit	3/8 - 7/8" (10-22 mm)	A5057
16766	16756	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A5057
16767	16768	Cut-to-fit	3/8 - 1 11/16" (10-43 mm)	A5057
16770	16760	Pre-cut	1 3/16" (21 mm)	A5057
16771	16761	Pre-cut	1" (25 mm)	A5057
16772	16762	Pre-cut	1 1/8" (28 mm)	A5057
16773	16757	Pre-cut	1 1/4" (31 mm)	A5057
NA	16763	Pre-cut	1 3/8" (35 mm)	A5057
NA	16764	Pre-cut	1 1/2" (38 mm)	A5057

SenSura® Mio Deep Convex MAXI Drainable Pouch

- Extended wear
- Belt tabs for use with 4237, 4247
- EasiClose™ WIDE Outlet
- FlexShell is designed for stability and flexibility
- **Without filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**



Transparent 11" (28 cm) / 655 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
16778	Cut-to-fit	3/8 - 1 11/16" (10-43 mm)	A4390

Indicates convexity

SenSura® Mio MAXI Drainable Pouch with Soft Outlet

- Extended wear
- **Special full-circle filter designed specifically for ileostomy output**
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or most standard night bags
- **Not made with natural rubber latex**



11" (28 cm) / 630 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

18665	18660	Cut-to-fit	$\frac{3}{8}$ - 2 $\frac{3}{4}$ " (10-70 mm)	A5056
18666	18661	Cut-to-fit	$\frac{3}{8}$ - 4" (10-100 mm)	A5056

SenSura® Mio Soft Convex MAXI Drainable Pouch with Soft Outlet

- Extended wear
- Belt tabs for use with 4237, 4247
- **Special full-circle filter designed specifically for ileostomy output**
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or any standard night bag
- **Not made with natural rubber latex**



11" (28 cm) / 630 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

18671	18670	Cut-to-fit	$\frac{3}{8}$ - 2" (10-50 mm)	A5057
-------	-------	------------	-------------------------------	-------

SenSura® Mio Convex Light MAXI Drainable Pouch with Soft Outlet

- Extended wear
- Belt tabs for use with 4237, 4247
- **Special full-circle filter designed specifically for ileostomy output**
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or any standard night bag
- **Not made with natural rubber latex**



11" (28 cm) / 630 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

18676	18675	Cut-to-fit	$\frac{3}{8}$ - 1 $\frac{11}{16}$ " (10-43 mm)	A5057
-------	-------	------------	--	-------

Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

SenSura® Mio Convex Flip MAXI Drainable Pouch

- Extended wear
- Designed for curves, bulges and hernias
- Curved, star-shaped barrier conforms to outward areas with less creasing and folding
- “Flip to fit” design allows for easy application
- Opaque pouches feature a peek-a-boo inspection window for easy application
- **Full-circle filter**
- **Not made with natural rubber latex**



11" (28 cm) / 655 mL

Units: 10

Code	Code	Cut	Stoma Size	Max Cut	HCPCS*
------	------	-----	------------	---------	--------

Transparent Opaque

18305	NA	Cut-to-fit	$\frac{3}{8}$ - $1\frac{9}{16}$ " (10-40 mm)	$1\frac{9}{16}$ " (40 mm)	A5057
18321	18312	Cut-to-fit	$\frac{3}{8}$ - 2" (10-50 mm)	2" (50 mm)	A5057
18341	18330	Cut-to-fit	$\frac{3}{8}$ - $2\frac{3}{8}$ " (10-60 mm)	$2\frac{3}{8}$ " (60 mm)	A5057
NA	18301	Pre-cut	1" (25 mm)	$1\frac{9}{16}$ " (40 mm)	A5057
NA	18302	Pre-cut	$1\frac{3}{16}$ " (30 mm)	$1\frac{9}{16}$ " (40 mm)	A5057
18322	18313	Pre-cut	1" (25 mm)	2" (50 mm)	A5057
18323	18314	Pre-cut	$1\frac{3}{16}$ " (30 mm)	2" (50 mm)	A5057
18324	18315	Pre-cut	$1\frac{3}{8}$ " (35 mm)	2" (50 mm)	A5057
18342	18331	Pre-cut	$1\frac{9}{16}$ " (40 mm)	$2\frac{3}{8}$ " (60 mm)	A5057
18343	18332	Pre-cut	$1\frac{3}{4}$ " (45 mm)	$2\frac{3}{8}$ " (60 mm)	A5057
18344	18333	Pre-cut	2" (50 mm)	$2\frac{3}{8}$ " (60 mm)	A5057

SenSura® Mio Convex Flip MAXI Drainable Pouch

- Extended wear
- **Without filter**
- Designed for curves, bulges and hernias
- Curved, star-shaped barrier conforms to outward areas with less creasing and folding
- “Flip to fit” design allows for easy application
- **Not made with natural rubber latex**



Transparent 11" (28 cm) / 655 mL

Units: 10

Code	Cut	Stoma Size	Max Cut	HCPCS*
------	-----	------------	---------	--------

18320	Cut-to-fit	$\frac{3}{8}$ - 2" (10-50 mm)	2" (50 mm)	A4390
18340	Cut-to-fit	$\frac{3}{8}$ - $2\frac{3}{8}$ " (10-60 mm)	$2\frac{3}{8}$ " (60 mm)	A4390

 Indicates Convex Flip

SenSura® Mio Flat MAXI Closed Pouch

- Full-circle filter
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window allows easy application
- Not made with natural rubber latex

**Opaque** 8 ¼" (21 cm) / 550 mL

Units: 30

Code	Cut	Stoma Size	HCPCS*
10881	Cut-to-fit	5/8 - 1 3/4" (15-45 mm)	A4416
10885***	Cut-to-fit	1 9/16 - 2 1/8" (40-55 mm)	A4416
10883	Pre-cut	1 3/16" (30 mm)	A4416
10884	Pre-cut	1 3/8" (35 mm)	A4416
10885***	Pre-cut	1 9/16" (40 mm)	A4416

***Because of its larger starter hole, 10885 functions as both a pre-cut and cut-to-fit solution.

SenSura® Mio Flat MIDI Closed Pouch

- Full-circle filter
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window allows easy application
- Not made with natural rubber latex

**Opaque** 7 ¼" (19 cm) / 340 mL

Units: 30

Code	Cut	Stoma Size	HCPCS*
10871	Cut-to-fit	5/8 - 1 3/4" (15-45 mm)	A4416
10872	Pre-cut	1" (25 mm)	A4416
10873	Pre-cut	1 3/16" (30 mm)	A4416
10874	Pre-cut	1 3/8" (35 mm)	A4416

SenSura® Mio MINI Closed Pouch

- Full-circle filter
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window allows easy application
- Comfortable and discreet stoma cap alternative
- Not made with natural rubber latex

**Opaque** 5 ¼" (13 cm) / 110 mL

Units: 30

Code	Cut	Stoma Size	HCPCS*
10870	Cut-to-fit	5/8 - 1 3/4" (15-45 mm)	A4416

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

SenSura® Mio Soft Convex MAXI Closed Pouch 

- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**



8¼" (21 cm) / 550 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
NA	16305	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4417
16315	16306	Cut-to-fit	3/8 - 2" (10-50 mm)	A4417
NA	16310	Pre-cut	1" (25 mm)	A4417
NA	16311	Pre-cut	1 1/8" (28 mm)	A4417
NA	16312	Pre-cut	1 1/4" (31 mm)	A4417
NA	16313	Pre-cut	1 3/8" (35 mm)	A4417

SenSura® Mio Soft Convex MIDI Closed Pouch 

- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**

**Opaque** 7¼" (19 cm) / 340 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
16301	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4417

SenSura® Mio Convex Light MAXI Closed Pouch 

- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**



8¼" (21 cm) / 550 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
NA	16325	Cut-to-fit	3/8 - 7/8" (10-22 mm)	A4417
16335	NA	Cut-to-fit	3/8 - 1 11/16" (10-43 mm)	A4417
NA	16326	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4417
NA	16330	Pre-cut	1" (25 mm)	A4417
NA	16331	Pre-cut	1 1/8" (28 mm)	A4417
NA	16327	Pre-cut	1 1/4" (31 mm)	A4417
NA	16333	Pre-cut	1 3/8" (35 mm)	A4417
NA	16334	Pre-cut	1 1/2" (38 mm)	A4417

 Indicates convexity

SenSura® Mio Convex Light MIDI Closed Pouch

- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**

**Opaque** 7¼" (19 cm) / 340 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
16321	Cut-to-fit	¾ - 7/8" (10-22 mm)	A4417
16322	Cut-to-fit	5/8 - 1 ⁵ / ₁₆ " (15-33 mm)	A4417
16323**	Cut-to-fit / Pre-cut	1¼ - 1 ¹¹ / ₁₆ " (31-43 mm)	A4417

**Because of its larger starter hole, 16323 functions as both a pre-cut and cut-to-fit solution.

SenSura® Mio Deep Convex MAXI Closed Pouch

- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**



8¼" (21 cm) / 550 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
NA	16345	Cut-to-fit	¾ - 7/8" (10-22 mm)	A4417
16355	NA	Cut-to-fit	¾ - 1 ¹¹ / ₁₆ " (10-43 mm)	A4417
NA	16346	Cut-to-fit	5/8 - 1 ⁵ / ₁₆ " (15-33 mm)	A4417
NA	16350	Pre-cut	1" (25 mm)	A4417
NA	16351	Pre-cut	1 ¹ / ₈ " (28 mm)	A4417
NA	16347	Pre-cut	1¼" (31 mm)	A4417
NA	16353	Pre-cut	1 ³ / ₈ " (35 mm)	A4417
NA	16354	Pre-cut	1½" (38 mm)	A4417

SenSura® Mio Deep Convex MIDI Closed Pouch

- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**

**Opaque** 7¼" (19 cm) / 340 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
16341	Cut-to-fit	¾ - 7/8" (10-22 mm)	A4417
16342	Cut-to-fit	5/8 - 1 ⁵ / ₁₆ " (15-33 mm)	A4417
16343**	Cut-to-fit / Pre-cut	1¼ - 1 ¹¹ / ₁₆ " (31-43 mm)	A4417

**Because of its larger starter hole, 16343 functions as both a pre-cut and cut-to-fit solution.

Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

SenSura® Mio Convex Flip MAXI Closed Pouch 

- Designed for curves, bulges and hernias
- Curved, star-shaped barrier conforms to outward areas with less creasing and folding
- "Flip to fit" design allows for easy application
- Peek-a-boo inspection window for easy application
- **Not made with natural rubber latex**



Opaque 8¼" (21 cm) / 550 mL Units: 10

Code	Cut	Stoma Size	Max Cut	HCPCS*
18152	Cut-to-fit	¾ - 1⅞" (10-40 mm)	1⅞" (40 mm)	A4417
18170	Cut-to-fit	¾ - 2" (10-50 mm)	2" (50 mm)	A4417
18171	Pre-cut	1" (25 mm)	2" (50 mm)	A4417
18172	Pre-cut	1⅜" (30 mm)	2" (50 mm)	A4417
18173	Pre-cut	1⅜" (35 mm)	2" (50 mm)	A4417
18180	Cut-to-fit	¾ - 2⅜" (10-60 mm)	2⅜" (60 mm)	A4417
18181	Pre-cut	1⅞" (40 mm)	2⅜" (60 mm)	A4417
18182	Pre-cut	1¾" (45 mm)	2⅜" (60 mm)	A4417

Urostomy

SenSura® Mio Flat MAXI Urostomy Pouch

- Extended wear
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or most standard night bags
- Neutral gray, water-resistant textile on back
- **Not made with natural rubber latex**



Transparent 10½" (27 cm) / 400 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
10586	Cut-to-fit	¾ - 1¾" (10-45 mm)	A4391

 Indicates Convex Flip

SenSura® Mio Soft Convex MAXI Urostomy Pouch 

- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**



10½" (27 cm) / 400 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
Transparent Opaque			
13679	16805	Cut-to-fit 3/8 - 1 5/16" (10-33 mm)	A4393
16810	NA	Cut-to-fit 3/8 - 2" (10-50 mm)	A4393
13731	NA	Pre-cut 7/8" (21 mm)	A4393
13732	NA	Pre-cut 1" (25 mm)	A4393
13733	NA	Pre-cut 1 1/8" (28 mm)	A4393
13734	NA	Pre-cut 1 1/4" (31 mm)	A4393

SenSura® Mio Soft Convex MIDI Urostomy Pouch 

- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**



Opaque 9½" (24 cm) / 330 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
16801	Cut-to-fit	3/8 - 1 5/16" (10-33 mm)	A4393

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

SenSura® Mio Convex Light MAXI Urostomy Pouch 

- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**



10½" (27 cm) / 400 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
16835	16825	Cut-to-fit	¾ - 7/8" (10-22 mm)	A4393
16836	16826	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4393
16837	NA	Cut-to-fit	¾ - 1 11/16" (10-43 mm)	A4393
NA	16830	Pre-cut	1 1/16" (18 mm)	A4393
16840	16831	Pre-cut	1 3/16" (21 mm)	A4393
16841	16832	Pre-cut	1" (25 mm)	A4393
16842	16833	Pre-cut	1 1/8" (28 mm)	A4393
NA	16827	Pre-cut	1 1/4" (31 mm)	A4393

SenSura® Mio Convex Light MIDI Urostomy Pouch 

- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**



Opaque 9½" (24 cm) / 330 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
16821	Cut-to-fit	¾ - 7/8" (10-22 mm)	A4393
16822	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4393
16823	Pre-cut	1 1/4" (31 mm)	A4393

 Indicates convexity

SenSura® Mio Deep Convex MAXI Urostomy Pouch 

- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**



10½" (27 cm) / 400 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

16865	16855	Cut-to-fit	3/8 - 7/8" (10-22 mm)	A4393
16866	16856	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4393
16867	NA	Cut-to-fit	3/8 - 1 11/16" (10-43 mm)	A4393
NA	16860	Pre-cut	1 1/16" (18 mm)	A4393
16870	16861	Pre-cut	1 3/16" (21 mm)	A4393
16871	16862	Pre-cut	1" (25 mm)	A4393
16872	16863	Pre-cut	1 1/8" (28 mm)	A4393
NA	16857	Pre-cut	1 1/4" (31 mm)	A4393

SenSura® Mio Convex Flip MAXI Urostomy Pouch 

- Extended wear
- Designed for curves, bulges and hernias
- Curved, star-shaped barrier conforms to outward areas with less creasing and folding
- "Flip to fit" design allows for easy application
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- **Not made with natural rubber latex**



10½" (27 cm) / 400 mL Units: 10

Code	Code	Cut	Stoma Size	Max Cut	HCPCS*
------	------	-----	------------	---------	--------

Transparent Opaque

18405	18401	Cut-to-fit	3/8 - 1 9/16" (10-40 mm)	1 9/16" (40 mm)	A4393
NA	18402	Pre-cut	3/4" (20 mm)	1 9/16" (40 mm)	A4393
NA	18403	Pre-cut	1" (25 mm)	1 9/16" (40 mm)	A4393
NA	18404	Pre-cut	1 3/16" (30 mm)	1 9/16" (40 mm)	A4393

 Indicates convexity

 Indicates Convex Flip

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

SenSura® Flat MAXI Drainable Pouch

- Standard wear
- EasiClose™ WIDE Outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11½" (29½ cm) / 655 mL

Units: 20

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
15521	15531	Cut-to-fit	¾ - 3" (10-76 mm)	A4424
15522	15532	Pre-cut	1" (25 mm)	A4424
15523	15533	Pre-cut	1⅜" (30 mm)	A4424
15524	15534	Pre-cut	1⅝" (35 mm)	A4424

SenSura® Flat MAXI Drainable Pouch

- Standard wear
- EasiClose™ WIDE Outlet
- **Without filter**
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 11½" (29½ cm) / 655 mL

Units: 20

Code	Cut	Stoma Size	HCPCS*
15821	Cut-to-fit	¾ - 3" (10-76 mm)	A5061

SenSura® Xpro Flat MAXI Drainable Pouch

- Extended wear
- EasiClose™ WIDE Outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11½" (29½ cm) / 655 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
15981	NA	Cut-to-fit	¾ - 3" (10-76 mm)	A5056
NA	15991	Pre-cut	1" (25 mm)	A5056
NA	15992	Pre-cut	1⅜" (30 mm)	A5056
NA	15993	Pre-cut	1¼" (32 mm)	A5056
NA	15994	Pre-cut	1½" (38 mm)	A5056

SenSura® Convex Light MAXI Drainable Pouch

- Standard wear
- EasiClose™ WIDE Outlet
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11½" (29½ cm) / 655 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent **Opaque**

15606	15608	Cut-to-fit	5/8 - 1 ¹¹ / ₁₆ " (15-43 mm)	A4389
15623	15633	Pre-cut	1" (25 mm)	A4389
15624	15634	Pre-cut	1 ¹ / ₈ " (28 mm)	A4389

SenSura® Xpro Convex Light MAXI Drainable Pouch

- Extended wear
- EasiClose™ WIDE Outlet
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11½" (29½ cm) / 655 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent **Opaque**

15695	NA	Cut-to-fit	5/8 - 1 ⁵ / ₁₆ " (15-33 mm)	A5057
15696	NA	Cut-to-fit	5/8 - 1 ¹¹ / ₁₆ " (15-43 mm)	A5057
NA	15697	Pre-cut	7/8" (22 mm)	A5057
NA	15698	Pre-cut	1" (25 mm)	A5057

SenSura® Flat MAXI Drainable Pouch

- Standard wear
- Soft outlet provides effective drainage; connects to a standard drain bag or Coloplast fistula drain bag #14010
- **Filter**
- Soft cloth on front and back
- Peek-a-boo inspection window allows easy application
- **Not made with natural rubber latex**



Opaque 12½" (32 cm) / 630 mL

Units: 30

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

19004	Cut-to-fit	3/8 - 3" (10-76 mm)	A4424
-------	------------	---------------------	-------

For bed drainage bag, see page 16.

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

**Product has not been through PDAC verification

SenSura® Flat MAXI Closed Pouch

- Standard wear
- **Filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**



Opaque 8½" (21 cm) / 550 mL Units: 30

Code	Cut	Stoma Size	HCPCS*
15480	Cut-to-fit	¾ - 3" (10-76 mm)	A4416
15484	Pre-cut	1¾" (35 mm)	A4416

Urostomy

SenSura® Flat MAXI Urostomy Pouch

- Standard wear
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 10¾" (26 cm) / 400 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
11804	Cut-to-fit	¾ - 3" (10-76 mm)	A5071

SenSura® Flat MIDI Urostomy Pouch

- Standard wear
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 9½" (24 cm) / 330 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
11803	Cut-to-fit	¾ - 2½" (10-66 mm)	A5071

Sensura® Xpro Flat MAXI Urostomy Pouch

- Extended wear
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



10³/₈" (26 cm) / 400 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

11809	NA	Cut-to-fit	$\frac{3}{8}$ - 3" (10-76 mm)	A4391**
NA	11805	Pre-cut	$\frac{7}{8}$ " (22 mm)	A4391**
NA	11806	Pre-cut	1" (25 mm)	A4391**
NA	11807	Pre-cut	1 ³ / ₁₆ " (30 mm)	A4391**
NA	11808	Pre-cut	1 ¹ / ₄ " (32 mm)	A4391**

Sensura® Convex Light MAXI Urostomy Pouch

- Standard wear
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 10³/₈" (26 cm) / 400 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

11814	Cut-to-fit	$\frac{5}{8}$ - 1 ⁵ / ₁₆ " (15-33 mm)	A4392
--------------	------------	---	-------

Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

**Product has not been through PDAC verification

SenSura® Xpro Convex Light MAXI Urostomy Pouch

- Extended wear
- Belt tabs for use with 4215, 4220
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



10³/₈" (26 cm) / 400 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

11819	NA	Cut-to-fit	5/8 - 1 ⁵ / ₁₆ " (15-33 mm)	A4393**
NA	11816	Pre-cut	7/8" (22 mm)	A4393**
NA	11817	Pre-cut	1" (25 mm)	A4393**
NA	11818	Pre-cut	1 ¹ / ₄ " (32 mm)	A4393**

Urostomy Night Bag and Bag Hanger

See Accessories page 16 for more information.



EasiClose™ WIDE outlets offer a 1/2" wider integrated closure for easy, fast emptying on the go.

Assura® Flat MAXI Drainable Pouch

- Standard wear
- EasiClose™ WIDE Outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11¹/₂" (29 cm) / 530 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

15860	15870	Cut-to-fit	3/8 - 3" (10-75 mm)	A4424
NA	15874	Pre-cut	1" (25 mm)	A4424
NA	15875	Pre-cut	1 ³ / ₁₆ " (30 mm)	A4424
NA	15876	Pre-cut	1 ³ / ₈ " (35 mm)	A4424

 Indicates convexity

Assura® Flat MAXI Drainable Pouch

- Standard wear
- EasiClose™ WIDE Outlet
- **Without filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**



Opaque 11½" (29 cm) / 530 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
15871	Cut-to-fit	¾ - 3" (10-75 mm)	A5061

Assura® Flat MAXI Drainable Pouch

- Extended wear
- EasiClose™ WIDE Outlet
- **Without filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11¼" (28½ cm) / 530 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
15869	15879	Cut-to-fit	¾ - 3" (10-75 mm)	A4388

Assura® Convex Light MAXI Drainable Pouch

- Standard wear
- EasiClose™ WIDE Outlet
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 11½" (29 cm) / 530 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
14511	Cut-to-fit	⅝ - 1⅝" (15-33 mm)	A4389
14512	Cut-to-fit	⅝ - 1⅞" (15-43 mm)	A4389

Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

**Product has not been through PDAC verification

Assura® Flat MAXI Drainable Pouch

- Standard wear
- EasiClose™ Outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11¼" (28½ cm) / 530 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

13860	13870	Cut-to-fit	¾ - 3" (10-75 mm)	A4424
NA	13874	Pre-cut	1" (25 mm)	A4424
NA	13875	Pre-cut	1⅜" (30 mm)	A4424
NA	13876	Pre-cut	1⅜" (35 mm)	A4424

Assura® Flat MIDI Drainable Pouch

- Standard wear
- EasiClose™ Outlet
- **Filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**

**Opaque** 10¼" (26 cm) / 390 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

13840	Cut-to-fit	¾ - 2½" (10-65 mm)	A4424
13844	Pre-cut	1" (25 mm)	A4424
13845	Pre-cut	1⅜" (30 mm)	A4424
13846	Pre-cut	1⅜" (35 mm)	A4424

Assura® Convex Light MAXI Drainable Pouch

- Standard wear
- EasiClose™ Outlet
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11¼" (28½ cm) / 530 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

Transparent Opaque

14411	14414	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4389
14412	14415	Cut-to-fit	5/8 - 1 11/16" (15-43 mm)	A4389
NA	14421	Pre-cut	7/8" (22 mm)	A4389
NA	14422	Pre-cut	1" (25 mm)	A4389
NA	14423	Pre-cut	1 1/8" (28 mm)	A4389
NA	14424	Pre-cut	1 1/4" (31 mm)	A4389

Assura® Convex Light MIDI Drainable Pouch

- Standard wear
- EasiClose™ Outlet
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on front and back
- **Not made with natural rubber latex**



Opaque 10¼" (26 cm) / 390 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

14404	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4389
14405	Cut-to-fit	5/8 - 1 11/16" (15-43 mm)	A4389

Assura® Deep Convex MAXI Drainable Pouch

- Standard wear
- EasiClose™ Outlet
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11¼" (28½ cm) / 530 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

14103	14106	Cut-to-fit	5/8 - 1 11/16" (15-43 mm)	A4389
14163	14173	Pre-cut	7/8" (22 mm)	A4389
14164	14174	Pre-cut	1" (25 mm)	A4389
14165	14175	Pre-cut	1 1/8" (28 mm)	A4389
14166	14176	Pre-cut	1 1/4" (31 mm)	A4389

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

Assura® Deep Convex MIDI Drainable Pouch

- Standard wear
- EasiClose™ Outlet
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on front and back
- **Not made with natural rubber latex**



Opaque 10¼" (26 cm) / 390 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
14196	Cut-to-fit	5/8 - 1 ¹¹ / ₁₆ " (15-43 mm)	A4389

Assura® Convex Light MAXI Drainable Pouch

- Extended wear
- EasiClose™ Outlet
- **Without filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11¼" (28½ cm) / 530 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
14413	NA	Cut-to-fit	5/8 - 1 ⁵ / ₁₆ " (15-33 mm)	A4390
14416	14419	Cut-to-fit	5/8 - 1 ¹¹ / ₁₆ " (15-43 mm)	A4390

Assura® Flat MAXI Drainable Pouch

- Standard wear
- Clamp closure
- **Without filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11¼" (28½ cm) / 530 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
12660	12670	Cut-to-fit	3/8 - 2 ³ / ₄ " (10-70 mm)	A5061

Assura® Flat MINI Drainable Pouch

- Standard wear
- Clamp closure
- **Without filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



9¾" (24½ cm) / 230 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
Transparent Opaque			
12630	12620	Cut-to-fit	3/8 - 2¼" (10-55 mm)
			A5061

 Indicates convexity

Assura® Deep Convex MAXI Drainable Pouch 

- Standard wear
- Clamp closure
- **Without filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11¼" (28½ cm) / 530 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent **Opaque**

13706	12533	Cut-to-fit	5/8 - 1 ¹¹ / ₁₆ " (15-43 mm)	A4389
-------	-------	------------	--	-------

Assura® Deep Convex MINI Drainable Pouch 

- Standard wear
- Clamp closure
- **Without filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- **Not made with natural rubber latex**

**Transparent** 9¾" (24½ cm) / 230 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

13703	Cut-to-fit	5/8 - 1 ¹¹ / ₁₆ " (15-43 mm)	A4389
-------	------------	--	-------

Assura® Deep Convex MAXI Drainable Pouch 

- Extended wear
- Clamp closure
- **Without filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11¼" (28½ cm) / 530 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent **Opaque**

12416	12419	Cut-to-fit	5/8 - 1 ¹¹ / ₁₆ " (15-43 mm)	A4390
-------	-------	------------	--	-------

Assura® Flat MAXI Closed Pouch

- Standard wear
- **Filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**

**Opaque** 8½" (22 cm) / 530 mL

Units: 30

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

12170	Cut-to-fit	1 ³ / ₁₆ - 3" (20-75 mm)	A4416
-------	------------	--	-------

12174	Pre-cut	1" (25 mm)	A4416
-------	---------	------------	-------

12175	Pre-cut	1 ³ / ₁₆ " (30 mm)	A4416
-------	---------	--	-------

12176	Pre-cut	1 ³ / ₈ " (35 mm)	A4416
-------	---------	---	-------

12177	Pre-cut	1 ⁹ / ₁₆ " (40 mm)	A4416
-------	---------	--	-------

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

Assura® Flat MIDI Closed Pouch

- Standard wear
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



7" (18 cm) / 350 mL Units: 30

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
12130	12140	Cut-to-fit	1 ³ / ₁₆ - 2 ¹ / ₂ " (20-65 mm)	A4416
12134	12144	Pre-cut	1" (25 mm)	A4416
12135	12145	Pre-cut	1 ³ / ₁₆ " (30 mm)	A4416
12136	12146	Pre-cut	1 ³ / ₈ " (35 mm)	A4416
NA	12147	Pre-cut	1 ⁹ / ₁₆ " (40 mm)	A4416
NA	12149	Pre-cut	2" (50 mm)	A4416

Assura® Flat MINI Closed Pouch

- Standard wear
- **Filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**



Opaque 5" (12 cm) / 170 mL Units: 30

Code	Cut	Stoma Size	HCPCS*
12110	Cut-to-fit	1 ³ / ₁₆ - 2 ¹ / ₂ " (20-65 mm)	A4416

Assura® Convex Light MAXI Closed Pouch

- Standard wear
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



8 1/2" (22 cm) / 560 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
14441	14444	Cut-to-fit	5/8 - 1 ⁵ / ₁₆ " (15-33 mm)	A4417
14442	14445	Cut-to-fit	5/8 - 1 ¹¹ / ₁₆ " (15-43 mm)	A4417

Indicates convexity

Assura® Convex Light MIDI Closed Pouch

- Standard wear
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



7" (18 cm) / 350 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

Transparent Opaque

14431	14434	Cut-to-fit	5/8" - 1 5/16" (15-33 mm)	A4417
14432	14435	Cut-to-fit	5/8" - 1 11/16" (15-43 mm)	A4417
NA	14451	Pre-cut	7/8" (22 mm)	A4417
NA	14452	Pre-cut	1" (25 mm)	A4417
NA	14453	Pre-cut	1 1/8" (28 mm)	A4417
NA	14454	Pre-cut	1 1/4" (31 mm)	A4417

Assura® Flat Stoma Cap

- Standard Wear
- **Filter**
- Belt tabs
- Soft cloth on back
- Small pouch size for intimate moments or colostomates who irrigate
- **Not made with natural rubber latex**



Opaque

Units: 30

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

2501	Cut-to-fit	1 3/16" - 2 1/4" (20-55 mm)	A5055
------	------------	-----------------------------	-------

Urostomy

Assura® Flat MAXI Urostomy Pouch

- Standard wear
- Multi-chambers minimize sloshing
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 10 3/4" (27 1/2 cm) / 470 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

14222	Cut-to-fit	3/8" - 2 1/4" (10-55 mm)	A5071
-------	------------	--------------------------	-------

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

Assura® Convex Light MAXI Urostomy Pouch 

- Standard wear
- Multi-chambers minimize sloshing
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 10¾" (27½ cm) / 470 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
14717	Cut-to-fit	5/8 - 1 ¹¹ / ₁₆ " (15-43 mm)	A4392

Assura® Flat MAXI Urostomy Pouch

- Standard wear
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



10¾" (27 cm) / 465 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
5585	5580	Cut-to-fit	3/8 - 2 ¹ / ₄ " (10-55 mm)	A5071

Assura® Flat MIDI Urostomy Pouch

- Standard wear
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



7" (18 cm) / 265 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
5575	5570	Cut-to-fit	3/8 - 2 ¹ / ₄ " (10-55 mm)	A5071

Assura® Flat MAXI Urostomy Pouch

- Extended wear
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 9½" (24 cm) / 375 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
12474	Cut-to-fit	3/8 - 2 ¹ / ₄ " (10-55 mm)	A4391

 Indicates convexity

Assura® Flat MIDI Urostomy Pouch

- Extended wear
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



9 ½" (24 cm) / 265 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

12472	12471	Cut-to-fit	3/8 - 2 1/4" (10-55 mm)	A4391
-------	-------	------------	-------------------------	-------

Assura® Deep Convex MAXI Urostomy Pouch

- Standard wear
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- **Not made with natural rubber latex**

**Transparent** 10 ¾" (27 cm) / 465 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
12595	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4392
12596	Cut-to-fit	5/8 - 1 11/16" (15-43 mm)	A4392
12991	Pre-cut	5/8" (15 mm)	A4392
12992	Pre-cut	3/4" (18 mm)	A4392
12993	Pre-cut	7/8" (22 mm)	A4392
12994	Pre-cut	1" (25 mm)	A4392
12995	Pre-cut	1 1/8" (28 mm)	A4392
12996	Pre-cut	1 1/4" (31 mm)	A4392
12997	Pre-cut	1 3/8" (35 mm)	A4392

Assura® Deep Convex MAXI Urostomy Pouch

- Extended wear
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- **Not made with natural rubber latex**

**Transparent** 10 ¾" (27 cm) / 465 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
12432	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4393
12433	Cut-to-fit	5/8 - 1 11/16" (15-43 mm)	A4393

Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 120


Download the Coloplast Stoma Care App



Designed with nurses in mind

Features:

1. Easy access to Coloplast® Care
2. BodyCheck Assessment
3. Step-by-step guides with how-to-use videos
4. Common problems and recommended actions
5. General information on urostomies, ileostomies and colostomies



Download the Coloplast Stoma Care app now!

Search for “**Stoma Care**” in the App Store and Google Play to download for free today.





Soft
outlet



Tap
outlet

Hospital Assortment



Elastic adhesive fits over folds and contours, yet remains securely in place



Two outlet types for different options:

- Soft (for liquid stool; fits standard drain bags)
- Tap (for thicker output; fits fecal drain bag with wide-lumen tubing)



Printed measurement scale on post op pouches helps visually track output volume



Security and comfort from hospital to home.

Your patients deserve all the benefits of SenSura[®] Mio from day one.

Getting used to a stoma is challenging enough without the added complexity of having to learn a different pouching system. Why not give your patients the security and comfort of Mio - right from the start? Options are available in different sizes, outlet types, and coupling systems to meet patient needs.



2-piece Click
Soft | Tap
**Also available
in smaller size!*



2-piece Flex
Soft | Tap



1-piece
Soft
*Available in flat,
soft convex, and
light convex!*



Post Op with window
Soft
Sterile | Non-Sterile

High Output Pouches




SenSura® Mio Click High Output Pouch with Soft Outlet

- **Special full-circle filter designed specifically for ileostomy output**
- Neutral gray, water-resistant textile
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Soft outlet provides effective drainage; connects to a standard drain bag or Coloplast fistula drain bag #14010
- **Not made with natural rubber latex**



12½" (32 cm) / 1080 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
18640	18630		A4413
18641	18631		A4413
18642	18632		A4413





SenSura® Mio Click High Output Pouch with Tap Outlet

- **Special full-circle filter designed specifically for ileostomy output**
- Neutral gray, water-resistant textile
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- **Not made with natural rubber latex**



12½" (32 cm) / 1080 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
18610	18600		A4413
18611	18601		A4413
18612	18602		A4413
18613	18603		A4413

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

SenSura® Mio Click Drainable Pouch with Tap Outlet

- Special full-circle filter designed specifically for ileostomy output
- Neutral gray, water-resistant textile
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- **Not made with natural rubber latex**



9" (23 cm) / 650 mL Units: 10

Code	Color	HCPCS*
Opaque		
18620		A4427
18621		A4427
18622		A4427

Assura® High Output Drainable Pouch

- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



12½" (32 cm) / 750 mL Units: 10

Code	Code	Color	HCPCS*
Transparent			
8114	NA		A4413
8115	2846		A4413
8116	2847		A4413

For bed drainage bag and bag hanger, see page 16.

SenSura® Mio Flex High Output Pouch with Tap Outlet

- Special full-circle filter designed specifically for ileostomy output
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- **Not made with natural rubber latex**



12½" (32 cm) / 1080 mL Units: 10



Code	Code	Color	HCPCS*
Transparent			
Opaque			
18655	18650		A4413
18656	18651		A4413

SenSura® Mio Flex High Output Pouch with Soft Outlet

- Special full-circle filter designed specifically for ileostomy output
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to a standard drain bag or Coloplast fistula drain bag #14010
- **Not made with natural rubber latex**



12½" (32 cm) / 1080 mL Units: 10

Code	Code	Color	HCPCS*
Transparent		Opaque	
18657	18652		A4413
18658	18653		A4413

SenSura® Mio MAXI Drainable Pouch with Soft Outlet

- Extended wear
- Special full-circle filter designed specifically for ileostomy output
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or most standard night bags
- **Not made with natural rubber latex**



11" (28 cm) / 630 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent		Opaque		
18665	18660	Cut-to-fit	⅜ - 2 ¾" (10-70 mm)	A5056
18666	18661	Cut-to-fit	⅜ - 4" (10-100 mm)	A5056

SenSura® Mio Soft Convex MAXI Drainable Pouch with Soft Outlet

- Extended wear
- Belt tabs for use with 4237, 4247
- Special full-circle filter designed specifically for ileostomy output
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or any standard night bag
- **Not made with natural rubber latex**



11" (28 cm) / 630 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent		Opaque		
18671	18670	Cut-to-fit	⅜ - 2" (10-50 mm)	A5057

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

SenSura® Mio Convex Light MAXI Drainable Pouch with Soft Outlet

- Extended wear
- Belt tabs for use with 4237, 4247
- **Special full-circle filter designed specifically for ileostomy output**
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or any standard night bag
- **Not made with natural rubber latex**



11" (28 cm) / 630 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent **Opaque**

18676	18675	Cut-to-fit	$\frac{3}{8}$ - $1\frac{11}{16}$ " (10-43 mm)	A5057
-------	-------	------------	---	-------

SenSura® Flat MAXI Drainable Pouch

- Standard wear
- Soft outlet provides effective drainage; connects to a standard drain bag or Coloplast fistula drain bag #14010
- **Filter**
- Soft cloth on front and back
- Peek-a-boo inspection window allows easy application
- **Not made with natural rubber latex**



Opaque 12 $\frac{1}{2}$ " (32 cm) / 630 mL Units: 30

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

19004	Cut-to-fit	$\frac{3}{8}$ - 3" (10-76 mm)	A4424
-------	------------	-------------------------------	-------

For bed drainage bag, see page 16.

SenSura® Mio 1-piece Post Op Pouch with Window

- Neutral gray, water-resistant textile
- Printed scale provides visual indication of output volume
- Easy-to-open and close post-op window attached to pouch
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- **Not made with natural rubber latex**



Transparent Units: 5

Code	Code	Cut	Volume	Stoma Size
------	------	-----	--------	------------

Non-sterile **Sterile**

18690	18680	Cut-to-fit	680mL	$\frac{3}{8}$ - $2\frac{3}{4}$ " (10-70 mm)
-------	-------	------------	-------	---

18691	18681	Cut-to-fit	725mL	$\frac{3}{8}$ - 4" (10-100 mm)
-------	-------	------------	-------	--------------------------------

 Indicates convexity

SenSura® Standard Wear 1-piece Post-Op Pouch with Window

- Flat
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- Easy-to-open and -close access window attached to pouch
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 13" (33 cm) / 830 mL Units: 5

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Non-sterile **Sterile**

19020	19010	Cut-to-fit	$\frac{3}{8}$ - 3" (10-76 mm)	A6154
--------------	--------------	------------	-------------------------------	-------

19021	19011	Cut-to-fit	$\frac{3}{8}$ - 4½" (10-115 mm)	A6154
--------------	--------------	------------	---------------------------------	-------

SenSura® Standard Wear 1-piece Post-Op Pouch

- Flat
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- Effective pouching solution for small fistulas
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 12¼" (31 cm) / 690 mL Units: 5

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Non-sterile **Sterile**

19030	19001	Cut-to-fit	$\frac{3}{8}$ - 3" (10-76 mm)	A6154
--------------	--------------	------------	-------------------------------	-------

For SenSura® Magnum high output pouches, see page 30.

Assura® Standard Wear 1-piece Post-Op Pouch with Window

- Flat
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- Easy-to-open and -close access window
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 13" (33 cm) / 830 mL Units: 5

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Non-sterile **Sterile**

12800	12808	Cut-to-fit	$\frac{1}{2}$ - 2¾" (13-70 mm)	A6154
--------------	--------------	------------	--------------------------------	-------

12802	12810	Cut-to-fit	$\frac{1}{2}$ - 4" (13-100 mm)	A6154
--------------	--------------	------------	--------------------------------	-------

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

Assura® Standard Wear 1-piece Post-Op Pouch

- Flat
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- Soft cloth on back
- **Not made with natural rubber latex**

Transparent 12¾" (31 cm) / 690 mL Units: 5

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Non-sterile **Sterile**

12805	12818	Cut-to-fit	½ - 2¾" (13-70 mm)	A6154
-------	-------	------------	--------------------	-------

Assura® Ileo Night Bag


- 2 liter capacity for very high output stomas
- Clamp closure
- Comes with bag hanger
- **Filter**
- **Not made with natural rubber latex**



2 L Units: 5

Code	Color	HCPCS*
------	-------	--------

Opaque

2836		A4413
------	---	-------

2837		A4413
------	---	-------



Coloplast® Fistula and Wound Management System

Each system includes:

- 1 pouch and tracing guide
- 1 drain port (extras sold separately on next page)
- 1 lid (extras sold separately on next page)
- 1 clear drape
- 1 clamp
- **Not made with natural rubber latex**

Each box includes one pump.



14050 Mini



14060 Midi



14070 Maxi

Units: 3

Units for 14050: 6

Code	Cut	Stoma Size	HCPCS*
14050	Mini	4 $\frac{1}{8}$ x 6 $\frac{1}{4}$ " (104-159 mm)	A6154
14060	Midi	6 $\frac{1}{8}$ x 9" (156-228 mm)	A6154
14070	Maxi	8 $\frac{1}{8}$ x 11 $\frac{3}{4}$ " (208-297 mm)	A6154

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

Fistula / Wound Pouch Accessories

Coloplast® Fistula Drain Port

Can be easily applied to the lid to allow for catheter insertion. For catheters ranging from 18FR - 22FR, the drain port is ready for use. For catheters ranging from 24FR - 30FR, pull off the top of the drain port - DO NOT CUT.



- Not made with natural rubber latex

Code	Size	Units
14015	One size	12

Coloplast® Bed Drainage Bag

Bag tubing specifically designed for fecal output. Holds up to 2000 mL. Long tubing can be cut easily, allowing nurses to find the right length for bed drainage.



- Not made with natural rubber latex

Code	Size	Units: 6	HCPCS*
14010	2000 mL		A4357

For matching bag hanger, see page 16.

Coloplast® Transparent Fistula Lids

Replacement lids match the transparent lid included with each fistula pouch.



- Not made with natural rubber latex

Code	Size	Units
14021	Mini (matches 14050)	10
14031	Midi (matches 14060)	10
14041	Maxi (matches 14070)	10

Coloplast® Pouch Clamp

Simple closure clamp is easy to manipulate, even for those with low dexterity.



- Not made with natural rubber latex

Code	Size	Units: 20	HCPCS*
9500	One size		A4363

2-Piece Click Single Use Kits

SenSura® Mio Click Barrier + Drainable Pouch


- Extended wear
- Barrier #10522 Flat with belt tabs
- Pouch #11483 MAXI / EasiClose™ WIDE Outlet
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit / Transparent

21008	¾ - 2¼" (10-55 mm)		A4410
			A4427

SenSura® Click Standard Wear Barrier + Drainable Pouch


- Barrier #10031 Flat with belt tabs
- Pouch #11196 MAXI / EasiClose™ WIDE Outlet
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit / Transparent

11009	¾ - 2¼" (10-55 mm)		A4414
			A4427

SenSura® Click Standard Wear Barrier + Drainable Pouch


- Barrier #11031 Convex Light with belt tabs
- Pouch #11196 MAXI / EasiClose™ WIDE Outlet
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit / Transparent

11010	⅝ - 1 ¹¹ / ₁₆ " (15-43 mm)		A4373
			A4427

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

SenSura® Click Xpro Extended Wear Barrier + Urostomy Pouch


- Barrier #10035 Flat with belt tabs
- Pouch #11856 MAXI / soft outlet
- Urostomy adapter
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit / Transparent

11004	⅜ - 2¼" (10-55 mm)		A4409
			A4433

2-Piece Flex Single Use Kits

SenSura® Mio Flex Barrier + MAXI Drainable Pouch


- Extended wear
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**




Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------


Cut-to-fit / Transparent

21009	⅜ - 1⅞" (10-48 mm)		A4410
			A4425


Barrier #10561 Flat, Pouch #12273 EasiClose™ WIDE Outlet / filter

21010	⅜ - 2¾" (10-68 mm)		A4410
			A4425

Barrier #10571 Flat, Pouch #12283 EasiClose™ WIDE Outlet / filter

21012	⅜ - 1⅞" (10-48 mm)		A4410
			A5063

Barrier #10561 Flat, Pouch #12274 EasiClose™ WIDE Outlet / without filter

21013	⅜ - 2¾" (10-68 mm)		A4410
			A5063

Barrier #10571 Flat, Pouch #12284 EasiClose™ WIDE Outlet / without filter

SenSura® Mio Flex + MAXI Urostomy Pouch


- Extended wear
- Urostomy adapter
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit / Transparent

110072	3/8 - 1 7/8" (10-48 mm)		A4410 A5073**
Barrier #10561 Flat, Pouch #12298 Soft outlet / multi-chamber			

SenSura® Mio Flex Convex Barrier + MAXI Drainable Pouch


- Extended wear
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**




Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit / Transparent

13660	5/8 - 1 9/16" (15-40 mm)		A4407 A4425
Barrier #16483 Convex Light, Pouch #12273 EasiClose™ WIDE Outlet / filter			

13661	5/8 - 2 1/16" (15-53 mm)		A4408 A4425
Barrier #16491 Convex Light, Pouch #12283 EasiClose™ WIDE Outlet / filter			

SenSura® Mio Flex Convex + MAXI Urostomy Pouch


- Extended wear
- Urostomy adapter
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit / Transparent

13662	5/8 - 1 9/16" (15-40 mm)		A4407 A5073**
Barrier #16483 Convex, Pouch #12298 Soft outlet / multi-chamber			

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

**Product has not been through PDAC verification

SenSura® Flex Standard Wear Barrier + MAXI Drainable Pouch

- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit / Transparent

11011K	⅜ - 1⅞" (10-48 mm)		A4414
			A4425
Barrier #10102 Flat, Pouch #11512 EasiClose™ WIDE Outlet / filter			
21001	⅜ - 2¾" (10-68 mm)		A4415
			A4425
Barrier #10103 Flat, Pouch #11517 EasiClose™ WIDE Outlet / filter			
21011	⅜ - 3½" (10-88 mm)		A4415
			A5063
Barrier #10108 Flat, Pouch #11690 EasiClose™ WIDE Outlet / without filter			

SenSura® Flex Xpro Extended Wear Barrier + MAXI Drainable Pouch



- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit / Transparent

21002	⅜ - 1⅞" (10-48 mm)		A4409
			A4425
Barrier #10105 Flat, Pouch #11512 EasiClose™ WIDE Outlet / filter			
21003	⅜ - 2¾" (10-68 mm)		A4410
			A4425
Barrier #10106 Flat, Pouch #11517 EasiClose™ WIDE Outlet / filter			

SenSura® Flex Xpro Extended Wear Barrier + MAXI Drainable Pouch


- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit / Transparent

21004	⅝ - 2⅛" (15-53 mm)		A4408
			A4425
Barrier #11308 Convex Light, Pouch #11517 EasiClose™ WIDE Outlet / filter			

 Indicates convexity

Assura® Standard Wear Barrier + Drainable Pouch


- Barrier #2883 Flat with belt tabs
- Pouch #15976 MAXI / EasiClose™ WIDE Outlet
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit / Transparent

11002	3/8 - 2 1/4" (10-55 mm)		A4414
			A4427

Assura® Standard Wear Barrier + Drainable Pouch


- Barrier #14263 Convex Light with belt tabs
- Pouch #15976 MAXI / EasiClose™ WIDE Outlet
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit / Transparent

11001	5/8 - 1 11/16" (15-43 mm)		A4373
			A4427

1-Piece Single Use Kits

SenSura® Mio 1-piece Single Use Kit

- Extended wear
- Pouch #10481 MAXI / EasiClose™ WIDE Outlet
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Transparent Units: 5

Code	Cut	Stoma Size	HCPCS*
21007	Cut-to-fit	3/8 - 2 1/8" (10-55 mm)	A5056

SenSura® Mio Convex Soft Single Use Kit

- Extended wear
- Pouch #16716 MAXI / EasiClose™ WIDE Outlet / Filter
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Transparent Units: 5

Code	Cut	Stoma Size	HCPCS*
13665	Cut-to-fit	3/8 - 2" (10-50 mm)	A5057

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

SenSura® Mio Convex Light Single Use Kit

- Extended wear
- Pouch #16739 MAXI / EasiClose™ WIDE Outlet / Filter
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Transparent Units: 5

Code	Cut	Stoma Size	HCPCS*
13664	Cut-to-fit	$\frac{3}{8}$ - $1\frac{11}{16}$ " (10-43 mm)	A5057

SenSura® Mio Soft Convex Urostomy Single Use Kit

- Extended wear
- Pouch #16837 MAXI / EasiClose™ WIDE Outlet / Filter
- Instructions for use
- Stoma measuring guide
- Urostomy adapter
- **Not made with natural rubber latex**



Transparent Units: 5

Code	Cut	Stoma Size	HCPCS*
13663	Cut-to-fit	$\frac{3}{8}$ - $1\frac{11}{16}$ " (10-43 mm)	A4393

SenSura® 1-piece Single Use Kit

- Pouch #15521 MAXI / EasiClose™ WIDE Outlet
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Transparent Units: 5

Code	Cut	Stoma Size	HCPCS*
11008	Cut-to-fit	$\frac{3}{8}$ - 3" (10-76 mm)	A4424

Assura® 1-piece Single Use Kit

- Pouch #15860 MAXI / EasiClose™ WIDE Outlet
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Transparent Units: 5

Code	Cut	Stoma Size	HCPCS*
11005	Cut-to-fit	$\frac{3}{8}$ - 3" (10-76 mm)	A4424

SenSura® Xpro 1-piece Urostomy Single Use Kit

- Extended wear
- Pouch #11809 MAXI / soft outlet
- Instructions for use
- Stoma measuring guide
- Urostomy adapter
- **Not made with natural rubber latex**



Transparent $10\frac{3}{8}$ " (26 cm) / 400 mL Units: 5

Code	Cut	Stoma Size	HCPCS*
21005	Cut-to-fit	$\frac{3}{8}$ - 3" (10-76 mm)	A4391**

 Indicates convexity

SenSura® Xpro 1-piece Convex Light Urostomy Single Use Kit

- Extended wear
- Pouch #11819 MAXI / soft outlet
- Instructions for use
- Stoma measuring guide
- Urostomy adapter
- **Not made with natural rubber latex**



Transparent 10³/₈" (26 cm) / 400 mL Units: 5

Code	Cut	Stoma Size	HCPCS*
21006	Cut-to-fit	5 ⁵ / ₈ - 1 ⁵ / ₁₆ " (15-33 mm)	A4393**

Extra Protection

Brava® Elastic Barrier Strips

Skin-friendly strips are used to "picture frame" around an ostomy barrier, providing extra security.

- Elastic (follow the contours of the body)
- Ensure the position of the barrier and minimize roll-up
- Skin-friendly and secure alternative to tape
- **Not made with natural rubber latex**



Brava® Elastic Barrier Strips- Curved

For round or oval barrier shapes

Code	Size	Units	HCPCS*
120700†	5½" strips	20 strips	A4362

†2 strips = 1 billable unit of A4362

Brava® Elastic Barrier Strips- Y shape

Designed for users with challenging body shapes, e.g. hernias or skin folds



Code	Size	Units	HCPCS*
120721†	5¾" strips	30 strips	A5121

†3 strips = 1 billable unit of A5121

Brava® Elastic Barrier Strips- Straight

For square or rectangular barrier shapes



Code	Size	Units	HCPCS*
120740†	5¾" strips	20 strips	A4362

†2 strips = 1 billable unit of A4362

Brava® Elastic Barrier Strip XL

Designed for increased security & support around the barrier. Wider shape accommodates the petals of Convex Flip.



Code	Size	Units	HCPCS*
12076**	6¾" strips	30 strips	A5121

**3 strips = 1 billable unit of A5121

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

**Product has not been through PDAC verification

Brava® Protective Seal

Moldable barrier ring designed to adhere securely and protect the skin around the stoma.

- Soft, flexible material adapts to body shape as the wearer bends and stretches
- Can be easily and cleanly removed
- Alcohol-free and sting-free
- Multiple sizes for a custom fit
- **Not made with natural rubber latex**



Code	Stoma size	Adhesive diameter	Thickness	Seal width	Units	HCPCS*
12035	¾"	1 7/8"	2.5 mm	Standard (5/8")	10	A4385
12036	¾"	2 ¼"	2.5 mm	Wide (¾")	10	A4385
12038	¾"	2 ½"	2.5 mm	Wide (7/8")	10	A4385
12037	1 1/8"	2 ¼"	2.5 mm	Standard (5/8")	10	A4385
12039	1 3/8"	2 ½"	2.5 mm	Standard (5/8")	10	A4385
12045	¾"	1 7/8"	4.2 mm	Standard (5/8")	10	A4385
12046	¾"	2 ¼"	4.2 mm	Wide (¾")	10	A4385
12048	¾"	2 ½"	4.2 mm	Wide (7/8")	10	A4385
12047	1 1/8"	2 ¼"	4.2 mm	Standard (5/8")	10	A4385
12049	1 3/8"	2 ½"	4.2 mm	Standard (5/8")	10	A4385
12032	¾"	3"	2.5 mm	Extra wide (1 1/8")	10	A4385
12033	1 1/8"	3"	2.5 mm	Extra wide (1 5/16")	10	A4385

Brava® Protective Seal Convex

Moldable barrier ring designed to provide convexity for flush or retracted stomas.

- High-tack adhesive that adheres quickly to the skin, filling small creases, folds, or scarring around the stoma
- Erosion resistant composition
- High flexibility to maintain seal during body movement
- Steep sides designed to provide convexity to peristomal area
- Thin edges fit snugly around peristomal area
- **Not made with natural rubber latex**



Code	Stoma size	Thickness	Units	HCPCS*
12090	¾"	8.2mm	10	A4411
12091	1"	8.2mm	10	A4411
12092	1 3/16"	8.2mm	10	A4411
12093	1 3/8"	8.2mm	10	A4411
12094	1 9/16"	8.2mm	10	A4411
12095	1 ¾"	8.2mm	10	A4411

Brava® Moldable Ring

Designed for dual-action leakage reduction - used to fill in uneven skin surfaces and protect the skin.

- Molds and adheres securely to dips, creases, and folds to help create a tight seal around the stoma
- Helps keep skin healthy by protecting it from output and absorbing moisture from the skin
- Easy to handle, mold and apply
- Alcohol-free and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
120307	2 mm (thin)	10	A4385
120427	4.2 mm (thick)	10	A4385

Brava® Strip Paste

Pectin-based product that protects peristomal skin.

- Protects peristomal skin from stoma output
- Used to fill uneven skin surfaces to create a flat pouching surface
- Alcohol-free and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
26555	2 oz. total package (.2 oz. per strip)	10 strips	A4406

Brava® Paste

A sting-free alternative to conventional alcohol-based paste.

- Creates a tight, effective seal between peristomal skin and the ostomy/fistula appliance
- No alcohol = no need to "set"; contains pectin
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
12050	2 oz.	1	A4406

Coloplast® Paste

Used to fill in uneven skin surfaces to create a flat pouching surface.

- Provides tight sealing around the stoma
- Protects against stoma output, minimizing risk of maceration
- Low alcohol content; no pectin
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
2650	2 oz.	12	A4405

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

Brava® Skin Cleanser Wipes

The wipe that cleans and cares. A convenient way to cleanse and maintain healthy skin when changing the ostomy pouching system.



- Contains ProVitamin B5 to provide a soothing and moisturizing effect
- Does not contain alcohol or perfumes
- No residue left on the skin that could affect barrier adhesion
- **Not made with natural rubber latex**

Code	Size	Units	HCPCS*
12080	1 wipe	15 wipes/pack	A9286**

Brava® Lubricating Deodorant

A pouch lubricant + deodorizer in one. Powder-fresh liquid reduces odor and allows output to slide to the bottom of the pouch for easier emptying.



- **Not made with natural rubber latex**

Code	Size	Units	HCPCS*
12060	.25 fl. oz. (7.5 mL) sachet	20 sachets	A4394
12061	8 oz. bottle	16	A4394

Brava® Adhesive Remover Spray

- Designed to lift the edges of the barrier for ease in removal
- Dry within seconds without leaving an oily residue behind
- Alcohol-free, silicone-based, and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
120105	1.7 oz. (50 mL)	1	A4455**

Brava® Adhesive Remover Wipes

- Wipe away adhesive residue
- Dry within seconds without leaving an oily residue behind
- Alcohol-free, silicone-based, and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
120115	1 wipe	30 wipes/box	A4456**

Brava® Skin Barrier Spray

- Protects against stoma output and adhesives
- Dry within seconds without building a thick or tacky layer
- Alcohol-free, silicone-based, and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
120205	1.7 oz. (50 mL)	1	A4369**

Brava® Skin Barrier Wipes

- Protects against stoma output and adhesives
- Dry within seconds without building a thick or tacky layer
- **Tip:** wring the wipe like a wet rag to expel the liquid, then wipe to disperse on skin
- Alcohol-free, silicone-based, and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
120215	1 wipe	30 wipes/box	A5120**

*For important information regarding the use of HCPCS codes listed and code definitions see page 120


**Product has not been through PDAC verification

BodyCheck

Identifying ostomy solutions based on individual needs

Together with 300 Wound, Ostomy, Contenance (WOC) nurses, we've developed BodyCheck body profile assessment which uses Body Profile Terminology, a common language to describe the peristomal body profile.

- 1 Assess your patient's overall peristomal area as a Regular, Inward, or Outward



- 2 Visit www.bodycheck.coloplast.us to complete the assessment



- 3 Receive personalized ostomy solutions *(may include accessory options)*



Visit www.bodycheck.coloplast.us


Henry

Making life easier_

For little fighters

Babies and kids present different pouching challenges than adults - so we've applied **BodyFit Technology** to develop specific solutions for different stages in a child's life.

Stepping up over time, as kids grow



Ask us about our pediatric education!



Pediatric FAQ

A quick tipsheet with answers to frequently asked questions, like "what do I do if friends and family - or other kids - have questions?"

and "what do I do if my child pulls the pouch off?"

Caregivers can request this via our Coloplast® Care program.



My Child Has An Ostomy: A Resource Sheet for Caregivers

A fill-in-the-blank sheet to help parents communicate about their child's needs to

day care providers, school nurses, and other caregivers.

Caregivers can request this via our Coloplast® Care program.




Life After Your Child's Ostomy

A high-level overview of anatomy, how to change a pouch, pouching system options, where to get

supplies, dietary notes, and other information in English or Spanish.

Caregivers can request this via our Coloplast® Care program.



When I Met Tipster: A coloring book for children with ostomies

A simple coloring book to help kids engage with their ostomy surgery and "new normal" - or learn more about a sibling's ostomy.

Available to hospitals via Coloplast sales representatives, or to caregivers via our Coloplast® Care program.



Jake Gets a Stoma

A children's book written by Jake's mom, who wanted a better way to explain the surgery and stoma management to her son and other kids going through the same thing.

Available to hospitals via Coloplast sales representatives, or to caregivers via our Coloplast® Care program.




SenSura® Mio Bunny

Coloplast offers a complimentary training bunny for product education use with both

SenSura® Mio Baby & Kids. Contact your sales representative for more information.

Limitations apply. The SenSura® Mio bunny is provided to healthcare providers to offer to families for teaching purposes only and should only be used with Coloplast products.



Teen Life with an Ostomy

Navigating school, what to wear, and hanging out with friends aren't covered by the average adult-oriented materials -

that's why we're so excited to feature tips from Paige, a teen who's been through it all!

Available to patients via our Coloplast® Care program.

SenSura® Mio Baby

2-piece pouching system designed with BodyFit Technology to tackle the challenges of preemies and babies from 30 weeks of gestation to 6 months of age.


SenSura® Mio Baby Flex 2-Piece Barrier

- Oval barrier accommodates one or multiple stomas
- No starter hole – custom cut
- Gentle adhesive sticks securely and removes gently with Brava® Baby Adhesive Remover
- Soft flange enables the pouch to be removed and replaced
- Sterilized surface, single-packed
- **Not made with natural rubber latex**



Cut-to-fit

Units: 5

Code	Stoma Size	Color	HCPCS*
18705	0 - 1 ⁹ / ₁₆ " (0-40 mm)		A4414

SenSura® Mio Baby Flex 2-Piece Pouch

- Adhesive coupling is easy to apply and remove
- Access window
- Soft outlet fits to enteral feeding systems
- Sterile, single-packed
- **Not made with natural rubber latex**



Transparent, No filter

Units: 10

Code	Color	HCPCS*
18700		A5063

Brava® Baby Adhesive Remover

- Oil-based adhesive remover designed specifically for babies
- Liquid sachet (no wipe) can be used with cotton swabs or gauze
- Gently removes adhesive / adhesive residues from SenSura® Mio Baby Barrier
- Not recommended for use with SenSura Mio Kids or adult products
- **Not made with natural rubber latex**



Units: 30

Code	Size	HCPCS*
12014	3 mL sachet	A4455

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

SenSura® Mio Kids

1-piece and 2-piece pouching systems designed with BodyFit Technology to fit kids' active movements and rounded bellies. Designed with options to step up as kids grow, the range covers kids 6 months to 4 years of age (larger options for older kids are found in the SenSura Mio adult line).

SenSura® Mio Kids Flex 2-Piece Barrier

- Star-shaped, elastic barrier can be rotated for flexibility in placement
- No starter hole – custom cut
- **Not made with natural rubber latex**



Cut-to-fit

Units: 5

Code	Stoma Size	Color	HCPCS*
18735	0 - 1 ¹ / ₁₆ " (0-27 mm)		A4414
18747	0 - 1 ³ / ₈ " (0-35 mm)		A4414

SenSura® Mio Kids Flex 2-Piece MINI Drainable Pouch

- EasiClose™ Outlet
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Small size enables kids to step up in pouch size over time
- Opaque includes inspection window
- **Not made with natural rubber latex**

6¹/₈" folded

Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
18721	18720		A4425



SenSura® Mio Kids Flex 2-Piece MIDI Drainable Pouch

- EasiClose™ Outlet
- Full-circle filter
- Neutral gray, water-resistant textile
- Opaque includes inspection window
- **Not made with natural rubber latex**



6¾" folded

Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
18723	18722		A4425
18725	18724		A4425



SenSura® Mio Kids Flex 2-piece MIDI Urostomy Pouch

- Neutral gray, water-resistant textile
- Soft outlet can be connected to a drain bag for nighttime use, or a leg bag for kids in a wheelchair
- **Not made with natural rubber latex**



Opaque 6¾" folded

Units: 10

Code	Color	HCPCS*
18730		A5073
18731		A5073

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

SenSura® Mio Kids 1-Piece MIDI Drainable Pouch

- Star-shaped, elastic barrier with no starter hole – custom cut
- EasiClose™ Outlet
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque includes inspection window
- **Not made with natural rubber latex**



6¾" folded | Cut-to-fit

Units: 10

Code	Code	Stoma size	HCPCS*
Transparent	Opaque		
18711	18710	0 – 1 ³ / ₈ " (0-35 mm)	A4424

SenSura® Mio Kids 1-Piece MIDI Urostomy Pouch

- Star-shaped, elastic barrier with starter hole
- Neutral gray, water-resistant textile
- Soft outlet can be connected to a drain bag for nighttime use, or a leg bag for kids in a wheelchair
- **Not made with natural rubber latex**



Units: 10

Code	Code	Stoma size	HCPCS*
Transparent	Neutral gray		
18712	18713	³ / ₈ - 1 ³ / ₈ " (10-35 mm)	A5071

Assura® AC (Adhesive Coupling) Pediatric

A simple, soft, secure system for neonates! This ostomy pouching system is based on adhesives that make it easy to apply the pouch to the barrier. The pouching system is easy to clean, and simple to handle with the EasiClose™ standard integrated outlet.

Assura® AC Pediatric Flat Barrier

- Standard wear
- No starter hole
- **Not made with natural rubber latex**



Cut-to-fit

Units: 10



Code	Stoma Size	Color	HCPCS*
14307	0 - 5/8" (0-17 mm)		A4414
14308	0 - 1" (0-27 mm)		A4414

Assura® AC Pediatric Drainable Pouch

- EasiClose™ Outlet
- Filter
- Soft cloth on back
- **Not made with natural rubber latex**



Units: 10

Code	Code	Length	Color	HCPCS*
	Transparent			
	Opaque (Teddy Bear)			
14681	14691	5½" (14 cm)		A4425
14682	14692	6½" (16 cm)		A4425

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

Assura® 2-piece Pediatric Pouching Systems

Specifically developed adhesive and a smaller barrier size are ideal for pediatrics, while a mix-and-match “green” flange size offers versatility. The SenSura® and Assura® pouches and barrier may be used interchangeably with “green”-sized adult SenSura® Click or Assura® products for a broad array of pouching options.

Assura® Pediatric Standard Wear Barrier

- Flat
- **Not made with natural rubber latex**



Cut-to-fit

Units: 5

Code		Color	HCPCS*
2182	3/8 - 1 3/8" (10-35 mm)		A4414

Assura® Pediatric Drainable Pouch

- Clamp closure
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 8 1/2" (21 cm) / 150 mL

Units: 10

Code		Color	HCPCS*
2156			A4426

Assura® Pediatric Closed Pouch

- Soft cloth on front and back
- **Not made with natural rubber latex**



Opaque 5 3/4" (14 1/2 cm)

Units: 30

Code		Color	HCPCS*
2161			A4423

Assura® Stoma Cap

- With filter
- Small size is optimal for swimming at the pool
- **Not made with natural rubber latex**



Opaque

Units: 30

Code		Color	HCPCS*
2801			A5055

Assura® Pediatric Urostomy Pouch

- Soft outlet can be connected to a drain bag for nighttime use, or a leg bag for kids in a wheelchair - see page 16 for options
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 8½" (21 cm) / 150 mL

Units: 10

Code	Color	HCPCS*
1789		A4433

Assura® Urostomy Micro-Pouch

- Miniature size prevents twisting
- Soft outlet can be connected to a drain bag for nighttime use, or a leg bag for people in a wheelchair
- Soft cloth on front
- **Not made with natural rubber latex**



Transparent 6" (15 cm)

Units: 10

Code	Color	HCPCS*
14204		A4433

Assura® Urostomy Minicap

- Ultra-absorbent material inside
- Small size is optimal for swimming at the pool
- **Not made with natural rubber latex**



Opaque 100 mL

Units: 30

Code	Color	HCPCS*
2807		A5055

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

ColoKids™ 1-piece Pediatric Pouching Systems

Specially developed Assura® “swiss roll” technology provides moisture absorption and balance while offering consistent tack and adhesion to the skin. Small pouch sizes provide a pediatric-friendly pouching option.

ColoKids™ Flat 1-piece Drainable Pouch

- Standard wear
- Clamp closure
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent

Units: 10

Code	Cut	Stoma Size	HCPCS*
2115	Cut-to-fit	3/8 - 1 3/8" (10-35 mm)	A5061

ColoKids™ Flat 1-piece Closed Pouch

- Standard wear
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent

Units: 30

Code	Cut	Stoma Size	HCPCS*
2125	Cut-to-fit	3/8 - 1 3/8" (10-35 mm)	A4416

ColoKids™ Flat 1-piece Urostomy Pouch

- Standard wear
- Soft outlet can be connected to a drain bag for nighttime use, or a leg bag for people in a wheelchair
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent

Units: 10

Code	Cut	Stoma Size	HCPCS*
8009	Cut-to-fit	5/8 - 1 3/8" (15-35 mm)	A5071

Brava® supporting products

Pediatric weight thresholds

How does Coloplast define pediatric age / weight parameters as relate to products?

- Premature: 30+ gestational weeks
- Neonate: 0-6 months
- Toddler: 6 months – 4 years

Product name	Caution	Conversion (kg to lbs.)	Approximate minimum age
Brava® Elastic Barrier Strips	Not recommended for children below 5.6 kg	12.34 lbs.	Toddler
Brava® Protective Seal	Not recommended for children below 5.6 kg	12.34 lbs.	Toddler
Brava® Moldable Ring	Not recommended for children below 7.2 kg	15.87 lbs.	Toddler
Brava® Strip Paste	Not recommended for children below 2.6 kg	5.73 lbs.	Neonate
Brava® Paste (tube paste)	Not recommended for children below 6 kg	13.22 lbs.	Toddler
Brava® Adhesive Remover Wipe	Not recommended for children below 2.4 kg	5.29 lbs.	Neonate
Brava® Adhesive Remover Spray	Not recommended for children below 5.6 kg	12.34 lbs.	Toddler
Brava® Baby Adhesive Remover	Not recommended for children below 2.4kg	5.29 lbs.	Neonate
Brava® Skin Barrier Wipe	Not recommended for children below 2.4 kg	5.29 lbs.	Neonate
Brava® Skin Barrier Spray	Not recommended for children below 5.6 kg	12.34 lbs.	Toddler
Brava® Skin Cleanser	Not recommended for children below 2.4 kg	5.29 lbs.	Neonate
Brava® Protective Sheet	Not recommended for children below 7 kg	15.43 lbs.	Toddler
Brava® Belt	Does not fit kids	-	-
Brava® Lubricating Deodorant	Not recommended for children	-	-
Brava® Powder	Not recommended for children below 5.6 kg	12.34 lbs.	Toddler

kg to lbs.

1 kg = 2.20 lbs.

2.4 kg = 5.29 lbs.

2.6 kg = 5.73 lbs.

5.6 kg = 12.34 lbs.

kg to lbs.

6.0 kg = 13.22 lbs.

7.0 kg = 15.43 lbs.

7.2 kg = 15.87 lbs.

*For important information regarding the use of HCPCS codes listed and code definitions see page 120

Night Drainage Bags

Urostomy Night Bag and Bag Hanger

See Accessories page 16 for more information.





Gabriela

Making life easier_

Coloplast® Care is the only manufacturer-provided patient support program with better patient outcomes supported by peer-reviewed data.



About the team

Coloplast® Care Advisors assist patients with their intimate healthcare needs. They are available Monday through Friday from 7 a.m. to 7 p.m. (CT). The team offers one-on-one assistance for the following:

Quality of life guidance

Our team can answer questions regarding product lines, routines that will fit their lifestyle, and encourage self-assessments for ongoing wellness support.

Connecting to the care network

Coloplast Care can help locate a medical supplies dealer in-network with your patient's insurance to ensure that they have access to the correct products.

Community resources

Coloplast Care Advisors also provide awareness of resources in communities where your patients live, such as support groups, events, or conferences. The Advisors can also offer free samples*, if needed, from the Coloplast ostomy and continence product lines.

For ostomy: visit www.ostomy.coloplastcare.us or call **1-877-858-2656** to enroll your patients today.

For continence: visit bladder.coloplastcare.us or call **1-866-226-6362** to enroll your patients today.

For Peristeen® inquiries, visit coloplastcare.us/bowel or call **1-855-605-7594** for more information.

*Limitations apply. SpeediCath®, Self-Cath® and SureCath® catheters are available by prescription only. Samples are not available for Peristeen.

Information from Coloplast Care is for educational purposes only. It is not intended to substitute for professional medical advice and should not be interpreted to contain treatment recommendations.

Coloplast® Care *Connect*

ColoplastCareConnect is the new easy online portal for enrolling your patients and access to team documentation. Scan below to enroll your patients or to view your dashboard of electronic submissions:



Contact your rep if you need access.

The following materials are presented for general purposes only. They do not constitute medical advice and should not be a substitute for consultation with a trained medical professional.

The information presented or discussed may not be representative of all patient outcomes. Each person's situation is unique and risks, outcomes, experience, and results may vary. Please see complete product instructions for use, including all product indications, contraindications, precautions and warnings, and adverse events.

Ostomy Self-Assessment Tools

- BodyCheck
- Ostomy Check
- Troubleshooter
- Apply-Remove-Check (ARC)

Diaries:

- Pouch change diary
- Exercise diary
- Hernia prevention action plan
- Ostomy Wellness Output Diary

Educational Guides

- Life After: Ileostomy
- Life After: Colostomy
- Life After: Urostomy
- Life After: Your Child's Ostomy
- Teen Life with an Ostomy

Product Pictorial Guide - Also available in Spanish

NOTE: These are abbreviated instructions designed to visually show how to use the product. Always refer to the product labeling (including the package insert) before using.

- SenSura® Mio 1-piece Pouching System
- SenSura® Mio Click 2-piece Pouching System
- SenSura® Mio Flex 2-piece Pouching System
- SenSura® 1-piece Pouching System
- SenSura® Click 2-piece Pouching System
- SenSura® Flex 2-piece Pouching System
- Fistula Wound Management System
- SenSura® Mio 1-piece Urostomy Pouching System
- SenSura® Mio Click 2-piece Urostomy Pouching System
- SenSura® Mio Flex 2-piece Urostomy Pouching System
- SenSura® Mio Baby 2-piece Flex Barrier, pouch, and Brava® Baby Adhesive Remover

Spanish Educational Tools

Educational Guides

- La Vida Después de la Colostomía Life After: Colostomy
- La Vida Después de la Ileostomía Life After: Ileostomy
- La Vida Después de la Urostomía Life After: Urostomy

Other languages also available in digital format only:
Arabic, Chinese, Russian and Vietnamese

For additional educational resources, visit
www.coloplast.us



**Online Learning with
CE Contact Hours**

Courses for Healthcare Professionals working in continence, ostomy and wound clinical spaces.



Expert Education Programs

Live or virtual education on a variety of clinical topics presented by content experts, providing interactive learning experiences.



Facility-led Education

Train the Trainer courses include slide presentations, training resources, toolkits and more to help you facilitate a meaningful training experience.



Education Resources

Videos, podcasts and printables designed to support your patients with product questions and needs in everyday situations and help you promote positive outcomes for your patients.



Register by using one of the following methods:

Website - www.coloplast.us/professional

Email - CProfessional_US@coloplast.us

Medicare Utilization - Ostomy

Usual maximum limits for some products from the Medicare Local Coverage Determinations (LCD) are shown below. The amount of ostomy supplies varies by patient based on the type of ostomy, its location, its construction, and the condition of the peristomal skin. The table below shows the maximum number of items that are usually medically necessary for patients with ostomies. It may be more or less depending on the factors that affect the frequency of barrier and pouch change. The medical necessity for a greater quantity of supplies must be clearly documented in the patient's medical record (physician's signed and dated order).

HCPCS*	Short Description	Usual Maximum
A4357	Bedside drainage bag	2 ea/month
A4361	Ostomy faceplate, each	3 ea/6 months
A4362	Skin barrier, solid, 4x4 or equivalent; each*	20 ea/month
A4363	Ostomy clamp replacement	none*
A4367	Ostomy belt, each	1 ea/month
A4368	Ostomy filter	none*
A4369	Ostomy skin barrier; liquid, per ounce	2 oz./month
A4371	Ostomy skin barrier powder, per ounce	10 oz/6 months
A4372	Ostomy skin barrier, 4 x 4, standard wear, CONVEX	none*
A4373	Ostomy skin barrier, with flange, CONVEX	none*
A4385	Ostomy skin barrier, 4 x 4, extended wear, FLAT	none*
A4388	One-piece ostomy drainable pouch and barrier, extended wear, FLAT	none*
A4389	One-piece ostomy drainable pouch and barrier, CONVEX	none*
A4390	One-piece ostomy drainable pouch, extended wear barrier, CONVEX	none*
A4391	One-piece urinary pouch, extended wear barrier, FLAT	none*
A4392	One-piece urinary pouch, standard wear barrier, FLAT	none*
A4394	Ostomy deodorant per ounce	none*
A4397	Ostomy irrigation supply; sleeve, each • HCPCS code A4397 will be discontinued effective 12/31/2021. Effective 1/1/2022 CMS has established two new HCPCS codes to distinguish between reusable and disposable sleeves: - A4436 Irrigation supply, sleeve, reusable - A4437 Irrigation supply, sleeve, disposable	4 ea/month TBD TBD
A4398	Ostomy irrigation supply; bag, each	2 ea/6 months
A4399	Ostomy irrigation supply; cone/catheter, with or without brush	2 ea/6 months
A4404	Ostomy ring, each	10 ea/month
A4405	Ostomy skin barrier, non-pectin based, paste, per ounce	4 oz/month
A4406	Ostomy skin barrier, pectin-based, per ounce	4 oz/month
A4407	Two-piece ostomy skin barrier, with flange, extended wear, CONVEX, (4 x 4 inches or smaller)	none*
A4408	Two-piece ostomy skin barrier, with flange extended wear, CONVEX, (larger than 4 x 4 inches)	none*
A4409	Two-piece ostomy skin barrier, with flange, extended wear, FLAT, (4 x 4 inches or smaller)	none*

* No published Medicare usual maximum. Monthly quantity subject to medical necessity documentation.

**Elastic Barrier Strips - Curved or Straight: 40 strips/month

www.coloplast.us

1-800-533-0464

HCPCS*	Short Description	Usual Maximum
A4410	Two-piece ostomy skin barrier, with flange, extended wear, FLAT, (larger than 4 x 4 inches)	none*
A4411	Two-piece ostomy solid skin barrier, extended wear, CONVEX	none*
A4412	Two-piece high output ostomy drainable pouch without filter	none*
A4413	Two-piece high output ostomy drainable pouch with filter	none*
A4414	Two-piece ostomy skin barrier with flange, FLAT (4x4 inches or smaller)	20 ea/month
A4415	Two-piece ostomy skin barrier, with flange, FLAT, greater than 4 x 4 inches	20 ea/month
A4416	One-piece ostomy closed pouch and barrier with filter, FLAT	60 ea/month
A4417	One-piece ostomy closed pouch and barrier with filter, CONVEX	60 ea/month
A4418	One-piece ostomy closed pouch with filter and barrier	60 ea/month
A4419	Two-piece ostomy closed pouch with filter for non-locking system	60 ea/month
A4423	Two-piece ostomy closed pouch with filter for locking system	60 ea/month
A4424	One-piece ostomy drainable pouch and barrier with filter, FLAT	20 ea/month
A4425	Two-piece ostomy drainable pouch with filter for non-locking system	20 ea/month
A4426	Two-piece ostomy drainable pouch without filter for locking system	20 ea/month
A4427	Two-piece ostomy drainable pouch with filter for locking system	20 ea/month
A4433	Two-piece ostomy urinary pouch for locking system	20 ea/month
A4436	Irrigation supply, sleeve, reusable, per month • Effective 1/1/2022	TBD
A4437	Irrigation supply, sleeve, reusable, per month • Effective 1/1/2022	TBD
A4455	Adhesive remover or solvent, per ounce	16 oz/6 months
A4456	Adhesive remover, wipes	none*
A5055	Stoma cap	31 ea/month
A5056	One-piece drainable ostomy pouch, extended wear with filter	40 ea/month
A5057	One-piece drainable ostomy pouch extended wear with filter, CONVEX	40 ea/month
A5061	One-piece ostomy drainable pouch and barrier, FLAT	20 ea/month
A5063	Two-piece ostomy drainable pouch for barrier with flange	20 ea/month
A5071	One-piece ostomy urinary pouch	20 ea/month
A5073	Two-piece ostomy urinary pouch	20 ea/month
A5081	Stoma plug or seal	31 ea/month
A5120	Skin barrier, wipes or swabs, each	150 ea/6 months
A5121	Solid skin barrier, 6x6 or equivalent, each**	20 ea/month
A5122	Solid skin barrier, 8x8 or equivalent, each	20 ea/month
A6154	Wound pouch, each	3/week

* No published Medicare usual maximum. Monthly quantity subject to medical necessity documentation.

**Elastic Barrier Strips - Y-shape: 60 strips/month

Reimbursement Disclaimer: Coloplast Corp. provides this information for convenience and your general reference only. It does not constitute legal advice or a recommendation regarding clinical practice. Reimbursement, coverage and payment policies can vary from one insurer and region to another and is subject to change without notice. The provider has the responsibility to determine medical necessity and to submit appropriate codes and charges for care provided. Coloplast does not guarantee coverage or payment of products and Coloplast makes no guarantee that the use of this information will prevent differences of opinion or disputes with Medicare, insurers, or other payers as to the correct form of billing or the amount that will be paid. This information is provided for your general information only and is not intended to replace any advice you receive from your own internal or external insurance coverage consultants, reimbursement specialists or legal counsel.

HCPCS Code Reference


A4357	Bedside drainage bag, day or night, with or without anti-reflux device, with or without tube, each
A4361	Ostomy faceplate, each
A4362	Solid skin barrier; solid, 4x4 or equivalent; each
A4363	Ostomy clamp, any type, replacement only, each
A4366	Ostomy vent, any type, each
A4367	Ostomy belt, each
A4368	Ostomy filter, any type, each
A4369	Ostomy skin barrier liquid (spray, brush, etc.), per ounce
A4371	Ostomy skin barrier powder, per ounce
A4372	Ostomy skin barrier, solid 4x4 or equivalent, standard wear, with built-in convexity, each
A4373	Ostomy skin barrier, with flange (solid, flexible or accordion), with built-in convexity, any size, each
A4385	Ostomy skin barrier, solid 4 x 4 or equivalent, extended wear, without built-in convexity, each
A4388	Ostomy pouch, drainable, with extended wear barrier attached, (1 piece), each
A4389	Ostomy pouch, drainable, with barrier attached, with built-in convexity (1-piece), each
A4390	Ostomy pouch, drainable, with extended wear barrier attached, with built-in convexity (1-piece), each
A4391	Ostomy pouch, urinary, with extended wear barrier attached (1-piece), each
A4392	Ostomy pouch, urinary with standard wear barrier attached, with built-in convexity (1-piece), each
A4393	Ostomy pouch, urinary, with extended wear barrier attached, with built-in convexity (1 pc), each
A4394	Ostomy deodorant, with or without lubricant, for use in ostomy pouch, per fluid ounce
A4396	Ostomy belt with peristomal hernia support
A4397	Irrigation supply; sleeve, each
A4398	Ostomy irrigation supply; bag, each
A4399	Ostomy Irrigation supply; cone /catheter, with or without brush
A4404	Ostomy ring, each
A4405	Ostomy skin barrier, non-pectin based, paste, per ounce
A4406	Ostomy skin barrier, pectin-based, paste, per ounce
A4407	Ostomy skin barrier, with flange (solid, flexible, or accordion), extended wear, with built-in convexity, 4x4 inches or smaller, each

HCPCS Code Reference

A4408	Ostomy skin barrier, with flange (solid, flexible or accordion), extended wear, with built-in convexity, larger than 4x4 inches, each
A4409	Ostomy skin barrier, with flange (solid, flexible or accordion), extended wear, without built in convexity, 4x4 inches or smaller, each
A4410	Ostomy skin barrier with flange (solid, flexible or accordion), extended wear, without built-in convexity, larger than 4x4 inches, each
A4411	Ostomy skin barrier, solid 4x4 or equivalent, extended wear, with built-in convexity, each
A4412	Ostomy pouch, drainable, high output, for use on a barrier with flange (2-piece system), without filter, each
A4413	Ostomy pouch, drainable, high output, for use on a barrier with flange (2 piece system) with filter, each
A4414	Ostomy skin barrier, with flange (solid, flexible or accordion), without built-in convexity, 4x4 inches or smaller, each
A4415	Ostomy skin barrier, with flange (solid, flexible or accordion), without built-in convexity, larger than 4x4 inches, each
A4416	Ostomy pouch, closed, with barrier attached, with filter (1-piece,) each
A4417	Ostomy pouch, closed, with barrier attached, with built-in convexity, with filter (1-piece), each
A4418	Ostomy pouch, drainable, high output, for use on a barrier with flange (2-piece system), without filter, each
A4419	Ostomy pouch, closed; for use on barrier with non-locking flange, with filter (2-piece), each
A4421	Ostomy supply; miscellaneous
A4422	Ostomy pouch absorbent material (sheet, pad, crystal packet) for use in ostomy pouch to thicken liquid stomal output, each
A4423	Ostomy pouch, closed; for use on barrier with locking flange, with filter (2-piece), each
A4424	Ostomy pouch, drainable, with barrier attached, with filter (1-piece), each
A4425	Ostomy pouch, drainable; for use on barrier with non-locking flange, with filter (2-piece system), each
A4426	Ostomy pouch, drainable; for use on barrier with locking flange (2-piece system), each
A4427	Ostomy pouch, drainable; for use on barrier with locking flange, with filter (2 piece system), each
A4428	Ostomy pouch, urinary, with extended wear barrier attached, with faucet-type tap with valve (1-piece), each
A4429	Ostomy pouch, urinary, with barrier attached, with built-in convexity, with faucet-type tap with valve (1-piece), each
A4430	Ostomy pouch, urinary, with extended wear barrier attached, with built-in convexity, with faucet-type tap with valve (1-piece), each
A4431	Ostomy pouch, urinary; with barrier attached, with faucet-type tap with valve (1- piece), each
A4432	Ostomy pouch, urinary; with barrier attached, with faucet-type tap with value (1- piece), each
A4433	Ostomy pouch, urinary; for use on barrier with locking flange (2-piece), each
A4434	Ostomy pouch, urinary; for use on barrier with locking flange, with faucet-type tap with valve (2-piece), each
A4435	Ostomy pouch, drainable, high output, with extended wear barrier (1-piece system), with or without filter, each

HCPCS Code Reference

A4455	Adhesive remover or solvent (for tape, cement or other adhesive), per ounce
A4456	Adhesive remover, wipes, any type, each
A5051	Ostomy pouch, closed; with barrier attached (1-piece), each
A5052	Ostomy pouch, closed; without barrier attached (1-piece), each
A5053	Ostomy pouch, closed; for use on faceplate, each
A5054	Ostomy pouch, closed; for use on barrier with flange (2-piece), each
A5055	Stoma cap
A5056	Ostomy pouch, drainable, with extended wear barrier attached, with filter (1-piece), each
A5057	Ostomy pouch, drainable, with extended wear barrier attached, with built in convexity, with filter, (1-piece), each
A5061	Ostomy pouch, drainable; with barrier attached, (1-piece), each
A5062	Ostomy pouch, drainable; without barrier attached (1-piece), each
A5063	Ostomy pouch, drainable; for use on barrier with flange (2-piece system), each
A5071	Ostomy pouch, urinary; with barrier attached (1-piece), each
A5072	Ostomy pouch, urinary; without barrier attached (1-piece), each
A5073	Ostomy pouch, urinary; for use on barrier with flange (2-piece), each
A5081	Stoma plug or seal, any type
A5093	Ostomy accessory; convex insert
A5120	Skin barrier, wipes or swabs, each
A5121	Skin barrier; solid, 6x6 or equivalent, each
A5122	Skin barrier; solid, 8x8 or equivalent, each



Mike

Making life easier_

Index

Making life *easier*

Code	Page #	Code	Page #	Code	Page #
0509	14	5580	78	10481	53
0925	14	5585	78	10482	53
1003	17	7583	13	10483	53
1110	18	8009	112	10484	53
1500	17	8114	37,84	10485	53
1500H	17	8115	37,84	10486	53
1751	40	8116	37,84	10487	53
1752	40	9500	14,90	10488	53
1754	40	10011	27	10489	54
1755	40	10012	27	10490	53
1756	40	10013	27	10491	53
1757	40	10015	27	10492	54
1758	40	10016	27	10493	54
1759	40	10017	27	10494	54
1789	111	10021	27	10495	54
2041	14	10022	27	10496	54
2115	112	10023	27	10502	21
2125	112	10025	27	10503	21
2156	110	10026	27	10504	21
2161	110	10027	27	10505	21
2182	110	10031	27	10506	21
2310	14	10032	27	10512	21
2315	14	10035	27	10513	21
2320	14	10036	27	10514	21
2325	14	10041	27	10515	21
2330	14	10045	27	10516	21
2340	14	10101	47	10522	21
2350	14	10102	47	10523	21
2501	77	10103	47	10524	21
2650	11,99	10104	47	10525	21
2801	39,110	10105	47	10531	21
2802	39	10106	47	10551	43
2803	39	10108	47	10561	43
2807	39,111	10154	30	10562	43
2808	39	10155	30	10563	43
2831	32	10156	30	10564	43
2832	32	10164	30	10565	43
2833	32	10165	30	10566	43
2836	38,88	10166	30	10567	43
2837	38,88	10167	30	10568	43
2846	37,84	10187	30	10571	43
2847	37,84	10461	53	10572	43
2881	32	10471	53	10573	43
2882	32	10472	53	10586	62
2883	32	10473	53	10870	59
4215	15	10474	53	10871	59
4220	15	10475	53	10872	59
4237	15	10476	53	10873	59
4247	15	10477	53	10874	59
5070	16	10478	53	10881	59
5570	78	10479	53	10883	59
5575	78	10480	53	10884	59

Code	Page #	Code	Page #	Code	Page #
10885	59	11307	47	11819	70
10911	49	11308	48	11841	31
10912	49	11312	48	11844	31
10914	49	11313	48	11845	31
10915	49	11314	48	11846	31
10917	49	11315	48	11851	31
10918	49	11401	25	11854	31
11001	95	11402	25	11855	31
11002	95	11412	25	11856	31
11004	92	11413	25	12003	15
11005	96	11421	25	12004	15
11008	96	11422	25	12005	15
11009	91	11425	25	12006	15
11010	91	11426	25	12007	15
11011	28	11461	23	12008	15
11011K	94	11462	23	12009	15
11012	28	11463	23	12014	105
11015	28	11471	23	12032	10, 98
11016	28	11472	23	12033	10, 98
11021	28	11473	23	12035	10, 98
11022	28	11474	23	12036	10, 98
11023	28	11482	23	12037	10, 98
11024	28	11483	23	12038	10, 98
11025	28	11484	23	12039	10, 98
11026	28	11486	23	12045	10, 98
11027	28	11487	23	12046	10, 98
11028	28	11488	23	12047	10, 98
11031	28	11491	26	12048	10, 98
11032	28	11492	26	12049	10, 98
11035	28	11497	26	12050	11, 99
11036	28	11498	26	12060	12, 100
11041	28	11499	26	12061	12, 100
11045	28	11500	26	12076	9, 97
11114	29	11512	48	12080	12, 100
11115	29	11514	48	12090	10, 98
11124	29	11515	48	12091	10, 98
11125	29	11517	48	12092	10, 98
11126	29	11518	48	12093	10, 98
11127	29	11590	48	12094	10, 98
11135	29	11591	48	12095	10, 98
11136	29	11612	49	12110	76
11185	29	11690	49	12130	76
11186	29	11803	68	12134	76
11194	29	11804	68	12135	76
11195	29	11805	69	12136	76
11196	29	11806	69	12140	76
11197	29	11807	69	12144	76
11301	47	11808	69	12145	76
11302	47	11809	69	12146	76
11303	47	11814	69	12147	76
11304	48	11816	70	12149	76
11305	48	11817	70	12170	75
11306	48	11818	70	12174	75

Code	Page #	Code	Page #	Code	Page #
12175	75	12670	74	13924	36
12176	75	12701	34	13925	36
12177	75	12702	34	13964	36
12202	46	12703	34	13965	36
12203	46	12704	34	13966	36
12213	46	12705	34	13974	36
12216	46	12706	34	13975	36
12222	46	12707	34	13976	36
12223	46	12708	34	13984	36
12262	44	12709	34	13985	36
12272	44	12716	34	13986	36
12273	44	12719	34	14010	16, 90
12274	44	12800	87	14015	16, 90
12281	44	12802	87	14021	16, 90
12282	44	12805	88	14031	16, 90
12283	44	12808	87	14041	16, 90
12284	44	12810	87	14050	89
12291	46	12818	88	14060	89
12292	46	12820	18	14070	89
12296	46	12830	17	14103	73
12297	46	12834	18	14106	73
12298	46	12835	18	14163	73
12344	38	12836	18	14164	73
12345	38	12841	32	14165	73
12346	38	12842	32	14166	73
12354	38	12844	32	14173	73
12355	38	12845	32	14174	73
12356	38	12847	32	14175	73
12374	38	12848	32	14176	73
12375	38	12991	79	14196	74
12376	38	12992	79	14204	39, 111
12384	38	12993	79	14205	39
12385	38	12994	79	14206	39
12386	38	12995	79	14222	77
12416	75	12996	79	14224	40
12419	75	12997	79	14225	40
12432	79	13660	93	14226	40
12433	79	13661	93	14227	40
12471	79	13662	93	14228	40
12472	79	13663	96	14229	40
12474	78	13664	96	14232	34
12527	37	13665	95	14233	34
12528	37	13703	75	14234	34
12529	37	13706	75	14235	34
12533	75	13840	72	14236	34
12577	37	13844	72	14237	34
12578	37	13845	72	14238	34
12579	37	13846	72	14239	34
12595	79	13860	72	14243	34
12596	79	13870	72	14246	34
12620	74	13874	72	14249	34
12630	74	13875	72	14250	32
12660	74	13876	72	14251	32

Code	Page #	Code	Page #	Code	Page #
14252	32	14454	77	16301	60
14254	32	14511	71	16305	60
14255	32	14512	71	16306	60
14257	32	14681	109	16310	60
14258	32	14682	109	16311	60
14261	33	14691	109	16312	60
14262	33	14692	109	16313	60
14263	33	14717	78	16315	60
14271	33	15480	68	16321	61
14272	33	15484	68	16322	61
14273	33	15521	66	16323	61
14274	33	15522	66	16325	60
14275	33	15523	66	16326	60
14276	33	15524	66	16327	60
14277	33	15531	66	16330	60
14278	33	15532	66	16331	60
14281	33	15533	66	16333	60
14282	33	15534	66	16334	60
14283	33	15606	67	16341	61
14291	33	15608	67	16342	61
14292	33	15623	67	16343	61
14293	33	15624	67	16345	61
14294	33	15633	67	16346	61
14295	33	15634	67	16347	61
14296	33	15695	67	16350	61
14297	33	15696	67	16351	61
14298	33	15697	67	16353	61
14307	109	15698	67	16354	61
14308	109	15821	66	16355	61
14404	73	15860	70	16471	43
14405	73	15869	71	16472	43
14411	73	15870	70	16481	43
14412	73	15871	71	16483	43
14413	74	15874	70	16484	43
14414	73	15875	70	16485	43
14415	73	15876	70	16486	43
14416	74	15879	71	16487	43
14419	74	15954	35	16488	43
14421	73	15955	35	16491	43
14422	73	15956	35	16493	43
14423	73	15974	35	16494	43
14424	73	15975	35	16705	54
14431	77	15976	35	16706	54
14432	77	15981	66	16715	54
14434	77	15982	35	16716	54
14435	77	15983	35	16725	55
14441	76	15984	35	16726	55
14442	76	15985	35	16727	55
14444	76	15986	35	16730	55
14445	76	15991	66	16731	55
14451	77	15992	66	16732	55
14452	77	15993	66	16733	55
14453	77	15994	66	16734	55

Code	Page #	Code	Page #	Code	Page #
16735	55	16867	65	18341	58
16737	55	16870	65	18342	58
16738	55	16871	65	18343	58
16739	55	16872	65	18344	58
16740	55	16901	21	18401	65
16741	55	16903	21	18402	65
16742	55	16911	21	18403	65
16748	55	16913	21	18404	65
16755	56	16914	21	18405	65
16756	56	16915	21	18500	22
16757	56	16921	21	18501	22
16760	56	16923	21	18510	22
16761	56	16924	21	18511	22
16762	56	16931	21	18512	22
16763	56	16933	21	18513	22
16764	56	16934	21	18520	22
16765	56	16941	22	18521	22
16766	56	16942	22	18522	22
16767	56	16951	22	18523	22
16768	56	16953	22	18524	22
16770	56	16954	22	18530	22
16771	56	16955	22	18531	22
16772	56	16961	22	18532	22
16773	56	16963	22	18600	24, 83
16778	56	16964	22	18601	24, 83
16801	63	16971	22	18602	24, 83
16805	63	16973	22	18603	24, 83
16810	63	16974	22	18610	24, 83
16821	64	18152	62	18611	24, 83
16822	64	18170	62	18612	24, 83
16823	64	18171	62	18613	24, 83
16825	64	18172	62	18620	25, 84
16826	64	18173	62	18621	25, 84
16827	64	18180	62	18622	25, 84
16830	64	18181	62	18630	24, 83
16831	64	18182	62	18631	24, 83
16832	64	18301	58	18632	24, 83
16833	64	18302	58	18640	24, 83
16835	64	18305	58	18641	24, 83
16836	64	18312	58	18642	24, 83
16837	64	18313	58	18650	45, 84
16840	64	18314	58	18651	45, 84
16841	64	18315	58	18652	45, 85
16842	64	18320	58	18653	45, 85
16855	65	18321	58	18655	45, 84
16856	65	18322	58	18656	45, 84
16857	65	18323	58	18657	45, 85
16860	65	18324	58	18658	45, 85
16861	65	18330	58	18660	57, 85
16862	65	18331	58	18661	57, 85
16863	65	18332	58	18665	57, 85
16865	65	18333	58	18666	57, 85
16866	65	18340	58	18670	57, 85

Index

<u>Code</u>	<u>Page #</u>	<u>Code</u>	<u>Page #</u>	<u>Code</u>	<u>Page #</u>
18671 57, 85	32155 13		
18675 57, 86	32205 13		
18676 57, 86	67589 13		
18680 86	95050 13		
18681 86	110072 93		
18690 86	120105	12, 100		
18691 86	120115	12, 100		
18700 105	120205	12, 101		
18705 105	120215	13, 101		
18710 108	120307	.. 11, 99		
18711 108	120427	.. 11, 99		
18712 108	120700 9, 97		
18713 108	120721 9, 97		
18720 106	120740 9, 97		
18721 106				
18722 107				
18723 107				
18724 107				
18725 107				
18730 107				
18731 107				
18735 106				
18747 106				
19001 87				
19004 67, 86				
19010 87				
19011 87				
19020 87				
19021 87				
19030 87				
19040 30				
19041 30				
19042 30				
19050 30				
19051 30				
19052 30				
19075 13				
21001 94				
21002 94				
21003 94				
21004 94				
21005 96				
21006 97				
21007 95				
21008 91				
21009 92				
21010 92				
21011 94				
21012 92				
21013 92				
21365 15				
26555 11, 99				
32105 13				