



Contenance, Ostomy,
Wound and Skin Care

Product Guide

Effective October 2021

Tabitha

Making life easier_



Making life easier—

About Tabitha

When Tabitha became paralyzed from the waist down in 2014, her life changed dramatically. From running around in high heels, going to work and out and about with friends, she was suddenly confined to a wheelchair, dependent on people around her for help with everything.

A fiercely independent person by nature and used to living a very active life, Tabitha was terrified by the prospect of losing her independence. “How am I going to get around? How am I even going to live?” This was the beginning of a massive learning curve for Tabitha. But as she went through rehabilitation, she realized that even though you don’t walk around on two legs, you can still do things. You just have to find new ways of doing them. “So I gave myself a stern talking to”, she recalls. **“I could either sit and cry about what I wasn’t able to do. Or I could get on with life.”** Tabitha chose the latter – and intermittent self-catheterization (ISC) helped to make that possible.

She’d been discharged from the hospital with an indwelling catheter. But being attached to tubes and having a urine bag strapped to her body felt invasive, making her embarrassed and self-conscious. It was yet another reminder of what she couldn’t do. The transition to a compact catheter was liberating. The fact that they are so discrete and easy to use, gave her the freedom and confidence to go out more. She regained control of her life in a way she hadn’t seen coming.

From having no interest in competitive sports when she could walk, Tabitha started practicing wheelchair tennis. She now travels around the world competing in international tournaments, with several trophies and medals to show for her talent. This achievement is proof of what Tabitha now knows to be true: **“that you can do what you want to do. You build your life.”** That’s the spirit she hopes to inspire in others, by sharing her experience and talking openly about catheterization.

Tabitha has received compensation from Coloplast to provide this information. Each person’s situation is unique so your experience may not be the same. Talk to your healthcare provider about whether this product is right for you.

Table of Contents

Continence Care.....	7-30
Straight Tip Intermittent Catheters	11-14
Coudé Tip Intermittent Catheters.....	15-16
Intermittent Catheters Closed Systems / Sets	17-20
Transanal Irrigation	21-22
Male External Catheters.....	23-25
Leg and Drain Bags	26-27
Indwelling Catheters	28-29
Accessories.....	30
Wound Care	31-44
Foam Dressings.....	35-37
Alginate Dressings.....	38-39
Hydrocolloid Dressings.....	40-41
Hydrophilic Dressing / Wound Cleansers / Hydrogels.....	42
Contact Layer	43
Super Absorbents.....	44
Skin Care	45-58
Cleansers.....	49-51
Moisturizers.....	52
Skin Protectants.....	53-54
Antifungals.....	55
Skin Fold Management.....	56
Hand Hygiene	57
Accessories.....	58
Ostomy Care	59-172
Accessories.....	63-72
Irrigation	72-74
2-Piece Click.....	75-96
2-Piece Flex.....	97-106
1-Piece.....	107-136
Hospital Assortment / Wound Management.....	137-158
Pediatric.....	159-170
Discontinued Items.....	171-172

Reference Tools 173-188

Education

Contenance Care 175

Ostomy Care 176-177

Contenance and Ostomy Education 178

Wound and Skin Care Education 179

Medicare Utilization

Contenance..... 180

Ostomy..... 181-182

Wound Care..... 183

HCPCS Codes 184-188

Index..... 189

Welcome to Coloplast

Coloplast is a leader in the market within intimate healthcare. For more than sixty years, we've worked closely with healthcare professionals and users to develop products and services that make life easier for people with intimate healthcare needs.

Our Mission

***Making life easier for
people with intimate
healthcare needs***

Making
life *easier*

Eva



Making life easier_

SpeediCath® catheters are prescribed for use for patients who require bladder drainage due to urine retention or post void residual volume (PVR). Before use, carefully read all of the instructions. Call your doctor if you think you have a UTI or can't pass the catheter into the bladder. For more information regarding risks, potential complications and product support, call Coloplast Corp. at 1-866-226-6362 and/ or consult the company website at www.coloplast.us. Caution: Federal law restricts this device to sale by or on the order of a physician.

Self-Cath® Plus Intermittent Urinary Catheters are intended for use in male or female patients requiring bladder drainage as determined by their physician. This device is indicated for those individuals unable to promote a natural urine flow or for those who have a significant volume of residual urine following a natural bladder-voiding episode. The device can be used by the patient once appropriate training has taken place or by a trained healthcare professional. Single-use only. See the device manual for detailed information regarding the procedure, contraindications, warnings, precautions, and potential complications/adverse events. For further information, call Coloplast Corp. at 1-866-226-6362 and/or consult the company website at www.coloplast.us. Caution: Federal law restricts this device to sale by or on the order of a physician.

Self-Cath® Intermittent Catheters are intended for use in male, female, and pediatric patients (neonates, infants, children, adolescents, and transitional adolescents) requiring bladder drainage as determined by their physician. This device is indicated for those individuals unable to promote a natural urine flow or for those individuals who have a significant volume of residual urine following a natural bladder-voiding episode. Single-use only. See the device manual for detailed information regarding the procedure, contraindications, warnings, precautions, and potential complications/adverse events. For further information, call Coloplast Corp. at 1-866-226-6362 and/ or consult the company website at www.coloplast.us. Caution: Federal law restricts this device to sale by or on the order of a physician.

SureCath® Catheters are indicated for use by patients for intermittent catheterization for the purpose of bladder drainage. See the device manual for detailed information regarding the instructions for use, contraindications, warnings, precautions, and potential complications/adverse events. For further information, call Coloplast Corp. at 1-866-226-6362 and/or consult the company website at www.coloplast.us. Caution: Federal law restricts this device to sale by or on the order of a physician.

Folysil® Indwelling Catheters empty urine from the bladder through the urethra via a small cut through the abdomen above the pubic bone (suprapubic) during a surgical procedure. The choice of the size, balloon volume, and catheter material is the responsibility of your healthcare professional based on your condition and the length of time to be used. Coloplast's indwelling catheters are available by prescription only. Coloplast's indwelling catheters are for single-use only; discard it after use. If you experience symptoms of a urinary tract infection, are unable to pass urine, or the balloon will not deflate, contact your healthcare professional. The risk information provided here is not comprehensive. To learn more, talk to your healthcare professional to understand the risks and benefits to determine if Folysil is right for you. Although many patients benefit from the use of this device, results may vary. For further questions, call Coloplast Corp. at 1-866-226-6362 and/or consult the company website at www.coloplast.us.

Indication: The **Peristeen® (Plus) Transanal Irrigation System** is intended to instill water into the colon through a rectal catheter, which incorporates an inflatable balloon, inserted into the rectum to promote evacuation of the contents of the lower colon. The Peristeen Transanal Irrigation System is indicated for use by children (2 years - <12 years old), adolescent (12 years - < 18 years old), transitional adolescent (18 - <21 years old) and adult patients with neurogenic bowel dysfunction who suffer from fecal incontinence, chronic constipation, and/or time-consuming bowel management procedures.

Contraindications: Peristeen (Plus) Transanal Irrigation must not be used in the following situations: known anal or colorectal stenosis, colorectal cancer, radiotherapy to the pelvis, and recent abdomino-perineal surgery, active inflammatory bowel disease, acute diverticulitis, severe diverticulitis, previous diverticulitis and diverticular abscess, chronic symptomatic diverticular disease, within 3 months of abdominal, anal or colorectal surgery, within 4 weeks of endoscopic polypectomy, recent colonic biopsy, recent endoscopic mucosal resection (EMR) and recent endoscopic sub-mucosal dissection (ESD), severe autonomic dysreflexia, ischemic colitis, during spinal cord shock phase, complex diverticular disease, in patients who are pregnant and have not used the system before*. Since the list is not exhaustive, the healthcare professional should always consider individual patient factors as well.

*If the patient is pregnant and has never used transanal irrigation before, they should not start the irrigation procedure during pregnancy.

Warning: Peristeen (Plus) Transanal Irrigation procedure should always be carried out with caution. Bowel perforation is an extremely rare, but serious and potentially lethal complication to transanal irrigation and will require immediate admission to a hospital, often requiring surgery. See the device manual for complete user instructions, contraindications, warnings, precautions, and potential complications/adverse events. For further information, call Coloplast Corp at 1-855-605-7594 or consult the company website at www.coloplast.us. Caution: Federal law restricts this device to sale by or on the order of a physician.

Ensure your patients get the product you recommend. Ask for SpeediCath® and Self-Cath® catheters by brand name.

Intermittent catheter users who have insurance can typically pay the same price when buying any catheter within the same category: straight tip,¹ coude tip,² or catheter sets.³



* Dispense as written; exact language to specify brand dispensed varies by state. Contact your Coloplast Sales Rep or National Association of Boards of Pharmacy for more information.

Advise your patients to find a supplier who accepts assignment for the intermittent catheter. If your patients need assistance in finding a supplier, they may contact Coloplast Consumer Care at 1-866-226-6362.

Steps to consider for patients to get catheter they want:

1. Include the following information on each prescription:

Patient Name _____
 Address _____

Type of catheter (straight, coude, kit / closed system)

Rx: SpeediCath Compact Set 28422
 Dx: Urine retention
 5 catheters / day x 30 days / month = 150 catheters / month

Brand of catheter

Frequency of catheterizations / day and total catheters needed / month

Do not refill _____
 Refill _____ Times NPI Number _____
 Date _____ Print Last Name _____

2. Enroll your patients in Coloplast® Care

- Call** - contact Coloplast Care at 1-866-226-6362
- Email** - send request to care-us@coloplast.com⁴
- Fax** - send patient enrollment form to 1-855-676-2594
- Website** - www.ColoplastCareEnroll.com

Information from Coloplast® Care is for educational purposes only. It is not intended to substitute for professional medical advice and should not be interpreted to contain treatment recommendations.

1. HCPCS A4351 (straight)**
2. HCPCS A4352 (coude)**
3. HCPCS A4353 (catheter set)**
4. We recommend encrypting emails and forms if sending over email.

⁴For important information regarding the use of HCPCS codes listed and code definitions see page 180 | For important safety information see page 8



All SpeediCath® catheters come with a unique hydrophilic coating, featuring our **Triple Action Coating Technology**. It consists of three strongly bonded elements, designed to protect the urethra in three ways:



Stays bonded

It stays bonded during insertion and withdrawal due to a proprietary base-coat that is covalently bonded to the catheter surface for an even coating without dry spots.



Stays smooth

It stays smooth because the polymer top-coat enables excellent water absorption for low friction, smooth insertion and reduced discomfort compared to uncoated catheters.



Stays hydrated

It stays hydrated with a special water-based solution that swells the coating and makes it instantly ready to use.

SpeediCath® Standard

The world's first ready-to-use hydrophilic-coated intermittent catheter

Easy to carry, store and use. The catheter is pre-hydrated and ready-to-use, eliminating wait time. The hydrophilic coating is smooth and uniform, designed to minimize urethral friction compared to uncoated catheters. The eyelets are polished for comfort. The uncoated Grip Zone™ area is created to offer control in handling. **This catheter is not made with natural rubber latex, DEHP or other phthalates.**

SpeediCath® Standard

SpeediCath catheter only. No insertion supplies.



Straight Tip

Code	Size	Units	Description	HCPCS*
28408	● 8 FR	30	14" length	A4351
28410	● 10 FR	30	14" length	A4351
28412	○ 12 FR	30	14" length	A4351
28414	● 14 FR	30	14" length	A4351
28416	● 16 FR	30	14" length	A4351
28418	● 18 FR	30	14" length	A4351

Female Straight Tip

28506	● 6 FR	30	6" length	A4351
28508	● 8 FR	30	6" length	A4351
28510	● 10 FR	30	6" length	A4351
28512	○ 12 FR	30	6" length	A4351
28514	● 14 FR	30	6" length	A4351
28516	● 16 FR	30	6" length	A4351

Pediatric Straight Tip

28706	● 6 FR	30	6" length	A4351
28708	● 8 FR	30	6" length	A4351
28710	● 10 FR	30	6" length	A4351

Boy Straight Tip

28606	● 6 FR	30	10" length	A4351
28608	● 8 FR	30	10" length	A4351
28610	● 10 FR	30	10" length	A4351
28612	○ 12 FR	30	10" length	A4351

*For important information regarding the use of HCPCS codes listed and code definitions see page 180 | For important safety information see page 8

SpeediCath® Soft

SpeediCath® Soft is designed to make catheterization simple and easy with fast preparation and a confident grip.

- Pre-lubricated intermittent catheter that is instantly ready-to-use with a unique hydrophilic coating, eliminating the need for additional lubricants or waiting for a catheter to hydrate.
- Insertion grip made from a soft, flexible material designed to ensure a confident and hygienic grip of the catheter, without touching the tube or tip.
- Soft catheter that is easy to use the first time and every time
- **This catheter is not made with natural rubber latex, DEHP or other phthalates.**



SpeediCath® Soft

SpeediCath Soft straight tip intermittent catheter.

Code	Size	Units	Description	HCPCS*
27010	● 10 FR	30	13" length	A4351
27012	○ 12 FR	30	13" length	A4351
27014	● 14 FR	30	13" length	A4351
27016	● 16 FR	30	13" length	A4351

SpeediCath® Compact Female and Compact Plus

The most compact, ready-to-use hydrophilic-coated intermittent catheters for women.

SpeediCath Compact Female and Compact Plus offer the benefits of our ready-to-use hydrophilic-coated intermittent catheters in a small, discreet package. Designed for women, the length of SpeediCath Compact Female matches the average length of the female urethra, and when unopened it is no longer than a tube of lipstick. **This catheter is not made with natural rubber latex, DEHP or other phthalates.**

SpeediCath® Compact Female

SpeediCath Compact straight tip intermittent catheter.

Code	Size	Units	Description	HCPCS*
28578	● 8 FR	30	2.75" length	A4351
28580	● 10 FR	30	2.75" length	A4351
28582	○ 12 FR	30	2.75" length	A4351
28584	● 14 FR	30	2.75" length	A4351

SpeediCath® Compact Plus

SpeediCath Compact Female with additional catheter length for women who want discretion but require a longer catheter.

Code	Size	Units	Description	HCPCS*
28810	● 10 FR	30	3.5" length	A4351
28812	○ 12 FR	30	3.5" length	A4351
28814	● 14 FR	30	3.5" length	A4351

Self-Cath® Plus

Self-Cath® Plus intermittent catheters have fire-polished eyelets designed for comfort and a siliconized surface created for smooth insertion. The catheter has a hydrophilic coating that activates immediately upon exposure to water. **This catheter is not made with natural rubber latex.**

Self-Cath® Plus

Hydrophilic-coated intermittent catheter.



Straight Tip

Code	Size	Units	Description	HCPCS*
4408	● 8 FR	30	16" length, funnel end	A4351
4410	● 10 FR	30	16" length, funnel end	A4351
4412	○ 12 FR	30	16" length, funnel end	A4351
4414	● 14 FR	30	16" length, funnel end	A4351
4416	● 16 FR	30	16" length, funnel end	A4351
4418	● 18 FR	30	16" length, funnel end	A4351

Female Straight Tip

4210	● 10 FR	30	6" length, funnel end	A4351
4212	○ 12 FR	30	6" length, funnel end	A4351
4214	● 14 FR	30	6" length, funnel end	A4351

Pediatric Straight Tip

4306	● 6 FR	30	10" length, funnel end	A4351
4308	● 8 FR	30	10" length, luer end	A4351
4310	● 10 FR	30	10" length, luer end	A4351

Self-Cath® Plus Soft

The Self-Cath® Plus Soft is the same as Self-Cath® Plus but designed for individuals that prefer a softer catheter.

Straight Tip

Code	Size	Units	Description	HCPCS*
4110	● 10 FR	30	16" length, funnel end	A4351
4112	○ 12 FR	30	16" length, funnel end	A4351
4114	● 14 FR	30	16" length, funnel end	A4351
4116	● 16 FR	30	16" length, funnel end	A4351

*For important information regarding the use of HCPCS codes listed and code definitions see page 180 | For important safety information see page 8

Self-Cath®

Self-Cath intermittent catheters have fire-polished eyelets designed for comfort, a siliconized surface created for smooth insertion, and are available in a variety of options. **This catheter is not made with natural rubber latex.**



Self-Cath®

Straight Tip

Code	Size	Units	Description	HCPCS*
408	● 8 FR	30	16" length, funnel end	A4351
410	● 10 FR	30	16" length, funnel end	A4351
412	○ 12 FR	30	16" length, funnel end	A4351
414	● 14 FR	30	16" length, funnel end	A4351
416	● 16 FR	30	16" length, funnel end	A4351
418	● 18 FR	30	16" length, funnel end	A4351
450	● 14 FR	50	16" length, funnel end, curved packaging	A4351
460	○ 12 FR	50	16" length, funnel end, curved packaging	A4351

Female Straight Tip

Code	Size	Units	Description	HCPCS*
208	● 8 FR	30	6" length, funnel end	A4351
210	● 10 FR	30	6" length, funnel end	A4351
212	○ 12 FR	30	6" length, funnel end	A4351
214	● 14 FR	30	6" length, funnel end	A4351
240	● 14 FR	30	6" length, luer end	A4351

Pediatric Straight Tip

Code	Size	Units	Description	HCPCS*
305	○ 5 FR	30	10" length, funnel end	A4351
306	● 6 FR	30	10" length, funnel end	A4351
308	● 8 FR	30	10" length, luer end	A4351
310	● 10 FR	30	10" length, luer end	A4351

Self-Cath® Soft

The Self-Cath® Soft is the same as Self-Cath® but designed for individuals that prefer a softer catheter.



Straight Tip

Code	Size	Units	Description	HCPCS*
110	● 10 FR	30	16" length, funnel end	A4351
112	○ 12 FR	30	16" length, funnel end	A4351
114	● 14 FR	30	16" length, funnel end	A4351
116	● 16 FR	30	16" length, funnel end	A4351

SpeediCath® Flex Coudé Pro

Every detail makes it easy. Hydrophilic-coated intermittent catheter that is instantly ready-to-use. Flexible, soft catheter that is easy to handle with a curved tip which can bend to follow the curves of the urethra. Hygienic with a protective dry sleeve and soft-squeeze grip. Discreet standard length catheter packaged in two compact sizes. **This catheter is not made with natural rubber latex, DEHP or other phthalates.**



Code	Size	Units	Description	HCPCS*
20010	● 10 FR	30	13" length	A4352
20012	○ 12 FR	30	13" length	A4352
20014	● 14 FR	30	13" length	A4352
20016	● 16 FR	30	13" length	A4352
20021	● 10 FR	30	13" length <i>Pocket</i>	A4352
20022	○ 12 FR	30	13" length <i>Pocket</i>	A4352
20024	● 14 FR	30	13" length <i>Pocket</i>	A4352

SpeediCath® Standard Coudé Tip

The world's first ready-to-use hydrophilic-coated intermittent catheter

This tapered coudé tip catheter is easy to carry, store and use. It is pre-hydrated and ready-to-use, eliminating wait time. The hydrophilic coating is smooth and uniform, designed to minimize urethral friction compared to uncoated catheters. The eyelets are polished and designed for comfort. The uncoated Grip Zone™ area is designed to offer control in handling. **This catheter is not made with natural rubber latex, DEHP or other phthalates.**



SpeediCath® Standard Tapered Coudé Tip

Code	Size	Units	Description	HCPCS*
28490	● 10 FR	30	14" length	A4352
28492	○ 12 FR	30	14" length	A4352
28494	● 14 FR	30	14" length	A4352
28496	● 16 FR	30	14" length	A4352

*For important information regarding the use of HCPCS codes listed and code definitions see page 180 | For important safety information see page 8

Self-Cath® Plus Coudé Tip

This coudé tip intermittent catheter has fire-polished eyelets designed for comfort and a siliconized surface created for smooth insertion. The catheter has a hydrophilic coating that activates immediately upon exposure to water. The blue Guide Stripe® along the top of the catheter allows monitoring of the coudé tip throughout the catheterization process. **This catheter is not made with natural rubber latex.**

**Self-Cath® Plus Tapered Coudé Tip with Guide Stripe®**

Code	Size	Units	Description	HCPCS*
4608	● 8 FR	30	16" length, funnel end	A4352
4610	● 10 FR	30	16" length, funnel end	A4352
4612	○ 12 FR	30	16" length, funnel end	A4352
4614	● 14 FR	30	16" length, funnel end	A4352

**Self-Cath® Plus Olive Coudé Tip with Guide Stripe®**

Code	Size	Units	Description	HCPCS*
4808	● 8 FR	30	16" length, funnel end	A4352
4810	● 10 FR	30	16" length, funnel end	A4352
4812	○ 12 FR	30	16" length, funnel end	A4352
4814	● 14 FR	30	16" length, funnel end	A4352
4816	● 16 FR	30	16" length, funnel end	A4352

Self-Cath® Coudé Tip

Self-Cath coudé tip intermittent catheters have fire-polished eyelets designed for comfort, a siliconized surface created for smooth insertion. The blue Guide Stripe® along the top of the catheter allows monitoring of the coudé tip throughout the catheterization process. They are available in a variety of options. **This catheter is not made with natural rubber latex.**

Self-Cath® Olive Coudé Tip with Guide Stripe®

Code	Size	Units	Description	HCPCS*
806	● 6 FR	30	16" length, funnel end	A4352
808	● 8 FR	30	16" length, funnel end	A4352
810	● 10 FR	30	16" length, funnel end	A4352
812	○ 12 FR	30	16" length, funnel end	A4352
814	● 14 FR	30	16" length, funnel end	A4352
816	● 16 FR	30	16" length, funnel end	A4352
818	● 18 FR	30	16" length, funnel end	A4352

**Self-Cath® Tapered Coudé Tip with Guide Stripe®**

Code	Size	Units	Description	HCPCS*
608	● 8 FR	30	16" length, funnel end	A4352
610	● 10 FR	30	16" length, funnel end	A4352
612	○ 12 FR	30	16" length, funnel end	A4352
614	● 14 FR	30	16" length, funnel end	A4352

SpeediCath® Compact Male with SpeediBag®

The most compact catheter for men.

Hydrophilic-coated catheter that is pre-hydrated and ready-to use. This catheter is hygienic with a no-touch insertion. Its telescoping mechanism allows expansion from its small modest tube into a full length intermittent catheter for males. A 750 mL SpeediBag is included. **This catheter is not made with natural rubber latex, DEHP or other phthalates.**

SpeediCath® Compact Male with SpeediBag®



Code	Size	Units	Description	HCPCS*
28702	○ 12/18 FR	20	13.2" length includes funnel	A4353

SpeediCath® Compact Set

The most compact set solution. An instantly ready-to-use hydrophilic-coated intermittent catheter in a compact design with an integrated bag available for both males and females. Provides ease of use and discretion and thanks to its award winning simple and compact design it is small enough to keep in a handbag or pocket. **This catheter is not made with natural rubber latex, DEHP or other phthalates.**



reddot design award
winner 2013

SpeediCath® Compact Set



Female

Code	Size	Units	Description	HCPCS*
28520	● 10 FR	20	3.5" length	A4353
28522	○ 12 FR	20	3.5" length	A4353
28524	● 14 FR	20	3.5" length	A4353

Male

28422	○ 12/18 FR	20	13.2" length includes funnel	A4353
-------	------------	----	------------------------------	-------

*For important information regarding the use of HCPCS codes listed and code definitions see page 180 | For important safety information see page 8

SpeediCath® Standard with Insertion Supplies

The world's first ready-to-use hydrophilic-coated intermittent catheter.

Easy to carry, store and use. The catheter is pre-hydrated and ready-to-use, eliminating wait time. The hydrophilic coating is smooth and uniform, designed to minimize urethral friction. The eyelets are polished for comfort. The uncoated Grip Zone™ area is designed to offer control in handling. **This catheter is not made with natural rubber latex, DEHP or other phthalates.**

SpeediCath® Standard with Insertion Supplies



SpeediCath Standard intermittent catheter with one pair of non-latex exam gloves, one waterproof drape, two iodine prep pads and urine collection bag.

Straight Tip

Code	Size	Units	Description	HCPCS*
28481	● 8 FR	100	14" length	A4353
28482	● 10 FR	100	14" length	A4353
28483	○ 12 FR	100	14" length	A4353
28484	● 14 FR	100	14" length	A4353
28485	● 16 FR	100	14" length	A4353
28486	● 18 FR	100	14" length	A4353

Female Straight Tip

28499	● 6 FR	100	6" length	A4353
28500	● 8 FR	100	6" length	A4353
28501	● 10 FR	100	6" length	A4353
28502	○ 12 FR	100	6" length	A4353
28503	● 14 FR	100	6" length	A4353
28504	● 16 FR	100	6" length	A4353

SureCath® Set

For intermittent self-catheterization.

Hydrophilic-coated intermittent catheter set with sterile solution enclosed in a single-use, collection bag. The catheter is available with and without insertion supplies. **This catheter is not made with natural rubber latex.**

SureCath® Set with Insertion Supplies

SureCath Set intermittent catheter and collection bag, one pair of exam gloves, one CSR wrap and two iodine prep pads.

Straight Tip 700mL

Code	Size	Units	Description	HCPCS*
28022S	● 8 FR	100	14" length	A4353
28023S	● 10 FR	100	14" length	A4353
28024S	○ 12 FR	100	14" length	A4353
28032	● 14 FR	20	14" length	A4353
28032S	● 14 FR	100	14" length	A4353
28033S	● 16 FR	100	14" length	A4353

Straight Tip 1200 mL

Code	Size	Units	Description	HCPCS*
28036	● 14 FR	20	14" length	A4353
28037	● 16 FR	20	14" length	A4353
28036S	● 14 FR	100	14" length	A4353
28037S	● 16 FR	100	14" length	A4353

SureCath® Set without Insertion Supplies



SureCath hydrophilic-coated intermittent catheter and collection bag only. No additional insertion supplies.

Straight Tip 700 mL

Code	Size	Units	Description	HCPCS*
28005S	● 8 FR	100	14" length	A4353
28006S	● 10 FR	100	14" length	A4353
28007S	○ 12 FR	100	14" length	A4353
28008S	● 14 FR	100	14" length	A4353
28009	● 16 FR	20	14" length	A4353
28009S	● 16 FR	100	14" length	A4353

Straight Tip 1200 mL

28026S	● 10 FR	100	14" length	A4353
28027	○ 12 FR	20	14" length	A4353
28027S	○ 12 FR	100	14" length	A4353
28028	● 14 FR	20	14" length	A4353
28028S	● 14 FR	100	14" length	A4353
28029S	● 16 FR	100	14" length	A4353

 **WARNING:** This product can expose you to chemicals including Di(2-ethylhexyl) phthalate (DEHP) and Methyl isobutyl ketone (MIBK), which are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov 

*For important information regarding the use of HCPCS codes listed and code definitions see page 180 | For important safety information see page 8

Self-Cath® Closed System

The Self-Cath Closed System has a Self-Cath pre-lubricated catheter enclosed in a single-use, collection bag. The catheter is available with and without insertion supplies. **This catheter is not made with natural rubber latex.**

Self-Cath® Closed Systems with Insertion Supplies

Self-Cath pre-lubricated intermittent catheter, gloves, 1100 mL seamless collection bag, povidone-iodine swabsticks, BZK wipe (towelette), underpad and instructions for use.



Straight Tip

Code	Size	Units	Description	HCPCS*
1008	● 8 FR	50	16" length	A4353
1010	● 10 FR	50	16" length	A4353
1012	○ 12 FR	50	16" length	A4353
1014	● 14 FR	50	16" length	A4353
1016	● 16 FR	50	16" length	A4353

Female Straight Tip

Code	Size	Units	Description	HCPCS*
3214	● 14 FR	50	6" length	A4353

Soft Straight Tip

3114	● 14 FR	50	16" length	A4353
------	---------	----	------------	-------

Olive Coudé Tip with Guide Stripe®

C3814	● 14 FR	50	16" length	A4353
C3816	● 16 FR	50	16" length	A4353

Tapered Coudé Tip with Guide Stripe®

3614	● 14 FR	50	16" length	A4353
------	---------	----	------------	-------

Self-Cath® Closed System Single Units

Self-Cath pre-lubricated intermittent catheter and 1100 mL collection bag only. No additional insertion supplies.



Straight Tip

Code	Size	Units	Description	HCPCS*
1108	● 8 FR	50	16" length	A4353
C1110	● 10 FR	50	16" length	A4353
1112	○ 12 FR	50	16" length	A4353
1114	● 14 FR	50	16" length	A4353
1116	● 16 FR	50	16" length	A4353

Female Straight Tip

2214	● 14 FR	50	6" length	A4353
------	---------	----	-----------	-------

Soft Straight Tip

2114	● 14 FR	50	16" length	A4353
------	---------	----	------------	-------

Olive Coudé Tip with Guide Stripe®

2814	● 14 FR	50	16" length	A4353
2816	● 16 FR	50	16" length	A4353

Tapered Coudé Tip with Guide Stripe®

2614	● 14 FR	50	16" length	A4353
------	---------	----	------------	-------

Peristeen® Plus Transanal Irrigation System

A better path to bowel control

Transanal Irrigation



Transanal irrigation with Peristeen is a well-documented bowel management technique that is effective at preventing fecal incontinence and chronic constipation for patients with neurogenic bowel dysfunction.

To learn more visit www.coloplastcare.us/bowel.

Introducing
**Peristeen
Plus**



Peristeen® Transanal Irrigation

Peristeen is a FDA Class II cleared transanal irrigation system intended to instill water into the colon through a rectal catheter to promote evacuation of the contents of the lower colon. It is indicated for use by children (2 and up), adolescents, transitional adolescents and adults with neurogenic bowel dysfunction.

Peristeen is specially designed for portability and ease of use. It offers independence and a sense of security to the user. In clinical trials with patients with neurogenic bowel dysfunction, Peristeen system reduced fecal incontinence and constipation and improved quality of life.^{1,2}



Peristeen® Transanal Irrigation Available through June 2022

Code

29129	System without catheters	1 control unit / 1 water bag (with lid) / 2 straps
29122	Accessory unit Regular	15 rectal balloon catheters Regular / 1 water bag
29127	Accessory unit Small	15 rectal balloon catheters Small / 1 water bag

HCPCS A4459

Manual pump enema system, includes balloon, catheter and all accessories, reusable, any type

A4459 would include 29129 (1 each) and either 29122 or 29127.

Available via prior authorization process for Medicaid and private insurance.

HCPCS A4453 (effective 10/1/2021)

Rectal catheter for use with the manual pump operated enema system, replacement only

Peristeen® Plus Transanal Irrigation Coming January 2022!



Code

29152	Peristeen Plus system without catheters	1 control unit / 1 water bag (with lid) / 2 straps / 1 tube
29142	Peristeen Plus Accessory unit Regular	15 rectal balloon catheters* Regular / 1 water bag
29149	Peristeen Plus Accessory unit Small	15 rectal balloon catheters* Small / 1 water bag

***Catheters and tubing from the original Peristeen system are not compatible with Peristeen Plus**

- Christensen P et al. A randomized, controlled trial of transanal irrigation versus conservative bowel management in spinal cord injured patients. Supported by Coloplast. *Gastroenterology* 2006; 131:738-747, n=87
- Mosiello G., Marshall D., Rolle U., Crétolle C., Santacruz B.G., Frischer J. and Benninga M. Consensus Review of Best Practice of Transanal Irrigation in Children. Supported by Coloplast. *JPGN*; 2017; 64(3):343-352

⚠ WARNING: This product can expose you to chemicals including Methyl isobutyl ketone (MIBK, which is known to the State of California to cause cancer and birth defects or other reproductive harm, and Di-n-butyl phthalate (DBP). For more information go to www.P65Warnings.ca.gov

This warning does not apply to catalog number 29129 but does apply to catalog numbers 29122 and 29127

Male External Catheters (Condom Catheters)

Coloplast offers a complete line of male external catheters. A variety of styles and sheath lengths are designed to suit end-user needs. A complete incontinence system can be achieved by combining a Coloplast male external catheter and a Coloplast leg bag.

Conveen® Optima

A self-adhering, 100% silicone clear male external catheter packaged in a discreet, recyclable smart pack.

- Easy on - the double strip grip provides an easy and accurate application
- Triple action system with anti-kink bellow, sure grip ribs and push ring provide increased security
- The PSX silicone is soft, flexible and breathable
- The GX adhesive is reliable and skin friendly
- **Not made with natural rubber latex**



Standard Length

Code	Size	Units	HCPCS*
22025	25 mm	30	A4349
22028	28 mm	30	A4349
22030	30 mm	30	A4349
22035	35 mm	30	A4349
22040	40 mm	30	A4349

Sport Length

22121	21 mm	30	A4349
22125	25 mm	30	A4349
22130	30 mm	30	A4349
22135	35 mm	30	A4349

Conveen® Security+ Self-Sealing Male External Catheter

A self-adhering male external catheter features a placement ring to assist with individual preference and an anti-kink design for security. Easy to apply with the applicator and roll-out strip. **Not made with natural rubber latex.**



Code	Size	Units	HCPCS*
5221	21 mm	35	A4349
5225	25 mm	35	A4349
5230	30 mm	35	A4349
5235	35 mm	35	A4349
5240	40 mm	35	A4349
5225H	25 mm	100	A4349
5230H	30 mm	100	A4349
5235H	35 mm	100	A4349
5240H	40 mm	100	A4349

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Freedom Clear Advantage® with Aloe

The Clear Advantage male external catheter is a self-adhering, silicone catheter. It features aloe in the adhesive and an anti-kink reservoir. **Not made with natural rubber latex.**



Code	Size	Units	HCPCS*
6100	23 mm	100	A4349
6200	28 mm	100	A4349
6300	31 mm	100	A4349
6400	35 mm	100	A4349
6500	40 mm	100	A4349
6130	23 mm	30	A4349
6230	28 mm	30	A4349
6330	31 mm	30	A4349
6430	35 mm	30	A4349
6530	40 mm	30	A4349

Freedom Clear®

The Freedom Clear male external catheter offers a silicone sheath with an anti-kink reservoir. **Not made with natural rubber latex.**



Code	Size	Units	HCPCS*
5100	23 mm	100	A4349
5200	28 mm	100	A4349
5300	31 mm	100	A4349
5400	35 mm	100	A4349
5500	40 mm	100	A4349

Freedom Clear® LS

The Freedom Clear LS (Long Seal) male external catheter is designed with a longer-length adhesive seal on a silicone sheath featuring an anti-kink reservoir. **Not made with natural rubber latex.**



Code	Size	Units	HCPCS*
5190	23 mm	100	A4349
5290	28 mm	100	A4349
5390	31 mm	100	A4349
5490	35 mm	100	A4349
5590	40 mm	100	A4349

Freedom Clear® SS

The Freedom Clear SS (Sport Sheath) male external catheter offers a shortened silicone sheath for less roll-up of the shaft of the catheter. This product features an anti-kink reservoir. **Not made with natural rubber latex.**



Code	Size	Units	HCPCS*
5110	23 mm	100	A4349
5210	28 mm	100	A4349
5310	31 mm	100	A4349
5410	35 mm	100	A4349

Freedom Cath®

Freedom Cath male external catheter is a self-adhering, one-piece latex catheter. Combines a latex sheath with a wide, water resistant adhesive seal featuring an anti-kink reservoir.

**Standard Length**

Code	Size	Units	HCPCS*
8000	23 mm	100	A4349
8200	28 mm	100	A4349
8205	31 mm	100	A4349
8400	35 mm	100	A4349
8030	23 mm	30	A4349
8230	28 mm	30	A4349
8235	31 mm	30	A4349
8430	35 mm	30	A4349

Freedom Active Cath®

Freedom Active Cath male external catheter is a self-adhering, extended wear, one-piece latex catheter. Combines a shortened latex sheath with a wide, water resistant adhesive seal. Features an anti-kink reservoir.

**Sport Length**

Code	Size	Units	HCPCS*
8100	23 mm	100	A4349
8300	28 mm	100	A4349
8305	31 mm	100	A4349
8500	35 mm	100	A4349

Manhood® Absorbent Pouch

Manhood Absorbent Pouch is designed for light to moderate incontinence. Using a super-absorbent polymer, the pouch can hold 216 mL of liquid.



Code	Size	Units	Description
27-4200-B	216 mL	30	One size fits most

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Conveen® Leg Bags

Features a discreet design and color scheme to be less visible under clothing. Two types are available, clamp or lever tap. Anti-kink tubing, soft cloth backing, straps and extended outlet tubing for urine emptying process.

Conveen® Security+ Leg Bag

Secure and reliable urine collection bag with anti-kink tubing reducing risk of accidental leakage.

Not made with natural rubber latex and is DEHP and phthalate free.



Code	Size	Units	Description	HCPCS*
Non-Sterile				
21026	17 oz. (500 mL)	10	20" adjustable tubing, clamp closure, flocked backing, 1" leg bag straps	A4358
21027	34 oz. (1000 mL)	10	20" adjustable tubing, clamp closure, flocked backing, 1" leg bag straps	A4358
Sterile				
21029	17 oz. (500 mL)	10	3" non-adjustable tubing, lever tap, 1" leg bag straps	A4358
21034	17 oz. (500 mL)	10	20" non-adjustable tubing, lever tap, 1" leg bag straps	A4358
21053	34 oz. (1000 mL)	10	3" non-adjustable tubing, lever tap, 1" leg bag straps	A4358
21054	34 oz. (1000 mL)	10	20" non-adjustable tubing, lever tap, 1" leg bag straps	A4358

Conveen® Contour Leg Bag*

Deluxe leg bag with pre-attached, adjustable, anti-kink tubing with fabric straps. Secure, user-friendly closure, soft cloth backing, and chambered panels for discretion and low profile. **Not made with natural rubber latex.**



Non-Sterile

Code	Size	Units	Description	HCPCS*
5170	21 oz. (600 mL)	10	Medium with 18" adjustable tubing, clamp closure, 1" leg bag straps	A4358
5174	28 oz. (800 mL)	10	Large with 18" adjustable tubing, clamp closure, 1" leg bag straps	A4358

Sterile

5171	21 oz. (600 mL)	10	Medium with 18" non-adjustable tubing, clamp closure, 1" leg bag straps	A4358
------	-----------------	----	---	-------

* **WARNING:** This product can expose you to Di(2-ethylhexyl)phthalate (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov

Conveen® Standard Drainage Bag

Large drainage bag with soft cloth backing. Secure, user-friendly closure. **Not made with natural rubber latex and is DEHP and phthalate free.**



Code	Size	Units	Description	HCPCS*
5062	51 oz. (1500 mL)	10	Non-sterile, extra large with 36" tubing, clamp closure, straps not included	A4357

Conveen® Security+ Drainage Bag (formerly Moveen®)

Large volume (2 liters) bedside drainage bag. Non-sterile package includes 30 bags and tubes and 3 hangers. Sterile package includes 10 bags and 2 hangers. **Not made with natural rubber latex and is DEHP and phthalate free.**



Non-Sterile

Code	Size	Units	Description	HCPCS*
21346	2 liters	30	52" adjustable tubing, slide closure	A4357

Sterile

21356	2 liters	10	Non-adjustable 52" tubing, slide closure	A4357
-------	----------	----	--	-------

Conveen® Active Leg Bag

The Conveen Active Leg bag provides a reliable and discreet solution for patients. It wraps around the thigh with the bag facing inward for minimal visibility, even when wearing shorts. Developed to fit a normal active lifestyle and easy to use. **Not made with natural rubber latex.**



Code	Size	Units	Description	HCPCS*
25502	8.5 oz. (250 mL)	30	Non-sterile, 3" adjustable tubing. Wrap 13" long and 4.5" wide	A4358

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Cysto-Care®

Folysil® 100% Silicone Indwelling Catheters.



Folysil® 2-way Pediatric

Code	Balloon Size	Size	Description	Units	HCPCS*
Straight Tip					
AA 6106	1.5 cc	6 FR	12" length	5	A4344
AA 6108	3 cc	8 FR	12" length	5	A4344
AA 6110	3 cc	10 FR	12" length	5	A4344
Coudé Tip					
AA 6308	3 cc	8 FR	12" length	5	A4340
AA 6310	3 cc	10 FR	12" length	5	A4340
Open-tip					
AA 6408	3 cc	8 FR	10" length	5	A4344
AA 6410	3 cc	10 FR	12" length	5	A4344



Folysil® 2-way Straight Tip

Code	Balloon Size	Size	Description	Units	HCPCS*
AA 6112	10 cc	12 FR	16" length	5	A4344
AA 6114	10 cc	14 FR	16" length	5	A4344
AA 6116	15 cc	16 FR	16" length	5	A4344
AA 6118	15 cc	18 FR	16" length	5	A4344
AA 6120	15 cc	20 FR	16" length	5	A4344
AA 6122	15 cc	22 FR	16" length	5	A4344
AA 6124	15 cc	24 FR	16" length	5	A4344
AA 6C18	30 cc	18 FR	16" length	5	A4344
AA 6C20	30 cc	20 FR	16" length	5	A4344
AA 6C22	30 cc	22 FR	16" length	5	A4344
AA 6C24	30 cc	24 FR	16" length	5	A4344



Folysil® 2-way Coudé Tip

Code	Balloon Size	Size	Description	Units	HCPCS*
AA 6312	10 cc	12 FR	16" length	5	A4340
AA 6314	10 cc	14 FR	16" length	5	A4340
AA 6316	15 cc	16 FR	16" length	5	A4340
AA 6318	15 cc	18 FR	16" length	5	A4340
AA 6320	15 cc	20 FR	16" length	5	A4340
AA 6322	15 cc	22 FR	16" length	5	A4340
AA 6324	15 cc	24 FR	16" length	5	A4340

Folysil® 2-way Open Tip

Code	Balloon Size	Size	Description	Units	HCPCS*
AA 6412	10 cc	12 FR	16" length	5	A4344
AA 6414	10 cc	14 FR	16" length	5	A4344
AA 6416	15 cc	16 FR	16" length	5	A4344
AA 6418	15 cc	18 FR	16" length	5	A4344
AA 6422	15 cc	22 FR	16" length	5	A4344
AA 6424	15 cc	24 FR	16" length	5	A4344

*For important information regarding the use of HCPCS codes listed and code definitions see page 180 | For important safety information see page 8

Intermittent Catheter Accessories

This tubing attaches to the end of an intermittent catheter. It is designed to allow the individual to catheterize when unable to reach the toilet from a slight distance.

Self-Cath® Intermittent Catheter Extension Tube



Code	Size	Units	HCPCS*
475	24"	12	A4331

⚠ WARNING: This product can expose you to Di(2-ethylhexyl)phthalate (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov

Male External Catheter/Bag Accessories

Bag Hanger

Reusable, sturdy, plastic drainage bag hanger secures tubing and ensures bag stability when attached to a bed or wheelchair.

Not made with natural rubber latex



Code	Size	Units	HCPCS*
5070	One size	1	Medicare A9900 Other A9999


Conveen® Fabric Leg Bag Straps

Soft, reusable, elastic fabric, button straps with hook and loop closure. Each pair includes a 1.5" width and 1" width strap.

Not made with natural rubber latex



Code	Size	Units	HCPCS*
50501	23" length	10 pair	A5114

A middle-aged man with a grey beard and balding head is shown in a maroon t-shirt and white work overalls. He is leaning forward, holding a long metal rod or tool horizontally across his chest. He is looking towards the camera with a slight smile. The background shows a construction site with a body of water, a stone wall, and some buildings in the distance under a cloudy sky. The overall tone is professional and focused.

Flemming

Wound

Making life easier_

Dressing selection chart

 = exudate level
 Ag = antimicrobial silver

PREPARATION

Cleanse

Sea-Cleans®
Wound
Cleanser



Rehydrate

Woun'Dres®
Hydrogel



Debride

Triad®
Wound
Dressing



EXUDATE

Superficial
Partial thickness

Biatain®
Silicone
Lite



<2cm depth
No filler needed

Biatain®
Silicone



or

Biatain®
Non-Adhesive
+ wrap or
stockinette






>2cm depth
Cavity or tunneling





Biatain®
Alginate
+ cover
dressing



BIOBURDEN

<p><2cm depth No filler needed</p>	<p>Biatain® Silicone Ag</p>  <p>or</p> <p>Biatain® Ag Non-Adhesive + wrap or stockinette</p>  <p>—</p>
<p>>2cm depth Cavity or tunneling</p>	<p>Biatain® Alginate Ag + cover dressing</p>  <p>—</p>

PROTECTION

<p>Granulation Or skin tears</p>	<p>Biatain® Contact + cover dressing</p> 
<p>Friction</p>	<p>Comfeel® Plus Transparent</p> 
<p>Stage 1 or 2 Pressure injuries</p>	<p>Comfeel® Plus</p> 
<p>Difficult-to-dress</p>	<p>Triad® Wound Dressing</p> 

Discover the difference of Biatain® dressings with 3DFit Technology™



Biatain foam dressings feature 3DFit Technology for optimal moist wound healing



Conforms to the wound bed

Upon contact with exudate, the unique foam conforms to the wound bed and reduces exudate pooling



Absorbs vertically

Microcapillaries go into action to absorb the exudate vertically and create optimal healing conditions



Retains exudate

Exudate is locked-away and retained even under compression, reducing risk of leakage and maceration

Visit www.biatainsiliconeus.com to learn more

Biatain® Silicone Foam Dressing

- Features 3DFit Technology™ for optimal moist wound healing
- Gentle and secure fit
- Soft and flexible
- Aseptic 3-piece non-touch opening



Code	Size	Units	HCPCS*
33434	3 x 3" (7.5 x 7.5 cm)	10	A6212
33435	4 x 4" (10 x 10 cm)	10	A6212
33436	5 x 5" (12.5 x 12.5 cm)	10	A6212
33437	6 x 6" (15 x 15 cm)	5	A6213
33438	7 x 7" (17.5 x 17.5 cm)	5	A6213
33400	4 x 8" (10 x 20 cm)	5	A6245
33401	4 x 12" (10 x 30 cm)	5	A6246
33406	Heel 7 x 7" (18 x 18 cm)	5	A6245
39000**	Sacral 9.8 x 9.8" (25 x 25 cm)	5	A6213
39001**	Sacral 6 x 7.5" (15 x 19 cm)	5	A6210
39010**	Multishape 5.5 x 7.6" (14 x 19.5 cm)	5	A6212

💧—💧💧 Exudate Level

**indicated for pressure injury prevention

Biatain® Silicone Lite Foam Dressing

- Features 3DFit Technology™ for optimal moist wound healing
- Gentle and secure fit
- Thin, soft, and ultra flexible
- Aseptic 3-piece non-touch opening



Code	Size	Units	HCPCS*
33452	2 x 2" (5 x 5 cm)	5	A6413
33453	2 x 5" (5 x 12.5 cm)	5	A6212
33444	3 x 3" (7.5 x 7.5 cm)	10	A6212
33445	4 x 4" (10 x 10 cm)	10	A6212
33446	5 x 5" (12.5 x 12.5 cm)	10	A6212

💧—💧 Exudate Level

Biatain® Non-Adhesive Foam Dressing

- Features 3DFit Technology™ for optimal moist wound healing
- Soft and flexible dressing with beveled edges
- Suitable for compression wraps and fragile skin



Code	Size	Units	HCPCS*
6105	2 x 2.75" (5 x 7 cm)	10	A6209
3410	4 x 4" (10 x 10 cm)	10	A6209
3413	6 x 6" (15 x 15 cm)	5	A6210
3416	8 x 8" (20 x 20 cm)	5	A6211

💧—💧💧 Exudate Level

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Biatain® Soft-Hold Foam Dressing

- Features 3DFit Technology™ for optimal moist wound healing
- Skin-friendly contact adherent provides a “third hand” during dressing changes
- Helps minimize adhesive trauma, easy to apply and remove



Code	Size	Units	HCPCS*
3473	2 x 2.75" (5 x 7 cm)	5	A6209
3470	4 x 4" (10 x 10 cm)	5	A6209
3475	6 x 6" (15 x 15 cm)	5	A6210

💧—💧💧 Exudate Level

Biatain® Adhesive Foam Dressing

- Features 3DFit Technology™ for optimal moist wound healing
- Absorptive and retentive foam with a hydrocolloid-based adhesive border
- Designed for wounds with healthy periwound skin that need extra dressing adhesion



Code	Size	Units	HCPCS*
3430	4 x 4" (10 x 10 cm) circle pad	10	A6212
3420	5 x 5" (12.5 x 12.5 cm)	10	A6212
3423	7 x 7" (18 x 18 cm)	5	A6213
3485	Sacral 9 x 9" (23 x 23 cm)	5	A6213
3488	Heel 7.5 x 8" (19 x 20 cm)	5	A6212

💧—💧💧 Exudate Level

Biatain® Ag Foam Dressings combine 3DFit Technology™ with controlled ionic silver-release. The dressings demonstrate in-vitro antibacterial activity in certain strains known to be detrimental to wound healing.¹ Clinical evidence shows a 45-56% ulcer size reduction within 4 weeks.²⁻⁵ **Not made with natural rubber latex; sterile.**

Biatain® Ag Non-Adhesive Foam Dressing

- Features 3DFit Technology™ for optimal moist wound healing
- Combines effective antimicrobial technology with exudate management
- Soft and flexible dressing with beveled edges
- Suitable for compression wraps and fragile skin



Code	Size	Units	HCPCS*
9622	4 x 4" (10 x 10 cm)	5	A6209
9625	6 x 6" (15 x 15 cm)	5	A6210

💧—💧💧 Exudate Level

Biatain® Ag Adhesive Foam Dressing

- Features 3DFit Technology™ for optimal moist wound healing
- Combines effective antimicrobial technology with exudate management
- Absorptive and retentive foam with a hydrocolloid-based adhesive border
- Designed for wounds with healthy periwound skin that need extra dressing adhesion



Code	Size	Units	HCPCS*
9632	5 x 5" (12.5 x 12.5 cm)	5	See PDAC
9635	7 x 7" (18 x 18 cm)	5	See PDAC
9641	Sacral 9 x 9" (23 x 23 cm)	5	See PDAC
9643	Heel 7.5 x 8" (19 x 20 cm)	5	See PDAC


💧—💧💧 Exudate Level

1. Data on file at Coloplast Corp.
2. Münter K-C et al. Effect of a sustained silver-releasing dressing on ulcers with delayed healing: the CONTOP study. *Journal of Wound Care* 2006, 15(5), 199-206.
3. Jørgensen B et al. The silver-releasing foam dressing, Contreet Foam, promotes faster healing of critically colonised venous leg ulcers: a randomised, controlled trial. *International Wound Journal* 2005, 2(1), 64-73.
4. Rayman et al. Sustained silver-releasing dressing in the treatment of diabetic foot ulcers. *British Journal of Nursing* 2005, 14(2), 109-114.
5. Karlsmark T. et al. Clinical performance of a new silver dressing, Contreet Foam, for chronic exuding venous leg ulcers. *Journal of Wound Care* 2003, 12(9), 351-354.

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Biatain® Alginate

Biatain® Alginate is a highly absorbent wound dressing composed of alginate and hydrocolloid (sodium carboxymethylcellulose). On contact with wound exudate it converts to a soft cohesive gel which provides an optimal moist wound healing environment. **Not made with natural rubber latex; sterile.**



Biatain® Alginate Dressing

- Fast-gelling alginate
- One-piece removal
- Recommended for cavity filling



Code	Size	Units	HCPCS*
3705	2 x 2" (5 x 5 cm)	30	A6196
3710	4 x 4" (10 x 10 cm)	10	A6196
3715	6 x 6" (15 x 15 cm)	10	A6197

💧 – 💧💧 Exudate Level

Biatain® Alginate Rope

- Fast-gelling alginate
- One-piece removal
- Recommended for cavity filling



Code	Size	Units	HCPCS*
3740	1 x 17.5" rope (3 x 44 cm)	6	A6199

💧 – 💧💧 Exudate Level

Biatain® Alginate Ag Dressing

- Fast-gelling alginate containing an ionic silver complex with broad antimicrobial profile
- One-piece removal
- Recommended for cavity filling



Code	Size	Units	HCPCS*
3755	2 x 2" (5 x 5 cm)	30	A6196
3760	4 x 4" (10 x 10 cm)	10	A6196
3765	6 x 6" (15 x 15 cm)	10	A6197

💧—💧💧 Exudate Level

Biatain® Alginate Ag Rope

- Fast-gelling alginate rope containing an ionic silver complex with broad antimicrobial profile
- One-piece removal
- Recommended for cavity filling



Code	Size	Units	HCPCS*
3780	1 x 17.5" rope (3 x 44 cm)	10	A6199

💧—💧💧 Exudate Level

The Comfeel® Plus portfolio offers a range of dressings to manage and protect non- to moderately-exuding wounds.

- Exudate management with fewer dressing changes and less disruption to wound bed
- Specialty shapes and sizes for difficult-to-dress areas
- Color guide indicates best time to change dressing
- Non-touch application system for ease of use
- **Not made with natural rubber latex; sterile**

Comfeel® Plus Transparent Dressing (Thin)

- Comfortable, flexible hydrocolloid adaptable to wound and body contours
- Clear dressing and mapping grid allow for visual inspection without disturbing wound
- Does not contain alginate



Code	Size	Units	HCPCS*
33530	2 x 2.75" (5 x 7 cm)	10	A6234
33533	4 x 4" (10 x 10 cm)	10	A6234
33536	3.5 x 5.5" (9 x 14 cm)	10	A6235
33539	6 x 6" (15 x 15 cm)	5	A6235

💧 Exudate Level

Comfeel® Plus Dressing

- Visual wound monitoring with mapping grid
- Contains alginate for enhanced absorption capacity



Code	Size	Units	HCPCS*
33146	1.5 x 2.5" (4 x 6 cm)	30	A6234
33110	4 x 4" (10 x 10 cm)	10	A6234
33115	6 x 6" (15 x 15 cm)	5	A6235
33120	8 x 8" (20 x 20 cm)	5	A6236

💧—💧 Exudate Level

Comfeel® Plus Contour Dressing

- Butterfly-shaped hydrocolloid dressing for better conformity and adhesion
- Specifically designed for difficult-to-dress areas
- Contains alginate for enhanced absorption



Code	Size	Units	HCPCS*
33280	24 sq. in. (6 x 8 cm)	5	A6235
33283	42 sq. in. (9 x 11 cm)	5	A6235

💧—💧 Exudate Level

Comfeel® Plus Sacral Dressing

- Triangular-shaped hydrocolloid dressing designed for the sacral area
- Contains alginate for enhanced absorption



Code	Size	Units	HCPCS*
33285	7 x 7" (17 x 17 cm)	5	A6235

💧—💧 Exudate Level

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Triad® Hydrophilic Wound Dressing

- Zinc oxide-based hydrophilic dressing absorbs moderate levels of wound exudate
- Facilitates autolytic debridement
- Ideal alternative for difficult-to-dress areas
- **Not made with natural rubber latex; sterile**



Code	Size	Units	HCPCS*
1964	2.5 oz. tube (71 g)	12	A6240
1967	6 oz. tube (170 g)	12	A6240

💧—💧 Exudate Level

Sea-Cleans® Wound Cleanser

- Saline-based solution for cleansing and irrigating acute and chronic wounds
- No-rinse, isotonic, pH-balanced formula
- Easy-to-use spray bottle with SPRAY and STREAM setting with irrigation pressure



Code	Size	Units	HCPCS*
1063	6 fl. oz. (178 mL)	12	A6260
1061	12 fl. oz. (355 mL)	12	A6260

💧—💧💧 Exudate Level

Woun'Dres® Collagen Hydrogel

- Rehydrates dry wounds and eschar
- Promotes moist wound environment and autolytic debridement
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
1166	1 oz. tube (28 g)	36	A6248
7690	3 oz. tube (84 g)	12	A6248

💧 Exudate Level

Purilon® Gel

- Sterile hydrogel with hydrating and absorbing properties to promote autolytic debridement
- Maintains a moist wound environment
- Contains alginate
- **Not made with natural rubber latex; sterile**



Code	Size	Units	HCPCS*
3906	0.28 oz. (8 g)	10	A6248
3900	0.5 oz. (15 g)	10	A6248
3903	0.88 oz. (25 g)	10	A6248

💧—💧 Exudate Level

Biatain® Contact

Biatain Contact is a one-sided silicone wound contact layer designed to protect wounds from secondary dressing changes.

- Easy application and removal
- Protects against in-growth
- Bigger pores for enhanced exudate transfer and reduced risk of pooling
- Improved transparency to facilitate wound assessment without removing the dressing



Code	Size	Units	HCPCS*
33560	2 x 3" (5 x 7.5 cm)	10	A6206
33561	3 x 4" (7.5 x 10 cm)	10	A6206
33562	4 x 7" (10 x 18 cm)	10	A6207
33563	6 x 10" (15 x 25 cm)	5	A6208
33564	7 x 12" (18 x 30 cm)	5	A6208

💧—💧💧 Exudate Level

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Biatain® Super, Adhesive

- Fast absorption, distribution and retention of exudate
- Suitable for highly exuding wounds surrounded by healthy skin



Code	Size	Units	HCPCS*
46100	4 x 4" (10 x 10 cm)	10	A6203
46120	5 x 5" (12.5 x 12.5 cm)	10	A6203
46250	5 x 8" (12.5 x 20 cm)	10	A6204
46150	6 x 6" (15 x 15 cm)	10	A6203
46200	8 x 8" (20 x 20 cm)	10	A6204

💧💧 — 💧💧💧 Exudate Level


Biatain® Super, Non-Adhesive

- Fast absorption, distribution and retention of exudate
- Suitable for highly exuding wounds surrounded by fragile skin



Code	Size	Units	HCPCS*
46300	4 x 4" (10 x 10 cm)	10	A6251
46320	5 x 5" (12.5 x 12.5 cm)	10	A6251
46450	5 x 8" (12.5 x 20 cm)	10	A6252
46350	6 x 6" (15 x 15 cm)	10	A6251
46390	8 x 8" (20 x 20 cm)	10	A6252

💧💧 — 💧💧💧 Exudate Level



Skin

Making life easier_

A scientifically simple approach to skin care

For over 35 years, Coloplast has been dedicated to creating and providing high quality skin care products formulated with scientifically simple ingredients.

When patients and caregivers asked for cleaner, simpler products free of harmful ingredients, Coloplast listened. Today, we are proud to offer one of the most simplified and clean skin care product lines on the market.

Skin

What does a scientifically simple skin care line mean?*

NO	Parabens	YES	Hypoallergenic
NO	Formaldehyde-donors	YES	Neonate suitable
NO	Scents	YES	CHG compatible
NO	Alcohol or latex	YES	pH-balanced
NO	Dyes	YES	Non-sensitizing
		YES	Non-irritating
		YES	Fewer ingredients than leading competitors



*Not all products in Coloplast's skin care line comply with these standards.

Three simple steps to promote healthy skin



At Coloplast, we believe that a simplified, structured routine decreases variation in the delivery of care among staff and promotes compliance for optimal outcomes.

Coloplast has developed a unique skin care routine that includes three elements of care: *Cleanse, Moisturize and Protect.*



Cleanse

To remove daily waste and environmental contaminants

Choose a cleanser that is mild, hypoallergenic and PH-balanced to match the acid mantle on the skin to minimize skin sensitivities, allergens and irritants.

Coloplast cleansers:

- EasiCleanse[®] Bath
- Bedside-Care[®] – Foam or Spray
- Bedside-Care[®] Sensitive Skin – Foam or Spray
- Gentle Rain[®] Extra Mild



Moisturize

To avoid skin complications

Choose a moisturizer without harmful ingredients like parabens and formaldehyde donors to avoid risk of skin complications.

Coloplast moisturizers:

- Sween[®] 24
- Sween Cream[®]
- Sween[®] Lotion
- Atrac-Tain[®] Cream



Protect

To avoid exposure to elements

Choose a skin protectant that has high quality ingredients to protect the skin from stool and urine, maintain skin hydration and help prevent skin maceration.

Coloplast protectants:

- Critic Aid[®] Clear
- Critic Aid[®] Paste
- Critic Aid[®] Clear Antifungal
- Baza[®] Clear Ointment
- Baza[®] Protect Cream
- Baza Cleanse and Protect[®] Lotion
- Baza[®] Antifungal Cream

Simple you can see

We make it simple for you to customize patient care with an array of easy-to-identify products.



brand name

product category

product description

uses

graphic symbol for quick identification

color-coded by category



Cleanse - Rinse



Cleanse - No-Rinse



Moisturize



Protect



Protect - Antifungal

Bedside-Care® Foam

No-rinse foam body wash, shampoo and incontinence cleanser

- Cleanses, moisturizes and conditions the skin
- pH-balanced to help maintain the natural acid mantle of the skin
- Alcohol-free
- CHG compatible



Product	Code	Size	Units
Scented	67143	4.1 fl. oz. (120mL)	36
	67145	8.1 fl. oz. (240mL)	12
Unscented	67147	4.1 fl. oz. (120mL)	36
	67146	8.1 fl. oz. (240mL)	12

Bedside-Care® Spray

No-rinse body wash, shampoo and incontinence cleanser

- Tear-free, hypoallergenic
- Cleanses and conditions the skin
- Alcohol-free
- pH-balanced to help maintain the natural acid mantle of the skin
- CHG compatible



Product	Code	Size	Units
Scented	61768	4.1 fl. oz. (120mL)	36
	61762	8.1 fl. oz. (240mL)	12
Unscented	61760	4.1 fl. oz. (120mL)	36
	61761	8.1 fl. oz. (240mL)	12

Bedside-Care® Sensitive Skin Foam

No-rinse foam body wash, shampoo and incontinence cleanser for all ages including neonates¹

- Tear-free, hypoallergenic
- Cleanses and conditions the skin
- Alcohol-free
- pH-balanced to help maintain the natural acid mantle of the skin
- CHG compatible



Code	Size	Units
7300	4 fl. oz. (118 mL)	36
7301	8 fl. oz. (237 mL)	12

Bedside-Care® Sensitive Skin Spray

No-rinse body wash, shampoo and incontinence cleanser for all ages including neonates¹

- Tear-free, hypoallergenic
- Cleanses and conditions the skin
- Alcohol-free
- pH-balanced to help maintain the natural acid mantle of the skin
- CHG compatible



Code	Size	Units
1300	4 fl. oz. (118 mL)	36

1. Association of Women's Health, Obstetrics and Neonatal Nurses. Neonatal Skin Care, Fourth Edition, Evidence-Based Clinical Practice Guideline, 2018

Smarter bedside bathing

EasiCleanse Bath is a no-rinse, self-sudsing, disposable washcloth offering a real bathing experience that is truly standardized.

EasiCleanse® Bath

No-rinse body wash, shampoo, meatal, perineal and incontinence cleanser in a self-sudsing disposable washcloth

- 3 ingredients with no preservatives
- No filmy residue left on skin
- No basin needed, to help avoid the risk of basin-associated contamination
- Suitable for all ages including neonates¹

Code	Size	Units
7055	30 washcloths/pk	30 pks/cs
7056	5 washcloths/pk	100 pks/cs
7058	8 washcloths/pk	80 pks/cs



Visit www.easicleansebath.com to learn more

1. Association of Women's Health, Obstetric and Neonatal Nurses. Neonatal Skin Care, Third Edition. Evidence-Based Clinical Practice Guideline, 2018.

Gentle Rain® Extra Mild

Moisturizing body wash, shampoo and hand wash

- Contains mild ingredients and added emollients for use on fragile skin
- Dermatologist tested, hypoallergenic and dye-free
- pH-balanced to help maintain the natural acid mantle of the skin



Cleaners



Code	Size	Units
7231	5 mL single application packet	300
7230	2 fl. oz. (59 mL)	36
7229	4 fl. oz. (118 mL)	36
7235	8 fl. oz. (237 mL)	36
7233	21 fl. oz. (621 mL) [†]	12
7234	34 fl. oz. (1000 mL) [†]	12

[†] See dispensers on page 58.

Sween® 24*Once-a-day moisturizing body cream*

- Once-a-day application
- Contains 6% dimethicone
- Fragrance, petrolatum, lanolin, alcohol and preservative-free
- CHG compatible
- Suitable for all ages including neonates¹



Code	Size	Units
7090	4 g single application packet	300
7091	2 oz. (57 g)	12
7092	5 oz. (142 g)	12
7095	9 oz. (256 g)	12

Sween Cream®*Moisturizing cream*

- Non-occlusive formula
- Will not interfere with tape adhesion
- CHG compatible



Code	Size	Units
0283	2 g single application packet	300
7067	3 oz. (85 g)	12
7068	6.5 oz. (184 g)	12
7069	12 oz. (340 g) jar	12

Sween® Lotion*Moisturizing lotion*

- Enriched with vitamin E
- Non-occlusive formula
- CHG compatible



Product	Code	Size	Units
Scented	0424	2 mL single application packet	300
	0407	2 fl. oz. (59 mL)	36
	0406	4 fl. oz. (118 mL)	36
	0402	8 fl. oz. (237 mL)	36
	0408	21 fl. oz. (621 mL) [†]	12
Unscented	7222	2 fl. oz. (59 mL)	36
	7223	4 fl. oz. (118 mL)	36

Atrac-Tain® Cream*Moisturizing cream with urea & alpha hydroxy acid (AHA)*

- 10% urea and 4% AHA
- Moisturizing cream for extremely dry, scaly cracked skin and fissured heels
- Fragrance, dye and preservative-free
- CHG compatible



Code	Size	Units
1843	2 g single application packet	300
1802	2 oz. (57 g)	12
1814	5 oz. (142 g)	12



Critic-Aid® Clear

Clear moisture barrier ointment

- Clear moisture barrier
- Minimal ingredients and preservative-free
- Adheres to intact and eroded skin
- Easy to apply and remove
- CHG compatible
- Suitable for all ages including neonates¹



Code	Size	Units
7565	4 g single application packet	300
7566	2.5 oz. (71 g)	12
7567	6 oz. (170 g)	12

Critic-Aid® Skin Paste

Thick moisture barrier paste

- Thick, zinc oxide-based skin protectant
- Adheres to intact and moist, denuded skin
- Contains only 4 ingredients to help minimize potential skin irritation
- CHG compatible
- Suitable for all ages including neonates¹



Code	Size	Units
1944	2.5 oz. (71 g)	12
1947	6 oz. (170 g)	12

Baza® Clear Ointment

Moisture barrier ointment

- Clear, petrolatum-based moisture barrier
- Enriched with skin conditioners
- Contains 96% petrolatum
- Clear formulation makes it easy to visually assess the skin
- CHG compatible



Code	Size	Units
7015	4 g single application packet	300
1005	1.75 oz. (50 g)	12
1006	5 oz. (142 g)	12

1. Association of Women's Health, Obstetric and Neonatal Nurses. Neonatal Skin Care, Third Edition. Evidence-Based Clinical Practice Guideline, 2018.

Baza® Protect Cream

Moisture barrier cream

- Thick, zinc oxide-based skin protectant
- Enriched with skin conditioners
- CHG compatible



Code	Size	Units
1873	4 g single application packet	300
1877	2 oz. (57 g)	12
1880	5 oz. (142 g)	12

Baza Cleanse & Protect® Lotion

All-in-one perineal lotion

- One step provides mild cleansing, a breathable barrier and moisturization
- Available with odor control
- pH-balanced to help maintain the natural acid mantle of the skin
- CHG compatible



Product	Code	Size	Units
Unscented	7710	4 g single application packet	300
	7712	8 fl. oz. (237 mL)	12
Scented	7725	8 fl. oz. (237 mL)	12

Sween® Body Powder

All-body powder

- Absorbs moisture and perspiration
- Non-talcum based formula
- Helps reduce friction and prevent shearing



Code	Size	Units
0506	3 oz. (85 g)	36
0505	8 oz. (227 g)	36

Critic-Aid® Clear Antifungal Ointment*Clear moisture barrier ointment with antifungal*

- 2% miconazole nitrate
- Contains only 5 ingredients to help minimize potential skin irritation
- Adheres to intact and moist, denuded skin
- Clear formulation makes it easy to visually inspect the skin
- CHG compatible



Code	Size	Units
7570	4 g single application packet	300
7571	2 oz. (57 g)	12
7572	5 oz. (142 g)	12

Baza® Antifungal Cream*Moisture barrier antifungal cream*

- 2% miconazole nitrate
- Enriched with skin conditioners
- CHG compatible



Code	Size	Units
1622	4 g single application packet	300
1611	2 oz. (57 g)	12
1607	5 oz. (142 g)	12

Micro-Guard® Powder*Antifungal powder*

- Non-talcum based formula with 2% miconazole nitrate
- Effective treatment of topical fungal infections
- Helps reduce the effects of friction and shear



Code	Size	Units
1337	3 oz. (85 g)	12

Micro-Guard® 4% CHG Solution

Combine with EasiCleanse Bath to deliver a real suds bath and the trusted protection of CHG

Code	Size	Units
67770	4 fl. oz. (118 mL)	48

InterDry®

Moisture-wicking fabric with antimicrobial silver



InterDry with FourFold Technology™ is an innovative moisture wicking fabric with antimicrobial silver.

- Manages all four causes of intertrigo: moisture, friction and bacterial or fungal organisms in skin folds¹ and other skin-to-skin and skin-to-device contact areas
- Manages the symptoms associated with intertrigo such as erythema, painful rash, denudement, satellite lesions and odor²
- Silver within the fabric provides effective antimicrobial action for up to 5 days

By targeting all four causes of Intertrigo, InterDry is clinically proven to provide complete, cost-effective relief within 5 days for as little as \$5³.

InterDry®

Code	Size	Units
7910	10 x 144"	10 rolls/cs
7912	10 x 36"	10 pks/box
7915	10 x 18"	40 pks/box



Visit www.interdry.com to learn more

1. Vorbeck, E. Evaluation of a Skin Fold Management Textile with Antimicrobial Silver Complex in a Variety of Case Studies
2. Freyberg J, et al. Moisture Management Challenges for the WOC Nurse presented at the WOCN Society 39th Annual Conference, June 2007
3. Kennedy-Evans KL, Viggiano B, Henn T, Smith D. Multi-site feasibility study using a new textile with silver for management of skin conditions located in skin folds. Poster presented at: 20th Annual Symposium on Advanced Wound Care; April 28 - May 1, 2007; Tampa, FL and 39th WOCN® Society Annual Conference; June 9-13, 2007; Salt Lake City, UT.

Isagel®

No-rinse, antiseptic hand cleansing gel

- Formulated with emollients for extra mildness
- Contains 60% ethanol alcohol
- CHG and latex glove compatible



Code	Size	Units
1644	4 fl. oz. (118 mL)	36
1645	21 fl. oz. (621 mL) [†]	12
7041	27 fl. oz. (800 mL) [†]	12



[†] See dispensers on page 58.

Bedside-Care® Gentle Touch Cloths

Disposable washcloths for use with all-body bathing and perineal cleansing products.

- Rayon, for softness and absorbency
- Polyester, for strength and tear resistance



Code	Size	Units
7057	25 washcloths/pack	24

Coloplast® Wall Dispenser

- Holds 27 fl. oz. and 34 fl. oz. (800 mL and 1000 mL) refill bags
- Completely closed to prevent possible contamination
- Metered 1 mL pump dispenses appropriate amount
- No-drip nozzle leaves facility clean and safe
- Easy-to-operate push pad meets Americans with Disabilities Act guidelines
- Keyless opening for simple cartridge replacement
- Installs quickly with choice of screws or tape



Code	Size	Units
7251	Holds 27 fl. oz. and 34 fl. oz. (800 mL and 1000 mL)	12
7262	Drip tray	60

Embrace™ Bracket

- Holds 21 fl. oz. (621 mL) Coloplast containers securely in place
- Improved tamper and theft-resistant design
- Optional locking collar conceals cap
- Installs quickly with choice of screws or tape



Code	Size	Units
0089	Holds 21 fl. oz. (621 mL)	12



Stina

Making life easier_

Ostomy

Through BodyFit Technology, we offer a secure fit to real bodies *to help prevent leakage and skin issues*

Leakage and skin issues are closely related challenges. The **right fit** helps **prevent** both.




Leakage

76%

of people with a stoma have experienced leakage within the last six months.¹

Leakage

is often caused by a poor fit to the body.²



Skin issues

Skin issues

There is a significant correlation between leakage frequency and the severity of the peristomal skin condition.²

61%

of users attributed their most recent skin issue to leakage under the barrier.³

1. Results from the 2014 Ostomy Life Study - an international survey aimed to better understand the challenges people living with an ostomy face (n=4138 ostomates). Published: Claessens I, Probert R, et. al. The Ostomy Life Study: the everyday challenges faced by people living with a stoma in a snapshot. Gastrointestinal Nursing. 2015;13(5), 18-25.
 2. Results from a multinational, randomized controlled cross-over trial comparing SenSura Mio Convex to the study participants' own appliance (n=129). Published: Walker H, Hopkins G, et. al. Raising the bar: New flexible convex ostomy appliance - a randomised controlled trial. World Council of Enterostomal Therapists Journal. 2016;36(1) S6-S11.
 3. Results from the 2018-2019 Ostomy Life Study - an international survey (n=3112 ostomates); Coloplast market research, Data on file.

The #1 ostomy brand in the world, **SenSura® Mio**, is now **Extended Wear***!



To learn more, visit us at www.sensuramio.us or call 1-888-726-7872.

SenSura Mio was designed to fit individual bodies – and now, a recent 2,359-patient survey demonstrated that SenSura Mio users report less leakage onto their clothes than other extended wear products, while experiencing similar or even longer wear time using SenSura Mio¹.

While SenSura Mio was initially launched into the standard wear category, recent user experience¹, clinician evaluation², and lab-based data³ reinforces SenSura Mio's performance as an extended wear product.

With this change, users who need or would prefer an extended wear product will now have access to SenSura Mio's unique technology, designed to reduce the worry about leakage & provide a comfortable fit:

- Full-circle filter to reduce ballooning
- Discreet neutral-gray, soft woven textile
- Wide outlet
- A flexible fit, thanks to BodyFit Technology, an elastic adhesive that is unique in all SenSura Mio barriers.

*Excluding 1-piece closed-end pouches

1. Ostomy pouch user survey, 2021, n=2,359. Dataset collected by Coloplast but includes all brands. Data on file.

2. Product evaluation of SenSura Mio Convex, 2016, n=98 clinicians.

3. Data from 48h bench testing performed by Coloplast; may not be indicative of clinical performance. Data on file.

Products of Distinction



SenSura® Mio Convex Flip

Designed to fit curves, bulges and hernias with less creasing and folding



SenSura® Mio Flex

- Floating flange
- Easy to apply
- Flexible coupling



Soft
6mm

Light
6mm

Deep
9mm

SenSura® Mio Convex

Three types and two depths to fit different bodies



SenSura® Mio High Output

The only 1-piece convex option with a spout outlet for fecal output. Also available in 1-piece flat, 2-piece Click, and 2-piece Flex

Brava[®]

Security to live *everyday life*

**“I was determined
to play sports
again. The
Brava[®] Protective
Seal gave me
confidence.”**

– Kevin*

*Kevin is a Brava user who has received compensation from Coloplast to provide this information.



Accessories

Brava® supporting products

*For customizing pouching system fit
and for issues such as leakage and
skin irritation.*

NEW! Brava® Protective Seal Convex

The new Brava Protective Seal Convex is designed to protect the skin against leakage for flush or retracted stomas and patients with folds, creases or scarring around their stoma.



Brava® Protective Seal

- Easily shaped
- Leaves minimal residue behind
- Follows body movements
- Significantly reduces leakage compared to competitor rings*

The **Brava® Ostomy Support Belt** brings relief by relieving the heavy sensation of hernias.

- Comfortable and breathable fabric
- Stretchable to fit body contours
- Easy-to-use pocket closure
- Silicone grip to prevent sliding and rolling



*Based on bench test data. May not be indicative of clinical performance.

Think again. Think Brava®

Designed to address leakage and care for your skin.

Brava® Elastic Barrier Strips

Skin-friendly strips are used to “picture frame” around an ostomy barrier, designed for extra security.

- Elastic (follow the contours of the body)
- Ensure the position of the barrier and designed to minimize roll-up
- Skin-friendly and secure alternative to tape
- **Not made with natural rubber latex**



Brava® Elastic Barrier Strips- Curved

For round or oval barrier shapes

Code	Size	Units	HCPCS*
120700**	5½" strips	20 strips	A4362

**2 strips = 1 billable unit of A4362

Brava® Elastic Barrier Strips- Y shape

Designed for users with challenging body shapes, e.g. hernias or skin folds



Code	Size	Units	HCPCS*
120721**	5¾" strips	30 strips	A5121

**3 strips = 1 billable unit of A5121

Brava® Elastic Barrier Strips- Straight

For square or rectangular barrier shapes



Code	Size	Units	HCPCS*
120740**	5¼" strips	20 strips	A4362

**2 strips = 1 billable unit of A4362

Brava® Elastic Barrier Strip XL

Designed for increased security & support around the barrier. Wider shape accommodates the petals of Convex Flip.



Code	Size	Units	HCPCS*
12076**	6¾" strips	30 strips	A5121

**3 strips = 1 billable unit of A5121

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Brava® Protective Seal

Moldable barrier ring designed to adhere securely and protect the skin around the stoma.

- Soft, flexible material adapts to body shape as the wearer bends and stretches
- Can be easily and cleanly removed
- Alcohol-free and sting-free
- Multiple sizes for a custom fit
- **Not made with natural rubber latex**



Code	Stoma size	Adhesive diameter	Thickness	Seal width	Units	HCPCS*
12035	¾"	1 7/8"	2.5 mm	Standard (5/8")	10	A4385
12036	¾"	2 ¼"	2.5 mm	Wide (¾")	10	A4385
12038	¾"	2 ½"	2.5 mm	Wide (7/8")	10	A4385
12037	1 1/8"	2 ¼"	2.5 mm	Standard (5/8")	10	A4385
12039	1 3/8"	2 ½"	2.5 mm	Standard (5/8")	10	A4385
12045	¾"	1 7/8"	4.2 mm	Standard (5/8")	10	A4385
12046	¾"	2 ¼"	4.2 mm	Wide (¾")	10	A4385
12048	¾"	2 ½"	4.2 mm	Wide (7/8")	10	A4385
12047	1 1/8"	2 ¼"	4.2 mm	Standard (5/8")	10	A4385
12049	1 3/8"	2 ½"	4.2 mm	Standard (5/8")	10	A4385
12032	¾"	3"	2.5 mm	Extra wide (1 1/8")	10	A4385
12033	1 1/8"	3"	2.5 mm	Extra wide (1 5/16")	10	A4385

Brava® Protective Seal Convex

Moldable barrier ring designed to provide convexity for flush or retracted stomas and patients with folds, creases or scarring around their stoma.

- High-tack adhesive that flows quickly into skin, filling small creases
- Erosion resistant composition
- High flexibility to maintain seal during body movement
- Steep sides designed to provide convexity to peristomal area
- Thin edges fit snugly around peristomal area
- **Not made with natural rubber latex**



Code	Stoma size	Thickness	Units	HCPCS*
12090	¾"	8.2mm	10	A4411
12091	1"	8.2mm	10	A4411
12092	1 3/16"	8.2mm	10	A4411
12093	1 3/8"	8.2mm	10	A4411
12094	1 9/16"	8.2mm	10	A4411
12095	1 3/4"	8.2mm	10	A4411

Brava® Moldable Ring

Designed for dual-action leakage reduction - used to fill in uneven skin surfaces and protect the skin.

- Molds and adheres securely to dips, creases, and folds to help create a tight seal around the stoma
- Helps keep skin healthy by protecting it from output and absorbing moisture from the skin
- Easy to handle, mold and apply
- Alcohol-free and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
120307	2 mm (thin)	10	A4385
120427	4.2 mm (thick)	10	A4385

Brava® Strip Paste

Pectin-based product that protects peristomal skin.

- Protects peristomal skin from stoma output
- Used to fill uneven skin surfaces to create a flat pouching surface
- Alcohol-free and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
26555	2 oz. total package (.2 oz. per strip)	10 strips	A4406

Brava® Paste

A sting-free alternative to conventional alcohol-based paste.

- Creates a tight, effective seal between peristomal skin and the ostomy/fistula appliance
- **No alcohol** = no need to "set"; **contains pectin**
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
12050	2 oz.	1	A4406

Coloplast® Paste

Used to fill in uneven skin surfaces to create a flat pouching surface.

- Provides tight sealing around the stoma
- Protects against stoma output, minimizing risk of maceration
- **Low alcohol content; no pectin**
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
2650	2 oz.	12	A4405

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Brava® Skin Cleanser Wipes

The wipe that cleans and cares. A convenient way to cleanse and maintain healthy skin when changing the ostomy pouching system.

- Contains ProVitamin B5 to provide a soothing and moisturizing effect
- Does not contain alcohol or perfumes
- No residue left on the skin that could affect barrier adhesion
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
12080	1 wipe	15 wipes/pack	A9286**

Brava® Lubricating Deodorant

A pouch lubricant + deodorizer in one. Powder-fresh liquid neutralizes odor and allows output to slide to the bottom of the pouch for easier emptying.

- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
12060	.25 fl. oz. (7.5 mL) sachet	20 sachets	A4394
12061	8 oz. bottle	16	A4394

Brava® Adhesive Remover Spray

- Designed to lift the edges of the barrier for ease in removal
- Dries within seconds without leaving an oily residue behind
- Alcohol-free, silicone-based, and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
120105	1.7 oz. (50 mL)	1	A4455**

Brava® Adhesive Remover Wipes

- Wipe away adhesive residue
- Dries within seconds without leaving an oily residue behind
- Alcohol-free, silicone-based, and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
120115	1 wipe	30 wipes/box	A4456**

Brava® Skin Barrier Spray

- Protects against stoma output and adhesives
- Dry within seconds without building a thick or tacky layer
- Alcohol-free, silicone-based, and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
120205	1.7 oz. (50 mL)	1	A4369**

Brava® Skin Barrier Wipes

- Protects against stoma output and adhesives
- Dry within seconds without building a thick or tacky layer
- **Tip:** wring the wipe like a wet rag to expel the liquid, then wipe to disperse on skin
- Alcohol-free, silicone-based, and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
120215	1 wipe	30 wipes/box	A5120**

Brava® Protective Sheet

Hydrocolloid skin barrier that protects the skin from stoma output.

- Offers effective protection for peristomal skin and "3-way stretch"
- Provides a ready surface for good adhesion
- Adheres to moist skin
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
32105	4 x 4" (10 x 10 cm)	10	A4362
32155	6 x 6" (15 x 15 cm)	5	A5121
32205	8 x 8" (20 x 20 cm)	5	A5122

Brava® Powder

Non-sterile protective powder for moisture absorption.

- May use in conjunction with Brava® Skin Barrier Spray or PREP™ to "crust", thereby protecting moist skin
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
19075	1 oz. (25 g)	16	A4371

Hex-On®

Odor Antagonist

- Helps control airborne odors
- Non-aerosol spray
- **Not made with natural rubber latex**



Product	Code	Size	Units
Fresh Linen	7583	2 fl. oz./59 mL	12
Fragrance Free	67589	2 fl. oz./59 mL	12

Ostomy Scissors

Curved ostomy scissors.

- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
95050	One size	30/bx	N/A

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

**Product has not been through PDAC verification

Coloplast® Skin Barrier Rings

Soft and flexible barrier rings provide skin protection for stomas, fistulas, and draining wounds.

- Non-pectin based
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
2310	⅜" (10 mm)	30	A4404
2315	⅔" (15 mm)	30	A4404
2320	¾" (20 mm)	30	A4404
2325	1" (25 mm)	30	A4404
2330	1⅛" (30 mm)	30	A4404
2340	1⅜" (40 mm)	30	A4404
2350	2" (50 mm)	30	A4404

Prep™ Protective Skin Barrier

Clear film forms a breathable barrier protecting against skin irritants such as adhesives, urine and feces. No need to remove between applications.

- **Not made with natural rubber latex**



Wipes

Code	Size	Units	HCPCS*
2041	54 wipes/pack	12 packs/box	A5120

Dabber

Code	Size	Units	HCPCS*
0925	2 fl. oz. (59 mL)	12	A4369

Pouch Filters (Filtrodor®)

Charcoal-activated, deodorizing filter with self-sealing foam center; adheres to pouch.

- **Not made with natural rubber latex**



Code	Units	HCPCS*
0509	50 <i>each box contains 5 packs of 10</i>	A4368

Drainable Pouch Clamp

Simple closure clamp is easy to manipulate, even for those with low dexterity.

- **Not made with natural rubber latex**



Code	Units	HCPCS*
9500	20	A4363

Brava® Ostomy Support Belt

This belt relieves the heavy sensation of hernias and may help prevent hernias post-operatively. Comfortable and breathable fabric, stretchable to fit body contours, intuitive pocket closure, and includes a silicone grip to prevent sliding and rolling.



- Not made with natural rubber latex
- All sizes are an 8" width

Code	Size	Units	HCPCS*
12003	S (29-32")	1	A4396
12004	M (31-35")	1	A4396
12005	L (34-39")	1	A4396
12006	XL (37-42")	1	A4396
12007	2XL (40-46")	1	A4396
12008	3XL (44-51")	1	A4396
12009	4XL (50-58")	1	A4396

Brava® Belt for SenSura® Mio

Neutral gray color was designed to blend in under clothing. Four-tab connection provides extra stability across the entire peristomal area. Designed specifically for use with SenSura® Mio products with belt tabs.



- Not made with natural rubber latex

Code	Size	Units	HCPCS*
4237	40" length	15	A4367
4247	61" length	15	A4367

Brava® Belt

Neutral gray color was designed to blend in under clothing. Comfortable, adjustable belt offers extra security for people with a flush or retracted stoma by supporting a convex barrier and allowing the stoma to protrude.



- Single tab attachment is designed to be compatible with most conventional ostomy products with belt tabs.
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
4215	43 $\frac{2}{3}$ " length	15	A4367
4220	49" length	15	A4367

Urostomy Night Bag

Includes a fixed connector that fits directly to SenSura® Mio, SenSura® and Assura® urostomy pouches. The bag holds up to 2 liters and has a 120 cm long flexible tube with anti-reflux valve. Bag hanger (#5070) not included. Order separately.



- Not made with natural rubber latex

Code	Units	HCPCS*
21365	10	A4357

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Bed Drainage Bag

Wide-lumen bag tubing specifically designed for fecal output. Holds up to 2000 mL. Long tubing can be cut easily and allows nurses to find the right length for bed drainage. Can be used with fistula / wound management systems or high-output pouches. Bag hanger (#5070) not included. Order separately.



- **Not made with natural rubber latex**

Code	Units	HCPCS*
14010	6	A4357

Bag Hanger

Reusable, sturdy, plastic drainage bag hanger secures tubing and ensures bag stability.



- **Not made with natural rubber latex**

Code	Size	Units	HCPCS*
5070	One Size	1	Medicare N/A Other A4421

Extra Drain Port

For use with Fistula Wound Management Systems (page 149). Drain port can be easily applied to the lid to allow for catheter insertion. For catheters ranging from 18FR - 22FR, the drain port is ready for use. For catheters ranging from 24FR-30FR, pull off the top of the drain port - DO NOT CUT.



- **Not made with natural rubber latex**

Code	Units
14015	12

Coloplast® Transparent Fistula Lids

Replacement lids match the transparent lid included with each fistula pouch.



- **Not made with natural rubber latex**

Code	Size	Units
14021	Mini (matches 14050)	10
14031	Midi (matches 14060)	10
14041	Maxi (matches 14070)	10

Irrigation Sets

Coloplast® Economy Version

Durable water bag with integrated thermometer, bold print, large opening. Can be used with 1- or 2-piece ostomy systems.

- **Not made with natural rubber latex**

Code	Units	Includes	HCPCS*
1500	1	1- Water Bag/Regulator	A4398
		1- Cone (#1110)	A4399

Assura® Hospital Version

Durable water bag with integrated thermometer, bold print, large opening. Can be used with 1- or 2-piece ostomy systems.

- **Not made with natural rubber latex**

Code	Units	Includes	HCPCS*
1500H	1	1- Water Bag/Regulator	A4398
		1- Cone (#1110)	A4399
		4- Adhesive Irrigation Sleeves (#1003)	A4397
		2- Closure Clips (#9500)	(effective 1/1/2022 this code will be A4437)



SenSura® or Assura® Deluxe Version

Components supplied in a travel bag.

- **Not made with natural rubber latex**

Code	Units	Includes	HCPCS*
12830	1	1- Water Bag/Regulator	A4398
		1- Cone (#1110)	A4399
		1- Irrigation Faceplate (#12820)	A4361
		1- Belt (#4215)	A4367
		2- Irrigation Sleeves (#12836)	A4397
			(effective 1/1/2022 this code will be A4436)

Irrigation Sleeves

Irrigation Sleeves (with adhesive)

Disposable irrigation sleeve (adhesive-backed).

- **Not made with natural rubber latex**

Transparent

Code	Units	HCPCS*
1003	10	A4397
		(effective 1/1/2022 this code will be A4437)






*For important information regarding the use of HCPCS codes listed and code definitions see page 180




**Irrigation Sleeves
(without adhesive)**



Irrigation sleeves for use with SenSura® Click or Assura® barriers or faceplate. Consult chart below to order the correct transparent sleeves.

• **Not made with natural rubber latex**

Code	Size	Units	HCPCS*
12834		5	A4397 (effective 1/1/2022 this code will be A4436)
12835		5	A4397 (effective 1/1/2022 this code will be A4436)
12836		5	A4397 (effective 1/1/2022 this code will be A4436)

Irrigation Sleeve #	Use with SenSura® Click Product Code	Units
12834	 - 10011, 10012, 10013, 10015-10017, 11011, 11012, 11015, 11016	5
12835	 - 10021, 10022, 10023, 10025-10027, 11021-11025	5
12836	 - 10031, 10032, 10035, 10036, 11031, 11032, 11035, 11036	5

Irrigation Accessories

Assura® Irrigation Faceplate

For use with ostomy belt 4215 or 4220.
(Use with #12836 Irrigation Sleeve) Note: not compatible with Brava® Belt for SenSura® Mio



• **Not made with natural rubber latex**

Code	Units	HCPCS*
12820	1	A4361

Coloplast Stoma Cone

Soft and flexible stoma cone with tubing.



• **Not made with natural rubber latex**

Code	Units	HCPCS*
1110	1	A4399



2-piece Click

(Mechanical coupling)



Pouch can be repositioned while it's still on the body



Tactile experience of pressing two pieces together can be reassuring



Lock ring gives an audible cue that it's securely locked

SenSura® Mio

The fit is
the difference
between sitting on the
edge and diving right in.

Mike is a SenSura® Mio Click user and has received compensation from Coloplast.

SenSura® Mio Click

BodyFit Technology is an elastic adhesive barrier found in all SenSura® Mio Barriers. This unique barrier quickly adheres to the skin, flows into creases and folds, is highly stretchable, and flexes with the body. This may help prevent channels and gaps and may result in reducing leakage and skin issues.



NOTE: SenSura® Mio Click is NOT compatible with SenSura® Click or Assura® systems.

SenSura® Mio Click Barrier

- Extended wear
- Flat
- Belt tabs for use with 4237, 4247
- **Not made with natural rubber latex**



Code	Stoma Size	Color	Units: 5 HCPCS*
Cut-to-fit			
10502	3/8-1 3/8" (10-35 mm)		A4409
10512	3/8-1 3/4" (10-45 mm)		A4409
10522	3/8-2 1/8" (10-55 mm)		A4410
10531	3/8-2 9/16" (10-65 mm)		A4410
Pre-cut			
10503	3/4" (20 mm)		A4409
10506	7/8" (22 mm)		A4409
10504	1" (25 mm)		A4409
10505	1 1/8" (28 mm)		A4409
10513	1 3/16" (30 mm)		A4409
10514	1 1/4" (32 mm)		A4409
10515	1 3/8" (35 mm)		A4409
10516	1 1/2" (38 mm)		A4409
10523	1 9/16" (40 mm)		A4409
10524	1 3/4" (45 mm)		A4409
10525	2" (50 mm)		A4409

SenSura® Mio Click Convex Light Barrier

- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Not made with natural rubber latex**



Code	Stoma Size	Color	Units: 5 HCPCS*
Cut-to-fit			
16901	3/8 - 3/4" (10-20 mm)		A4407
16911	5/8 - 1 3/16" (15-30 mm)		A4407
16921	5/8 - 1 9/16" (15-40 mm)		A4407
16931	5/8 - 2 1/16" (15-53 mm)		A4408
Pre-cut			
16903	1 1/16" (18 mm)		A4407
16913	1 3/16" (21 mm)		A4407
16914	1" (25 mm)		A4407
16915	1 1/8" (28 mm)		A4407
16923	1 1/4" (31 mm)		A4407
16924	1 3/8" (35 mm)		A4407
16933	1 1/2" (38 mm)		A4408
16934	1 5/8" (41 mm)		A4408













Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

SenSura® Mio Click Deep Convex Barrier 

- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Not made with natural rubber latex**































Code	Stoma Size	Color	Units: 5 HCPCS*
Cut-to-fit			
16941	3/8 - 3/4" (10-20 mm)		A4407
16951	5/8 - 1 1/16" (15-30 mm)		A4407
16961	5/8 - 1 9/16" (15-40 mm)		A4407
16971	5/8 - 2 1/16" (15-53 mm)		A4408
Pre-cut			
16942	1 1/16" (18 mm)		A4407
16953	1 3/16" (21 mm)		A4407
16954	1" (25 mm)		A4407
16955	1 1/8" (28 mm)		A4407
16963	1 1/4" (31 mm)		A4407
16964	1 3/8" (35 mm)		A4407
16973	1 1/2" (38 mm)		A4408
16974	1 5/8" (41 mm)		A4408

SenSura® Mio Click Convex Flip Barrier 

- Extended wear
- Designed for curves, bulges and hernias
- Curved, star-shaped barrier conforms to outward areas with less creasing and folding
- "Flip to fit" design allows for easy application
- Belt tabs for use with 4237, 4247
- **Not made with natural rubber latex**



Code	Stoma Size	Barrier Size	Color	Units: 5 HCPCS*	
Cut-to-fit					
18500	3/8 - 1 3/8" (10-35 mm)			A4407	
18510	3/8 - 1 3/4" (10-45 mm)			A4408	
18520	3/8 - 2 1/8" (10-55 mm)			A4408	
18530	3/8 - 2 9/16" (10-66 mm)			A4408	
Code	Stoma Size	Max Cut	Barrier Size	Color	HCPCS*
Pre-cut					
18501	3/4" (20 mm)	1 3/8" (35 mm)			A4407
18511	1" (25 mm)	1 3/4" (45 mm)			A4407
18512	1 3/16" (30 mm)	1 3/4" (45 mm)			A4407
18513	1 3/8" (35 mm)	1 3/4" (45 mm)			A4407
18521	1 3/16" (30 mm)	2 1/8" (55 mm)			A4407
18522	1 3/8" (35 mm)	2 1/8" (55 mm)			A4407
18523	1 9/16" (40 mm)	2 1/8" (55 mm)			A4407
18524	1 3/4" (45 mm)	2 1/8" (55 mm)			A4407
18531	1 3/4" (45 mm)	2 9/16" (66 mm)			A4408
18532	2" (50 mm)	2 9/16" (66 mm)			A4408

NOTE: SenSura® Mio Click is NOT compatible with SenSura® Click or Assura® systems.

SenSura® Mio Click MAXI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo inspection window for easy application
- **Not made with natural rubber latex**



11" (28 cm) / 655 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent			
Opaque			
11463	11462		A4427
11473	11472		A4427
11483	11482		A4427
11487	11486		A4427

SenSura® Mio Click MAXI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Without filter**
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**



11" (28 cm) / 655 mL

Units: 10

Code	Color	HCPCS*
Transparent		
11474		A4426
11484		A4426
11488		A4426

SenSura® Mio Click MIDI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window for easy application
- **Not made with natural rubber latex**



10¾" (26 cm) / 485 mL

Units: 10

Code	Color	HCPCS*
Opaque		
11461		A4427
11471		A4427

Indicates convexity Indicates Convex Flip





*For important information regarding the use of HCPCS codes listed and code definitions see page 180

SenSura® Mio Click High Output Pouch with Tap Outlet

- **Special full-circle filter designed specifically for ileostomy output**
- Neutral gray, water-resistant textile
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- **Not made with natural rubber latex**



12½" (32 cm) / 1080 mL Units: 10




Code	Code	Color	HCPCS*
Transparent			
18610	Opaque		A4413
18611	18601		A4413
18612	18602		A4413
18613	18603		A4413

SenSura® Mio Click High Output Pouch with Soft Outlet

- **Special full-circle filter designed specifically for ileostomy output**
- Neutral gray, water-resistant textile
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Soft outlet provides effective drainage; connects to a standard drain bag or Coloplast fistula drain bag #14010
- **Not made with natural rubber latex**



12 ½" (32 cm) / 1080 mL Units: 10

Code	Code	Color	HCPCS*
Transparent			
18640	Opaque		A4413
18641	18631		A4413
18642	18632		A4413




SenSura® Mio Click Drainable Pouch with Tap Outlet

- Special full-circle filter designed specifically for ileostomy output
- Neutral gray, water-resistant textile
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- **Not made with natural rubber latex**



9" (23 cm) / 650 mL

Units: 10

Code	Color	HCPCS*
Opaque		
18620		A4427
18621		A4427
18622		A4427





SenSura® Mio Click MAXI Closed Pouch

- Full-circle filter
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo inspection window for easy application
- **Not made with natural rubber latex**



8¼" (20½ cm) / 570 mL

Units: 30

Code	Code	Color	HCPCS*
Transparent			
Opaque			
N/A	11402		A4423
N/A	11413		A4423
N/A	11422		A4423
11426	11425		A4423




SenSura® Mio Click MIDI Closed Pouch

- Full-circle filter
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window for easy application
- **Not made with natural rubber latex**



7¼" (18½ cm) / 275 mL

Units: 30

Code	Color	HCPCS*
Opaque		
11401		A4423
11412		A4423
11421		A4423


*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Urostomy

SenSura® Mio Uro

Wide anti-reflux valves support varying ways of managing stents:

- Stents threaded into a pouch connected to a night drainage bag
- Rolling the stents in the top of the pouch, making a loop
- Stents that have been cut



Welded "multi-chambers" minimize ballooning




SenSura® Mio Click MAXI Urostomy Pouch

- Multi-chambers minimize sloshing
- Wider anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**



10½" (26½ cm) / 480 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
11492	11491		A4433
11498	11497		A4433
11500	11499		A4433

NOTE: SenSura® Mio Click is NOT compatible with SenSura® Click or Assura® systems. SenSura® Click and Assura® are compatible.

SenSura® Click Flat Barrier


- Standard wear
- Belt tabs for use with 4215, 4220
- **Not made with natural rubber latex**








Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit

10011	3/8 - 1 3/8" (10-35 mm)		A4414
10021	3/8 - 1 3/4" (10-45 mm)		A4414
10031	3/8 - 2 1/4" (10-55 mm)		A4414
10041	3/8 - 2 1/2" (10-65 mm)		A4415**

Pre-cut

10012	1 3/16" (20 mm)		A4414
10013	1" (25 mm)		A4414
10022	1 3/16" (30 mm)		A4414
10023	1 3/8" (35 mm)		A4414
10032	1 9/16" (40 mm)		A4414

SenSura® Click Xpro Flat Barrier


- Extended wear
- Belt tabs
- **Not made with natural rubber latex**






Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit

10015	3/8 - 1 3/8" (10-35 mm)		A4409
10025	3/8 - 1 3/4" (10-45 mm)		A4409
10035	3/8 - 2 1/4" (10-55 mm)		A4409
10045	3/8 - 2 1/2" (10-65 mm)		A4410

Pre-cut

10016	1 3/16" (20 mm)		A4409
10017	1" (25 mm)		A4409
10026	1 3/16" (30 mm)		A4409
10027	1 3/8" (35 mm)		A4409
10036	1 9/16" (40 mm)		A4409

*For important information regarding the use of HCPCS codes listed and code definitions see page 180










**Product has not been through PDAC verification

SenSura® Click Convex Light Barrier 

- Standard wear
- Belt tabs
- Not made with natural rubber latex



Units: 5










Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
11011	5/8 - 7/8" (15-22 mm)		A4373
11021	5/8 - 1 5/16" (15-33 mm)		A4373
11031	5/8 - 1 11/16" (15-43 mm)		A4373
11041	5/8 - 2 1/16" (15-53 mm)		A4373**
Pre-cut			
11012	7/8" (22 mm)		A4373
11022	1" (25 mm)		A4373
11023	1 1/8" (28 mm)		A4373
11024	1 1/4" (31 mm)		A4373
11032	1 3/8" (35 mm)		A4373

SenSura® Click Xpro Convex Light Barrier 

- Extended wear
- Belt tabs
- Not made with natural rubber latex



Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
11015	5/8 - 7/8" (15-22 mm)		A4407
11025	5/8 - 1 5/16" (15-33 mm)		A4407
11035	5/8 - 1 11/16" (15-43 mm)		A4407
11045	5/8 - 2 1/16" (15-53 mm)		A4408
Pre-cut			
11016	7/8" (22 mm)		A4407
11026	1" (25 mm)		A4407
11027	1 1/8" (28 mm)		A4407
11028	1 1/4" (31 mm)		A4407
11036	1 3/8" (35 mm)		A4407

For Assura® CONVEX barriers, see pages 89-90.





SenSura® Click MAXI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11½" (29½ cm) / 655 mL

Units: 20

Code	Code	Color	HCPCS*
Transparent	Opaque		
11194	11124		A4427
11195	11125		A4427
11196	11126		A4427
11197	11127		A4427**



SenSura® Click MIDI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**



10½" (26½ cm) / 485 mL

Units: 20

Code	Color	HCPCS*
Opaque		
11114		A4427
11115		A4427



SenSura® Click MAXI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Without filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11½" (29½ cm) / 655 mL

Units: 20

Code	Code	Color	HCPCS*
Transparent	Opaque		
11185	11135		A4426
11186	11136		A4426

*For important information regarding the use of HCPCS codes listed and code definitions see page 180




**Product has not been through PDAC verification

NOTE: SenSura® Click is compatible with Assura® but NOT with SenSura® Mio Click.

SenSura® Click Magnum High Output Drainable Pouch

- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**







12½" (32 cm) / 1160 mL			Units: 10
Code	Code	Color	HPCPS*
Transparent	Opaque		
19050	19040		A4413**
19051	19041		A4413**
19052	19042		A4413**

For Assura® Ileo Night Bag, see page 94.

SenSura® Click MAXI Closed Pouch

- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**






8½" (21 cm) / 570 mL			Units: 30
Code	Code	Color	HPCPS*
Transparent	Opaque		
NA	10164		A4423
NA	10165		A4423
NA	10166		A4423
10187	10167		A4423**

SenSura® Click MIDI Closed Pouch

- **Filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**



7" (18 cm) / 275 mL		Units: 30
Code	Color	HPCPS*
Opaque		
10154		A4423
10155		A4423
10156		A4423

For Assura® stoma caps, see page 95.

Urostomy



SenSura® Click Uro MAXI Urostomy Pouch

- Multi-chambers minimize sloshing
- Wider anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



10³/₈" (26 cm) / 480 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
11854	11844		A4433
11855	11845		A4433
11856	11846		A4433


SenSura® Click Uro MIDI Urostomy Pouch

- Multi-chambers minimize sloshing
- Wider anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



9¹/₂" (24 cm) / 430 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
11851	11841		A4433

For Assura® urostomy micro-pouches and minicaps, see page 95.

Urostomy Night Bag and Bag Hanger

See Accessories page 72 for more information.



NOTE: Assura® is compatible with SenSura® Click , but **NOT** with SenSura® Mio Click.

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

**Product has not been through PDAC verification

Assura® Flat Barrier

- Standard wear
- Belt tabs for use with 4215, 4220
- **Not made with natural rubber latex**



Code	Stoma Size	Color	Units: 5 HCPCS*
Cut-to-fit			
2881	3/8 - 1 3/8" (10-35 mm)		A4414
2882	3/8 - 1 3/4" (10-45 mm)		A4414
2883	3/8 - 2 1/4" (10-55 mm)		A4414
Pre-cut			
12841	1 3/16" (20 mm)		A4414
12842	1" (25 mm)		A4414
12844	1 3/16" (30 mm)		A4414
12845	1 3/8" (35 mm)		A4414
12847	1 9/16" (40 mm)		A4414
12848	1 3/4" (45 mm)		A4414

Assura® Flat Barrier

- Extended wear
- Belt tabs for use with 4215, 4220
- **Not made with natural rubber latex**



Code	Stoma Size	Color	Units: 5 HCPCS*
Cut-to-fit			
2831	3/8 - 1 3/8" (10-35 mm)		A4409
2832	3/8 - 1 3/4" (10-45 mm)		A4409
2833	3/8 - 2 1/4" (10-55 mm)		A4409
Pre-cut			
14250	5/8" (15 mm)		A4409
14251	1 3/16" (20 mm)		A4409
14252	1" (25 mm)		A4409
14254	1 3/16" (30 mm)		A4409
14255	1 3/8" (35 mm)		A4409
14257	1 9/16" (40 mm)		A4409
14258	1 3/4" (45 mm)		A4409












NOTE: Assura® is compatible with SenSura® Click, but NOT with SenSura® Mio Click.

Assura® Convex Light Barrier

- Standard wear
- Belt tabs for use with 4215, 4220
- **Not made with natural rubber latex**



Units: 5












Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
14261	5/8 - 15/16" (15-23 mm)		A4373
14262	5/8 - 15/16" (15-33 mm)		A4373
14263	5/8 - 111/16" (15-43 mm)		A4373
Pre-cut			
14271	3/4" (18 mm)		A4373
14272	7/8" (22 mm)		A4373
14273	1" (25 mm)		A4373
14274	1 1/8" (28 mm)		A4373
14275	1 1/4" (31 mm)		A4373
14276	1 3/8" (35 mm)		A4373
14277	1 1/2" (38 mm)		A4373
14278	1 5/8" (41 mm)		A4373

Assura® Convex Light Barrier

- Extended wear
- Belt tabs for use with 4215, 4220
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
14281	5/8 - 15/16" (15-23 mm)		A4407
14282	5/8 - 15/16" (15-33 mm)		A4407
14283	5/8 - 111/16" (15-43 mm)		A4407
Pre-cut			
14291	3/4" (18 mm)		A4407
14292	7/8" (22 mm)		A4407
14293	1" (25 mm)		A4407
14294	1 1/8" (28 mm)		A4407
14295	1 1/4" (31 mm)		A4407
14296	1 3/8" (35 mm)		A4407
14297	1 1/2" (38 mm)		A4407
14298	1 5/8" (41 mm)		A4407

 Indicates convexity












*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Assura® Deep Convex Barrier 

- Standard wear
- Belt tabs for use with 4215, 4220
- **Not made with natural rubber latex**



Units: 5












Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
12716	5/8" - 1 ⁵ / ₁₆ " (15-33 mm)		A4373
12719	5/8" - 1 ¹¹ / ₁₆ " (15-43 mm)		A4373
Pre-cut			
12701	5/8" (15 mm)		A4373
12702	3/4" (18 mm)		A4373
12703	7/8" (22 mm)		A4373
12704	1" (25 mm)		A4373
12705	1 ¹ / ₈ " (28 mm)		A4373
12706	1 ¹ / ₄ " (31 mm)		A4373
12707	1 ³ / ₈ " (35 mm)		A4373
12708	1 ¹ / ₂ " (38 mm)		A4373
12709	1 ⁵ / ₈ " (41 mm)		A4373

Assura® Deep Convex Barrier 

- Extended wear
- Belt tabs for use with 4215, 4220
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
14243	5/8" - 1 ⁵ / ₁₆ " (15-23 mm)		A4407
14246	5/8" - 1 ⁵ / ₁₆ " (15-33 mm)		A4407
14249	5/8" - 1 ¹¹ / ₁₆ " (15-43 mm)		A4407
Pre-cut			
14232	3/4" (18 mm)		A4407
14233	7/8" (22 mm)		A4407
14234	1" (25 mm)		A4407
14235	1 ¹ / ₈ " (28 mm)		A4407
14236	1 ¹ / ₄ " (31 mm)		A4407
14237	1 ³ / ₈ " (35 mm)		A4407
14238	1 ¹ / ₂ " (38 mm)		A4407
14239	1 ⁵ / ₈ " (41 mm)		A4407

 Indicates convexity




NOTE: Assura® is compatible with SenSura® Click, but NOT with SenSura® Mio Click.

Assura® MAXI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11½" (29½ cm) / 655 mL Units: 10



Code	Color	HCPCS*
Transparent	Opaque	
15974	15984 	A4427
15975	15985 	A4427
15976	15986 	A4427

Assura® MAXI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Without filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**



11½" (29½ cm) / 655 mL Units: 10




Code	Color	HCPCS*
Opaque		
15982		A4426
15983		A4426

Assura® MIDI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**



10½" (26½ cm) / 485 mL Units: 10

Code	Color	HCPCS*
Opaque		
15954		A4427
15955		A4427
15956		A4427




*For important information regarding the use of HCPCS codes listed and code definitions see page 182

NOTE: Assura® is compatible with SenSura® Click, but NOT with SenSura® Mio Click.

Assura® MAXI Drainable Pouch

- EasiClose™ Outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**






11½" (29½ cm) / 655 mL			Units: 10
Code	Code	Color	HCPCS*
Transparent	Opaque		
13974	13984		A4427
13975	13985		A4427
13976	13986		A4427

Assura® MIDI Drainable Pouch

- EasiClose™ Outlet
- **Filter**
- Soft cloth on front and back
- Graphic pattern
- **Not made with natural rubber latex**





10½" (26½ cm) / 485 mL			Units: 10
Code		Color	HCPCS*
Opaque			
13964			A4427
13965			A4427
13966			A4427

Assura® MINI Drainable Pouch

- EasiClose™ Outlet
- **Filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**



9¾" (23½ cm) / 280 mL			Units: 10
Code		Color	HCPCS*
Opaque			
13924			A4427
13925			A4427

NOTE: Assura® is compatible with SenSura® Click, but NOT with SenSura® Mio Click.




Assura® MAXI Drainable Pouch

- Clamp closure
- **Without filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**



11¼" (28½ cm) / 760 mL

Units: 10

Code	Color	HCPCS*
Opaque		
12577		A4426
12578		A4426
12579		A4426




Assura® MIDI Drainable Pouch

- Clamp closure
- **Without filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**



9¾" (24½ cm) / 400 mL

Units: 10

Code	Color	HCPCS*
Opaque		
12527		A4426
12528		A4426
12529		A4426




Assura® High Output Drainable Pouch

- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



12½" (32 cm) / 750 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent			
Opaque			
8114	NA		A4413
8115	2846		A4413
8116	2847		A4413



For bed drainage bag and bag hanger, see page 72.

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Assura® Ileo Night Bag

- 2 liter capacity for very high output stomas
- Clamp closure
- Comes with bag hanger
- **Filter**
- **Not made with natural rubber latex**






2 L			Units: 5
Code		Color	HCPCS*
Opaque			
2836			A4413
2837			A4413

Assura® MAXI Closed Pouch

- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**






8½" (21 cm) / 570 mL			Units: 30
Code	Code	Color	HCPCS*
Transparent			
Opaque			
12374	12384		A4423
12375	12385		A4423
12376	12386		A4423

Assura® MIDI Closed Pouch

- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



7" (18 cm) / 275 mL			Units: 30
Code	Code	Color	HCPCS*
Transparent			
Opaque			
12344	12354		A4423
12345	12355		A4423
12346	12356		A4423

NOTE: Assura® is compatible with SenSura® Click, but NOT with SenSura® Mio Click.




Ileostomy / Colostomy

Assura® Stoma Cap

- **Filter**
- Small size is optimal for intimate moments or colostomates who irrigate
- **Not made with natural rubber latex**



Units: 30

Code	Color	HCPCS*
Opaque		
2801		A5055
2802		A5055
2803		A5055**

Urostomy




Assura® Urostomy Micro-Pouch

- Miniature size prevents twisting
- Soft outlet can be connected to a drain bag for nighttime use, or a leg bag for patients in a wheelchair
- Soft cloth on front
- **Not made with natural rubber latex**



6" (15 cm)

Units: 10

Code	Color	HCPCS*
Transparent		
14204		A4433
14205		A4433
14206		A4433

Try a leg bag for discreet day wear — available on page 26-27. For nighttime, use a urostomy night bag found on page 72.



Assura® Urostomy Minicap

- Ultra-absorbent material inside
- Small size is optimal for intimate use
- **Not made with natural rubber latex**



100 mL

Units: 30

Code	Color	HCPCS*
Opaque		
2807		A5055
2808		A5055

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

**Product has not been through PDAC verification

NOTE: Assura® is compatible with SenSura® Click, but NOT with SenSura® Mio Click.




Assura® MAXI Urostomy Pouch

- Multi-chambers minimize sloshing
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



10½" (26½ cm) / 488 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent			
Opaque			
14227	14224		A4433
14228	14225		A4433
14229	14226		A4433




Assura® MAXI Urostomy Pouch

- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



10¾" (27¼ cm) / 375 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent			
Opaque			
1757	1754		A4433
1758	1755		A4433
1759	1756		A4433



Assura® MIDI Urostomy Pouch

- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back
- **Not made with natural rubber latex**



9½" (24 cm) / 150 mL

Units: 10

Code	Color	HCPCS*
Transparent		
1751		A4433
1752		A4433

Urostomy Night Bag and Bag Hanger

See Accessories page 72 for more information.





2-piece Flex

(Adhesive coupling)



Flexible system bends with body movements



The low profile of a 1-piece – the convenience of a 2-piece



Floating flange = minimal pressure to apply



Built-in security channel + strong adhesive = strong leakage reduction

SenSura® Mio

The fit is
the difference
between staying in
and working out

Deirdre is a SenSura® Mio user.

SenSura® Mio Flex

The BodyFit Technology in all SenSura Mio barriers combines an elastic adhesive, that bends and stretches with the body, with an innovative barrier designed to provide a secure and comfortable fit.



2-Piece Flex

NOTE: SenSura® Mio Flex and SenSura® Flex are compatible with one another

SenSura® Mio Flex Flat Barrier

- Extended wear
- **Not made with natural rubber latex**



Units: 5
18390 - 10 units

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
10551	3/8 - 1 5/16" (10-33 mm)		A4409
10561	3/8 - 1 7/8" (10-48 mm)		A4410
10571	3/8 - 2 11/16" (10-68 mm)		A4410

Coming soon!

18390	3/8 - 3 1/2" (10-88mm)		A4410**
-------	------------------------	--	---------

Pre-cut

10562	1" (25 mm)		A4409
10566	1 1/8" (28 mm)		A4409
10563	1 3/16" (30 mm)		A4409
10567	1 1/4" (32 mm)		A4409
10564	1 3/8" (35 mm)		A4409
10568	1 1/2" (38 mm)		A4409
10565	1 9/16" (40 mm)		A4409
10572	1 3/4" (45 mm)		A4410
10573	2" (50 mm)		A4410

SenSura® Mio Flex Convex Light Barrier

- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
16471	3/8 - 3/4" (10-20 mm)		A4407
16481	5/8 - 1 3/16" (15-30 mm)		A4407
16483	5/8 - 1 9/16" (15-40 mm)		A4407
16491	5/8 - 2 1/16" (15-53 mm)		A4408

Pre-cut

16472	1 1/16" (18 mm)		A4407
16484	1 3/16" (21 mm)		A4407
16485	1" (25 mm)		A4407
16486	1 1/8" (28 mm)		A4407
16487	1 1/4" (31 mm)		A4407
16488	1 3/8" (35 mm)		A4407
16493	1 1/2" (38 mm)		A4408
16494	1 5/8" (41 mm)		A4408

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

**Product has not been through PDAC verification

NOTE: SenSura® Mio pouches **are compatible** with SenSura Xpro extended wear barriers





SenSura® Mio Flex MAXI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo inspection window for easy application
- **Not made with natural rubber latex**



11" (28 cm) / 655 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent			
NA	Opaque		
12262	12262		A4425
12273	12272		A4425
12283	12282		A4425
Coming soon!			
18389	NA		A4425**



SenSura® Mio Flex MAXI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Without filter**
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**



11" (28 cm) / 655 mL

Units: 10

Code	Color	HCPCS*
Transparent		
12274		A5063
12284		A5063


SenSura® Mio Flex MIDI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window for easy application
- **Not made with natural rubber latex**



10¼" (26 cm) / 485 mL

Units: 10

Code	Color	HCPCS*
Opaque		
12281		A4425



SenSura® Mio Flex High Output Pouch with Tap Outlet

- **Special full-circle filter designed specifically for ileostomy output**
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- **Not made with natural rubber latex**



12½" (32 cm) / 1080 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
18655	18650		A4413
18656	18651		A4413



SenSura® Mio Flex High Output Pouch with Soft Outlet

- **Special full-circle filter designed specifically for ileostomy output**
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- **Not made with natural rubber latex**



12½" (32 cm) / 1080 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
18657	18652		A4413
18658	18653		A4413

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

**Product has not been through PDAC verification

NOTE: SenSura® Mio Flex and SenSura® Flex are compatible with one another




SenSura® Mio Flex MAXI Closed Pouch

- Full-circle filter
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo inspection window for easy application
- Not made with natural rubber latex



8¼" (20½ cm) / 570 mL

Units: 30

Code	Code	Color	HCPCS*
Transparent			
Opaque			
12203	12202		A4419
12216	12213		A4419
12223	12222		A4419



SenSura® Mio Flex MAXI Urostomy Pouch

- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- Not made with natural rubber latex



10½" (27 cm) / 400 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent			
Opaque			
12292	12291		A5073**
12298	12297		A5073**

SenSura® Mio Flex MIDI Urostomy Pouch

- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- Not made with natural rubber latex



9½" (24 cm) / 330 mL

Units: 10

Code	Color	HCPCS*
Opaque		
12296		A5073**





NOTE: SenSura® Mio Flex and SenSura® Flex are compatible with one another

SenSura® Flex Flat Barrier

- Standard wear
- **Not made with natural rubber latex**



Units: 10
10108** - 5 units




Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
10101	3/8 - 1 5/16" (10-33 mm)		A4414
10102	3/8 - 1 7/8" (10-48 mm)		A4414
10103	3/8 - 2 3/4" (10-68 mm)		A4415
10108**	3/8 - 3 1/2" (10-88 mm)		A4415

SenSura® Flex Xpro Flat Barrier

- Extended wear
- Compatible with SenSura Mio pouches
- **Not made with natural rubber latex**



Units: 10





Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
10104	3/8 - 1 5/16" (10-33 mm)		A4409
10105	3/8 - 1 7/8" (10-48 mm)		A4409
10106	3/8 - 2 3/4" (10-68 mm)		A4410

SenSura® Flex Convex Light Barrier

- Standard wear
- Belt tabs for use with 4215, 4220
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
11301	5/8 - 1 5/16" (15-23 mm)		A4373
11302	5/8 - 1 5/16" (15-33 mm)		A4373
11303	5/8 - 1 11/16" (15-43 mm)		A4373
11307	5/8 - 2 1/16" (15-53 mm)		A4373

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 180





**Product has not been through PDAC verification

SenSura® Flex Convex Light Barrier 

- Standard wear
- Belt tabs for use with 4215, 4220
- **Not made with natural rubber latex**



Units: 5





Code	Stoma Size	Color	HCPCS*
Pre-cut			
11312	7/8" (22 mm)		A4373
11313	1" (25 mm)		A4373
11314	1 1/8" (28 mm)		A4373
11315	1 1/4" (31 mm)		A4373

SenSura® Flex Xpro Convex Light Barrier 

- Extended wear
- Compatible with SenSura Mio pouches
- Belt tabs for use with 4215, 4220
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
11304	5/8 - 15/16" (15-23 mm)		A4407
11305	5/8 - 1 1/16" (15-33 mm)		A4407
11306	5/8 - 1 11/16" (15-43 mm)		A4407
11308	5/8 - 2 1/16" (15-53 mm)		A4408





SenSura® Flex MAXI Drainable Pouch

- EasiClose™ WIDE Outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11 1/2" (29 1/2 cm) / 655 mL

Units: 20

Code	Code	Color	HCPCS*
Transparent			
NA	11514		A4425
11512	11515		A4425
11517	11518		A4425
11 1/2" (29 1/2 cm) / 750 mL XXL pouch			
11590	11591		A4425

 Indicates convexity

NOTE: SenSura® Mio Flex and SenSura® Flex are compatible with one another


SenSura® Flex MAXI Drainable Pouch

- EasiClose™ WIDE Outlet
- Without filter
- Soft cloth on back
- Not made with natural rubber latex



11½" (29½ cm) / 655 mL

Units: 20

Code	Color	HCPCS*
Transparent		
11612		A5063
11½" (29½cm) / 750 mL XXL pouch		
11690		A5063




SenSura® Flex MAXI Closed Pouch

- Filter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex



8½" (21 cm) / 570 mL

Units: 30

Code	Code	Color	HCPCS*
Transparent		Opaque	
10911	10914		A4419
10912	10915		A4419
10917	10918		A4419

Note: Assura AC (excluding pediatrics) has been discontinued. For a full list of discontinued products and suggested replacement products please refer to the table starting on page 171.

*For important information regarding the use of HCPCS codes listed and code definitions see page 180



1-piece



The choice for flexibility:
bends with body movements



No coupling = pouch and
barrier are welded together
and apply as one piece



Low-profile

SenSura[®] Mio

The fit is
the difference
between asking for help
and doing it yourself

Walton is a SenSura[®] Mio Convex Flip user and has received compensation from Coloplast.

SenSura[®] Mio 1-piece Pouching System

The BodyFit Technology in all SenSura Mio barriers combines an elastic adhesive, that bends and stretches with the body, with an innovative barrier designed to provide a secure and comfortable fit.



1-piece

SenSura® Mio Flat MAXI Drainable Pouch

- Extended wear
- EasiClose™ WIDE Outlet
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo inspection window for easy application
- **Not made with natural rubber latex**



11" (28 cm) / 655 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

10481	10471	Cut-to-fit	$\frac{3}{8}$ - $2\frac{1}{8}$ " (10-55 mm)	A5056
10482	10472	Pre-cut	1" (25 mm)	A5056
10488	10473	Pre-cut	$1\frac{1}{8}$ " (28 mm)	A5056
10483	10474	Pre-cut	$1\frac{3}{16}$ " (30 mm)	A5056
10490	10475	Pre-cut	$1\frac{1}{4}$ " (32 mm)	A5056
10484	10476	Pre-cut	$1\frac{3}{8}$ " (35 mm)	A5056
10491	10480	Pre-cut	$1\frac{1}{2}$ " (38 mm)	A5056
10487	10479	Pre-cut	$1\frac{9}{16}$ " (40 mm)	A5056
10485	10477	Pre-cut	$1\frac{3}{4}$ " (45 mm)	A5056
10486	10478	Pre-cut	2" (50 mm)	A5056

SenSura® Mio Flat MIDI Drainable Pouch

- Extended wear
- EasiClose™ WIDE Outlet
- **Full-circle filter**
- Peek-a-boo inspection window allows easy application
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**

**Opaque** 10 $\frac{1}{4}$ " (26 cm) / 485 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

10461	Cut-to-fit	$\frac{3}{8}$ - $1\frac{3}{4}$ " (10-45 mm)	A5056
-------	------------	---	-------

SenSura® Mio Flat MAXI Drainable Pouch

- Extended wear
- EasiClose™ WIDE Outlet
- **Without filter**
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**

**Transparent** 11" (28 cm) / 655 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
10489	Cut-to-fit	$\frac{3}{8}$ - $2\frac{1}{8}$ " (10-55 mm)	A4388
10492	Pre-cut	1" (25 mm)	A4388
10493	Pre-cut	$1\frac{1}{8}$ " (28 mm)	A4388
10494	Pre-cut	$1\frac{1}{4}$ " (32 mm)	A4388
10495	Pre-cut	$1\frac{3}{8}$ " (35 mm)	A4388
10496	Pre-cut	$1\frac{1}{2}$ " (38 mm)	A4388

SenSura® Mio Soft Convex MAXI Drainable Pouch

- Extended wear
- Belt tabs for use with 4237, 4247
- EasiClose™ WIDE Outlet
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**



11" (28 cm) / 655 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent				
16715	16705	Cut-to-fit	$\frac{5}{8}$ - $1\frac{5}{16}$ " (15-33 mm)	A5057
16716	16706	Cut-to-fit	$\frac{3}{8}$ - 2" (10-50 mm)	A5057

Coming soon!

NA	13671	Pre-cut	$\frac{7}{8}$ " (21 mm)	A5057**
NA	13672	Pre-cut	1" (25 mm)	A5057**
NA	13673	Pre-cut	$1\frac{1}{8}$ " (28 mm)	A5057**
NA	13674	Pre-cut	$1\frac{1}{4}$ " (31 mm)	A5057**
NA	13675	Pre-cut	$1\frac{3}{8}$ " (35 mm)	A5057**
NA	13676	Pre-cut	$1\frac{1}{2}$ " (38 mm)	A5057**

 Indicates convexity

SenSura® Mio Convex Light MAXI Drainable Pouch

- Extended wear
- Belt tabs for use with 4237, 4247
- EasiClose™ WIDE Outlet
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**



11" (28 cm) / 655 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
16737	16725	Cut-to-fit	$\frac{3}{8}$ - $\frac{7}{8}$ " (10-22 mm)	A5057
16738	16726	Cut-to-fit	$\frac{5}{8}$ - $1\frac{5}{16}$ " (15-33 mm)	A5057
16739	16733	Cut-to-fit	$\frac{3}{8}$ - $1\frac{11}{16}$ " (10-43 mm)	A5057
16740	16730	Pre-cut	$1\frac{3}{16}$ " (21 mm)	A5057
16741	16731	Pre-cut	1" (25 mm)	A5057
16742	16732	Pre-cut	$1\frac{1}{8}$ " (28 mm)	A5057
NA	16727	Pre-cut	$1\frac{1}{4}$ " (31 mm)	A5057
NA	16734	Pre-cut	$1\frac{3}{8}$ " (35 mm)	A5057
NA	16735	Pre-cut	$1\frac{1}{2}$ " (38 mm)	A5057

SenSura® Mio Convex Light MAXI Drainable Pouch

- Extended wear
- Belt tabs for use with 4237, 4247
- EasiClose™ WIDE Outlet
- FlexShell is designed for stability and flexibility
- **Without filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**



Transparent 11" (28 cm) / 655 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
16748	Cut-to-fit	$\frac{3}{8}$ - $1\frac{11}{16}$ " (10-43 mm)	A4390

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

**Product has not been through PDAC verification

SenSura® Mio Deep Convex MAXI Drainable Pouch

- Extended wear
- Belt tabs for use with 4237, 4247
- EasiClose™ WIDE Outlet
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**



11" (28cm) / 655 mL Units: 10

Code	Code	Cut	Stoma Size	HPCS*
Transparent Opaque				
16765	16755	Cut-to-fit	3/8 - 7/8" (10-22 mm)	A5057
16766	16756	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A5057
16767	16768	Cut-to-fit	3/8 - 1 11/16" (10-43 mm)	A5057
16770	16760	Pre-cut	1 3/16" (21 mm)	A5057
16771	16761	Pre-cut	1" (25 mm)	A5057
16772	16762	Pre-cut	1 1/8" (28 mm)	A5057
16773	16757	Pre-cut	1 1/4" (31 mm)	A5057
NA	16763	Pre-cut	1 3/8" (35 mm)	A5057
NA	16764	Pre-cut	1 1/2" (38 mm)	A5057

SenSura® Mio Deep Convex MAXI Drainable Pouch

- Extended wear
- Belt tabs for use with 4237, 4247
- EasiClose™ WIDE Outlet
- FlexShell is designed for stability and flexibility
- **Without filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**



Transparent 11" (28 cm) / 655 mL Units: 10

Code	Cut	Stoma Size	HPCS*
16778	Cut-to-fit	3/8 - 1 11/16" (10-43 mm)	A4390

Indicates convexity

SenSura® Mio MAXI Drainable Pouch with Soft Outlet

- Extended wear
- **Special full-circle filter designed specifically for ileostomy output**
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or most standard night bags
- **Not made with natural rubber latex**



11" (28 cm) / 630 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

18665	18660	Cut-to-fit	$\frac{3}{8}$ - 2 $\frac{3}{4}$ " (10-70 mm)	A5056
18666	18661	Cut-to-fit	$\frac{3}{8}$ - 4" (10-100 mm)	A5056

SenSura® Mio Soft Convex MAXI Drainable Pouch with Soft Outlet

- Extended wear
- Belt tabs for use with 4237, 4247
- **Special full-circle filter designed specifically for ileostomy output**
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or any standard night bag
- **Not made with natural rubber latex**



11" (28 cm) / 630 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

18671	18670	Cut-to-fit	$\frac{3}{8}$ - 2" (10-50 mm)	A5057
-------	-------	------------	-------------------------------	-------

SenSura® Mio Convex Light MAXI Drainable Pouch with Soft Outlet

- Extended wear
- Belt tabs for use with 4237, 4247
- **Special full-circle filter designed specifically for ileostomy output**
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or any standard night bag
- **Not made with natural rubber latex**



11" (28 cm) / 630 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

18676	18675	Cut-to-fit	$\frac{3}{8}$ - 1 $\frac{11}{16}$ " (10-43 mm)	A5057
-------	-------	------------	--	-------

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 180












SenSura® Mio Convex Flip MAXI Drainable Pouch

- Extended wear
- Designed for curves, bulges and hernias
- Curved, star-shaped barrier conforms to outward areas with less creasing and folding
- "Flip to fit" design allows for easy application
- Opaque pouches feature a peek-a-boo inspection window for easy application
- **Full-circle filter**
- **Not made with natural rubber latex**



11" (28 cm) / 655 mL

Units: 10

Code	Code	Cut	Stoma Size	Max Cut	Barrier Size	HCPCS*
Transparent Opaque						
18305	NA	Cut-to-fit	3/8 - 1 9/16" (10-40 mm)	1 9/16" (40 mm)		A5057
18321	18312	Cut-to-fit	3/8 - 2" (10-50 mm)	2" (50 mm)		A5057
18341	18330	Cut-to-fit	3/8 - 2 3/8" (10-60 mm)	2 3/8" (60 mm)		A5057
NA	18301	Pre-cut	1" (25 mm)	1 9/16" (40 mm)		A5057
NA	18302	Pre-cut	1 3/16" (30 mm)	1 9/16" (40 mm)		A5057
18322	18313	Pre-cut	1" (25 mm)	2" (50 mm)		A5057
18323	18314	Pre-cut	1 3/16" (30 mm)	2" (50 mm)		A5057
18324	18315	Pre-cut	1 3/8" (35 mm)	2" (50 mm)		A5057
18342	18331	Pre-cut	1 9/16" (40 mm)	2 3/8" (60 mm)		A5057
18343	18332	Pre-cut	1 3/4" (45 mm)	2 3/8" (60 mm)		A5057
18344	18333	Pre-cut	2" (50 mm)	2 3/8" (60 mm)		A5057



SenSura® Mio Convex Flip MAXI Drainable Pouch

- Extended wear
- **Without filter**
- Designed for curves, bulges and hernias
- Curved, star-shaped barrier conforms to outward areas with less creasing and folding
- "Flip to fit" design allows for easy application
- **Not made with natural rubber latex**



Transparent 11" (28 cm) / 655 mL

Units: 10

Code	Cut	Stoma Size	Max Cut	Barrier Size	HCPCS*
18320	Cut-to-fit	3/8 - 2" (10-50 mm)	2" (50 mm)		A4390
18340	Cut-to-fit	3/8 - 2 3/8" (10-60 mm)	2 3/8" (60 mm)		A4390

 Indicates Convex Flip

SenSura® Mio Flat MAXI Closed Pouch

- Full-circle filter
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window allows easy application
- Not made with natural rubber latex

**Opaque** 8 ¼" (21 cm) / 550 mL

Units: 30

Code	Cut	Stoma Size	HCPCS*
10881	Cut-to-fit	5/8 - 1 3/4" (15-45 mm)	A4416
10885***	Cut-to-fit	1 9/16 - 2 1/8" (40-55 mm)	A4416
10883	Pre-cut	1 3/16" (30 mm)	A4416
10884	Pre-cut	1 3/8" (35 mm)	A4416
10885***	Pre-cut	1 9/16" (40 mm)	A4416

***Because of its larger starter hole, 10885 functions as both a pre-cut and cut-to-fit solution.

SenSura® Mio Flat MIDI Closed Pouch

- Full-circle filter
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window allows easy application
- Not made with natural rubber latex

**Opaque** 7 ¼" (19 cm) / 340 mL

Units: 30

Code	Cut	Stoma Size	HCPCS*
10871	Cut-to-fit	5/8 - 1 3/4" (15-45 mm)	A4416
10872	Pre-cut	1" (25 mm)	A4416
10873	Pre-cut	1 3/16" (30 mm)	A4416
10874	Pre-cut	1 3/8" (35 mm)	A4416

SenSura® Mio MINI Closed Pouch

- Full-circle filter
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window allows easy application
- Comfortable and discreet stoma cap alternative
- Not made with natural rubber latex

**Opaque** 5 ¼" (13 cm) / 110 mL

Units: 30

Code	Cut	Stoma Size	HCPCS*
10870	Cut-to-fit	5/8 - 1 3/4" (15-45 mm)	A4416

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

SenSura® Mio Soft Convex MAXI Closed Pouch 

- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**



8¾" (21 cm) / 550 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
NA	16305	Cut-to-fit	5/8 - 1 ⁵ / ₁₆ " (15-33 mm)	A4417
16315	16306	Cut-to-fit	3/8 - 2" (10-50 mm)	A4417
NA	16310	Pre-cut	1" (25 mm)	A4417
NA	16311	Pre-cut	1 ¹ / ₈ " (28 mm)	A4417
NA	16312	Pre-cut	1 ¹ / ₄ " (31 mm)	A4417
NA	16313	Pre-cut	1 ³ / ₈ " (35 mm)	A4417

SenSura® Mio Soft Convex MIDI Closed Pouch 

- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**

**Opaque** 7¼" (19 cm) / 340 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
16301	Cut-to-fit	5/8 - 1 ⁵ / ₁₆ " (15-33 mm)	A4417

SenSura® Mio Convex Light MAXI Closed Pouch 

- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**



8¾" (21 cm) / 550 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
NA	16325	Cut-to-fit	3/8 - 7/8" (10-22 mm)	A4417
16335	NA	Cut-to-fit	3/8 - 1 ¹¹ / ₁₆ " (10-43 mm)	A4417
NA	16326	Cut-to-fit	5/8 - 1 ⁵ / ₁₆ " (15-33 mm)	A4417
NA	16330	Pre-cut	1" (25 mm)	A4417
NA	16331	Pre-cut	1 ¹ / ₈ " (28 mm)	A4417
NA	16327	Pre-cut	1 ¹ / ₄ " (31 mm)	A4417
NA	16333	Pre-cut	1 ³ / ₈ " (35 mm)	A4417
NA	16334	Pre-cut	1 ¹ / ₂ " (38 mm)	A4417

 Indicates convexity

SenSura® Mio Convex Light MIDI Closed Pouch 

- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**

**Opaque** 7¼" (19 cm) / 340 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
16321	Cut-to-fit	¾ - 7/8" (10-22 mm)	A4417
16322	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4417
16323**	Cut-to-fit / Pre-cut	1 ¼ - 1 11/16" (31-43 mm)	A4417

**Because of its larger starter hole, 16323 functions as both a pre-cut and cut-to-fit solution.

SenSura® Mio Deep Convex MAXI Closed Pouch 

- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**



8¼" (21 cm) / 550 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
NA	16345	Cut-to-fit	¾ - 7/8" (10-22 mm)	A4417
16355	NA	Cut-to-fit	¾ - 1 11/16" (10-43 mm)	A4417
NA	16346	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4417
NA	16350	Pre-cut	1" (25 mm)	A4417
NA	16351	Pre-cut	1 1/8" (28 mm)	A4417
NA	16347	Pre-cut	1 ¼" (31 mm)	A4417
NA	16353	Pre-cut	1 3/8" (35 mm)	A4417
NA	16354	Pre-cut	1 ½" (38 mm)	A4417

SenSura® Mio Deep Convex MIDI Closed Pouch 

- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**

**Opaque** 7¼" (19 cm) / 340 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
16341	Cut-to-fit	¾ - 7/8" (10-22 mm)	A4417
16342	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4417
16343**	Cut-to-fit / Pre-cut	1 ¼ - 1 11/16" (31-43 mm)	A4417

**Because of its larger starter hole, 16343 functions as both a pre-cut and cut-to-fit solution.

 Indicates convexity









*For important information regarding the use of HCPCS codes listed and code definitions see page 180

SenSura® Mio Convex Flip MAXI Closed Pouch 

- Designed for curves, bulges and hernias
- Curved, star-shaped barrier conforms to outward areas with less creasing and folding
- "Flip to fit" design allows for easy application
- Peek-a-boo inspection window for easy application
- **Not made with natural rubber latex**

**Opaque** 8¼" (21 cm) / 550 mL

Units: 10

Code	Cut	Stoma Size	Max Cut	Barrier Size	HCPCS*
18152	Cut-to-fit	¾ - 1⅞" (10-40 mm)	1⅞" (40 mm)	 S	A4417
18170	Cut-to-fit	¾ - 2" (10-50 mm)	2" (50 mm)	 M	A4417
18171	Pre-cut	1" (25 mm)	2" (50 mm)	 M	A4417
18172	Pre-cut	1⅜" (30 mm)	2" (50 mm)	 M	A4417
18173	Pre-cut	1⅝" (35 mm)	2" (50 mm)	 M	A4417
18180	Cut-to-fit	¾ - 2⅜" (10-60 mm)	2⅜" (60 mm)	 L	A4417
18181	Pre-cut	1⅞" (40 mm)	2⅜" (60 mm)	 L	A4417
18182	Pre-cut	1¾" (45 mm)	2⅜" (60 mm)	 L	A4417

Urostomy**SenSura® Mio Flat MAXI Urostomy Pouch**

- Extended wear
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or most standard night bags
- Neutral gray, water-resistant textile on back
- **Not made with natural rubber latex**

**Transparent** 10½" (27 cm) / 400 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
10586	Cut-to-fit	¾ - 1¾" (10-45 mm)	A4391

 Indicates Convex Flip

SenSura® Mio Soft Convex MAXI Urostomy Pouch 

- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**



10½" (27 cm) / 400 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

Transparent **Opaque**

13679***	16805	Cut-to-fit	$\frac{3}{8}$ - $1\frac{5}{16}$ " (10-33 mm)	A4393
16810	NA	Cut-to-fit	$\frac{3}{8}$ - 2" (10-50 mm)	A4393

Coming soon!

13731	NA	Pre-cut	$\frac{7}{8}$ " (21 mm)	A4393**
13732	NA	Pre-cut	1" (25 mm)	A4393**
13733	NA	Pre-cut	$1\frac{1}{8}$ " (28 mm)	A4393**
13734	NA	Pre-cut	$1\frac{3}{4}$ " (31 mm)	A4393**

SenSura® Mio Soft Convex MIDI Urostomy Pouch 

- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**



Opaque 9½" (24 cm) / 330 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

16801	Cut-to-fit	$\frac{3}{8}$ - $1\frac{5}{16}$ " (10-33 mm)	A4393
--------------	------------	--	-------

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

**Product has not been through PDAC verification

***This cut-to-fit size is coming soon

SenSura® Mio Convex Light MAXI Urostomy Pouch

- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**



10½" (27 cm) / 400 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
16835	16825	Cut-to-fit	⅜ - 7⁄8" (10-22 mm)	A4393
16836	16826	Cut-to-fit	5⁄8 - 1 5⁄16" (15-33 mm)	A4393
16837	NA	Cut-to-fit	⅜ - 1 11⁄16" (10-43 mm)	A4393
NA	16830	Pre-cut	1 1⁄16" (18 mm)	A4393
16840	16831	Pre-cut	1 3⁄16" (21 mm)	A4393
16841	16832	Pre-cut	1" (25 mm)	A4393
16842	16833	Pre-cut	1 1⁄8" (28 mm)	A4393
NA	16827	Pre-cut	1 3⁄4" (31 mm)	A4393

SenSura® Mio Convex Light MIDI Urostomy Pouch

- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**



Opaque 9½" (24 cm) / 330 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
16821	Cut-to-fit	⅜ - 7⁄8" (10-22 mm)	A4393
16822	Cut-to-fit	5⁄8 - 1 5⁄16" (15-33 mm)	A4393
16823	Pre-cut	1 3⁄4" (31 mm)	A4393

 Indicates convexity

SenSura® Mio Deep Convex MAXI Urostomy Pouch

- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**



10½" (27 cm) / 400 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

16865	16855	Cut-to-fit	⅜ - ⅞" (10-22 mm)	A4393
16866	16856	Cut-to-fit	⅝ - 1⅝" (15-33 mm)	A4393
16867	NA	Cut-to-fit	⅜ - 1⅞" (10-43 mm)	A4393
NA	16860	Pre-cut	1⅞" (18 mm)	A4393
16870	16861	Pre-cut	1⅜" (21 mm)	A4393
16871	16862	Pre-cut	1" (25 mm)	A4393
16872	16863	Pre-cut	1⅛" (28 mm)	A4393
NA	16857	Pre-cut	1¼" (31 mm)	A4393

SenSura® Mio Convex Flip MAXI Urostomy Pouch

- Extended wear
- Designed for curves, bulges and hernias
- Curved, star-shaped barrier conforms to outward areas with less creasing and folding
- "Flip to fit" design allows for easy application
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- **Not made with natural rubber latex**







10½" (27 cm) / 400 mL

Units: 10

Code	Code	Cut	Stoma Size	Max Cut	Barrier Size	HCPCS*
------	------	-----	------------	---------	--------------	--------

Transparent Opaque

18405	18401	Cut-to-fit	⅜ - 1⅞" (10-40 mm)	1⅞" (40 mm)		A4393
NA	18402	Pre-cut	¾" (20 mm)	1⅞" (40 mm)		A4393
NA	18403	Pre-cut	1" (25 mm)	1⅞" (40 mm)		A4393
NA	18404	Pre-cut	1⅜" (30 mm)	1⅞" (40 mm)		A4393

 Indicates convexity

 Indicates Convex Flip

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

SenSura® Flat MAXI Drainable Pouch

- Standard wear
- EasiClose™ WIDE Outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11½" (29½ cm) / 655 mL

Units: 20

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

15521	15531	Cut-to-fit	¾ - 3" (10-76 mm)	A4424
15522	15532	Pre-cut	1" (25 mm)	A4424
15523	15533	Pre-cut	1⅜" (30 mm)	A4424
15524	15534	Pre-cut	1⅝" (35 mm)	A4424

SenSura® Flat MAXI Drainable Pouch

- Standard wear
- EasiClose™ WIDE Outlet
- **Without filter**
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 11½" (29½ cm) / 655 mL

Units: 20

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

15821	Cut-to-fit	¾ - 3" (10-76 mm)	A5061
-------	------------	-------------------	-------

SenSura® Xpro Flat MAXI Drainable Pouch

- Extended wear
- EasiClose™ WIDE Outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11½" (29½ cm) / 655 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

15981	NA	Cut-to-fit	¾ - 3" (10-76 mm)	A5056**
NA	15991	Pre-cut	1" (25 mm)	A5056**
NA	15992	Pre-cut	1⅜" (30 mm)	A5056**
NA	15993	Pre-cut	1¼" (32 mm)	A5056**
NA	15994	Pre-cut	1½" (38 mm)	A5056**

SenSura® Convex Light MAXI Drainable Pouch

- Standard wear
- EasiClose™ WIDE Outlet
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11½" (29½ cm) / 655 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
15606	15608	Cut-to-fit	5/8 - 1 ¹¹ / ₁₆ " (15-43 mm)	A4389
15623	15633	Pre-cut	1" (25 mm)	A4389
15624	15634	Pre-cut	1 ¹ / ₈ " (28 mm)	A4389

SenSura® Xpro Convex Light MAXI Drainable Pouch

- Extended wear
- EasiClose™ WIDE Outlet
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11½" (29½ cm) / 655 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
15695	15697	Cut-to-fit	5/8 - 1 ⁵ / ₁₆ " (15-33 mm)	A5057
15696	15698	Cut-to-fit	5/8 - 1 ¹¹ / ₁₆ " (15-43 mm)	A5057

SenSura® Flat MAXI Drainable Pouch

- Standard wear
- Soft outlet provides effective drainage; connects to a standard drain bag or Coloplast fistula drain bag #14010
- **Filter**
- Soft cloth on front and back
- Peek-a-boo inspection window allows easy application
- **Not made with natural rubber latex**



Opaque 12½" (32 cm) / 630 mL

Units: 30

Code	Cut	Stoma Size	HCPCS*
19004	Cut-to-fit	3/8 - 3" (10-76 mm)	A4424

For bed drainage bag, see page 72.

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

**Product has not been through PDAC verification

SenSura® Flat MAXI Closed Pouch

- Standard wear
- **Filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**



Opaque 8½" (21 cm) / 550 mL Units: 30

Code	Cut	Stoma Size	HCPCS*
15480	Cut-to-fit	¾ - 3" (10-76 mm)	A4416
15484	Pre-cut	1¾" (35 mm)	A4416

Urostomy

SenSura® Flat MAXI Urostomy Pouch

- Standard wear
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 10¾" (26 cm) / 400 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
11804	Cut-to-fit	¾ - 3" (10-76 mm)	A5071

SenSura® Flat MIDI Urostomy Pouch

- Standard wear
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 9½" (24 cm) / 330 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
11803	Cut-to-fit	¾ - 2½" (10-66 mm)	A5071

SenSura® Xpro Flat MAXI Urostomy Pouch

- Extended wear
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



10³/₈" (26 cm) / 400 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

11809	NA	Cut-to-fit	3/8 - 3" (10-76 mm)	A4391**
NA	11805	Pre-cut	7/8" (22 mm)	A4391**
NA	11806	Pre-cut	1" (25 mm)	A4391**
NA	11807	Pre-cut	1 3/16" (30 mm)	A4391**
NA	11808	Pre-cut	1 1/4" (32 mm)	A4391**

SenSura® Convex Light MAXI Urostomy Pouch



- Standard wear
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- **Not made with natural rubber latex**

Transparent 10³/₈" (26 cm) / 400 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

11814	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4392
-------	------------	--------------------------	-------

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

**Product has not been through PDAC verification

SenSura® Xpro Convex Light MAXI Urostomy Pouch

- Extended wear
- Belt tabs for use with 4215, 4220
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



10³/₈" (26 cm) / 400 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

11819	NA	Cut-to-fit	5/8 - 1 ⁵ / ₁₆ " (15-33 mm)	A4393**
NA	11816	Pre-cut	7/8" (22 mm)	A4393**
NA	11817	Pre-cut	1" (25 mm)	A4393**
NA	11818	Pre-cut	1 ¹ / ₄ " (32 mm)	A4393**

Urostomy Night Bag and Bag Hanger

See Accessories page 72 for more information.



EasiClose™ WIDE outlets offer a 1/2" wider integrated closure for easy, fast emptying on the go.

Assura® Flat MAXI Drainable Pouch

- Standard wear
- EasiClose™ WIDE Outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11¹/₂" (29 cm) / 530 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

15860	15870	Cut-to-fit	3/8 - 3" (10-75 mm)	A4424
NA	15874	Pre-cut	1" (25 mm)	A4424
NA	15875	Pre-cut	1 ³ / ₁₆ " (30 mm)	A4424
NA	15876	Pre-cut	1 ³ / ₈ " (35 mm)	A4424

 Indicates convexity

Assura® Flat MAXI Drainable Pouch

- Standard wear
- EasiClose™ WIDE Outlet
- **Without filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**



Opaque 11½" (29 cm) / 530 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
15871	Cut-to-fit	¾ - 3" (10-75 mm)	A5061

Assura® Flat MAXI Drainable Pouch

- Extended wear
- EasiClose™ WIDE Outlet
- **Without filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11¼" (28½ cm) / 530 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
15869	15879	Cut-to-fit	¾ - 3" (10-75 mm)	A4388

Assura® Convex Light MAXI Drainable Pouch

- Standard wear
- EasiClose™ WIDE Outlet
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 11½" (29 cm) / 530 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
14511	Cut-to-fit	⅝ - 1⅝" (15-33 mm)	A4389
14512	Cut-to-fit	⅝ - 1⅞" (15-43 mm)	A4389

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

**Product has not been through PDAC verification

Assura® Flat MAXI Drainable Pouch

- Standard wear
- EasiClose™ Outlet
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11¼" (28½ cm) / 530 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

13860	13870	Cut-to-fit	¾ - 3" (10-75 mm)	A4424
NA	13874	Pre-cut	1" (25 mm)	A4424
NA	13875	Pre-cut	1⅜" (30 mm)	A4424
NA	13876	Pre-cut	1⅜" (35 mm)	A4424

Assura® Flat MIDI Drainable Pouch

- Standard wear
- EasiClose™ Outlet
- **Filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**

**Opaque** 10¼" (26 cm) / 390 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

13840	Cut-to-fit	¾ - 2½" (10-65 mm)	A4424
13844	Pre-cut	1" (25 mm)	A4424
13845	Pre-cut	1⅜" (30 mm)	A4424
13846	Pre-cut	1⅜" (35 mm)	A4424

Assura® Convex Light MAXI Drainable Pouch

- Standard wear
- EasiClose™ Outlet
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11¼" (28½ cm) / 530 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

Transparent Opaque

14411	14414	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4389
14412	14415	Cut-to-fit	5/8 - 1 11/16" (15-43 mm)	A4389
NA	14421	Pre-cut	7/8" (22 mm)	A4389
NA	14422	Pre-cut	1" (25 mm)	A4389
NA	14423	Pre-cut	1 1/8" (28 mm)	A4389
NA	14424	Pre-cut	1 1/4" (31 mm)	A4389

Assura® Convex Light MIDI Drainable Pouch

- Standard wear
- EasiClose™ Outlet
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on front and back
- **Not made with natural rubber latex**



Opaque 10¼" (26 cm) / 390 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

14404	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4389
14405	Cut-to-fit	5/8 - 1 11/16" (15-43 mm)	A4389

Assura® Deep Convex MAXI Drainable Pouch

- Standard wear
- EasiClose™ Outlet
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11¼" (28½ cm) / 530 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

14103	14106	Cut-to-fit	5/8 - 1 11/16" (15-43 mm)	A4389
14163	14173	Pre-cut	7/8" (22 mm)	A4389
14164	14174	Pre-cut	1" (25 mm)	A4389
14165	14175	Pre-cut	1 1/8" (28 mm)	A4389
14166	14176	Pre-cut	1 1/4" (31 mm)	A4389

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Assura® Deep Convex MIDI Drainable Pouch

- Standard wear
- EasiClose™ Outlet
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on front and back
- **Not made with natural rubber latex**



Opaque 10¼" (26 cm) / 390 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
14196	Cut-to-fit	5/8 - 1 ¹¹ / ₁₆ " (15-43 mm)	A4389

Assura® Convex Light MAXI Drainable Pouch

- Extended wear
- EasiClose™ Outlet
- **Without filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11¼" (28½ cm) / 530 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

14413	NA	Cut-to-fit	5/8 - 1 ⁵ / ₁₆ " (15-33 mm)	A4390
14416	14419	Cut-to-fit	5/8 - 1 ¹¹ / ₁₆ " (15-43 mm)	A4390

Assura® Flat MAXI Drainable Pouch

- Standard wear
- Clamp closure
- **Without filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11¼" (28½ cm) / 530 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

12660	12670	Cut-to-fit	3/8 - 2 ³ / ₄ " (10-70 mm)	A5061
-------	-------	------------	--	-------

Assura® Flat MINI Drainable Pouch

- Standard wear
- Clamp closure
- **Without filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



9¾" (24½ cm) / 230 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

Transparent Opaque

12630	12620	Cut-to-fit	3/8 - 2¼" (10-55 mm)	A5061
-------	-------	------------	----------------------	-------

 Indicates convexity

Assura® Deep Convex MAXI Drainable Pouch 

- Standard wear
- Clamp closure
- **Without filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11¼" (28½ cm) / 530 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

13706	12533	Cut-to-fit	5/8" - 1 ¹¹ / ₁₆ " (15-43 mm)	A4389
-------	-------	------------	---	-------

Assura® Deep Convex MINI Drainable Pouch 

- Standard wear
- Clamp closure
- **Without filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 9¾" (24½ cm) / 230 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

13703	Cut-to-fit	5/8" - 1 ¹¹ / ₁₆ " (15-43 mm)	A4389
-------	------------	---	-------

Assura® Deep Convex MAXI Drainable Pouch 

- Extended wear
- Clamp closure
- **Without filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



11¼" (28½ cm) / 530 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

12416	12419	Cut-to-fit	5/8" - 1 ¹¹ / ₁₆ " (15-43 mm)	A4390
-------	-------	------------	---	-------

Assura® Flat MAXI Closed Pouch 

- Standard wear
- **Filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**



Opaque 8½" (22 cm) / 530 mL

Units: 30

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

12170	Cut-to-fit	1 ³ / ₁₆ " - 3" (20-75 mm)	A4416
-------	------------	--	-------

12174	Pre-cut	1" (25 mm)	A4416
-------	---------	------------	-------

12175	Pre-cut	1 ³ / ₁₆ " (30 mm)	A4416
-------	---------	--	-------

12176	Pre-cut	1 ³ / ₈ " (35 mm)	A4416
-------	---------	---	-------

12177	Pre-cut	1 ⁹ / ₁₆ " (40 mm)	A4416
-------	---------	--	-------

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Assura® Flat MIDI Closed Pouch

- Standard wear
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



7" (18 cm) / 350 mL Units: 30

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
12130	12140	Cut-to-fit	1 ³ / ₁₆ - 2 ¹ / ₂ " (20-65 mm)	A4416
12134	12144	Pre-cut	1" (25 mm)	A4416
12135	12145	Pre-cut	1 ³ / ₁₆ " (30 mm)	A4416
12136	12146	Pre-cut	1 ³ / ₈ " (35 mm)	A4416
NA	12147	Pre-cut	1 ⁹ / ₁₆ " (40 mm)	A4416
NA	12149	Pre-cut	2" (50 mm)	A4416

Assura® Flat MINI Closed Pouch

- Standard wear
- **Filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**



Opaque 5" (12 cm) / 170 mL Units: 30

Code	Cut	Stoma Size	HCPCS*
12110	Cut-to-fit	1 ³ / ₁₆ - 2 ¹ / ₂ " (20-65 mm)	A4416

Assura® Convex Light MAXI Closed Pouch

- Standard wear
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



8 1/2" (22 cm) / 560 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent Opaque				
14441	14444	Cut-to-fit	5/8 - 1 ⁵ / ₁₆ " (15-33 mm)	A4417
14442	14445	Cut-to-fit	5/8 - 1 ¹¹ / ₁₆ " (15-43 mm)	A4417

Indicates convexity

Assura® Convex Light MIDI Closed Pouch 

- Standard wear
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



7" (18 cm) / 350 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

Transparent Opaque

14431	14434	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4417
14432	14435	Cut-to-fit	5/8 - 1 11/16" (15-43 mm)	A4417
NA	14451	Pre-cut	7/8" (22 mm)	A4417
NA	14452	Pre-cut	1" (25 mm)	A4417
NA	14453	Pre-cut	1 1/8" (28 mm)	A4417
NA	14454	Pre-cut	1 1/4" (31 mm)	A4417

Assura® Flat Stoma Cap

- Standard Wear
- **Filter**
- Belt tabs
- Soft cloth on back
- Small pouch size for intimate moments or colostomates who irrigate
- **Not made with natural rubber latex**



Opaque

Units: 30

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

2501	Cut-to-fit	1 3/16 - 2 1/4" (20-55 mm)	A5055
------	------------	----------------------------	-------

Urostomy

Assura® Flat MAXI Urostomy Pouch

- Standard wear
- Multi-chambers minimize sloshing
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 10 3/4" (27 1/2 cm) / 470 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

14222	Cut-to-fit	3/8 - 2 1/4" (10-55 mm)	A5071
-------	------------	-------------------------	-------

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Assura® Convex Light MAXI Urostomy Pouch 

- Standard wear
- Multi-chambers minimize sloshing
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 10¾" (27½ cm) / 470 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
14717	Cut-to-fit	5/8 - 1 ¹¹ / ₁₆ " (15-43 mm)	A4392

Assura® Flat MAXI Urostomy Pouch

- Standard wear
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



10¾" (27 cm) / 465 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

5585	5580	Cut-to-fit	3/8 - 2 ¹ / ₄ " (10-55 mm)	A5071
------	------	------------	--	-------

Assura® Flat MIDI Urostomy Pouch

- Standard wear
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



7" (18 cm) / 265 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

5575	5570	Cut-to-fit	3/8 - 2 ¹ / ₄ " (10-55 mm)	A5071
------	------	------------	--	-------

Assura® Flat MAXI Urostomy Pouch

- Extended wear
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 9½" (24 cm) / 375 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
12474	Cut-to-fit	3/8 - 2 ¹ / ₄ " (10-55 mm)	A4391

 Indicates convexity

Assura® Flat MIDI Urostomy Pouch

- Extended wear
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



9 ½" (24 cm) / 265 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent Opaque

12472	12471	Cut-to-fit	3/8" - 2 1/4" (10-55 mm)	A4391
-------	-------	------------	--------------------------	-------

Assura® Deep Convex MAXI Urostomy Pouch

- Standard wear
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- **Not made with natural rubber latex**

**Transparent** 10 ¾" (27 cm) / 465 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
12595	Cut-to-fit	5/8" - 1 5/16" (15-33 mm)	A4392
12596	Cut-to-fit	5/8" - 1 11/16" (15-43 mm)	A4392
12991	Pre-cut	5/8" (15 mm)	A4392
12992	Pre-cut	3/4" (18 mm)	A4392
12993	Pre-cut	7/8" (22 mm)	A4392
12994	Pre-cut	1" (25 mm)	A4392
12995	Pre-cut	1 1/8" (28 mm)	A4392
12996	Pre-cut	1 1/4" (31 mm)	A4392
12997	Pre-cut	1 3/8" (35 mm)	A4392

Assura® Deep Convex MAXI Urostomy Pouch

- Extended wear
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- **Not made with natural rubber latex**

**Transparent** 10 ¾" (27 cm) / 465 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
12432	Cut-to-fit	5/8" - 1 5/16" (15-33 mm)	A4393
12433	Cut-to-fit	5/8" - 1 11/16" (15-43 mm)	A4393

Note: SenSura Mio Original has been discontinued. For a full list of discontinued products and suggested replacement products please refer to the table starting on page 171.

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 180


Download the Coloplast Stoma Care App



Designed with nurses in mind

Features:

1. Easy access to Coloplast® Care
2. BodyCheck™ Assessment
3. Step-by-step guides with how-to-use videos
4. Common problems and recommended actions
5. General information on urostomies, ileostomies and colostomies



Download the Coloplast Stoma Care app now!

Search for “**Stoma Care**” in the App Store and Google Play to download for free today.





Soft
outlet



Tap
outlet

Hospital Assortment



Elastic adhesive fits over folds and contours, yet remains securely in place



Two outlet types for different options:

- Soft (for liquid stool; fits standard drain bags)
- Tap (for thicker output; fits fecal drain bag with wide-lumen tubing)



Printed measurement scale on post op pouches helps visually track output volume



Security and comfort from hospital to home.

Your patients deserve all the benefits of SenSura[®] Mio from day one.

Getting used to a stoma is challenging enough without the added complexity of having to learn a different pouching system. Why not give your patients the security and comfort of Mio - right from the start? Options are available in different sizes, outlet types, and coupling systems to meet patient needs.



2-piece Click
Soft | Tap
**Also available
in smaller size!*



2-piece Flex
Soft | Tap



1-piece
Soft
*Available in flat,
soft convex, and
light convex!*



Post Op with window
Soft
Sterile | Non-Sterile

High Output Pouches




SenSura® Mio Click High Output Pouch with Soft Outlet

- Special full-circle filter designed specifically for ileostomy output
- Neutral gray, water-resistant textile
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Soft outlet provides effective drainage; connects to a standard drain bag or Coloplast fistula drain bag #14010
- **Not made with natural rubber latex**



12½" (32 cm) / 1080 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
18640	18630		A4413
18641	18631		A4413
18642	18632		A4413





SenSura® Mio Click High Output Pouch with Tap Outlet

- Special full-circle filter designed specifically for ileostomy output
- Neutral gray, water-resistant textile
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- **Not made with natural rubber latex**



12½" (32 cm) / 1080 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
18610	18600		A4413
18611	18601		A4413
18612	18602		A4413
18613	18603		A4413

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

SenSura® Mio Click Drainable Pouch with Tap Outlet


- Special full-circle filter designed specifically for ileostomy output
- Neutral gray, water-resistant textile
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- **Not made with natural rubber latex**



9" (23 cm) / 650 mL Units: 10

Code	Color	HCPCS*
------	-------	--------

Opaque

18620		A4427
18621		A4427
18622		A4427

Assura® High Output Drainable Pouch




- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



12½" (32 cm) / 750 mL Units: 10

Code	Code	Color	HCPCS*
------	------	-------	--------

Transparent Opaque

8114	NA		A4413
8115	2846		A4413
8116	2847		A4413

For bed drainage bag and bag hanger, see page 72.

SenSura® Mio Flex High Output Pouch with Tap Outlet



- Special full-circle filter designed specifically for ileostomy output
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- **Not made with natural rubber latex**



12½" (32 cm) / 1080 mL Units: 10

Code	Code	Color	HCPCS*
------	------	-------	--------

Transparent Opaque



18655	18650		A4413
18656	18651		A4413

SenSura® Mio Flex High Output Pouch with Soft Outlet



- **Special full-circle filter designed specifically for ileostomy output**
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to a standard drain bag or Coloplast fistula drain bag #14010
- **Not made with natural rubber latex**

12½" (32 cm) / 1080 mL Units: 10

Code	Code	Color	HCPCS*
Transparent		Opaque	
18657	18652		A4413
18658	18653		A4413

SenSura® Mio MAXI Drainable Pouch with Soft Outlet



- Extended wear
- **Special full-circle filter designed specifically for ileostomy output**
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or most standard night bags
- **Not made with natural rubber latex**

11" (28 cm) / 630 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent		Opaque		
18665	18660	Cut-to-fit	⅜ - 2 ¾" (10-70 mm)	A5056
18666	18661	Cut-to-fit	⅜ - 4" (10-100 mm)	A5056

SenSura® Mio Soft Convex MAXI Drainable Pouch with Soft Outlet



- Extended wear
- Belt tabs for use with 4237, 4247
- **Special full-circle filter designed specifically for ileostomy output**
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or any standard night bag
- **Not made with natural rubber latex**

11" (28 cm) / 630 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent		Opaque		
18671	18670	Cut-to-fit	⅜ - 2" (10-50 mm)	A5057

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

SenSura® Mio Convex Light MAXI Drainable Pouch with Soft Outlet

- Extended wear
- Belt tabs for use with 4237, 4247
- **Special full-circle filter designed specifically for ileostomy output**
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or any standard night bag
- **Not made with natural rubber latex**



11" (28 cm) / 630 mL Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Transparent **Opaque**

18676	18675	Cut-to-fit	$\frac{3}{8}$ - $1\frac{11}{16}$ " (10-43 mm)	A5057
-------	-------	------------	---	-------

SenSura® Flat MAXI Drainable Pouch

- Standard wear
- Soft outlet provides effective drainage; connects to a standard drain bag or Coloplast fistula drain bag #14010
- **Filter**
- Soft cloth on front and back
- Peek-a-boo inspection window allows easy application
- **Not made with natural rubber latex**



Opaque 12 $\frac{1}{2}$ " (32 cm) / 630 mL Units: 30

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

19004	Cut-to-fit	$\frac{3}{8}$ - 3" (10-76 mm)	A4424
-------	------------	-------------------------------	-------

For bed drainage bag, see page 72.

SenSura® Mio 1-piece Post Op Pouch with Window

- Neutral gray, water-resistant textile
- Printed scale provides visual indication of output volume
- Easy-to-open and close post-op window attached to pouch
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- **Not made with natural rubber latex**



Transparent Units: 5

Code	Code	Cut	Volume	Stoma Size
------	------	-----	--------	------------

Non-sterile **Sterile**

18690	18680	Cut-to-fit	680mL	$\frac{3}{8}$ - $2\frac{3}{4}$ " (10-70 mm)
-------	-------	------------	-------	---

18691	18681	Cut-to-fit	725mL	$\frac{3}{8}$ - 4" (10-100 mm)
-------	-------	------------	-------	--------------------------------

 Indicates convexity

SenSura® Standard Wear 1-piece Post-Op Pouch with Window

- Flat
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- Easy-to-open and -close access window attached to pouch
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 13" (33 cm) / 830 mL Units: 5

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Non-sterile **Sterile**

19020	19010	Cut-to-fit	3/8 - 3" (10-76 mm)	A6154
19021	19011	Cut-to-fit	3/8 - 4 1/2" (10-115 mm)	A6154

SenSura® Standard Wear 1-piece Post-Op Pouch

- Flat
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- Effective pouching solution for small fistulas
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 12 1/4" (31 cm) / 690 mL Units: 5

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Non-sterile **Sterile**

19030	19001	Cut-to-fit	3/8 - 3" (10-76 mm)	A6154
--------------	--------------	------------	---------------------	-------

For SenSura® Magnum high output pouches, see page 86.

Assura® Standard Wear 1-piece Post-Op Pouch with Window

- Flat
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- Easy-to-open and -close access window
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 13" (33 cm) / 830 mL Units: 5

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Non-sterile **Sterile**

12800	12808	Cut-to-fit	1/2 - 2 3/4" (13-70 mm)	A6154
12802	12810	Cut-to-fit	1/2 - 4" (13-100 mm)	A6154

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Assura® Standard Wear 1-piece Post-Op Pouch

- Flat
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- Soft cloth on back
- **Not made with natural rubber latex**

Transparent 12¾" (31 cm) / 690 mL Units: 5

Code	Code	Cut	Stoma Size	HCPCS*
------	------	-----	------------	--------

Non-sterile **Sterile**

12805	12818	Cut-to-fit	½ - 2¾" (13-70 mm)	A6154
-------	-------	------------	--------------------	-------

Assura® Ileo Night Bag

- 2 liter capacity for very high output stomas
- Clamp closure
- Comes with bag hanger
- **Filter**
- **Not made with natural rubber latex**



2 L Units: 5

Code	Color	HCPCS*
------	-------	--------

Opaque

2836		A4413
------	---	-------

2837		A4413
------	---	-------



Coloplast® Fistula and Wound Management System

Each system includes:

- 1 pouch and tracing guide
- 1 drain port (extras sold separately on next page)
- 1 lid (extras sold separately on next page)
- 1 clear drape
- 1 clamp
- **Not made with natural rubber latex**

Each box includes one pump.



14050 Mini



14060 Midi



14070 Maxi

Units: 3
Units for 14050: 6

Code	Cut	Stoma Size	HCPCS*
14050	Mini	4 $\frac{1}{8}$ x 6 $\frac{1}{4}$ " (104-159 mm)	A6154
14060	Midi	6 $\frac{1}{8}$ x 9" (156-228 mm)	A6154
14070	Maxi	8 $\frac{1}{8}$ x 11 $\frac{3}{4}$ " (208-297 mm)	A6154

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Fistula / Wound Pouch Accessories

Coloplast® Fistula Drain Port

Can be easily applied to the lid to allow for catheter insertion. For catheters ranging from 18FR - 22FR, the drain port is ready for use. For catheters ranging from 24FR - 30FR, pull off the top of the drain port - DO NOT CUT.



- Not made with natural rubber latex

Code	Size	Units
14015	One size	12

Coloplast® Bed Drainage Bag

Bag tubing specifically designed for fecal output. Holds up to 2000 mL. Long tubing can be cut easily, allowing nurses to find the right length for bed drainage.



- Not made with natural rubber latex

Code	Size	Units: 6	HCPCS*
14010	2000 mL		A4357

For matching bag hanger, see page 72.

Coloplast® Transparent Fistula Lids

Replacement lids match the transparent lid included with each fistula pouch.



- Not made with natural rubber latex

Code	Size	Units
14021	Mini (matches 14050)	10
14031	Midi (matches 14060)	10
14041	Maxi (matches 14070)	10

Coloplast® Pouch Clamp

Simple closure clamp is easy to manipulate, even for those with low dexterity.



- Not made with natural rubber latex

Code	Size	Units: 20	HCPCS*
9500	One size		A4363

2-Piece Click Single Use Kits

SenSura® Mio Click Barrier + Drainable Pouch


- Extended wear
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit / Transparent

21008	3/8 - 2 1/8" (10-55 mm)		A4410
			A4427

Barrier #10522 Flat with belt tabs, Pouch #11483 MAXI / EasiClose™ WIDE Outlet

SenSura® Click Standard Wear Barrier + Drainable Pouch


- Barrier #10031 Flat with belt tabs
- Pouch #11196 MAXI / EasiClose™ WIDE Outlet
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit / Transparent

11009	3/8 - 2 1/4" (10-55 mm)		A4414
			A4427

SenSura® Click Standard Wear Barrier + Drainable Pouch


- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit / Transparent

11010	5/8 - 1 11/16" (15-43 mm)		A4373
			A4427

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

SenSura® Click Xpro Extended Wear Barrier + Urostomy Pouch


- Barrier #10035 Flat with belt tabs
- Pouch #11856 MAXI / soft outlet
- Urostomy adapter
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit / Transparent

11004	⅜ - 2¼" (10-55 mm)		A4409
			A4433

2-Piece Flex Single Use Kits

SenSura® Mio Flex Barrier + MAXI Drainable Pouch


- Extended wear
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**




Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------


Cut-to-fit / Transparent

21009	⅜ - 1⅞" (10-48 mm)		A4410
			A4425


Barrier #10561 Flat, Pouch #12273 EasiClose™ WIDE Outlet / filter

21010	⅜ - 2¾" (10-68 mm)		A4410
			A4425

Barrier #10571 Flat, Pouch #12283 EasiClose™ WIDE Outlet / filter

21012	⅜ - 1⅞" (10-48 mm)		A4410
			A5063

Barrier #10561 Flat, Pouch #12274 EasiClose™ WIDE Outlet / without filter

21013	⅜ - 2¾" (10-68 mm)		A4410
			A5063


Barrier #10571 Flat, Pouch #12284 EasiClose™ WIDE Outlet / without filter

SenSura® Mio Flex + MAXI Urostomy Pouch

- Extended wear
- Urostomy adapter
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5



Code	Stoma Size	Color	HCPCS*
Cut-to-fit / Transparent			
110072	3/8 - 1 1/8" (10-48 mm)		A4410 A5073**
Barrier #10561 Flat, Pouch #12298 Soft outlet / multi-chamber			

SenSura® Mio Flex Convex Barrier + MAXI Drainable Pouch 

- Extended wear
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit / Transparent			
13660	5/8 - 1 9/16" (15-40 mm)		A4407 A4425
Barrier #16483 Convex Light, Pouch #12273 EasiClose™ WIDE Outlet / filter			
13661	5/8 - 2 1/16" (15-53 mm)		A4408 A4425
Barrier #16491 Convex Light, Pouch #12283 EasiClose™ WIDE Outlet / filter			

SenSura® Mio Flex Convex + MAXI Urostomy Pouch

- Extended wear
- Urostomy adapter
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit / Transparent			
13662	5/8 - 1 9/16" (15-40 mm)		A4407 A5073**
Barrier #16483 Convex, Pouch #12298 Soft outlet / multi-chamber			

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

**Product has not been through PDAC verification

SenSura® Flex Standard Wear Barrier + MAXI Drainable Pouch

- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit / Transparent

11011K	⅜ - 1⅞" (10-48 mm)		A4414
			A4425
Barrier #10102 Flat, Pouch #11512 EasiClose™ WIDE Outlet / filter			
21001	⅜ - 2¾" (10-68 mm)		A4415
			A4425
Barrier #10103 Flat, Pouch #11517 EasiClose™ WIDE Outlet / filter			
21011	⅜ - 3½" (10-88 mm)		A4415
			A5063
Barrier #10108 Flat, Pouch #11690 EasiClose™ WIDE Outlet / without filter			

SenSura® Flex Xpro Extended Wear Barrier + MAXI Drainable Pouch



- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit / Transparent

21002	⅜ - 1⅞" (10-48 mm)		A4409
			A4425
Barrier #10105 Flat, Pouch #11512 EasiClose™ WIDE Outlet / filter			
21003	⅜ - 2¾" (10-68 mm)		A4410
			A4425
Barrier #10106 Flat, Pouch #11517 EasiClose™ WIDE Outlet / filter			

SenSura® Flex Xpro Extended Wear Barrier + MAXI Drainable Pouch


- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit / Transparent

21004	⅝ - 2⅛" (15-53 mm)		A4408
			A4425
Barrier #11308 Convex Light, Pouch #11517 EasiClose™ WIDE Outlet / filter			

 Indicates convexity

Assura® Standard Wear Barrier + Drainable Pouch


- Barrier #2883 Flat with belt tabs
- Pouch #15976 MAXI / EasiClose™ WIDE Outlet
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit / Transparent

11002	3/8 - 2 1/4" (10-55 mm)		A4414
			A4427

Assura® Standard Wear Barrier + Drainable Pouch


- Barrier #14263 Convex Light with belt tabs
- Pouch #15976 MAXI / EasiClose™ WIDE Outlet
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Units: 5

Code	Stoma Size	Color	HCPCS*
------	------------	-------	--------

Cut-to-fit / Transparent

11001	5/8 - 1 11/16" (15-43 mm)		A4373
			A4427

1-Piece Single Use Kits

SenSura® Mio 1-piece Single Use Kit

- Extended wear
- Pouch #10481 MAXI / EasiClose™ WIDE Outlet
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Transparent Units: 5

Code	Cut	Stoma Size	HCPCS*
21007	Cut-to-fit	3/8 - 2 1/8" (10-55 mm)	A5056

SenSura® Mio Convex Soft Single Use Kit

- Extended wear
- Pouch #16716 MAXI / EasiClose™ WIDE Outlet / Filter
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Transparent Units: 5

Code	Cut	Stoma Size	HCPCS*
13665	Cut-to-fit	3/8 - 2" (10-50 mm)	A5057

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

SenSura® Mio Convex Light Single Use Kit 

- Extended wear
- Pouch #16739 MAXI / EasiClose™ WIDE Outlet / Filter
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Transparent Units: 5

Code	Cut	Stoma Size	HCPCS*
13664	Cut-to-fit	$\frac{3}{8}$ - $1\frac{11}{16}$ " (10-43 mm)	A5057

SenSura® Mio Soft Convex Urostomy Single Use Kit 

- Extended wear
- Pouch #16837 MAXI / EasiClose™ WIDE Outlet / Filter
- Instructions for use
- Stoma measuring guide
- Urostomy adapter
- **Not made with natural rubber latex**



Transparent Units: 5

Code	Cut	Stoma Size	HCPCS*
13663	Cut-to-fit	$\frac{3}{8}$ - $1\frac{11}{16}$ " (10-43 mm)	A4393

SenSura® 1-piece Single Use Kit

- Pouch #15521 MAXI / EasiClose™ WIDE Outlet
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Transparent Units: 5

Code	Cut	Stoma Size	HCPCS*
11008	Cut-to-fit	$\frac{3}{8}$ - 3" (10-76 mm)	A4424

Assura® 1-piece Single Use Kit

- Pouch #15860 MAXI / EasiClose™ WIDE Outlet
- Instructions for use
- Stoma measuring guide
- **Not made with natural rubber latex**



Transparent Units: 5

Code	Cut	Stoma Size	HCPCS*
11005	Cut-to-fit	$\frac{3}{8}$ - 3" (10-76 mm)	A4424

SenSura® Xpro 1-piece Urostomy Single Use Kit

- Extended wear
- Pouch #11809 MAXI / soft outlet
- Instructions for use
- Stoma measuring guide
- Urostomy adapter
- **Not made with natural rubber latex**



Transparent $10\frac{3}{8}$ " (26 cm) / 400 mL Units: 5

Code	Cut	Stoma Size	HCPCS*
21005	Cut-to-fit	$\frac{3}{8}$ - 3" (10-76 mm)	A4391**

 Indicates convexity

SenSura® Xpro 1-piece Convex Light Urostomy Single Use Kit

- Extended wear
- Pouch #11819 MAXI / soft outlet
- Instructions for use
- Stoma measuring guide
- Urostomy adapter
- **Not made with natural rubber latex**



Transparent 10³/₈" (26 cm) / 400 mL Units: 5

Code	Cut	Stoma Size	HCPCS*
21006	Cut-to-fit	5 ⁵ / ₈ - 1 ⁵ / ₁₆ " (15-33 mm)	A4393**

Extra Protection

Brava® Elastic Barrier Strips

Skin-friendly strips are used to "picture frame" around an ostomy barrier, providing extra security.

- Elastic (follow the contours of the body)
- Ensure the position of the barrier and minimize roll-up
- Skin-friendly and secure alternative to tape
- **Not made with natural rubber latex**



Brava® Elastic Barrier Strips- Curved

For round or oval barrier shapes

Code	Size	Units	HCPCS*
120700**	5 ¹ / ₂ " strips	20 strips	A4362

**2 strips = 1 billable unit of A4362

Brava® Elastic Barrier Strips- Y shape

Designed for users with challenging body shapes, e.g. hernias or skin folds



Code	Size	Units	HCPCS*
120721**	5 ³ / ₄ " strips	30 strips	A5121

**3 strips = 1 billable unit of A5121

Brava® Elastic Barrier Strips- Straight

For square or rectangular barrier shapes



Code	Size	Units	HCPCS*
120740**	5 ¹ / ₄ " strips	20 strips	A4362

**2 strips = 1 billable unit of A4362

Brava® Elastic Barrier Strip XL

Designed for increased security & support around the barrier. Wider shape accommodates the petals of Convex Flip.



Code	Size	Units	HCPCS*
12076**	6 ³ / ₄ " strips	30 strips	A5121

**3 strips = 1 billable unit of A5121

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

**Product has not been through PDAC verification

Brava® Protective Seal

Moldable barrier ring designed to adhere securely and protect the skin around the stoma.

- Soft, flexible material adapts to body shape as the wearer bends and stretches
- Can be easily and cleanly removed
- Alcohol-free and sting-free
- Multiple sizes for a custom fit
- **Not made with natural rubber latex**



Code	Stoma size	Adhesive diameter	Thickness	Seal width	Units	HCPCS*
12035	¾"	1 7/8"	2.5 mm	Standard (5/8")	10	A4385
12036	¾"	2 ¼"	2.5 mm	Wide (¾")	10	A4385
12038	¾"	2 ½"	2.5 mm	Wide (7/8")	10	A4385
12037	1 1/8"	2 ¼"	2.5 mm	Standard (5/8")	10	A4385
12039	1 3/8"	2 ½"	2.5 mm	Standard (5/8")	10	A4385
12045	¾"	1 7/8"	4.2 mm	Standard (5/8")	10	A4385
12046	¾"	2 ¼"	4.2 mm	Wide (¾")	10	A4385
12048	¾"	2 ½"	4.2 mm	Wide (7/8")	10	A4385
12047	1 1/8"	2 ¼"	4.2 mm	Standard (5/8")	10	A4385
12049	1 3/8"	2 ½"	4.2 mm	Standard (5/8")	10	A4385
12032	¾"	3"	2.5 mm	Extra wide (1 1/8")	10	A4385
12033	1 1/8"	3"	2.5 mm	Extra wide (1 5/16")	10	A4385

Brava® Protective Seal Convex

Moldable barrier ring designed to provide convexity for flush or retracted stomas.

- High-tack adhesive that adheres quickly to the skin, filling small creases, folds, or scarring around the stoma
- Erosion resistant composition
- High flexibility to maintain seal during body movement
- Steep sides designed to provide convexity to peristomal area
- Thin edges fit snugly around peristomal area
- **Not made with natural rubber latex**



Code	Stoma size	Thickness	Units	HCPCS*
12090	¾"	8.2mm	10	A4411
12091	1"	8.2mm	10	A4411
12092	1 3/16"	8.2mm	10	A4411
12093	1 3/8"	8.2mm	10	A4411
12094	1 9/16"	8.2mm	10	A4411
12095	1 3/4"	8.2mm	10	A4411

Brava® Moldable Ring

Designed for dual-action leakage reduction - used to fill in uneven skin surfaces and protect the skin.

- Molds and adheres securely to dips, creases, and folds to help create a tight seal around the stoma
- Helps keep skin healthy by protecting it from output and absorbing moisture from the skin
- Easy to handle, mold and apply
- Alcohol-free and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
120307	2 mm (thin)	10	A4385
120427	4.2 mm (thick)	10	A4385

Brava® Strip Paste

Pectin-based product that protects peristomal skin.

- Protects peristomal skin from stoma output
- Used to fill uneven skin surfaces to create a flat pouching surface
- Alcohol-free and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
26555	2 oz. total package (.2 oz. per strip)	10 strips	A4406

Brava® Paste

A sting-free alternative to conventional alcohol-based paste.

- Creates a tight, effective seal between peristomal skin and the ostomy/fistula appliance
- No alcohol = no need to "set"; contains pectin
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
12050	2 oz.	1	A4406

Coloplast® Paste

Used to fill in uneven skin surfaces to create a flat pouching surface.

- Provides tight sealing around the stoma
- Protects against stoma output, minimizing risk of maceration
- Low alcohol content; no pectin
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
2650	2 oz.	12	A4405

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Brava® Skin Cleanser Wipes

The wipe that cleans and cares. A convenient way to cleanse and maintain healthy skin when changing the ostomy pouching system.



- Contains ProVitamin B5 to provide a soothing and moisturizing effect
- Does not contain alcohol or perfumes
- No residue left on the skin that could affect barrier adhesion
- **Not made with natural rubber latex**

Code	Size	Units	HCPCS*
12080	1 wipe	15 wipes/pack	A9286**

Brava® Lubricating Deodorant

A pouch lubricant + deodorizer in one. Powder-fresh liquid reduces odor and allows output to slide to the bottom of the pouch for easier emptying.



- **Not made with natural rubber latex**

Code	Size	Units	HCPCS*
12060	.25 fl. oz. (7.5 mL) sachet	20 sachets	A4394
12061	8 oz. bottle	16	A4394

Brava® Adhesive Remover Spray

- Designed to lift the edges of the barrier for ease in removal
- Dry within seconds without leaving an oily residue behind
- Alcohol-free, silicone-based, and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
120105	1.7 oz. (50 mL)	1	A4455**

Brava® Adhesive Remover Wipes

- Wipe away adhesive residue
- Dry within seconds without leaving an oily residue behind
- Alcohol-free, silicone-based, and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
120115	1 wipe	30 wipes/box	A4456**

Brava® Skin Barrier Spray

- Protects against stoma output and adhesives
- Dry within seconds without building a thick or tacky layer
- Alcohol-free, silicone-based, and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
120205	1.7 oz. (50 mL)	1	A4369**

Brava® Skin Barrier Wipes

- Protects against stoma output and adhesives
- Dry within seconds without building a thick or tacky layer
- **Tip:** wring the wipe like a wet rag to expel the liquid, then wipe to disperse on skin
- Alcohol-free, silicone-based, and sting-free
- **Not made with natural rubber latex**



Code	Size	Units	HCPCS*
120215	1 wipe	30 wipes/box	A5120**

*For important information regarding the use of HCPCS codes listed and code definitions see page 180


**Product has not been through PDAC verification

BodyCheck™

Identifying ostomy solutions based on individual needs

Together with 300 Wound, Ostomy, Contenance (WOC) nurses, we've developed BodyCheck™ body profile assessment which uses Body Profile Terminology, a common language to describe the peristomal body profile.

- 1 Assess your patient's overall peristomal area as a Regular, Inward, or Outward



- 2 Visit www.bodycheck.coloplast.us to complete the assessment



- 3 Receive personalized ostomy solutions *(may include accessory options)*



Visit www.bodycheck.coloplast.us


Henry

Making life easier_

For little fighters

Babies and kids present different pouching challenges than adults - so we've applied **BodyFit Technology** to develop specific solutions for different stages in a child's life.

Stepping up over time, as kids grow



Ask us about our pediatric education!



Pediatric FAQ

A quick tipsheet with answers to frequently asked questions, like "what do I do if friends and family - or other kids - have questions?"

and "what do I do if my child pulls the pouch off?"

Caregivers can request this via our Coloplast® Care program.




My Child Has An Ostomy: A Resource Sheet for Caregivers

A fill-in-the-blank sheet to help parents communicate about their child's needs to

day care providers, school nurses, and other caregivers.

Caregivers can request this via our Coloplast® Care program.




Life After Your Child's Ostomy

A high-level overview of anatomy, how to change a pouch, pouching system options, where to get

supplies, dietary notes, and other information in English or Spanish.

Caregivers can request this via our Coloplast® Care program.



When I Met Tipster: A coloring book for children with ostomies

A simple coloring book to help kids engage with their ostomy surgery and "new normal" - or learn more about a sibling's ostomy.

Available to hospitals via Coloplast sales representatives, or to caregivers via our Coloplast® Care program.



Jake Gets a Stoma

A children's book written by Jake's mom, who wanted a better way to explain the surgery and stoma management to her son and other kids going through the same thing.

Available to hospitals via Coloplast sales representatives, or to caregivers via our Coloplast® Care program.




SenSura® Mio Bunny

Coloplast offers a complimentary training bunny for product education use with both

SenSura® Mio Baby & Kids. Contact your sales representative for more information.

Limitations apply. The SenSura® Mio bunny is provided to healthcare providers to offer to families for teaching purposes only and should only be used with Coloplast products.



Teen Life with an Ostomy

Navigating school, what to wear, and hanging out with friends aren't covered by the average adult-oriented materials -

that's why we're so excited to feature tips from Paige, a teen who's been through it all!

Available to patients via our Coloplast® Care program.

SenSura® Mio Baby

2-piece pouching system designed with BodyFit Technology to tackle the challenges of preemies and babies from 30 weeks of gestation to 6 months of age.


SenSura® Mio Baby Flex 2-Piece Barrier

- Oval barrier accommodates one or multiple stomas
- No starter hole – custom cut
- Gentle adhesive sticks securely and removes gently with Brava® Baby Adhesive Remover
- Soft flange enables the pouch to be removed and replaced
- Sterilized surface, single-packed
- **Not made with natural rubber latex**



Cut-to-fit

Units: 5

Code	Stoma Size	Color	HCPCS*
18705	0 - 1 ⁹ / ₁₆ " (0-40 mm)		A4414

SenSura® Mio Baby Flex 2-Piece Pouch

- Adhesive coupling is easy to apply and remove
- Access window
- Soft outlet fits to enteral feeding systems
- Sterile, single-packed
- **Not made with natural rubber latex**



Transparent, No filter

Units: 10

Code	Color	HCPCS*
18700		A5063

Brava® Baby Adhesive Remover

- Oil-based adhesive remover designed specifically for babies
- Liquid sachet (no wipe) can be used with cotton swabs or gauze
- Gently removes adhesive / adhesive residues from SenSura® Mio Baby Barrier
- Not recommended for use with SenSura Mio Kids or adult products
- **Not made with natural rubber latex**



Units: 30

Code	Size	HCPCS*
12014	3 mL sachet	A4455

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

SenSura® Mio Kids

1-piece and 2-piece pouching systems designed with BodyFit Technology to fit kids' active movements and rounded bellies. Designed with options to step up as kids grow, the range covers kids 6 months to 4 years of age (larger options for older kids are found in the SenSura Mio adult line).

SenSura® Mio Kids Flex 2-Piece Barrier

- Star-shaped, elastic barrier can be rotated for flexibility in placement
- No starter hole – custom cut
- **Not made with natural rubber latex**



Cut-to-fit

Units: 5

Code	Stoma Size	Color	HCPCS*
18735	0 - 1 ¹ / ₁₆ " (0-27 mm)		A4414
18747	0 - 1 ³ / ₈ " (0-35 mm)		A4414

SenSura® Mio Kids Flex 2-Piece MINI Drainable Pouch

- EasiClose™ Outlet
- Full-circle filter
- Neutral gray, water-resistant textile
- Small size enables kids to step up in pouch size over time
- Opaque includes inspection window
- **Not made with natural rubber latex**

6¹/₈" folded

Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque		
18721	18720		A4425



SenSura® Mio Kids Flex 2-Piece MIDI Drainable Pouch

- EasiClose™ Outlet
- Full-circle filter
- Neutral gray, water-resistant textile
- Opaque includes inspection window
- **Not made with natural rubber latex**



6¾" folded

Units: 10

Code	Code	Color	HCPCS*
Transparent Opaque			
18723	18722		A4425
18725	18724		A4425



SenSura® Mio Kids Flex 2-piece MIDI Urostomy Pouch

- Neutral gray, water-resistant textile
- Soft outlet can be connected to a drain bag for nighttime use, or a leg bag for kids in a wheelchair
- **Not made with natural rubber latex**



Opaque 6¾" folded

Units: 10

Code	Color	HCPCS*
18730		A5073
18731		A5073

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

SenSura® Mio Kids 1-Piece MIDI Drainable Pouch

- Star-shaped, elastic barrier with no starter hole – custom cut
- EasiClose™ Outlet
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque includes inspection window
- **Not made with natural rubber latex**



6¾" folded | Cut-to-fit Units: 10

Code	Code	Stoma size	HCPCS*
Transparent	Opaque		
18711	18710	0 – 1 ³ / ₈ " (0-35 mm)	A4424

SenSura® Mio Kids 1-Piece MIDI Urostomy Pouch

- Star-shaped, elastic barrier with starter hole
- Neutral gray, water-resistant textile
- Soft outlet can be connected to a drain bag for nighttime use, or a leg bag for kids in a wheelchair
- **Not made with natural rubber latex**



Units: 10

Code	Code	Stoma size	HCPCS*
Transparent	Neutral gray		
18712	18713	³ / ₈ - 1 ³ / ₈ " (10-35 mm)	A5071

Assura® AC (Adhesive Coupling) Pediatric

A simple, soft, secure system for neonates! This ostomy pouching system is based on adhesives that make it easy to apply the pouch to the barrier. The pouching system is easy to clean, and simple to handle with the EasiClose™ standard integrated outlet.

Assura® AC Pediatric Flat Barrier

- Standard wear
- No starter hole
- **Not made with natural rubber latex**



Cut-to-fit

Units: 10



Code	Stoma Size	Color	HCPCS*
14307	0 - 5/8" (0-17 mm)		A4414
14308	0 - 1" (0-27 mm)		A4414

Assura® AC Pediatric Drainable Pouch

- EasiClose™ Outlet
- Filter
- Soft cloth on back
- **Not made with natural rubber latex**



Units: 10

Code	Code	Length	Color	HCPCS*
Transparent Opaque (Teddy Bear)				
14681	14691	5½" (14 cm)		A4425
14682	14692	6½" (16 cm)		A4425

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Assura® 2-piece Pediatric Pouching Systems

Specifically developed adhesive and a smaller barrier size are ideal for pediatrics, while a mix-and-match “green” flange size offers versatility. The SenSura® and Assura® pouches and barrier may be used interchangeably with “green”-sized adult SenSura® Click or Assura® products for a broad array of pouching options.


Assura® Pediatric Standard Wear Barrier

- Flat
- **Not made with natural rubber latex**



Cut-to-fit

Units: 5

Code		Color	HCPCS*
2182	3/8 - 1 3/8" (10-35 mm)		A4414

Assura® Pediatric Drainable Pouch

- Clamp closure
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 8 1/2" (21 cm) / 150 mL

Units: 10

Code		Color	HCPCS*
2156			A4426


Assura® Pediatric Closed Pouch

- Soft cloth on front and back
- **Not made with natural rubber latex**



Opaque 5 3/4" (14 1/2 cm)

Units: 30

Code		Color	HCPCS*
2161			A4423

Assura® Stoma Cap

- With filter
- Small size is optimal for swimming at the pool
- **Not made with natural rubber latex**



Opaque

Units: 30

Code		Color	HCPCS*
2801			A5055

Assura® Pediatric Urostomy Pouch

- Soft outlet can be connected to a drain bag for nighttime use, or a leg bag for kids in a wheelchair - see page 72 for options
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent 8½" (21 cm) / 150 mL

Units: 10

Code	Color	HCPCS*
1789		A4433

Assura® Urostomy Micro-Pouch

- Miniature size prevents twisting
- Soft outlet can be connected to a drain bag for nighttime use, or a leg bag for people in a wheelchair
- Soft cloth on front
- **Not made with natural rubber latex**



Transparent 6" (15 cm)

Units: 10

Code	Color	HCPCS*
14204		A4433

Assura® Urostomy Minicap

- Ultra-absorbent material inside
- Small size is optimal for swimming at the pool
- **Not made with natural rubber latex**



Opaque 100 mL

Units: 30

Code	Color	HCPCS*
2807		A5055

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

ColoKids™ 1-piece Pediatric Pouching Systems

Specially developed Assura® “swiss roll” technology provides moisture absorption and balance while offering consistent tack and adhesion to the skin. Small pouch sizes provide a pediatric-friendly pouching option.

ColoKids™ Flat 1-piece Drainable Pouch

- Standard wear
- Clamp closure
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent

Units: 10

Code	Cut	Stoma Size	HCPCS*
2115	Cut-to-fit	$\frac{3}{8}$ - $1\frac{3}{8}$ " (10-35 mm)	A5061

ColoKids™ Flat 1-piece Closed Pouch

- Standard wear
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent

Units: 30

Code	Cut	Stoma Size	HCPCS*
2125	Cut-to-fit	$\frac{3}{8}$ - $1\frac{3}{8}$ " (10-35 mm)	A4416

ColoKids™ Flat 1-piece Urostomy Pouch

- Standard wear
- Soft outlet can be connected to a drain bag for nighttime use, or a leg bag for people in a wheelchair
- Soft cloth on back
- **Not made with natural rubber latex**



Transparent

Units: 10

Code	Cut	Stoma Size	HCPCS*
8009	Cut-to-fit	$\frac{5}{8}$ - $1\frac{3}{8}$ " (15-35 mm)	A5071

Brava® supporting products

Pediatric weight thresholds

How does Coloplast define pediatric age / weight parameters as relate to products?

- Premature: 30+ gestational weeks
- Neonate: 0-6 months
- Toddler: 6 months – 4 years

Product name	Caution	Conversion (kg to lbs.)	Approximate minimum age
Brava® Elastic Barrier Strips	Not recommended for children below 5.6 kg	12.34 lbs.	Toddler
Brava® Protective Seal	Not recommended for children below 5.6 kg	12.34 lbs.	Toddler
Brava® Moldable Ring	Not recommended for children below 7.2 kg	15.87 lbs.	Toddler
Brava® Strip Paste	Not recommended for children below 2.6 kg	5.73 lbs.	Neonate
Brava® Paste (tube paste)	Not recommended for children below 6 kg	13.22 lbs.	Toddler
Brava® Adhesive Remover Wipe	Not recommended for children below 2.4 kg	5.29 lbs.	Neonate
Brava® Adhesive Remover Spray	Not recommended for children below 5.6 kg	12.34 lbs.	Toddler
Brava® Baby Adhesive Remover	Not recommended for children below 2.4kg	5.29 lbs.	Neonate
Brava® Skin Barrier Wipe	Not recommended for children below 2.4 kg	5.29 lbs.	Neonate
Brava® Skin Barrier Spray	Not recommended for children below 5.6 kg	12.34 lbs.	Toddler
Brava® Skin Cleanser	Not recommended for children below 2.4 kg	5.29 lbs.	Neonate
Brava® Protective Sheet	Not recommended for children below 7 kg	15.43 lbs.	Toddler
Brava® Belt	Does not fit kids	-	-
Brava® Lubricating Deodorant	Not recommended for children	-	-
Brava® Powder	Not recommended for children below 5.6 kg	12.34 lbs.	Toddler

kg to lbs.

1 kg = 2.20 lbs.

2.4 kg = 5.29 lbs.

2.6 kg = 5.73 lbs.

5.6 kg = 12.34 lbs.

kg to lbs.

6.0 kg = 13.22 lbs.

7.0 kg = 15.43 lbs.

7.2 kg = 15.87 lbs.

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Night Drainage Bags

Urostomy Night Bag and Bag Hanger

See Accessories page 71-72 for more information.



Bathing

Bedside-Care® EasiCleanse® (page 51)

Moisturizers

Sween Cream® (page 53)

Sween® 24 (page 53)

Sween® Lotion (page 53)

Skin Protectants

Critic-Aid® Clear (page 54)

Critic-Aid® Skin Paste (page 54)

Wound Dressings

Triad® Hydrophilic Wound Dressing (page 42)

All product recommendations on this page are based on:
Association of Women's Health, Obstetric and Neonatal Nurses.
Neonatal Skin Care, Fourth Edition. Evidence-Based Clinical
Practice Guideline, 2018.

Previous item code	New item code	Difference	New item page #
10811	10870	Textile grey pouch with pre-filter	115
10821	10871	Textile grey pouch with pre-filter	115
10822	10872	Textile grey pouch with pre-filter	115
10823	10873	Textile grey pouch with pre-filter	115
10824	10874	Textile grey pouch with pre-filter	115
10831	10881	Textile grey pouch with pre-filter	115
10833	10883	Textile grey pouch with pre-filter	115
10834	10884	Textile grey pouch with pre-filter	115
10835	10885	Textile grey pouch with pre-filter	115
14301	10551	Double-layer adhesive technology	99
14302	10561	Double-layer adhesive technology	99
14303	10571	Double-layer adhesive technology	99
14304	10104	Double-layer adhesive technology	103
14305	10105	Double-layer adhesive technology	103
14306	10106	Double-layer adhesive technology	103
14309	10108	Double-layer adhesive technology	103
14316	12202	Textile grey pouch with pre-filter	102
14317	12213	Textile grey pouch with pre-filter	102
14319	12222	Textile grey pouch with pre-filter	102
14326	12203	Textile grey pouch with pre-filter	102
14327	12216	Textile grey pouch with pre-filter	102
14328	12216	Textile grey pouch with pre-filter	102
14329	12223	Textile grey pouch with pre-filter	102
14331	12202	Textile grey pouch with pre-filter	102
14332	12213	Textile grey pouch with pre-filter	102
14333	12213	Textile grey pouch with pre-filter	102
14334	12222	Textile grey pouch with pre-filter	102
14342	12273	Textile grey pouch with pre-filter	100
14344	12283	Textile grey pouch with pre-filter	100
14346	12262	Textile grey pouch with pre-filter	100
14347	12272	Textile grey pouch with pre-filter	100
14349	12281	Textile grey pouch with pre-filter	100
14355	12222	Textile grey pouch with pre-filter	102
14356	12262	Textile grey pouch with pre-filter	100
14357	12273	Textile grey pouch with pre-filter	100
14359	12283	Textile grey pouch with pre-filter	100
14361	12262	Textile grey pouch with pre-filter	100
14362	12272	Textile grey pouch with pre-filter	100
14364	12282	Textile grey pouch with pre-filter	100
14379	11591	EasiClose™ WIDE Outlet	104

*For important information regarding the use of HCPCS codes listed and code definitions see page 180

Discontinued Items

Previous item code	New item code	Difference	New item page #
14401	16471	Double-layer adhesive technology	99
14403	16483	Double-layer adhesive technology	99
14406	16491	Double-layer adhesive technology	99
14408	11308	Double-layer adhesive technology	104
14521	12274	Textile grey pouch with pre-filter	100
14525	12274	Textile grey pouch with pre-filter	100
14527	11690	EasiClose™ WIDE Outlet	105
14552	12292	Textile grey pouch with anti-reflux	102
14553	12291	Textile grey pouch with anti-reflux	102
14554	12298	Textile grey pouch with anti-reflux	102
14555	12297	Textile grey pouch with anti-reflux	102
14601	11304	Double-layer adhesive technology	104
14602	11305	Double-layer adhesive technology	104
14603	11306	Double-layer adhesive technology	104
14642	16472	Double-layer adhesive technology	99
14643	16485	Double-layer adhesive technology	99
14644	16486	Double-layer adhesive technology	99
14645	16487	Double-layer adhesive technology	99
14646	16488	Double-layer adhesive technology	99



Gabriela

Making life easier_

Coloplast® Care is the only manufacturer-provided patient support program with better patient outcomes supported by peer-reviewed data.



About the team

Coloplast® Care Advisors assist patients with their intimate healthcare needs. They are available Monday through Friday from 7 a.m. to 7 p.m. (CT). The team offers one-on-one assistance for the following:

Quality of life guidance

Our team can answer questions regarding product lines, routines that will fit their lifestyle, and encourage self-assessments for ongoing wellness support.

Connecting to the care network

Coloplast Care can help locate a medical supplies dealer in-network with your patient's insurance to ensure that they have access to the correct products.

Community resources

Coloplast Care Advisors also provide awareness of resources in communities where your patients live, such as support groups, events, or conferences. The Advisors can also offer free samples*, if needed, from the Coloplast ostomy and continence product lines.

For ostomy: visit www.ostomy.coloplastcare.us or call **1-877-858-2656** to enroll your patients today.

For continence: visit bladder.coloplastcare.us or call **1-866-226-6362** to enroll your patients today.

For Peristeen® inquiries, visit coloplastcare.us/bowel or call **1-855-605-7594** for more information.

*Limitations apply. SpeediCath®, Self-Cath® and SureCath® catheters are available by prescription only. Samples are not available for Peristeen.

Information from Coloplast Care is for educational purposes only. It is not intended to substitute for professional medical advice and should not be interpreted to contain treatment recommendations.

Coloplast® Care **Connect**

ColoplastCareConnect is the new easy online portal for enrolling your patients and access to team documentation. Scan below to enroll your patients or to view your dashboard of electronic submissions:



Contact your rep if you need access.

The following materials are presented for general purposes only. They do not constitute medical advice and should not be a substitute for consultation with a trained medical professional.

The information presented or discussed may not be representative of all patient outcomes. Each person's situation is unique and risks, outcomes, experience, and results may vary. Please see complete product instructions for use, including all product indications, contraindications, precautions and warnings, and adverse events.

Videos

Intermittent Self-Catheterization and Bowel Management

- Visit www.coloplastcare.us/bowel to view interactive how-to-use videos for Peristeen.
- Visit bladderguides.coloplastcare.us for interactive how-to-use videos for intermittent catheters.

Measuring Guides

- Conveen® Optima Male External Catheter Selection and Fitting Guide
- Freedom Male External Catheter Selection and Fitting Guide

Educational Guides

- Making Peristeen® Work for Your Child, Pediatric Patient Information Kit, English & Spanish
- Welcome to Peristeen®, Patient Information Kit - Adult, English & Spanish
- Peristeen® Training Guide for Healthcare Professionals
- Male IC Care Kit, Patient Information Kit- Adult, English & Spanish
- Female IC Care Kit, Patient Information Kit- Adult, English & Spanish

Wellness Education

- Socializing and Peristeen®, English & Spanish
- Traveling and Peristeen®, English & Spanish
- Hand Hygiene, English & Spanish
- Basics of Intermittent Self Catheterization
- Bladder Hygiene, Better Bladder Health
- Medicare Reimbursement Guidelines
- Understanding UTIs
- Catheterization Reference Guide
- Basics of MECs
- Why Intermittent Self-Catheterization
- Bladder Diary

Coloplast® Care **Connect**

ColoplastCareConnect is the new easy online portal for enrolling your patients and access to team documentation. Scan below to enroll your patients or to view your dashboard of electronic submissions:



Contact your rep if you need access.

The following materials are presented for general purposes only. They do not constitute medical advice and should not be a substitute for consultation with a trained medical professional.

The information presented or discussed may not be representative of all patient outcomes. Each person's situation is unique and risks, outcomes, experience, and results may vary. Please see complete product instructions for use, including all product indications, contraindications, precautions and warnings, and adverse events.

Ostomy Self-Assessment Tools

- BodyCheck™
- Ostomy Check
- Troubleshooter
- Apply-Remove-Check (ARC)

Diaries:

- Pouch change diary
- Exercise diary
- Hernia prevention action plan
- Ostomy Wellness Output Diary

Educational Guides

- Life After: Ileostomy
- Life After: Colostomy
- Life After: Urostomy
- Life After: Your Child's Ostomy
- Teen Life with an Ostomy

Product Pictorial Guide - Also available in Spanish

NOTE: These are abbreviated instructions designed to visually show how to use the product. Always refer to the product labeling (including the package insert) before using.

- SenSura® Mio 1-piece Pouching System
- SenSura® Mio Click 2-piece Pouching System
- SenSura® Mio Flex 2-piece Pouching System
- SenSura® 1-piece Pouching System
- SenSura® Click 2-piece Pouching System
- SenSura® Flex 2-piece Pouching System
- Fistula Wound Management System
- SenSura® Mio 1-piece Urostomy Pouching System
- SenSura® Mio Click 2-piece Urostomy Pouching System
- SenSura® Mio Flex 2-piece Urostomy Pouching System
- SenSura® Mio Baby 2-piece Flex Barrier, pouch, and Brava® Baby Adhesive Remover

Spanish Educational Tools

Educational Guides

- La Vida Después de la Colostomía Life After: Colostomy
- La Vida Después de la Ileostomía Life After: Ileostomy
- La Vida Después de la Urostomía Life After: Urostomy

Other languages also available in digital format only:
Arabic, Chinese, Russian and Vietnamese

For additional educational resources, visit
www.coloplast.us

Coloplast® Professional



Focus on
YOUR
education

Online Learning with CE Contact Hours

Courses for Healthcare Professionals working in continence, ostomy and wound clinical spaces.



Expert Education Programs

Live or virtual education on a variety of clinical topics presented by content experts, providing interactive learning experiences.



Facility-led Education

Train the Trainer courses include slide presentations, training resources, toolkits and more to help you facilitate a meaningful training experience.



Education Resources

Videos, podcasts and printables designed to support your patients with product questions and needs in everyday situations and help you promote positive outcomes for your patients.



Register by using one of the following methods:

Website - www.coloplast.us/professional

Email - CProfessional_US@coloplast.us

Coloplast® Professional

Coloplast

Professional is a comprehensive program that suits a variety of learning needs and provides ongoing support to deliver best-in-class education training, and resources.



Product Training

This hands-on training program is tailored to your staff and your schedule to establish proper product usage and compliance to protocols. Coloplast provides product training at implementation and as needed to meet your needs.



Expert Education Programs

A Coloplast Clinical Consultant or Key Opinion Leader in the field of wound and skin care will provide a live learning experience at your facility, off-site, or via webinar. Learn best practices while sharing experiences with colleagues.



Online Learning

Coloplast HEAL is an online education program that offers CE credits for healthcare professionals. The program aims to expand knowledge of modern wound healing principles and improve the standard of care for patients. To learn more and register, visit www.coloplast.us/professional



Facility-led Education

This train-the-trainer program will provide expertise and materials to key individuals in your organization so they can drive education across the organization. Materials provided include presentations, training tools, and more to help you facilitate a meaningful training experience.



To learn more about the education programs offered by Coloplast, contact Customer Service at **1-800-533-0464**.

The information in this document is informational only, general in nature, and does not cover all payers' rules or policies. This information was obtained from third party sources and is subject to change without notice as a result of changes in reimbursement regulations and payer policies. This document represents no promise or guarantee by Coloplast Corp. regarding coverage or payment for products or procedures by CMS or other payers. Providers are responsible for reporting the codes that most accurately describe the patient's medical condition, procedures performed and products used. Providers should check Medicare bulletins, manuals, program memoranda and Medicare guidelines to ensure compliance with Medicare requirements. Inquires should be directed to the appropriate other payer for non-Medicare coverage situations.

Medicare Utilization - Continence Care

Usual maximum limits for some products from the Medicare Local Coverage Determinations (LCD) are shown below. Suppliers should examine their applicable Medicare contractor policies for full examination of policy requirements, including medical necessity and documentation.

HCPCS*	Long Description	Usual Maximum
A4310	Insert tray without drainage bag and without catheter (accessories only)	1 insertion tray per episode of indwelling catheter insertion
A4332	Lubricant, individual sterile packet, each	200/month
A4349	Male external catheter, with or without adhesive, disposable, each	35/month
A4351	Intermittent urinary catheter; straight tip, with or without coating (Teflon, silicone, silicone elastomer, or hydrophilic, etc.), each	200/month
A4352	Intermittent urinary catheter; Coudé (curved) tip, with or without coating (Teflon, silicone, silicone elastomer, or hydrophilic, etc.) each	200/month
A4353	Intermittent urinary catheter, with insertion supplies	200/month
A4357	Bedside drainage bag, day or night, with or without anti-reflux device, with or without tube, each	2/month
A4358	Urinary drainage bag, leg or abdomen, vinyl, with or without tube, with straps, each	2/month

Medicare Utilization - Ostomy

Usual maximum limits for some products from the Medicare Local Coverage Determinations (LCD) are shown below. The amount of ostomy supplies varies by patient based on the type of ostomy, its location, its construction, and the condition of the peristomal skin. The table below shows the maximum number of items that are usually medically necessary for patients with ostomies. It may be more or less depending on the factors that affect the frequency of barrier and pouch change. The medical necessity for a greater quantity of supplies must be clearly documented in the patient's medical record (physician's signed and dated order).

HCPCS*	Short Description	Usual Maximum
A4357	Bedside drainage bag	2 ea/month
A4361	Ostomy faceplate, each	3 ea/6 months
A4362	Skin barrier, solid, 4x4 or equivalent; each*	20 ea/month
A4363	Ostomy clamp replacement	none*
A4367	Ostomy belt, each	1 ea/month
A4368	Ostomy filter	none*
A4369	Ostomy skin barrier; liquid, per ounce	2 oz./month
A4371	Ostomy skin barrier powder, per ounce	10 oz/6 months
A4372	Ostomy skin barrier, 4 x 4, standard wear, CONVEX	none*
A4373	Ostomy skin barrier, with flange, CONVEX	none*
A4385	Ostomy skin barrier, 4 x 4, extended wear, FLAT	none*
A4388	One-piece ostomy drainable pouch and barrier, extended wear, FLAT	none*
A4389	One-piece ostomy drainable pouch and barrier, CONVEX	none*
A4390	One-piece ostomy drainable pouch, extended wear barrier, CONVEX	none*
A4391	One-piece urinary pouch, extended wear barrier, FLAT	none*
A4392	One-piece urinary pouch, standard wear barrier, FLAT	none*
A4394	Ostomy deodorant per ounce	none*
A4397	Ostomy irrigation supply; sleeve, each • HCPCS code A4397 will be discontinued effective 12/31/2021. Effective 1/1/2022 CMS has established two new HCPCS codes to distinguish between reusable and disposable sleeves: - A4436 Irrigation supply, sleeve, reusable - A4437 Irrigation supply, sleeve, disposable	4 ea/month TBD TBD
A4398	Ostomy irrigation supply; bag, each	2 ea/6 months
A4399	Ostomy irrigation supply; cone/catheter, with or without brush	2 ea/6 months
A4404	Ostomy ring, each	10 ea/month
A4405	Ostomy skin barrier, non-pectin based, paste, per ounce	4 oz/month
A4406	Ostomy skin barrier, pectin-based, per ounce	4 oz/month
A4407	Two-piece ostomy skin barrier, with flange, extended wear, CONVEX, (4 x 4 inches or smaller)	none*
A4408	Two-piece ostomy skin barrier, with flange extended wear, CONVEX, (larger than 4 x 4 inches)	none*
A4409	Two-piece ostomy skin barrier, with flange, extended wear, FLAT, (4 x 4 inches or smaller)	none*

* No published Medicare usual maximum. Monthly quantity subject to medical necessity documentation.

**Elastic Barrier Strips - Curved or Straight: 40 strips/month

HCPCS*	Short Description	Usual Maximum
A4410	Two-piece ostomy skin barrier, with flange, extended wear, FLAT, (larger than 4 x 4 inches)	none*
A4411	Two-piece ostomy solid skin barrier, extended wear, CONVEX	none*
A4412	Two-piece high output ostomy drainable pouch without filter	none*
A4413	Two-piece high output ostomy drainable pouch with filter	none*
A4414	Two-piece ostomy skin barrier with flange, FLAT (4x4 inches or smaller)	20 ea/month
A4415	Two-piece ostomy skin barrier, with flange, FLAT, greater than 4 x 4 inches	20 ea/month
A4416	One-piece ostomy closed pouch and barrier with filter, FLAT	60 ea/month
A4417	One-piece ostomy closed pouch and barrier with filter, CONVEX	60 ea/month
A4418	One-piece ostomy closed pouch with filter and barrier	60 ea/month
A4419	Two-piece ostomy closed pouch with filter for non-locking system	60 ea/month
A4423	Two-piece ostomy closed pouch with filter for locking system	60 ea/month
A4424	One-piece ostomy drainable pouch and barrier with filter, FLAT	20 ea/month
A4425	Two-piece ostomy drainable pouch with filter for non-locking system	20 ea/month
A4426	Two-piece ostomy drainable pouch without filter for locking system	20 ea/month
A4427	Two-piece ostomy drainable pouch with filter for locking system	20 ea/month
A4433	Two-piece ostomy urinary pouch for locking system	20 ea/month
A4436	Irrigation supply, sleeve, reusable, per month • Effective 1/1/2022	TBD
A4437	Irrigation supply, sleeve, reusable, per month • Effective 1/1/2022	TBD
A4455	Adhesive remover or solvent, per ounce	16 oz/6 months
A4456	Adhesive remover, wipes	none*
A5055	Stoma cap	31 ea/month
A5056	One-piece drainable ostomy pouch, extended wear with filter	40 ea/month
A5057	One-piece drainable ostomy pouch extended wear with filter, CONVEX	40 ea/month
A5061	One-piece ostomy drainable pouch and barrier, FLAT	20 ea/month
A5063	Two-piece ostomy drainable pouch for barrier with flange	20 ea/month
A5071	One-piece ostomy urinary pouch	20 ea/month
A5073	Two-piece ostomy urinary pouch	20 ea/month
A5081	Stoma plug or seal	31 ea/month
A5120	Skin barrier, wipes or swabs, each	150 ea/6 months
A5121	Solid skin barrier, 6x6 or equivalent, each**	20 ea/month
A5122	Solid skin barrier, 8x8 or equivalent, each	20 ea/month
A6154	Wound pouch, each	3/week

* No published Medicare usual maximum. Monthly quantity subject to medical necessity documentation.

**Elastic Barrier Strips - Y-shape: 60 strips/month

Medicare Utilization - Wound Care

Usual maximum limits from the Medicare Local Coverage Determinations (LCD) are shown below. Suppliers should examine their applicable Medicare contractor policies for full examination of policy requirements, including medical necessity and documentation.

HCPCS*	Short Description	Usual Maximum
A6196	Alginate dressing ≤ 16 sq. in.	1/day
A6197	Alginate dressing > 16 ≤ 48 sq. in.	1/day
A6199	Alginate dressing wound filler, per 6 inches	1/day
A6203	Sterile composite dressing, ≤ 16 sq. in. with border	3/week
A6204	Sterile composite dressing >16 ≤ 48 sq. in. with border	3/week
A6206	Sterile contact layer, ≤ 16 sq. in.	1/week
A6207	Sterile contact layer > 16 ≤ 48 sq. in.	1/week
A6209	Sterile foam dressing ≤ 16 sq. in. without border	3/week
A6210	Sterile foam dressing > 16 ≤ 48 sq. in. without border	3/week
A6211	Sterile foam dressing > 48 sq. in. without border	3/week
A6212	Sterile foam dressing ≤ 16 sq. in with border	3/week
A6213	Sterile foam dressing > 16 ≤ 48 sq. in. with border	3/week
A6234	Sterile hydrocolloid dressing ≤ 6 sq. in. without border	3/week
A6235	Sterile hydrocolloid dressing > 16 ≤ 48 sq. in. without border	3/week
A6236	Sterile hydrocolloid dressing > 48 sq. in. without border	3/week
A6237	Sterile hydrocolloid dressing ≤ 16 sq. in. with border	3/week
A6238	Sterile hydrocolloid dressing > 16 ≤ 48 sq. in. with border	3/week
A6240	Sterile hydrocolloid dressing filler, paste, per ounce	3/week
A6248	Hydrogel dressing, wound filler, gel, per fluid ounce	3 fl. oz. per wound/30 days
A6251	Sterile specialty absorptive dressing, wound cover, ≤ 16 sq. in. without border	1/day
A6252	Sterile specialty absorptive dressing, wound cover, > 16 ≤ 48 without border	1/day
A6260	Sterile wound cleanser, any type/size	Not covered under the surgical dressing benefit

Reimbursement Disclaimer: Coloplast Corp. provides this information for convenience and your general reference only. It does not constitute legal advice or a recommendation regarding clinical practice. Reimbursement, coverage and payment policies can vary from one insurer and region to another and is subject to change without notice. The provider has the responsibility to determine medical necessity and to submit appropriate codes and charges for care provided. Coloplast does not guarantee coverage or payment of products and Coloplast makes no guarantee that the use of this information will prevent differences of opinion or disputes with Medicare, insurers, or other payers as to the correct form of billing or the amount that will be paid. This information is provided for your general information only and is not intended to replace any advice you receive from your own internal or external insurance coverage consultants, reimbursement specialists or legal counsel.

HCPCS Code Reference

A4322	Irrigation syringe, bulb or piston, each
A4331	Extension drainage tubing, any type, any length, with connector/adapter, for use with urinary leg bag or urostomy pouch, each
A4332	Lubricant, individual sterile packet, each
A4334	Urinary catheter anchoring device, leg strap, each
A4335	Incontinence supply; miscellaneous
A4338	Indwelling catheter; Foley type, two-way latex with coating (Teflon, silicone, silicone elastomer, or hydrophilic, etc.), each
A4340	Indwelling catheter, specialty type, (e.g., Coudé, mushroom, wing, etc.), each
A4344	Indwelling catheter; Foley type, two-way, all silicone, each
A4346	Indwelling catheter; Foley type, three-way for continuous irrigation, each
A4349	Male external catheter, with or without adhesive, disposable, each
A4351	Intermittent urinary catheter; straight tip, with or without coating (Teflon, silicone, silicone elastomer, or hydrophilic, etc.), each
A4352	Intermittent urinary catheter; Coudé (curved) tip, with or without coating (Teflon, silicone, silicone elastomer, or hydrophilic, etc.), each
A4353	Intermittent urinary catheter, with insertion supplies
A4357	Bedside drainage bag, day or night, with or without anti-reflux device, with or without tube, each
A4358	Urinary drainage bag, leg or abdomen, vinyl, with or without tube, with straps, each
A4361	Ostomy faceplate, each
A4362	Solid skin barrier; solid, 4x4 or equivalent; each
A4363	Ostomy clamp, any type, replacement only, each
A4366	Ostomy vent, any type, each
A4367	Ostomy belt, each
A4368	Ostomy filter, any type, each
A4369	Ostomy skin barrier liquid (spray, brush, etc.), per ounce
A4371	Ostomy skin barrier powder, per ounce
A4372	Ostomy skin barrier, solid 4x4 or equivalent, standard wear, with built-in convexity, each
A4373	Ostomy skin barrier, with flange (solid, flexible or accordion), with built-in convexity, any size, each
A4385	Ostomy skin barrier, solid 4 x 4 or equivalent, extended wear, without built-in convexity, each

HCPCS Code Reference

A4388	Ostomy pouch, drainable, with extended wear barrier attached, (1 piece), each
A4389	Ostomy pouch, drainable, with barrier attached, with built-in convexity (1-piece), each
A4390	Ostomy pouch, drainable, with extended wear barrier attached, with built-in convexity (1-piece), each
A4391	Ostomy pouch, urinary, with extended wear barrier attached (1-piece), each
A4392	Ostomy pouch, urinary with standard wear barrier attached, with built-in convexity (1-piece), each
A4393	Ostomy pouch, urinary, with extended wear barrier attached, with built-in convexity (1 pc), each
A4394	Ostomy deodorant, with or without lubricant, for use in ostomy pouch, per fluid ounce
A4396	Ostomy belt with peristomal hernia support
A4397	Irrigation supply; sleeve, each
A4398	Ostomy irrigation supply; bag, each
A4399	Ostomy Irrigation supply; cone /catheter, with or without brush
A4402	Lubricant, per ounce
A4404	Ostomy ring, each
A4405	Ostomy skin barrier, non-pectin based, paste, per ounce
A4406	Ostomy skin barrier, pectin-based, paste, per ounce
A4407	Ostomy skin barrier, with flange (solid, flexible, or accordion), extended wear, with built-in convexity, 4x4 inches or smaller, each
A4408	Ostomy skin barrier, with flange (solid, flexible or accordion), extended wear, with built-in convexity, larger than 4x4 inches, each
A4409	Ostomy skin barrier, with flange (solid, flexible or accordion), extended wear, without built in convexity, 4x4 inches or smaller, each
A4410	Ostomy skin barrier with flange (solid, flexible or accordion), extended wear, without built-in convexity, larger than 4x4 inches, each
A4411	Ostomy skin barrier, solid 4x4 or equivalent, extended wear, with built-in convexity, each
A4412	Ostomy pouch, drainable, high output, for use on a barrier with flange (2-piece system), without filter, each
A4413	Ostomy pouch, drainable, high output, for use on a barrier with flange (2 piece system) with filter, each
A4414	Ostomy skin barrier, with flange (solid, flexible or accordion), without built-in convexity, 4x4 inches or smaller, each
A4415	Ostomy skin barrier, with flange (solid, flexible or accordion), without built-in convexity, larger than 4x4 inches, each
A4416	Ostomy pouch, closed, with barrier attached, with filter (1-piece,) each
A4417	Ostomy pouch, closed, with barrier attached, with built-in convexity, with filter (1-piece), each
A4418	Ostomy pouch, drainable, high output, for use on a barrier with flange (2-piece system), without filter, each
A4419	Ostomy pouch, closed; for use on barrier with non-locking flange, with filter (2-piece), each
A4421	Ostomy supply; miscellaneous
A4422	Ostomy pouch absorbent material (sheet, pad, crystal packet) for use in ostomy pouch to thicken liquid stomal output, each


HCPCS Code Reference	
A4423	Ostomy pouch, closed; for use on barrier with locking flange, with filter (2-piece), each
A4424	Ostomy pouch, drainable, with barrier attached, with filter (1-piece), each
A4425	Ostomy pouch, drainable; for use on barrier with non-locking flange, with filter (2-piece system), each
A4426	Ostomy pouch, drainable; for use on barrier with locking flange (2-piece system), each
A4427	Ostomy pouch, drainable; for use on barrier with locking flange, with filter (2 piece system), each
A4428	Ostomy pouch, urinary, with extended wear barrier attached, with faucet-type tap with valve (1-piece), each
A4429	Ostomy pouch, urinary, with barrier attached, with built-in convexity, with faucet-type tap with value (1-piece), each
A4430	Ostomy pouch, urinary, with extended wear barrier attached, with built-in convexity, with faucet-type tap with valve (1-piece), each
A4431	Ostomy pouch, urinary; with barrier attached, with faucet-type tap with valve (1- piece), each
A4432	Ostomy pouch, urinary; with barrier attached, with faucet-type tap with value (1- piece), each
A4433	Ostomy pouch, urinary; for use on barrier with locking flange (2-piece), each
A4434	Ostomy pouch, urinary; for use on barrier with locking flange, with faucet-type tap with valve (2-piece), each
A4435	Ostomy pouch, drainable, high output, with extended wear barrier (1-piece system), with or without filter, each
A4453	Rectal catheter for use with manual pump-operated enema system, replacement only
A4455	Adhesive remover or solvent (for tape, cement or other adhesive), per ounce
A4456	Adhesive remover, wipes, any type, each
A4459	Manual pump-operated enema system, includes balloon, catheter and all accessories, reusable, any type
A5051	Ostomy pouch, closed; with barrier attached (1-piece), each
A5052	Ostomy pouch, closed; without barrier attached (1-piece), each
A5053	Ostomy pouch, closed; for use on faceplate, each
A5054	Ostomy pouch, closed; for use on barrier with flange (2-piece), each
A5055	Stoma cap
A5056	Ostomy pouch, drainable, with extended wear barrier attached, with filter (1-piece), each
A5057	Ostomy pouch, drainable, with extended wear barrier attached, with built in convexity, with filter, (1-piece), each
A5061	Ostomy pouch, drainable; with barrier attached, (1-piece), each
A5062	Ostomy pouch, drainable; without barrier attached (1-piece), each
A5063	Ostomy pouch, drainable; for use on barrier with flange (2-piece system), each
A5071	Ostomy pouch, urinary; with barrier attached (1-piece), each
A5072	Ostomy pouch, urinary; without barrier attached (1-piece), each

HCPCS Code Reference

A5073	Ostomy pouch, urinary; for use on barrier with flange (2-piece), each
A5081	Stoma plug or seal, any type
A5093	Ostomy accessory; convex insert
A5113	Leg strap; latex, replacement only, per set
A5114	Leg strap; foam or fabric, replacement only, per set
A5120	Skin barrier, wipes or swabs, each
A5121	Skin barrier; solid, 6x6 or equivalent, each
A5122	Skin barrier; solid, 8x8 or equivalent, each
A5131	Appliance cleaner, incontinence and ostomy appliances, per 16 ounces
A6154	Wound pouch, each
A6196	Alginate or other fiber gelling dressing, wound cover, sterile, pad size 16 sq. in. or less, each dressing
A6197	Alginate or other fiber gelling dressing wound cover, sterile, pad size more than 16 sq. in. but less than or equal to 48 sq. in, each dressing
A6199	Alginate or other fiber gelling dressing, wound filler, sterile, per 6 inches
A6203	Composite dressing, sterile, pad size 16 sq. in. or less, with any size adhesive border, each dressing
A6204	Composite dressing, sterile, pad size more than 16 sq. in. but less than or equal to 48 sq. in., with any size adhesive border, each dressing
A6206	Contact layer, sterile, 16 sq. in. or less, each dressing
A6207	Contact layer, sterile, more than 16 sq. in. but less than or equal to 48 sq. in., each dressing
A6209	Foam dressing, wound cover, sterile, pad size 16 sq. in. or less, without adhesive border, each dressing
A6210	Foam dressing, wound cover, sterile, pad size more than 16 sq. in but less than or equal to 48 sq. in. without adhesive border, each dressing
A6211	Foam dressing, wound cover, sterile, pad size more than 48 sq. in. without adhesive border, each dressing
A6212	Foam dressing, wound cover, sterile, pad size 16 sq. in. or less, with any size adhesive border, each dressing
A6213	Foam dressing, wound cover, sterile, pad size more than 16 sq. in but less than or equal to 48 sq. in. with any size adhesive border, each dressing
A6234	Hydrocolloid dressing, wound cover, sterile, pad size 16 sq. in. or less without adhesive border, each dressing
A6235	Hydrocolloid dressing, wound cover, sterile, pad size more than 16 but less than or equal to 48 sq. in. without adhesive border, each dressing
A6236	Hydrocolloid dressing, wound cover, sterile, pad size more than 48 sq. in. without adhesive border, each dressing
A6237	Hydrocolloid dressing, wound cover, sterile, pad size 16 sq. in. or less, with any size adhesive border, each dressing
A6238	Hydrocolloid dressing, wound cover, sterile, pad size more than 16 sq. in. but less than or equal to 48 sq. in., with any size adhesive border, each dressing
A6239	Hydrocolloid dressing, wound cover, sterile, pad size more than 48 sq. in., with any size adhesive border, each dressing
A6240	Hydrocolloid dressing, wound filler, paste, sterile, per ounce
A6248	Hydrogel dressing, wound filler, gel, per fluid ounce

HCPCS Code Reference

A6251	Specialty absorptive dressing, wound cover, sterile, pad size 16 sq. in. or less, without border adhesive border, each dressing
A6252	Specialty absorptive dressing, wound cover, sterile, pad size more than 16 sq. in but less than or equal to 48 sq. in without adhesive border, each dressing
A6260	Wound cleansers, any type, any size
A6413	Adhesive bandage, first-aid type, any size, each Not covered under the surgical dressing benefit



Mike

Making life easier_

Index

Making life *easier*

Code	Page #	Code	Page #	Code	Page #
Continenence					
27-4200-B ...	25	4210	13	6430	24
110	14	4212	13	6500	24
112	14	4214	13	6530	24
114	14	4306	13	8000	25
116	14	4308	13	8030	25
208	14	4310	13	8100	25
210	14	4408	13	8200	25
212	14	4410	13	8205	25
214	14	4412	13	8230	25
240	14	4414	13	8235	25
305	14	4416	13	8300	25
306	14	4418	13	8305	25
308	14	4608	16	8400	25
310	14	4610	16	8430	25
408	14	4612	16	8500	25
410	14	4614	16	20010	15
412	14	4808	16	20012	15
414	14	4810	16	20014	15
416	14	4812	16	20016	15
418	14	4814	16	20021	15
450	14	4816	15, 16	20022	15
460	14	5062	27	20024	15
475	30	5070	30	21026	26
608	16	5100	24	21027	26
610	16	5110	24	21029	26
612	16	5170	26	21034	26
614	16	5171	26	21053	26
806	16	5174	26	21054	26
808	16	5190	24	21346	27
810	16	5200	24	21356	27
812	16	5210	24	22025	23
814	16	5221	23	22028	23
816	16	5225	23	22030	23
818	16	5225H	23	22035	23
1008	20	5230	23	22040	23
1010	20	5230H	23	22121	23
1012	20	5235	23	22125	23
1014	20	5235H	23	22130	23
1016	20	5240	23	22135	23
1108	20	5240H	23	25502	27
1112	20	5290	24	27010	12
1114	20	5300	24	27012	12
1116	20	5310	24	27014	12
2114	20	5390	24	27016	12
2214	20	5400	24	28005S	19
2614	20	5410	24	28006S	19
2814	20	5490	24	28007S	19
2816	20	5500	24	28008S	19
3114	20	5590	24	28009	19
3214	20	6100	24	28009S	19
3614	20	6130	24	28022S	19
4110	13	6200	24	28023S	19
4112	13	6230	24	28024S	19
4114	13	6300	24	28026S	19
4116	13	6330	24	28027	19
		6400	24	28027S	19

Code	Page #	Code	Page #	Code	Page #
28028	19	28814	12	1758	96
28028S	19	29122	22	1759	96
28029S	19	29127	22	1789	167
28032	19	29129	22	2041	70
28032S	19	29142	22	2115	168
28033S	19	29149	22	2125	168
28036	19	29152	22	2156	166
28036S	19	50501	30	2161	166
28037	19	AA 6C18	28	2182	166
28037S	19	AA 6C20	28	2310	70
28408	11	AA 6C22	28	2315	70
28410	11	AA 6C24	28	2320	70
28412	11	AA 6106	28	2325	70
28414	11	AA 6108	28	2330	70
28416	11	AA 6110	28	2340	70
28418	11	AA 6112	28	2350	70
28422	17	AA 6114	28	2501	133
28481	18	AA 6116	28	2650	67, 155
28482	18	AA 6118	28	2801	95, 166
28483	18	AA 6120	28	2802	95
28484	18	AA 6122	28	2803	95
28485	18	AA 6124	28	2807	95, 167
28486	18	AA 6308	28	2808	95
28490	15	AA 6310	28	2831	88
28492	15	AA 6312	29	2832	88
28494	15	AA 6314	29	2833	88
28496	15	AA 6316	29	2836	94, 144
28499	18	AA 6318	29	2837	94, 144
28500	18	AA 6320	29	2846	93, 140
28501	18	AA 6322	29	2847	93, 140
28502	18	AA 6324	29	2881	88
28503	18	AA 6408	28	2882	88
28504	18	AA 6410	28	2883	88
28506	11	AA 6412	29	4215	71
28508	11	AA 6414	29	4220	71
28510	11	AA 6416	29	4237	71
28512	11	AA 6418	29	4247	71
28514	11	AA 6422	29	5070	72
28516	11	AA 6424	29	5570	134
28520	17	C1110	20	5575	134
28522	17	C3814	20	5580	134
28524	17	C3816	20	5585	134
28578	12			7583	69
28580	12			8009	168
28582	12	Ostomy		8114	93, 140
28584	12	0509	70	8115	93, 140
28606	11	0925	70	8116	93, 140
28608	11	1003	73	9500	70, 146
28610	11	1110	74	10011	83
28612	11	1500	73	10012	83
28702	17	1500H	73	10013	83
28706	11	1751	96	10015	83
28708	11	1752	96	10016	83
28710	11	1754	96	10017	83
28810	12	1755	96	10021	83
28812	12	1756	96	10022	83
		1757	96		

Code	Page #	Code	Page #	Code	Page #
10023	83	10506	77	11011K	150
10025	83	10512	77	11012	86
10026	83	10513	77	11015	86
10027	83	10514	77	11016	86
10031	83	10515	77	11021	86
10032	83	10516	77	11022	86
10035	83	10522	77	11023	86
10036	83	10523	77	11024	86
10041	83	10524	77	11025	86
10045	83	10525	77	11026	86
10101	103	10531	77	11027	86
10102	103	10551	99	11028	86
10103	103	10561	99	11031	86
10104	103	10562	99	11032	86
10105	103	10563	99	11035	86
10106	103	10564	99	11036	86
10108	103	10565	99	11041	86
10154	86	10566	99	11045	86
10155	86	10567	99	11114	85
10156	86	10568	99	11115	85
10164	86	10571	99	11124	85
10165	86	10572	99	11125	85
10166	86	10573	99	11126	85
10167	86	10586	118	11127	85
10187	86	10811	171	11135	85
10461	109	10821	171	11136	85
10471	109	10822	171	11185	85
10472	109	10823	171	11186	85
10473	109	10824	171	11194	85
10474	109	10831	171	11195	85
10475	109	10833	171	11196	85
10476	109	10834	171	11197	85
10477	109	10835	171	11301	103
10478	109	10870	115	11302	103
10479	109	10871	115	11303	103
10480	109	10872	115	11304	104
10481	109	10873	115	11305	104
10482	109	10874	115	11306	104
10483	109	10881	115	11307	103
10484	109	10883	115	11308	104
10485	109	10884	115	11312	104
10486	109	10885	115	11313	104
10487	109	10911	105	11314	104
10488	109	10912	105	11315	104
10489	110	10914	105	11401	81
10490	109	10915	105	11402	81
10491	109	10917	105	11412	81
10492	110	10918	105	11413	81
10493	110	11001	151	11421	81
10494	110	11002	151	11422	81
10495	110	11004	148	11425	81
10496	110	11005	152	11426	81
10502	77	11008	152	11461	79
10503	77	11009	147	11462	79
10504	77	11010	147	11463	79
10505	77	11011	86	11471	79

Code	Page #	Code	Page #	Code	Page #
11472	79	12037	66, 154	12346	94
11473	79	12038	66, 154	12354	94
11474	79	12039	66, 154	12355	94
11482	79	12045	66, 154	12356	94
11483	79	12046	66, 154	12374	94
11484	79	12047	66, 154	12375	94
11486	79	12048	66, 154	12376	94
11487	79	12049	66, 154	12384	94
11488	79	12050	67, 155	12385	94
11491	82	12060	68, 156	12386	94
11492	82	12061	68, 156	12416	131
11497	82	12076	65, 153	12419	131
11498	82	12080	68, 156	12432	135
11499	82	12090	66, 154	12433	135
11500	82	12091	66, 154	12471	135
11512	104	12092	66, 154	12472	135
11514	104	12093	66, 154	12474	134
11515	104	12094	66, 154	12527	93
11517	104	12095	66, 154	12528	93
11518	104	12110	132	12529	93
11590	104	12130	132	12533	131
11591	104	12134	132	12577	93
11612	105	12135	132	12578	93
11690	105	12136	132	12579	93
11803	124	12140	132	12595	135
11804	124	12144	132	12596	135
11805	125	12145	132	12620	130
11806	125	12146	132	12630	130
11807	125	12147	132	12660	130
11808	125	12149	132	12670	130
11809	125	12170	131	12701	90
11814	125	12174	131	12702	90
11816	126	12175	131	12703	90
11817	126	12176	131	12704	90
11818	126	12177	131	12705	90
11819	126	12202	102	12706	90
11841	87	12203	102	12707	90
11844	87	12213	102	12708	90
11845	87	12216	102	12709	90
11846	87	12222	102	12716	90
11851	87	12223	102	12719	90
11854	87	12262	100	12800	143
11855	87	12272	100	12802	143
11856	87	12273	100	12805	144
12003	71	12274	100	12808	143
12004	71	12281	100	12810	143
12005	71	12282	100	12818	144
12006	71	12283	100	12820	74
12007	71	12284	100	12830	73
12008	71	12291	102	12834	74
12009	71	12292	102	12835	74
12014	161	12296	102	12836	74
12032	66, 154	12297	102	12841	88
12033	66, 154	12298	102	12842	88
12035	66, 154	12344	94	12844	88
12036	66, 154	12345	94	12845	88

Code	Page #	Code	Page #	Code	Page #
12847	88	14205	95	14308	165
12848	88	14206	95	14309	171
12991	135	14222	133	14316	171
12992	135	14224	96	14317	171
12993	135	14225	96	14319	171
12994	135	14226	96	14326	171
12995	135	14227	96	14327	171
12996	135	14228	96	14328	171
12997	135	14229	96	14329	171
13660	149	14232	90	14331	171
13661	149	14233	90	14332	171
13663	152	14234	90	14333	171
13664	152	14235	90	14334	171
13665	151	14236	90	14342	171
13703	131	14237	90	14344	171
13706	131	14238	90	14346	171
13840	128	14239	90	14347	171
13844	128	14243	90	14349	171
13845	128	14246	90	14355	171
13846	128	14249	90	14356	171
13860	128	14250	88	14357	171
13870	128	14251	88	14359	171
13874	128	14252	88	14361	171
13875	128	14254	88	14362	171
13876	128	14255	88	14364	171
13924	92	14257	88	14379	171
13925	92	14258	88	14401	172
13964	92	14261	89	14403	172
13965	92	14262	89	14404	129
13966	92	14263	89	14405	129
13974	92	14271	89	14406	172
13975	92	14272	89	14408	172
13976	92	14273	89	14411	129
13984	92	14274	89	14412	129
13985	92	14275	89	14413	130
13986	92	14276	89	14414	129
14010	.. 72, 146	14277	89	14415	129
14015	.. 72, 146	14278	89	14416	130
14021	.. 72, 146	14281	89	14419	130
14031	.. 72, 146	14282	89	14421	129
14041	.. 72, 146	14283	89	14422	129
14050	145	14291	89	14423	129
14060	145	14292	89	14424	129
14070	145	14293	89	14431	133
14103	129	14294	89	14432	133
14106	129	14295	89	14434	133
14163	129	14296	89	14435	133
14164	129	14297	89	14441	132
14165	129	14298	89	14442	132
14166	129	14301	171	14444	132
14173	129	14302	171	14445	132
14174	129	14303	171	14451	133
14175	129	14304	171	14452	133
14176	129	14305	171	14453	133
14196	130	14306	171	14454	133
14204	.. 95, 167	14307	165	14511	127

Code	Page #	Code	Page #	Code	Page #
14512	127	15981	122	16726	111
14521	172	15982	91	16727	111
14525	172	15983	91	16730	111
14527	172	15984	91	16731	111
14552	172	15985	91	16732	111
14553	172	15986	91	16733	111
14554	172	15991	122	16734	111
14555	172	15992	122	16735	111
14601	172	15993	122	16737	111
14602	172	15994	122	16738	111
14603	172	16301	116	16739	111
14642	172	16305	116	16740	111
14643	172	16306	116	16741	111
14644	172	16310	116	16742	111
14645	172	16311	116	16748	111
14646	172	16312	116	16755	112
14681	165	16313	116	16756	112
14682	165	16315	116	16757	112
14691	165	16321	117	16760	112
14692	165	16322	117	16761	112
14717	134	16323	117	16762	112
15480	124	16325	116	16763	112
15484	124	16326	116	16764	112
15521	122	16327	116	16765	112
15522	122	16330	116	16766	112
15523	122	16331	116	16767	112
15524	122	16333	116	16768	112
15531	122	16334	116	16770	112
15532	122	16341	117	16771	112
15533	122	16342	117	16772	112
15534	122	16343	117	16773	112
15606	123	16345	117	16778	112
15608	123	16346	117	16801	119
15623	123	16347	117	16805	119
15624	123	16350	117	16810	119
15633	123	16351	117	16821	120
15634	123	16353	117	16822	120
15695	123	16354	117	16823	120
15696	123	16355	117	16825	120
15697	123	16471	99	16826	120
15698	123	16472	99	16827	120
15821	122	16481	99	16830	120
15860	126	16483	99	16831	120
15869	127	16484	99	16832	120
15870	126	16485	99	16833	120
15871	127	16486	99	16835	120
15874	126	16487	99	16836	120
15875	126	16488	99	16837	120
15876	126	16491	99	16840	120
15879	127	16493	99	16841	120
15954	91	16494	99	16842	120
15955	91	16705	110	16855	121
15956	91	16706	110	16856	121
15974	91	16715	110	16857	121
15975	91	16716	110	16860	121
15976	91	16725	111	16861	121

Code	Page #	Code	Page #	Code	Page #
16862	121	18340	114	18676	113, 142
16863	121	18341	114	18680	142
16865	121	18342	114	18681	142
16866	121	18343	114	18690	142
16867	121	18344	114	18691	142
16870	121	18401	121	18700	161
16871	121	18402	121	18705	161
16872	121	18403	121	18710	164
16901	77	18404	121	18711	164
16903	77	18405	121	18712	164
16911	77	18500	78	18713	164
16913	77	18501	78	18720	162
16914	77	18510	78	18721	162
16915	77	18511	78	18722	163
16921	77	18512	78	18723	163
16923	77	18513	78	18724	163
16924	77	18520	78	18725	163
16931	77	18521	78	18730	163
16933	77	18522	78	18731	163
16934	77	18523	78	18735	162
16941	78	18524	78	18747	162
16942	78	18530	78	19001	143
16951	78	18531	78	19004	123, 142
16953	78	18532	78	19010	143
16954	78	18600	80, 139	19011	143
16955	78	18601	80, 139	19020	143
16961	78	18602	80, 139	19021	143
16963	78	18603	80, 139	19030	143
16964	78	18610	80, 139	19040	86
16971	78	18611	80, 139	19041	86
16973	78	18612	80, 139	19042	86
16974	78	18613	80, 139	19050	86
18152	118	18620	81, 140	19051	86
18170	118	18621	81, 140	19052	86
18171	118	18622	81, 140	19075	69
18172	118	18630	80, 139	21001	150
18173	118	18631	80, 139	21002	150
18180	118	18632	80, 139	21003	150
18181	118	18640	80, 139	21004	150
18182	118	18641	80, 139	21005	152
18301	114	18642	80, 139	21006	153
18302	114	18650	101, 140	21007	151
18305	114	18651	101, 140	21008	147
18312	114	18652	101, 141	21009	148
18313	114	18653	101, 141	21010	148
18314	114	18655	101, 140	21011	150
18315	114	18656	101, 140	21012	148
18320	114	18657	101, 141	21013	148
18321	114	18658	101, 141	21365	71
18322	114	18660	113, 141	26555	67, 155
18323	114	18661	113, 141	32105	69
18324	114	18665	113, 141	32155	69
18330	114	18666	113, 141	32205	69
18331	114	18670	113, 141	67589	69
18332	114	18671	113, 141	95050	69
18333	114	18675	113, 142	110072	149

Code	Page #	Code	Page #	Code	Page #
120105	68, 156	7234	51	3906	42
120115	68, 156	7235	51	6105	35
120205	68, 157	7251	58	7690	42
120215	69, 157	7262	58	9622	37
120307	67, 155	7300	49	9625	37
120427	67, 155	7301	49	9632	37
120700	65, 153	7565	53	9635	37
120721	65, 153	7566	53	9641	37
120740	65, 153	7567	53	9643	37
		7570	55	33110	40
		7571	55	33115	40
		7572	55	33120	40
		7710	54	33146	40
		7712	54	33280	41
		7725	54	33283	41
		7910	56	33285	41
		7912	56	33400	35
		7915	56	33401	35
		61760	49	33406	35
		61761	49	33434	35
		61762	49	33435	35
		61768	49	33436	35
		67143	49	33437	35
		67145	49	33438	35
		67146	49	33444	35
		67147	49	33445	35
		67770	55	33446	35
				33452	35
				33453	35
				33530	40
				33533	40
				33536	40
				33539	40
				33560	43
				33561	43
				33562	43
				33563	43
				33564	43
				39000	35
				39001	35
				39010	35
				46100	44
				46120	44
				46150	44
				46200	44
				46250	44
				46300	44
				46320	44
				46350	44
				46390	44
				46450	44

Skin Care		Wound Care	
0089	58	1061	42
0283	52	1063	42
0402	52	1166	42
0406	52	1964	42
0407	52	1967	42
0408	52	3410	35
0424	52	3413	35
0505	54	3416	35
0506	54	3420	36
1005	53	3423	36
1006	53	3430	36
1300	49	3444	35
1337	55	3470	36
1607	55	3473	36
1611	55	3475	36
1622	55	3485	36
1644	57	3488	36
1645	57	3705	38
1802	52	3710	38
1814	52	3715	38
1843	52	3740	38
1873	54	3755	39
1877	54	3760	39
1880	54	3765	39
1944	53	3780	39
1947	53	3900	42
7015	53	3903	42
7041	57		
7055	50		
7056	50		
7057	58		
7058	50		
7067	52		
7068	52		
7069	52		
7090	52		
7091	52		
7092	52		
7095	52		
7222	52		
7223	52		
7229	51		
7230	51		
7231	51		
7233	51		

Making life *easier*

Patient Education Resources

Scan here to enroll your patients into
the Coloplast® Care program



HCP Login (coloplastcareconnect.us)

Scan here for ostomy patient education
resources including How-To videos.



From going to the grocery store to biking up a mountain, users' hopes and aspirations fuel our passion of finding new ways to help make their lives easier.

We constantly question the status quo across every aspect of intimate healthcare and seek to understand our users, to better anticipate and advocate for their needs. Ever since the founding of Coloplast in 1957, we have worked closely with healthcare professionals and the users of our products to discover new solutions as well as new ways to create solutions that are sensitive to their special needs. Together, we move care forward.

Our business includes Ostomy Care, Continence Care, Wound & Skin Care and Interventional Urology. We operate globally and employ about 13,000 employees.

For additional information,
contact us at:

1-800-533-0464
info-us@coloplast.com
www.coloplast.us

