



Peristeen rectal balloon catheter – product utilization

Consideration / Request:

- 15 rectal catheters and 1 water bag are packaged per box
- Small and regular sizes available
- Prescription options, based on clinical evaluation and medical necessity:
 - 15 units per month – treatment provided every other day
 - 30 units per month – treatment provided daily



Peristeen rectal balloon catheter – size options, based on patient population

SKU	Product description	Catheter length (to blue insertion line)	Catheter diameter
29122	 Accessory unit 15 rectal catheters Regular / 1 water bag	9.5 cm	10 mm
29127	 Accessory unit 15 rectal catheters Small / 1 water bag	8 cm	10 mm



Indication: The Peristeen® Transanal Irrigation System is intended to instill water into the colon through a rectal catheter, which incorporates an inflatable balloon, inserted into the rectum to promote evacuation of the contents of the lower colon. The Peristeen Transanal Irrigation System is indicated for use by children (2 years - <12 years old), adolescent (12 years - < 18 years old), transitional adolescent (18 - <21 years old) and adult patients with neurogenic bowel dysfunction who suffer from fecal incontinence, chronic constipation, and/or time-consuming bowel management procedures.

Contraindications: Peristeen Transanal Irrigation must not be used in the following situations: known anal or colorectal stenosis, colorectal cancer, radiotherapy to the pelvis, and recent abdomino-perineal surgery, active inflammatory bowel disease, acute diverticulitis, severe diverticulitis, previous diverticulitis and diverticular abscess, chronic symptomatic diverticular disease, within 3 months of abdominal, anal or colorectal surgery, within 4 weeks of endoscopic polypectomy, recent colonic biopsy, recent endoscopic mucosal resection (EMR) and recent endoscopic sub-mucosal dissection (ESD), severe autonomic dysreflexia, ischemic colitis, during Spinal Cord Shock Phase, complex diverticular disease, in patients who are pregnant and have not used the system before*. Since the list is not exhaustive, the healthcare professional should always consider individual patient factors as well.

*If the patient is pregnant and has never used transanal irrigation before, they should not start the irrigation procedure during pregnancy.

Warning: Transanal irrigation procedure should always be carried out with care. Bowel perforation is an extremely rare, but serious and potentially lethal complication to transanal irrigation and will require immediate admission to a hospital, often requiring surgery. See the device manual for complete user instructions, contraindications, warnings, precautions, and potential complications/adverse events. For further information, call Coloplast Corp at 1-855-605-7594 or consult the company website at www.coloplast.us. **Caution: Federal law restricts this device to sale by or on the order of a physician.**