

Biatain Ag - Adhesive Dressing (Hydrocolloid border)– Code 33459



2009-01-19 Short Material Safety Data Sheet

Based on template version 1.0

Identification of the substance/mixture and of the company/undertaking

Product name: **Biatain Ag - Adhesive Dressing (Hydrocolloid border)**
Product code: **33459, 33541, 33543, 33562, 36434, 39631, 39632, 39635, 39641, 39643**
Product information: **Wound Care**
Manufacturer: **Coloplast A/S
Holtedam 1
DK-3050 Humlebaek
Denmark
Telephone +45 49111111
msds@coloplast.com**

USA

Coloplast Corp.
1601 West River Road.
Minneapolis, MN 55411
Telephone: +1-800-533-0464
Fax: +1-612-337-7803
www.us.coloplast.com

Canada

Coloplast Canada Corporation
3300 Ridgeway Drive, Unit 12
Mississauga, Ont. L5L 5Z9
Telephone: +1-800-533-0464
Fax: +1-877-820-8206

Europe

Coloplast A/S
Holtedam 1
DK-3050 Humlebaek
Telephone: +45 49 11 11 11
Fax: +45 49 11 19 32
www.coloplast.com

Hazards identification

This product consists primarily of polymer materials. The products pose no immediate hazard.

Composition/information on ingredients

This product contains a silver complex, which should be regarded as hazardous to the environment. However, the concentration is <1% and according to EC classification criteria this product should not be classified as dangerous to the environment. Main ingredients and packaging materials are listed below.

Biatain Ag - Adhesive Dressing (Hydrocolloid border)– Code 33459

Chemical name	CAS no.	EC no.
Product:		
PU (polyurethane) foam	-	
PU (polyurethane) film	-	
Silver complex	-	
TPE (thermoplastic elastomer)	-	
Adhesive	-	
DOA	103-23-1	
PET (polyester)	25038-59-9	
Packaging:		
Aluminum	-	
PE (polyethylene)	9002-88-4	
PET (polyester)	25038-59-9	

In addition all Substances of Very High Concern (SVHC) present in the product at a concentration above 0.1% should be stated in this short MSDS in order to comply with EC 1907/200, article IV, §33.

Disposal considerations

The recommended disposal technology at any approved facility. The disposal should always be in compliance with national, federal, state and local regulations. The product should not be discharged to the environment.

US

This product does not meet the criteria for hazardous waste as defined under the Resource Conservation and Recovery Act (RCRA) 40 CFR 261. Under normal private use the product may be disposed of together with other household waste per RCRA 40 RFT 261.4.B1.

European Union

Per The European Waste Catalogue (EWC), in accordance with EC Directive 75/422/EEC, the following Waste Code can be used: 18 01 04 00 wastes whose collection and disposal is not subject to special requirements in view of the prevention of infection (e.g. dressings, plaster casts, linen, disposable clothing, diapers). However, if the waste in view of the prevention of infection needs special requirements, other Waste Codes should be used. It is the responsibility of the disposed of together with other household waste.

Handling and storage

Handling: See instruction for use
Storage: Store until use as supplied and at room temperature unless other information is stated on the packaging or on the leaflet.

Other information

This MSDS is supplied as an additional service to the customer. The product is a medical device, which is regulated under the Council Directive 93/42/EEC, Medical Device Directive. Medical devices do not require a MSDS for hazard communication and therefore this MSDS is not a 16

Biatain Ag - Adhesive Dressing (Hydrocolloid border)– Code 33459

point data sheet. The product has been evaluated according to the requirements of medical devices. According to current knowledge this product is considered non-toxic. For further information please contact Coloplast A/S.